# AESCULAP®

**en Instructions for use/Technical description**

A-Space SIBD and Arcadius<sup>XP</sup> L Instruments

**USA**

Note for U.S. users

This Instructions for Use is NOT intended for United States users. Please discard. The Instructions for Use for United States users can be obtained by visiting our website at [www.aesculapimplant-systemsifus.com](http://www.aesculapimplant-systemsifus.com). If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-866-229-3002. A paper copy will be provided to you upon request at no additional cost.

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


\* Instrument offered upon customer request only, not included in the instrument set.

### Contents

1.	About this document . . . . .	3
1.1	Scope . . . . .	3
1.2	Safety messages . . . . .	3
2.	Clinical use . . . . .	3
2.1	Areas of use and limitations of use . . . . .	3
2.1.1	Intended use . . . . .	3
2.1.2	Indications . . . . .	4
2.1.3	Contraindications . . . . .	4
2.2	Safety information . . . . .	4
2.2.1	Clinical user . . . . .	4
2.2.2	Product-specific safety information . . . . .	4
2.2.3	Sterility . . . . .	4
2.3	Application . . . . .	4
2.3.1	Distraction . . . . .	5
2.3.2	Selection of Implant Size Using Trial Implants . . . . .	5
2.3.3	Filling the Implant . . . . .	5
2.3.4	Implant Placement and Positioning . . . . .	6
2.3.5	Bone Screw Placement . . . . .	6
2.3.6	Revision . . . . .	8
2.3.7	General . . . . .	8
3.	Validated reprocessing procedure . . . . .	9
3.1	General safety instructions . . . . .	9
3.2	General information . . . . .	9
3.3	Reusable products . . . . .	10
3.4	Preparations at the place of use . . . . .	10
3.5	Preparing for cleaning . . . . .	10
3.6	Disassembly . . . . .	10
3.6.1	Straight Ball Hex Screwdriver (ME013R) and Straight Hex Screwdriver (ME016R) . . . . .	10
3.6.2	U-Joint screwdriver (ME014R), U-Joint Bone Awl (ME190R) and Implant Extraction Instrument (ME018R) . . . . .	10
3.6.3	Implant Inserter Manipulator (ME015R) and Implant Introducer (SJ605R) . . . . .	11
3.6.4	Trial Insertion Instrument (ME020R) . . . . .	11
3.6.5	U-Joint Drill (ME189R) . . . . .	12
3.7	Cleaning/Disinfection . . . . .	12
3.7.1	Product-specific safety information on the reprocessing method . . . . .	12

3.7.2	Validated cleaning and disinfection procedure .....	13
3.8	Manual cleaning/disinfection .....	15
3.8.1	Manual cleaning with immersion disinfection .....	15
3.8.2	Manual cleaning with ultrasound and immersion disinfection .....	16
3.9	Mechanical cleaning/disinfection .....	18
3.9.1	Mechanical alkaline cleaning and thermal disinfecting .....	18
3.10	Mechanical cleaning/disinfection with manual pre-cleaning .....	18
3.10.1	Manual pre-cleaning with a brush .....	19
3.10.2	Manual pre-cleaning with ultrasound and brush .....	19
3.10.3	Mechanical alkaline cleaning and thermal disinfecting .....	20
3.11	Inspection .....	21
3.11.1	Visual inspection .....	21
3.11.2	Functional test .....	21
3.12	Assembly .....	21
3.12.1	Straight Ball Hex Screwdriver (ME013R) and Straight Hex Screwdriver (ME016R) .....	21
3.12.2	U-Joint screwdriver (ME014R) and Implant Extraction Instrument (ME018R) .....	22
3.12.3	Implant Inserter Manipulator (ME015R) and Implant Introducer (SJ605R) .....	22
3.12.4	Trial Insertion Instrument (ME020R) .....	23
3.12.5	U-Joint Drill (ME189R) .....	23
3.12.6	U-Joint Bone Awl (ME190R) .....	24
3.13	Packaging .....	24
3.14	Steam sterilization .....	25
3.15	Storage .....	25
4.	Technical service .....	25
5.	Disposal .....	25

## 1. About this document

### Note

General risk factors associated with surgical procedures are not described in these instructions for use.

### 1.1 Scope

These instructions for use apply for the A-Space SIBD and Arcadius<sup>XP</sup> L instruments.

### Note

The applicable CE mark for the product can be seen on the label or packaging of the product.

### Note

Instructions for use and further information about B. Braun / AESCULAP products can be found on the B. Braun eIFU website at [efu.bbraun.com](http://efu.bbraun.com)

### 1.2 Safety messages

Safety messages make clear the dangers to patient, user and/or product that could arise during the use of the product. Safety messages are labeled as follows:

#### WARNING

Indicates a possible threat of danger. If not avoided, minor or moderate injury may result.

#### CAUTION

Indicates a possible threat of material damage. If not avoided, the product may be damaged.

## 2. Clinical use

### 2.1 Areas of use and limitations of use

#### 2.1.1 Intended use

The A-Space SIBD and Arcadius<sup>XP</sup> L Instruments are used for the implantation of A-Space SIBD and Arcadius<sup>XP</sup> L devices into the lumbar spine.

### 2.1.2 Indications

#### Note

*The manufacturer is not responsible for any use of the product against the specified indications and/or the described applications.*

For indications, see Intended use.

### 2.1.3 Contraindications

No known contraindications.

## 2.2 Safety information

### 2.2.1 Clinical user

#### General safety information

To prevent damage caused by improper setup or operation, and to not compromise the manufacturer warranty and liability:

- ▶ Use the product only according to these instructions for use.
- ▶ Follow the safety and maintenance instructions.
- ▶ Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge and experience.
- ▶ Store any new or unused products in a dry, clean, and safe place.
- ▶ Prior to use, check that the product is in good working order.
- ▶ Keep the instructions for use accessible for the user.

#### Note

*The user is obligated to report all severe events in connection with the product to the manufacturer and the responsible authorities of the state in which the user is located.*

#### Notes on surgical procedures

It is the user's responsibility to ensure that the surgical procedure is performed correctly.

Appropriate clinical training as well as a theoretical and practical proficiency of all the required operating techniques, including the use of this product, are prerequisites for the successful use of this product.

The user is required to obtain information from the manufacturer if there is an unclear preoperative situation regarding the use of the product.

### 2.2.2 Product-specific safety information

- ▶ Always follow the instructions for use of the A-Space SIBD and Arcadius<sup>XP</sup> L implants (TA013175 and TA015555) and the respective surgical techniques.
- ▶ Remove the transport packaging and clean the new product, either manually or mechanically, prior to its initial sterilization.
- ▶ Prior to each use, inspect the product for loose, bent, broken, cracked, worn, or fractured components.
- ▶ Do not use the product if it is damaged or defective. Set aside the product if it is damaged.
- ▶ Replace any damaged components immediately with original spare parts.

#### Products ME014R, ME018R, ME189R, ME190R

#### Note

*The product is partially made of NiCoCrMo (MP35N). The material contains cobalt.*

### 2.2.3 Sterility

The product is delivered in an unsterile condition.

- ▶ Clean the new product after removing its transport packaging and prior to its initial sterilization.

## 2.3 Application

### WARNING

**Risk of injury and/or malfunction!**

- ▶ **Prior to each use, inspect the product for loose, bent, broken, cracked, worn, or fractured components.**
- ▶ **Always carry out a function test prior to each use of the product.**

### 2.3.1 Distraction

#### WARNING

Risk of damage to the spinal cord, spinal nerve roots, adjacent intervertebral bodies or soft tissue when inserting and using the Distractor!

- ▶ Do not over insert the Distractor into the disc space.
- ▶ Do not over distract the disc space.
- ▶ Desired distraction of the vertebral bodies can be achieved by using the Distractors (SJ020R–SJ030R).
- ▶ Desired working height can be gradually achieved by coupling a Distractor to a T-Handle (SJ033R), inserting horizontally into the disc space and rotating by 90°.

### 2.3.2 Selection of Implant Size Using Trial Implants

- ▶ Attach selected Trial Implant to the Trial Insertion Instrument (ME020R) by inserting the Instrument into the anterior face of the Trial Implant. Turn the wheel clockwise to engage the Trial Implant onto the Trial Insertion Instrument.

#### WARNING

Risk of over distracting, damaging the endplates, or adjacent vertebral bodies in the course of incorrectly applying the Trial Insertion Instrument!

- ▶ During the Trial Implant insertion, do not apply an excessive bending/tilting/levering force to the Trial Insertion Instrument (ME020R).
- ▶ Keep the instrument directly in the plane of the disc.
- ▶ To advance the Trial Implant into the disc space use either a mallet or a Slap Hammer Extension (SJ709R) assembled onto a Slap Hammer Handle (SJ708T).

#### WARNING

If the Trial implant is inserted too deep there is a risk of damage to the spinal cord, spinal nerve roots, posterior structures!

- ▶ Always use the Trial Insertion Instrument (ME020R) with a depth stop.
- ▶ Do not apply excessive hammering force.
- ▶ Under lateral and anterior-posterior fluoroscopy confirm the:
  - Implant depth, height, and lordosis
  - Endplate coverage (anterior-posterior/medial-lateral)
  - Midline placement
  - Implant rotation

#### WARNING

Risk of selecting incorrect implant size!

- ▶ X-ray is mandatory to confirm the final Trial Implant size.
- ▶ Trial Implant must provide good contact with the inferior and superior endplates. Therefore the proper footprint, height, and lordotic angle should be chosen for the disc space.
- ▶ Select the implant size that corresponds with the final Trial Implant size chosen.

### 2.3.3 Filling the Implant

- ▶ Fill the implant with bone or bone replacement material by using the Packing Block (SJ604R) and Tamp (SJ608R).

### 2.3.4 Implant Placement and Positioning

- ▶ Securely attach the implant (corresponding to the final Trial Implant size) to the distal end of the insertion instrument by turning the wheel clockwise.

The laser markings (CRANIAL/CAUDAL) on the implant insertion instrument indicate the correct orientation for implant engagement:

- The Implant Inserter/Manipulator (ME015R) will grip the implant utilizing the two lateral implant screw holes.
- The Implant Introducer (SJ605R) will grip the implant utilizing the two medial implant screw holes.

#### WARNING

**Inaccurate marking of the midline may result in incorrect positioning of the implant!**

- ▶ Always mark the midline under X-ray visualization.
- ▶ Determine the center of the vertebral disc using the midline marker, under X-Ray visualization.

#### WARNING

**If the implant is inserted too deep, the spinal canal and other posterior elements may be compressed!**

- ▶ Always use the Implant Inserter/Manipulator (ME015R) with a depth stop.
- ▶ Excessive hammering may over-insert implant.
- ▶ Use appropriate care when inserting the implant.
  - Drive the Implant into the disc space using either a mallet or the Slap Hammer Extension (SJ709R) assembled onto the Slap Hammer Handle (SJ708T).

- ▶ Confirm an anatomically suitable position and orientation of the implant. It is recommended to confirm implant position prior to removing the implant inserter.
  - Obtain an AP fluoroscopic image to confirm midline placement of the device.
  - Obtain a lateral fluoroscopic image to confirm that the anterior edge of the implant is seated flush with the anterior border of the vertebral body.
  - Observe the X-Ray markers in both the AP and lateral views to ensure that the implant is not rotated within the disc space.
- ▶ Use Impactor (SJ606R) to adjust the position (in AP direction) of the implant as necessary.

### 2.3.5 Bone Screw Placement

- ▶ Use the bone awls (Straight Bone Awl (ME017R) or U-Joint Bone Awl (ME190R)) and the drills (Straight Drill (SJ725R) or U-Joint Drill (ME189R)) to prepare a pilot hole for the bone screws at the intended screw placement site.

#### WARNING

**Risk of damaging biological structures (spinal cord, spinal nerve roots, ligaments, vessels, soft tissues, etc.) with the drills and the bone awls, especially if the PEEK cage is placed off midline!**

- ▶ The instruments are sharp, always use the bone awls/drills under X-Ray control.

#### WARNING

**Risk of damaging biological structures (spinal cord, spinal nerve roots, ligaments, vessels, soft tissues, etc.) due to Straight Bone Awl tip or Straight Drill tip breaking!**

- ▶ Do not apply an excessive bending/tilting/levering force to the Straight Bone Awl (ME017R) or Straight Drill (SJ725R).

**⚠ WARNING**

Risk of damaging PEEK cage with incorrectly applied Straight Drill Guide (SJ724R)!

- ▶ Do not apply an excessive bending/tilting/levering force to the instrument.

**⚠ WARNING**

Applying a mallet to advance the U-Joint Bone Awl (ME190R) instrument into the bone may result in damage to the A-Space SIBD or the Arcadius<sup>XP</sup> L cages!

- ▶ Do not use a mallet with the U-Joint Bone Awl.
- ▶ The U-Joint Bone Awl must be advanced into the bone by hand only.
- ▶ Always use the Drills and Straight Bone Awl with the drill guides (Straight Drill Guide (SJ724R) or Angled Drill Guide (SJ722R)).

**⚠ WARNING**

Engaging the screwdriver incorrectly when turning the bone screw into the PEEK cage may result in damage to the bone screws!

- ▶ Fully insert the tip of the screwdriver into the bone screw.
- ▶ Attach bone screw of appropriate length to the selected self-retaining screwdriver (Straight Hex Screwdriver (ME016R), Straight Ball Hex Screwdriver (ME013R), or U-Joint Screwdriver (ME014R)).

**⚠ WARNING**

Risk of damaging biological structures (spinal cord, spinal nerve roots, ligaments, vessels, soft tissues, etc.) while advancing incorrectly sized bone screws, especially if the PEEK cage is placed off midline!

- ▶ Always advance the bone screws, especially the diverging lateral screws, under X-Ray control.
- ▶ Follow Surgical technique, or training by Aesculap.
- ▶ Place a bone screw through selected screw hole on the anterior face of the implant.

**⚠ WARNING**

Applying too much torque may result in damage to the bone screws and the PEEK cages!

- ▶ Do not apply an excessive screw tightening torque.

**⚠ WARNING**

Risk of damage to the PEEK cage, or the bone screw head in the course of incorrectly driving the bone screw!

- ▶ Do not apply an excessive bending/tilting/levering force to the Straight Hex Screwdriver (ME016R) or Straight Ball Hex Screwdriver (ME013R) engaged in the bone screw.
- ▶ Keep the screwdrivers directly aligned with the axis of the bone screws.
- ▶ Turn the screwdriver in a clockwise motion to advance the bone screw into the vertebral body.

**⚠ WARNING**

The screws could back out or loosen if they are not fully inserted into the cage!

- ▶ Insert the bone screw until it is fully engaged.
- ▶ Follow Surgical technique, or training by Aesculap.
- ▶ Insert the bone screws until they reach the final seated position, ensuring full engagement of the two locking mechanisms.

**⚠ WARNING**

Risk of insufficient stability or implant failure if fewer than four screws are used!

- ▶ Apply all four screws or use an additional supplemental spinal fixation system such as the Aesculap S<sup>4</sup> Spinal System.
- ▶ Repeat the bone screw placement steps outlined above to insert the remaining bone screws through the implant.

### 2.3.6 Revision

#### WARNING

Risk of injury to the patient!

- ▶ Remove implant in the direction it was inserted.
- ▶ Excessive force during extraction may cause implant to dislodge suddenly and impact the surrounding tissues.
- ▶ Locate implant, which is to be revised, and engage selected self-retaining screwdriver (Straight Hex Screwdriver (ME016R), Straight Ball Hex Screwdriver (ME013R), or U-Joint Screwdriver (ME014R)) in the bone screw head.
- ▶ Retract the bone screw from the vertebral body by turning the screw in a counterclockwise motion.
- ▶ Repeat the bone screw removal process with the remaining bone screws in the PEEK cage.
- ▶ Engage Implant Extraction Instrument (ME018R) or the Implant Inserter/Manipulator (ME015R) to the implant and remove it from the disc space.

#### Note

If a fully seated bone screw is removed from the implant, a small piece of PEEK debris from the locking rim in the locking mechanism may be present.

### 2.3.7 General

#### WARNING

Activating the "PRESS" button on the handle while actively using the instruments can result in the instrument shaft detaching from the handle!

- ▶ Do not cover the "PRESS" button with hand when utilizing fully assembled instruments.
- ▶ Caution should be taken while operating the instruments with quick connect modular handles: T-Handle (SJ033R), Handle (FW440R), Ratchet Handle (SJ705R) and Slap Hammer Handle (SJ708T).

#### WARNING

Risk of not being able to securely engage implant or damaging the instrument with incorrectly assembled Implant Inserter/Manipulator (ME015R) or Introducer (SJ605R)!

- ▶ Ensure that the instruments are correctly assembled prior to the surgery.

#### WARNING

Danger from improperly assembled U-Joint Drill (ME189R): the drill bit may disengage from the instrument shaft!

- ▶ Prior to the surgery, ensure that the drill bit is secured to the instrument shaft by firmly tightening the nut until it is fully threaded.

#### WARNING

Danger from improperly assembled Straight Hex Screwdriver (ME016R) or Straight Ball Hex Screwdriver (ME013R): the screwdriver bit may disengage from the instrument shaft!

- ▶ Prior to surgery, ensure that the screwdriver bit is secured to the instrument shaft, by firmly tightening the nut until it is fully threaded.
- ▶ Caution should be taken for assembling the following modular instruments: Implant Inserter/Manipulator (ME015R), Implant Introducer (SJ605R), U-Joint Drill (ME189R), Straight Hex Screwdriver (ME016R) and Straight Ball Hex Screwdriver (ME013R).

### 3. Validated reprocessing procedure

#### 3.1 General safety instructions

##### Note

*Adhere to national statutory regulations, national and international standards and directives, and local, clinical hygiene instructions for sterile processing.*

##### Note

*For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of products.*

##### Note

*Mechanical reprocessing should be favored over manual cleaning as it gives better and more reliable results.*

##### Note

*Successful processing of this medical device can only be ensured if the processing method is first validated. The operator/sterile processing technician is responsible for this.*

##### Note

*If there is no final sterilization, then a virucidal disinfectant must be used.*

##### Note

*For up-to-date information about reprocessing and material compatibility, see B. Braun eFU at [eifu.bbraun.com](http://eifu.bbraun.com)*

*The validated steam sterilization procedure was carried out in the Aesculap sterile container system.*

#### 3.2 General information

Dried or affixed surgical residues can make cleaning more difficult or ineffective and lead to corrosion. Therefore the time interval between application and processing should not exceed 6 h; also, neither fixating pre-cleaning temperatures >45 °C nor fixating disinfecting agents (active ingredient: aldehydes/alcohols) should be used.

Excessive measures of neutralizing agents or basic cleaners may result in a chemical attack and/or to fading and the laser marking becoming unreadable visually or by machine for stainless steel.

Residues containing chlorine or chlorides e.g. in surgical residues, medicines, saline solutions and in the service water used for cleaning, disinfection and sterilization will cause corrosion damage (pitting, stress corrosion) and result in the destruction of stainless steel products. These must be removed by rinsing thoroughly with demineralized water and then drying.

Additional drying, if necessary.

Only process chemicals that have been tested and approved (e.g. VAH or FDA approval or CE mark) and which are compatible with the product's materials according to the chemical manufacturers' recommendations may be used for processing the product. All the chemical manufacturer's application specifications must be strictly observed. Failure to do so can result in the following problems:

- Optical changes of materials, e.g. fading or discoloration of titanium or aluminum. For aluminum, the application/process solution only needs to be of pH >8 to cause visible surface changes.
- Material damage such as corrosion, cracks, fracturing, premature aging or swelling.
- ▶ Do not use metal cleaning brushes or other abrasives that would damage the product surfaces and could cause corrosion.
- ▶ Further detailed advice on hygienically safe and material-/value-preserving reprocessing can be found at [www.a-k-i.org](http://www.a-k-i.org), link to "AKI-Brochures", "Red brochure".

### 3.3 Reusable products

Influences of the reprocessing which lead to damage to the product are not known.

A careful visual and functional inspection before the next use is the best opportunity to recognize a product that is no longer functional, see Inspection.

### 3.4 Preparations at the place of use

- ▶ If applicable, rinse non-visible surfaces preferably with deionized water, with a disposable syringe for example.
- ▶ Remove any visible surgical residues to the extent possible with a damp, lint-free cloth.
- ▶ Transport the dry product in a sealed waste container for cleaning and disinfection within 6 hours.

### 3.5 Preparing for cleaning

- ▶ Disassemble the product prior to cleaning, see Disassembly.

### 3.6 Disassembly

#### 3.6.1 Straight Ball Hex Screwdriver (ME013R) and Straight Hex Screwdriver (ME016R)

- ▶ Unthread the nut **24** from the distal end of the shaft.
- ▶ Remove the screwdriver bit **25** from the shaft.
- ▶ Remove the soft tissue protection sleeve **26** from the shaft, see Fig. 1.


Fig. 1 Disassembling ME013R and ME016R

#### 3.6.2 U-Joint screwdriver (ME014R), U-Joint Bone Awl (ME190R) and Implant Extraction Instrument (ME018R)

##### ⚠ WARNING

Risk of injury and/or damage to the gloves and/or spring when removing the spring!

- ▶ Do not use too much force to remove the spring.
- ▶ Apply both rotational and translational motion to the spring.
- ▶ Wear two pairs of gloves while removing the spring.
- ▶ Remove the spring **27** covering the u-joint from the distal end of the shaft, see Fig. 2, see Fig. 3 and see Fig. 4.


Fig. 2 Disassembling ME014R


Fig. 3 Disassembling ME018R


Fig. 4 Disassembling ME190R

### 3.6.3 Implant Inserter Manipulator (ME015R) and Implant Introducer (SJ605R)

- ▶ While holding the proximal end (mechanism with thumb wheel) upright, unthread the nut **28** counterclockwise towards the proximal end.
- ▶ With the nut disengaged from its threading, turn the thumb wheel **29** clockwise ("LOCK" - direction) until the sub-assembly **30** disengages from the outer sheath **31**.
- ▶ Pull out the sub-assembly **30** from the outer sheath **31**.
- ▶ Remove the depth stopper **32** by pushing it distally parallel along the shaft until it disengages from its rail, see Fig. 5 and see Fig. 6.


Fig. 5 Disassembling ME015R


Fig. 6 Disassembling SJ605R

### 3.6.4 Trial Insertion Instrument (ME020R)

- ▶ Press the locking mechanism button **33** near the thumb wheel and pull out the outer sheath **34**.
- ▶ Remove the inner shaft **35** from the thumb wheel.
- ▶ Remove the depth stopper **36** by pushing it distally parallel along the shaft until it disengages from its rail.


Fig. 7 Disassembling ME020R

### 3.6.5 U-Joint Drill (ME189R)

#### **⚠ WARNING**

Risk of injury and/or damage to the gloves and/or spring when removing the spring!

- ▶ Do not use too much force to remove the spring.
- ▶ Apply both rotational and translational motion to the spring.
- ▶ Wear two pairs of gloves while removing the spring.
- ▶ Unthread the nut **37** from the distal end of the shaft.
- ▶ Remove the drill bit **38** from the shaft end and slide the spring **39** off the u-joint, see Fig. 8.


Fig. 8 Disassembling ME189R

## 3.7 Cleaning/Disinfection

### 3.7.1 Product-specific safety information on the reprocessing method

Damage to or destruction of the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!

- ▶ Following the manufacturer's instructions, use cleaning and disinfecting agents that are approved for high-grade steel.
- ▶ Observe specifications regarding concentration, temperature and exposure time.
- ▶ Manually pre-clean the product (using a brush) if bone or tissue residues or residues of additives (e.g. flasking plaster) are present.

### 3.7.2 Validated cleaning and disinfection procedure

Validated procedure	Specific requirements	Reference
Manual cleaning with immersion disinfection <ul style="list-style-type: none"> <li>■ ME013R</li> <li>■ ME014R</li> <li>■ ME016R</li> <li>■ ME017R</li> <li>■ ME018R</li> <li>■ ME189R</li> <li>■ ME190R</li> <li>■ SJ033R</li> <li>■ SJ604R</li> <li>■ SJ606R</li> <li>■ SJ608R</li> <li>■ SJ722R</li> <li>■ SJ724R</li> <li>■ SJ725R</li> <li>■ SJ020R – SJ030R (6 sizes)</li> <li>■ SJ664T – SJ798T (36 sizes)</li> </ul>	<ul style="list-style-type: none"> <li>■ Cleaning brush: e.g. TA011944, TE654202, GK469200</li> <li>■ Disposable syringe 20 ml</li> <li>■ When cleaning instruments with movable hinges, ensure that these are in an open position and, if applicable, move the hinge while cleaning.</li> <li>■ Drying phase: Use a lint-free cloth or medical compressed air</li> </ul>	Chapter Manual cleaning/disinfection and subsection: <ul style="list-style-type: none"> <li>■ Chapter Manual cleaning with immersion disinfection</li> </ul>
Manual cleaning with ultrasound and immersion disinfection <ul style="list-style-type: none"> <li>■ ME015R</li> <li>■ ME020R</li> <li>■ SJ605R</li> <li>■ SJ705R</li> <li>■ SJ708T</li> <li>■ SJ709R</li> <li>■ SJ724R</li> </ul>	<ul style="list-style-type: none"> <li>■ Cleaning brush: e.g. TA011944, TE654202, GK469200</li> <li>■ Disposable syringe 20 ml</li> <li>■ When cleaning instruments with movable hinges, ensure that these are in an open position and, if applicable, move the hinge while cleaning.</li> <li>■ Drying phase: Use a lint-free cloth or medical compressed air</li> </ul>	Chapter Manual cleaning/disinfection and subsection: <ul style="list-style-type: none"> <li>■ Chapter Manual cleaning with ultrasound and immersion disinfection</li> </ul>

Validated procedure	Specific requirements	Reference
Mechanical alkaline cleaning and thermal disinfection <ul style="list-style-type: none"> <li>■ ME014R</li> <li>■ ME017R</li> <li>■ ME018R</li> <li>■ ME189R</li> <li>■ ME190R</li> <li>■ SJ033R</li> <li>■ SJ604R</li> <li>■ SJ606R</li> <li>■ SJ608R</li> <li>■ SJ725R</li> <li>■ SJ020R – SJ030R (6 sizes)</li> </ul>	<ul style="list-style-type: none"> <li>■ Place the product on a tray that is suitable for cleaning (avoid rinsing blind spots).</li> <li>■ To flush the product: Use a flushing nozzle or flushing sleeve.</li> <li>■ Place the product on the tray with all product links and joints open.</li> </ul>	Chapter Mechanical cleaning/disinfection and subsection: <ul style="list-style-type: none"> <li>■ Chapter Mechanical alkaline cleaning and thermal disinfecting</li> </ul>
Manual pre-cleaning with brush and subsequent mechanical alkaline cleaning and thermal disinfection <ul style="list-style-type: none"> <li>■ ME013R</li> <li>■ ME016R</li> <li>■ SJ705R</li> <li>■ SJ709R</li> <li>■ SJ722R</li> <li>■ SJ724R</li> <li>■ SJ708T</li> <li>■ SJ664T – SJ798T (36 sizes)</li> </ul>	<ul style="list-style-type: none"> <li>■ Cleaning brush: e.g. TA011944, TE654202, GK469200</li> <li>■ Place the product on a tray that is suitable for cleaning (avoid rinsing blind spots).</li> <li>■ To flush the product: Use a flushing nozzle or flushing sleeve.</li> <li>■ Place the product on the tray with all product links and joints open.</li> </ul>	Chapter Mechanical cleaning/disinfection with manual pre-cleaning and subsection: <ul style="list-style-type: none"> <li>■ Chapter Manual pre-cleaning with a brush</li> <li>■ Chapter Mechanical alkaline cleaning and thermal disinfecting</li> </ul>
Manual pre-cleaning with ultrasound and brush, and subsequent mechanical alkaline cleaning and thermal disinfection <ul style="list-style-type: none"> <li>■ FW440R</li> <li>■ ME015R</li> <li>■ ME020R</li> <li>■ SJ605R</li> </ul>	<ul style="list-style-type: none"> <li>■ Cleaning brush: e.g. TA011944, TE654202, GK469200</li> <li>■ Disposable syringe 20 ml</li> <li>■ Place the product on a tray that is suitable for cleaning (avoid rinsing blind spots).</li> <li>■ To flush the product: Use a flushing nozzle or flushing sleeve.</li> <li>■ Place the product on the tray with all product links and joints open.</li> </ul>	Chapter Mechanical cleaning/disinfection with manual pre-cleaning and subsection: <ul style="list-style-type: none"> <li>■ Chapter Manual pre-cleaning with ultrasound and brush</li> <li>■ Chapter Mechanical alkaline cleaning and thermal disinfecting</li> </ul>

### 3.8 Manual cleaning/disinfection

- ▶ Prior to manual disinfecting, allow water to drip off for a sufficient length of time to prevent dilution of the disinfecting solution.
- ▶ After manual cleaning/disinfection, check visible surfaces visually for residues.
- ▶ Repeat the cleaning/disinfection process if necessary.

#### 3.8.1 Manual cleaning with immersion disinfection

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	<b>Disinfecting cleaning</b>	RT (cold)	>15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
II	<b>Intermediate rinse</b>	RT (cold)	1	-	D-W	-
III	<b>Disinfection</b>	RT (cold)	5	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
IV	<b>Final rinse</b>	RT (cold)	1	-	FD-W	-
V	<b>Drying</b>	RT	-	-	-	-

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

RT: Room temperature

\*Recommended: BBraun Stabimed fresh

- ▶ Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

#### Phase I

- ▶ Fully immerse the product in the cleaning/disinfectant for at least 15 min. Ensure that all accessible surfaces are moistened.
- ▶ Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- ▶ If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
- ▶ Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- ▶ Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

#### Phase II

- ▶ Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- ▶ Drain any remaining water fully.

#### Phase III

- ▶ Fully immerse the product in the disinfectant solution.
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- ▶ Rinse lumens at least 5 times at the beginning of the exposure time using an appropriate disposable syringe. Ensure that all accessible surfaces are moistened.

**Phase IV**

- ▶ Rinse/flush the product thoroughly (all accessible surfaces).
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during final rinse.
- ▶ Rinse lumens with an appropriate disposable syringe at least five times.
- ▶ Drain any remaining water fully.

**Phase V**

- ▶ Dry the product in the drying phase with suitable equipment (e.g. cloth, compressed air), see Validated cleaning and disinfection procedure.

**3.8.2 Manual cleaning with ultrasound and immersion disinfection**

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Ultrasonic cleaning	RT (cold)	>15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
II	Intermediate rinse	RT (cold)	1	-	D-W	-
III	Disinfection	RT (cold)	5	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
IV	Final rinse	RT (cold)	1	-	FD-W	-
V	Drying	RT	-	-	-	-

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

RT: Room temperature

\*Recommended: BBraun Stabimed fresh

- ▶ Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

**Phase I**

- ▶ Clean the product in an ultrasonic cleaning bath (frequency 35 kHz) for at least 15 min. Ensure that all accessible surfaces are immersed and acoustic shadows are avoided.
- ▶ Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- ▶ If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
- ▶ Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- ▶ Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

**Phase II**

- ▶ Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- ▶ Drain any remaining water fully.

**Phase III**

- ▶ Fully immerse the product in the disinfectant solution.
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- ▶ Rinse lumens at least five times at the beginning of the exposure time with an appropriate disposable syringe. Ensure that all accessible surfaces are moistened.

**Phase IV**

- ▶ Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during final rinse.
- ▶ Rinse lumens with an appropriate disposable syringe at least five times.
- ▶ Drain any remaining water fully.

**Phase V**

- ▶ Dry the product in the drying phase with suitable equipment (e.g. cloth, compressed air), see Validated cleaning and disinfection procedure.

### 3.9 Mechanical cleaning/disinfection

#### Note

The cleaning and disinfection device must be of tested and approved effectiveness (e.g. FDA approval or CE mark according to DIN EN ISO 15883).

#### Note

The cleaning and disinfection device used for processing must be serviced and checked at regular intervals.

#### 3.9.1 Mechanical alkaline cleaning and thermal disinfecting

Machine type: single-chamber cleaning/disinfection device without ultrasound

Phase	Step	T [°C/°F]	t [min]	Water quality	Chemical/Note
I	Prerinse	<25/77	3	D-W	-
II	Cleaning	55/131	10	FD-W	<ul style="list-style-type: none"> <li>■ Concentrate, alkaline:               <ul style="list-style-type: none"> <li>- pH = 13</li> <li>- &lt;5 % anionic surfactant</li> </ul> </li> <li>■ working solution 0.5%               <ul style="list-style-type: none"> <li>- pH = 11*</li> </ul> </li> </ul>
III	Intermediate rinse	> 10/50	1	FD-W	-
IV	Thermal disinfecting	90/194	5	FD-W	-
V	Drying	-	-	-	In accordance with the program for the cleaning and disinfecting machine

DW: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

\*Recommended: BBraun Helimatic Cleaner alkaline

- ▶ Check visible surfaces for residues after mechanical cleaning/disinfecting.

### 3.10 Mechanical cleaning/disinfection with manual pre-cleaning

#### Note

The cleaning and disinfection device must be of tested and approved effectiveness (e.g. FDA approval or CE mark according to DIN EN ISO 15883).

#### Note

The cleaning and disinfection device used for processing must be serviced and checked at regular intervals.

### 3.10.1 Manual pre-cleaning with a brush

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	<b>Disinfecting cleaning</b>	RT (cold)	>15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
II	<b>Rinsing</b>	RT (cold)	1	-	D-W	-

D-W: Drinking water

RT: Room temperature

\*Recommended: BBraun Stabimed fresh

- ▶ Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

#### Phase I

- ▶ Fully immerse the product in the cleaning/disinfectant for at least 15 min. Ensure that all accessible surfaces are moistened.
- ▶ Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- ▶ If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
- ▶ Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- ▶ Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

#### Phase II

- ▶ Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.

### 3.10.2 Manual pre-cleaning with ultrasound and brush

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	<b>Ultrasonic cleaning</b>	RT (cold)	>15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
II	<b>Rinsing</b>	RT (cold)	1	-	D-W	-

D-W: Drinking water

RT: Room temperature

\*Recommended: BBraun Stabimed fresh

- ▶ Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

**Phase I**

- ▶ Clean the product in an ultrasonic cleaning bath (frequency 35 kHz) for at least 15 min. Ensure that all accessible surfaces are immersed and acoustic shadows are avoided.
- ▶ Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- ▶ If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
- ▶ Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- ▶ Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

**Phase II**

- ▶ Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.

**3.10.3 Mechanical alkaline cleaning and thermal disinfecting**

Machine type: single-chamber cleaning/disinfection device without ultrasound

Phase	Step	T [°C/°F]	t [min]	Water quality	Chemical
I	Prerinse	<25/77	3	D-W	-
II	Cleaning	55/131	10	FD-W	<ul style="list-style-type: none"> <li>■ Concentrate, alkaline:                             <ul style="list-style-type: none"> <li>- pH = 13</li> <li>- &lt;5 % anionic surfactant</li> </ul> </li> <li>■ 0.5 % working solution                             <ul style="list-style-type: none"> <li>- pH ~ 11*</li> </ul> </li> </ul>
III	Intermediate rinse	>10/50	1	FD-W	-
IV	Thermal disinfection	90/194	5	FD-W	-
V	Drying	-	-	-	According to the program for cleaning and disinfection device

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

\*Recommended: BBraun Helimatic Cleaner alkaline

- ▶ Check visible surfaces for residues after mechanical cleaning/disinfecting.

### 3.11 Inspection

- ▶ Allow the product to cool down to room temperature.
- ▶ Dry the product if it is wet or damp.

#### 3.11.1 Visual inspection

- ▶ Ensure that all soiling has been removed. In particular, pay attention to mating surfaces, hinges, shafts, recessed areas, drill grooves and the sides of the teeth on rasps.
- ▶ If the product is dirty: repeat the cleaning and disinfection process.
- ▶ Check the product for damage, e.g. insulation or corroded, loose, bent, broken, cracked, worn or severely scratched and fractured components.
- ▶ Check the product for missing or faded labels.
- ▶ Check the products with long, slim shapes (in particular rotating instruments) for deformities.
- ▶ Check the product for damage to the spiral element.
- ▶ Check the cutting edges for continuity, sharpness, nicks and other damage.
- ▶ Check the surfaces for rough spots.
- ▶ Check the product for burrs that could damage tissue or surgical gloves.
- ▶ Check the product for loose or missing parts.
- ▶ Immediately put aside damaged or inoperative products and send them to Aesculap Technical Service, see Technical service.

### 3.11.2 Functional test

#### CAUTION

**Damage (metal cold welding/friction corrosion) to the product caused by insufficient lubrication!**

- ▶ Prior to function checks, lubricate moving parts (e.g. joints, pusher components and threaded rods) at the marked points with maintenance oil suitable for the respective sterilization process (e.g. for steam sterilization: STERILIT® I oil spray JG600 or STERILIT® I drip lubricator JG598).
- ▶ Assemble disassembled products, see Assembly.
- ▶ Check that the product functions correctly.
- ▶ Check that all moving parts are working properly (e.g. hinges, locks/latches, sliding parts etc.).
- ▶ Check rotating products (e.g. reusable drills and cutters) for bends and deformities. To do this, roll the product on an even surface.
- ▶ Check for compatibility with associated products.
- ▶ Immediately put aside inoperative products and send them to Aesculap Technical Service, see Technical service.

### 3.12 Assembly

#### 3.12.1 Straight Ball Hex Screwdriver (ME013R) and Straight Hex Screwdriver (ME016R)

#### WARNING

**Danger from improperly assembled Straight Hex Screwdriver (ME016R) or Straight Ball Hex Screwdriver (ME013R): the screwdriver bit may disengage from the instrument shaft during shipment or operation!**

- ▶ Ensure the screwdriver bit **25** is secured to the shaft, by firmly tightening the nut **24** until it is fully threaded.
- ▶ Place the soft tissue protection sleeve **26** over the inner shaft.
- ▶ Insert the screwdriver bit **25** into the internal hexagon of the distal end of the shaft.

- ▶ Firmly thread the nut **24** proximally over the screwdriver bit until it is fully threaded, see Fig. 9.


Fig. 9 Assembling ME013R and ME016R

### 3.12.2 U-Joint screwdriver (ME014R) and Implant Extraction Instrument (ME018R)

- ▶ Slide the spring **27** over the distal end until it is completely covering the u-joint, see Fig. 10 and see Fig. 11. The ends of the spring should sit within the flange on the shaft and the flange on the distal head, see Fig. 12.


Fig. 10 Assembling ME014R


Fig. 11 Assembling ME018R


Fig. 12 Correct spring position

### 3.12.3 Implant Inserter Manipulator (ME015R) and Implant Introducer (SJ605R)

#### ⚠ WARNING

Risk of not being able to securely engage the implant or damaging the instrument with incorrectly assembled Implant Inserter/Manipulator (ME015R) or Introducer (SJ605R)!

- ▶ Ensure that there is no gap between components **30** and **31**.
- ▶ Firmly tighten the nut **28** until it is fully threaded.
- ▶ Place the depth stopper **32** on its rail and press it proximally along the shaft until it is fully engaged.
- ▶ Insert the sub-assembly **30** all the way into the outer sheath **31**.
- ▶ Start turning the thumb wheel **29** counterclockwise ("RELEASE"-direction) while pressing the sub-assembly towards the thumb wheel. The threads between the end of the sub-assembly shaft and the inner thread in the thumb wheel will start to engage.
- ▶ Continue turning the thumb wheel until the protrusions of the sub-assembly are flush with the slots in the outer sheath.
- ▶ Visually verify that both components (**30** and **31**) are fully engaged, see Fig. 15.
- ▶ Firmly thread the nut **28** towards the distal end until it is fully threaded, see Fig. 13 and see Fig. 14.


Fig. 13 Assembling ME015R


Fig. 14 Assembling SJ605R


Fig. 15 Correct engagement of sub-assembly and outer sheath

### 3.12.4 Trial Insertion Instrument (ME020R)

- ▶ Place the depth stop **36** on its rail and press it proximally along the shaft until it is fully engaged.
- ▶ Insert the hexagonal end of the inner shaft **35** into the inner hexagon of the thumb wheel **37**.
- ▶ Place the outer sheath **34** over the inner shaft **35**.
- ▶ Press the locking mechanism button **33** and fully insert the outer sheath into the locking mechanism (until it reaches the hard stop).
- ▶ Release the button **33**, see Fig. 16.


Fig. 16 Assembling ME020R

### 3.12.5 U-Joint Drill (ME189R)

#### ⚠ WARNING

Risk of using wrong spring!

- ▶ Check markings on the spring before assembling the instrument. The spring is marked with the item numbers ME189R and ME190R.

#### ⚠ WARNING

Danger from improperly assembled U-Joint Drill (ME189R): the drill bit may disengage from the instrument shaft during shipment or operation!

- ▶ Ensure that the drill bit **38** is secured to instrument shaft by firmly tightening the nut **37** until it is fully threaded.
- ▶ Slide the spring **39** over the distal end until it is completely covering the u-joint. The end of the spring must sit in the undercut **41**.

- ▶ Insert the drill bit **38** into the internal hexagon of the distal end of the u-joint shaft.
- ▶ Firmly thread the nut **37** over the drill bit towards the proximal end until it is fully threaded, see Fig. 17 and see Fig. 18.


Fig. 17 Assembling ME189R


Fig. 18 Correct spring position

### 3.12.6 U-Joint Bone Awl (ME190R)

**⚠ WARNING**

Risk of using wrong spring!

- ▶ Check markings on the spring before assembling the instrument. The spring is marked with the item numbers ME189R and ME190R.

**⚠ WARNING**

Risk of injury during spring assembly due to sharp bone awl tip.

- ▶ Use extra caution when sliding the spring over the sharp tip.

**⚠ WARNING**

Danger from improperly assembled U-Joint Bone Awl (ME190R): the spring may disengage from the instrument shaft during shipment or operation!

- ▶ Ensure the end of the spring is seated in the undercut.

- ▶ Slide the spring **27** over the distal end until it is completely covering the u-joint. The end of the spring must sit in the undercut **42**, see Fig. 19 and see Fig. 20.


Fig. 19 Assembling ME190R


Fig. 20 Correct spring position

### 3.13 Packaging

- ▶ Appropriately protect products with fine working tips for the following instruments: ME190R, ME017R, ME189R, SJ725R, ME013R, ME014R and ME016R.
- ▶ Place the product in its holder or on a suitable tray. Ensure that sharp edges are covered.
- ▶ Package trays appropriately for the sterilization process (e.g. in Aesculap sterile containers).
- ▶ Ensure that the packaging provides sufficient protection against contamination of the product during storage.

### 3.14 Steam sterilization

- ▶ Check to ensure that the sterilizing agent will come into contact with all external and internal surfaces (e.g., by opening any valves and faucets).
- ▶ Validated sterilization process
  - Steam sterilization using fractional vacuum process
  - Steam sterilizer according to DIN EN 285 and validated according to DIN EN ISO 17665
  - Sterilization using fractionated vacuum process at 134 °C/holding time 5 min
- ▶ If several devices are sterilized at the same time in the same steam sterilizer: Ensure that the maximum permitted load according to the manufacturers' specifications is not exceeded.

### 3.15 Storage

- ▶ Store sterile products in germ-proof packaging, protected from dust, in a dry, dark, temperature-controlled area.

## 4. Technical service

### CAUTION

**Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.**

- ▶ **Do not modify the product.**
- ▶ **For service and repairs, please contact your national B. Braun/Aesculap agency.**

#### Service addresses

Aesculap Technischer Service

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78532 Tuttlingen / Germany

Phone: +49 7461 95-1601

Fax: +49 7461 16-2887

E-Mail: [ats@aesculap.de](mailto:ats@aesculap.de)

Other service addresses can be obtained from the address indicated above.

## 5. Disposal

### WARNING

**Risk of infection due to contaminated products!**

- ▶ **Adhere to national regulations when disposing of or recycling the product, its components and its packaging.**

### WARNING

**Risk of injury due to sharp-edged and/or pointed products!**

- ▶ **When disposing of or recycling the product, ensure that the packaging prevents injury by the product.**

#### Note

*The user institution is obliged to reprocess the product before its disposal, see Validated reprocessing procedure.*

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