# AESCULAP®

**en Instructions for use/Technical description**  
Instruments for CASPAR anterior cervical fusion

**USA Note for U.S. users**  
This Instructions for Use is NOT intended for United States users. Please discard. The Instructions for Use for United States users can be obtained by visiting our website at [www.aesculapimplantsystemsifus.com](http://www.aesculapimplantsystemsifus.com). If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-866-229-3002. A paper copy will be provided to you upon request at no additional cost.

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




## 1. About this document

**Note**  
General risk factors associated with surgical procedures are not described in these instructions for use.

### 1.1 Scope

These instructions for use apply for instruments for CASPAR anterior cervical fusion.  
► For article-specific instructions for use as well as information on material compatibility and lifetime see B. Braun eIFU at eifu.bbraun.com

### 1.2 Safety messages

Safety messages make clear the dangers to patient, user and/or product that could arise during the use of the product. Safety messages are labeled as follows:

#### ⚠ WARNING

Indicates a possible threat of danger. If not avoided, minor or moderate injury may result.

#### ⚠ CAUTION

Indicates a possible threat of material damage. If not avoided, the product may be damaged.

## 2. Clinical use

### 2.1 Areas of use and limitations of use

#### 2.1.1 Intended use

The instruments for anterior cervical fusion have been designed for cervical spine plating with CASPAR plates and screws for anterior cervical fusion (ACF).

They are used for positioning, drilling, tapping, and measuring during application of the CASPAR plates and screws.

#### 2.1.2 Indications

**Note**  
The manufacturer is not responsible for any use of the product against the specified indications and/or the described applications.

For indications, see Intended use.

#### 2.1.3 Contraindications

No known contraindications.

### 2.2 Safety information

#### 2.2.1 Clinical user

##### General safety information

To prevent damage caused by improper setup or operation, and to not compromise the manufacturer warranty and liability:

- Use the product only according to these instructions for use.
- Follow the safety and maintenance instructions.
- Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge and experience.
- Store any new or unused products in a dry, clean, and safe place.
- Prior to use, check that the product is in good working order.
- Keep the instructions for use accessible for the user.

**Note**  
The user is obligated to report all severe events in connection with the product to the manufacturer and the responsible authorities of the state in which the user is located.

##### Notes on surgical procedures

It is the user's responsibility to ensure that the surgical procedure is performed correctly.

Appropriate clinical training as well as a theoretical and practical proficiency of all the required operating techniques, including the use of this product, are prerequisites for the successful use of this product.

The user is required to obtain information from the manufacturer if there is an unclear preoperative situation regarding the use of the product.

#### 2.2.2 Product

##### Product-specific safety information

Risk of injury and/or damage to the instruments due to incorrect application!

- Use the instruments for anterior cervical fusion only in connection with CASPAR cervical spine plates and screws for anterior cervical fusion.

Application of the instruments requires detailed knowledge of spine surgery and spine biomechanics. The operating surgeon must be familiar with all theoretical and practical aspects of the operating technique for the CASPAR evolution system.

- The patient must be informed about the advantages and disadvantages of anterior cervical fusion.
- Read and follow the operating manual.

**Note**  
The operating manual can be ordered under article number O43102 from the address given below.

#### 2.2.3 Sterility

The product is delivered in an unsterile condition.

- Clean the new product after removing its transport packaging and prior to its initial sterilization.

## 2.3 Application

#### ⚠ WARNING

Risk of injury and/or malfunction!

- Prior to each use, inspect the product for loose, bent, broken, cracked, worn, or fractured components.
- Always carry out a function test prior to each use of the product.

#### ⚠ WARNING

Risk of injury when using the product beyond the field of view!

- Apply the product only under visual control.

#### 2.3.1 Single drill guide FF885R (Fig. 1)

► Adjust drill sleeve 1 to the required drilling depth (e.g., 25 mm) by turning in +/- direction.  
The drill C is limited by the stop 2 on the drill sleeve 1 once the set drilling depth 6 is reached. The set drilling depth is reached in conjunction with a CASPAR plate 5.

**Note**  
Always check the set drilling depth in advance using a caliper. To do this, insert drill C into drill sleeve 1.

#### 2.3.2 Double drill guide FF886R (Fig. 2)

- Adjust drill sleeve 1 to the required drilling depth (e.g., 24 mm) by turning in +/- direction.  
The drill C is limited by the stop 2 on the drill sleeve 1 once the set drilling depth 6 is reached. The set drilling depth is reached in conjunction with a CASPAR plate 5.

**Note**  
Always check the set drilling depth in advance using a caliper. To do this, insert drill C into drill sleeve 1.

#### 2.3.3 Drill FG412R


Fig. 1 Drill FG412R

#### Legend

- 2 Stop
- 3 Coupling stem  $\varnothing$  2,35 mm
- 7 Effective length

#### ⚠ WARNING

Risk of injury and/or damage to instruments due to changes in the drilling depth/effective length 7 by repointing drill C! Agreement between the set drilling depth and the actual drilling depth will not be ensured anymore.

- Do not repoint drill C.

**Note**  
Only use the drill in conjunction with drill guides A and B.

**Note**  
The drill can be used in conjunction with the INTRA handpiece GD456 and the Aesculap microsPEED® uni or ELAN EC motor system.

#### 2.3.4 Tap FF975/FG413R (Fig. 3)

- Pull spacer sleeve 1 in the direction of the arrow against spring pressure.  
The form fit between pin 2 and flute 3 is released.
- Turn spacer sleeve 1 to the required cutting depth 4 (e.g., 22 mm) and let it slide back into the form fit.

#### 2.3.5 Depth gauge FF965R (Fig. 4)

Interaction of the depth gauge with a CASPAR plate 5 and a CASPAR screw 8 for measuring example 9 of a 20 mm CASPAR screw:

- Loosen knurled nut 2.
- Slide caliper 1 in the required direction (e.g., 20 mm).
- Tighten knurled nut 2.  
The gage is fixed at the measured depth of 20 mm.

## 3. Validated reprocessing procedure

### 3.1 General safety instructions

**Note**  
Adhere to national statutory regulations, national and international standards and directives, and local, clinical hygiene instructions for sterile processing.

**Note**  
For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of products.

**Note**  
Mechanical reprocessing should be favored over manual cleaning as it gives better and more reliable results.

**Note**  
Successful processing of this medical device can only be ensured if the processing method is first validated. The operator/sterile processing technician is responsible for this.

**Note**  
If there is no final sterilization, then a virucidal disinfectant must be used.

**Note**  
For up-to-date information about reprocessing and material compatibility, see B. Braun eIFU at eifu.bbraun.com. The validated steam sterilization procedure was carried out in the Aesculap sterile container system.

### 3.2 General information

Dried or affixed surgical residues can make cleaning more difficult or ineffective and lead to corrosion. Therefore the time interval between application and processing should not exceed 6 h; also, neither fixating pre-cleaning temperatures >45 °C nor fixating disinfecting agents (active ingredient: aldehydes/alcohols) should be used.

Excessive measures of neutralizing agents or basic cleaners may result in a chemical attack and/or to fading and the laser marking becoming unreadable visually or by machine for stainless steel.

Residues containing chlorine or chlorides e.g. in surgical residues, medicines, saline solutions and in the service water used for cleaning, disinfection and sterilization will cause corrosion damage (pitting, stress corrosion) and result in the destruction of stainless steel products. These must be removed by rinsing thoroughly with demineralized water and then drying.

Additional drying, if necessary.

Only process chemicals that have been tested and approved (e.g. VAH or FDA approval or CE mark) and which are compatible with the product's materials according to the chemical manufacturers' recommendations may be used for processing the product. All the chemical manufacturer's application specifications must be strictly observed. Failure to do so can result in the following problems:

- Optical changes of materials, e.g. fading or discoloration of titanium or aluminum. For aluminum, the application/process solution only needs to be of pH >8 to cause visible surface changes.
- Material damage such as corrosion, cracks, fracturing, premature aging or swelling.
- Do not use metal cleaning brushes or other abrasives that would damage the product surfaces and could cause corrosion.
- Further detailed advice on hygienically safe and material-/value-preserving reprocessing can be found at [www.a-k-i.org](http://www.a-k-i.org), link to "AKI-Brochures", "Red brochure".

### 3.3 Reusable products

Influences of the reprocessing which lead to damage to the product are not known. A careful visual and functional inspection before the next use is the best opportunity to recognize a product that is no longer functional, see Inspection.

### 3.4 Preparations at the place of use

- If applicable, rinse non-visible surfaces preferably with deionized water, with a disposable syringe for example.
- Remove any visible surgical residues to the extent possible with a damp, lint-free cloth.
- Transport the dry product in a sealed waste container for cleaning and disinfection within 6 hours.

### 3.5 Preparing for cleaning

- Disassemble the product immediately after use, as described in the respective instructions for use.
- Disassemble the product prior to cleaning, see Disassembly.
- Open up product with hinges.

### 3.6 Disassembly

#### 3.6.1 Single drill guide FF885R (Fig. 5)

- ▶ Unscrew drill sleeve 1 clockwise.

##### Note

Drill sleeve 1 audibly engages at every half turn.

##### Note

The drill sleeves 1 of drill guides A and B are interchangeable.

#### 3.6.2 Double drill guide FF886R (Fig. 6)

- ▶ Unscrew drill sleeve 1 clockwise.

##### Note

Drill sleeve 1 audibly engages at every half turn.

##### Note

The drill sleeves 1 of drill guides A and B are interchangeable.

#### 3.6.3 Tap FF975R/FG413R (Fig. 7)

- ▶ Loosen spacer sleeve 1 at the screw coupling and remove it in the direction of the arrow.

#### 3.6.4 Depth gauge FF965R (Fig. 8)

- ▶ Loosen knurled nut 2.
- ▶ Pull caliper 1 out in the direction of the arrow.
- ▶ Unscrew knurled nut 2.

### 3.7 Cleaning/Disinfection

#### 3.7.1 Product-specific safety information on the reprocessing method

Damage to or destruction of the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!

- ▶ Following the manufacturer's instructions, use cleaning and disinfecting agents
  - that are approved for (e.g. aluminum, plastics, high-grade steel),
  - that do not attack softeners (e.g., in silicone).
- ▶ Observe specifications regarding concentration, temperature and exposure time.
- ▶ Do not exceed the maximum allowable disinfection temperature of 95 °C.

For products with plasma layers (e.g. Noir instruments), the layer is attacked or removed when special cleaning procedures with the addition of oxidizing chemicals (e.g. hydrogen peroxide H<sub>2</sub>O<sub>2</sub>) are used.

- ▶ Do not use oxidizing chemicals for cleaning.
- ▶ Use suitable cleaning/disinfecting agents if the product is disposed of in a wet condition. To prevent foaming and degradation of the efficacy of the process chemicals: prior to mechanical cleaning and disinfection, rinse the product thoroughly with running water
- ▶ If the microsurgical products can be securely fixed in machines or storage devices in such a way that they will be cleaned thoroughly, clean and disinfect them mechanically.

#### 3.7.2 Validated cleaning and disinfection procedure

For mechanical cleaning:

Validated procedure	Specific requirements	Reference
Manual cleaning with immersion disinfection ■ FF885R ■ FF886R ■ FF975R ■ FG413R ■ FG412R ■ FF965R	<ul style="list-style-type: none"> <li>■ Keep working ends open for cleaning.</li> <li>■ When cleaning instruments with movable hinges, ensure that these are in an open position and, if applicable, move the hinge while cleaning.</li> <li>■ Drying phase: Use a lint-free cloth or medical compressed air</li> </ul>	Chapter Manual cleaning/disinfection and subsection: ■ Chapter Manual cleaning with immersion disinfection
Mechanical alkaline cleaning and thermal disinfection ■ FG412R	<ul style="list-style-type: none"> <li>■ Place the product on a tray that is suitable for cleaning (avoid rinsing blind spots).</li> </ul>	Chapter Mechanical cleaning/disinfection and subsection: ■ Chapter Mechanical alkaline cleaning and thermal disinfecting
Manual pre-cleaning with brush and subsequent mechanical alkaline cleaning and thermal disinfection ■ FF885R ■ FF886R ■ FF975R ■ FG413R ■ FF965R	<ul style="list-style-type: none"> <li>■ Place the product on a tray that is suitable for cleaning (avoid rinsing blind spots).</li> </ul>	Chapter Mechanical cleaning/disinfection with manual pre-cleaning and subsection: ■ Chapter Manual pre-cleaning with a brush ■ Chapter Mechanical alkaline cleaning and thermal disinfecting

### 3.8 Manual cleaning/disinfection

- ▶ Prior to manual disinfecting, allow water to drip off for a sufficient length of time to prevent dilution of the disinfecting solution.
- ▶ After manual cleaning/disinfection, check visible surfaces visually for residues.
- ▶ Repeat the cleaning/disinfection process if necessary.

#### 3.8.1 Manual cleaning with immersion disinfection

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Disinfecting cleaning	RT (cold)	>15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
II	Intermediate rinse	RT (cold)	1	-	D-W	-
III	Disinfection	RT (cold)	5	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
IV	Final rinse	RT (cold)	1	-	FD-W	-
V	Drying	RT	-	-	-	-

D-W: Drinking water  
 FD-W: Fully desalinated water (demineralized, microbiological, at least of drinking water quality)  
 RT: Room temperature  
 \*Recommended: BBraun Stabimed fresh

- ▶ Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

#### Phase I

- ▶ Fully immerse the product in the cleaning/disinfectant for at least 15 min. Ensure that all accessible surfaces are moistened.
- ▶ Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- ▶ If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
- ▶ Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- ▶ Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

#### Phase II

- ▶ Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- ▶ Drain any remaining water fully.

#### Phase III

- ▶ Fully immerse the product in the disinfectant solution.
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- ▶ Rinse lumens at least 5 times at the beginning of the exposure time using an appropriate disposable syringe. Ensure that all accessible surfaces are moistened.

#### Phase IV

- ▶ Rinse/flush the product thoroughly (all accessible surfaces).
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during final rinse.
- ▶ Rinse lumens with an appropriate disposable syringe at least five times.
- ▶ Drain any remaining water fully.

#### Phase V

- ▶ Dry the product in the drying phase with suitable equipment (e.g. cloth, compressed air), see Validated cleaning and disinfection procedure.

### 3.9 Mechanical cleaning/disinfection

##### Note

The cleaning and disinfection device must be of tested and approved effectiveness (e.g. FDA approval or CE mark according to DIN EN ISO 15883).

##### Note

The cleaning and disinfection device used for processing must be serviced and checked at regular intervals.

#### 3.9.1 Mechanical alkaline cleaning and thermal disinfecting

Machine type: single-chamber cleaning/disinfecting machine without ultrasound

Phase	Step	D [°C/°F]	t [min]	Water quality	Chemical/Note
I	Prerinse	<25/77	3	D-W	-
II	Cleaning	55/131	10	FD-W	<ul style="list-style-type: none"> <li>■ Concentrate, alkaline:                             <ul style="list-style-type: none"> <li>– pH ~ 13</li> <li>– &lt;5 % anionic surfactant</li> </ul> </li> <li>■ working solution 0.5%                             <ul style="list-style-type: none"> <li>– pH = 11*</li> </ul> </li> </ul>
III	Intermediate rinse	>10/50	1	FD-W	-
IV	Thermal disinfecting	90/194	5	FD-W	-
V	Drying	-	-	-	According to the program for cleaning and disinfection device

DW: Drinking water  
 FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

\*Recommended: BBraun Helimatic Cleaner alkaline

- ▶ Check visible surfaces for residues after mechanical cleaning/disinfecting.

### 3.10 Mechanical cleaning/disinfection with manual pre-cleaning

##### Note

The cleaning and disinfection device must be of tested and approved effectiveness (e.g. FDA approval or CE mark according to DIN EN ISO 15883).

##### Note

The cleaning and disinfection device used for processing must be serviced and checked at regular intervals.

#### 3.10.1 Manual pre-cleaning with a brush

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Disinfecting cleaning	RT (cold)	>15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
II	Rinsing	RT (cold)	1	-	D-W	-

D-W: Drinking water  
 RT: Room temperature  
 \*Recommended: BBraun Stabimed fresh

- ▶ Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

#### Phase I

- ▶ Fully immerse the product in the cleaning/disinfectant for at least 15 min. Ensure that all accessible surfaces are moistened.
- ▶ Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- ▶ If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
- ▶ Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- ▶ Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

#### Phase II

- ▶ Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.

### 3.10.2 Mechanical alkaline cleaning and thermal disinfecting

Machine type: single-chamber cleaning/disinfection device without ultrasound

Phase	Step	T [°C/°F]	t [min]	Water quality	Chemical
I	Prerinse	<25/77	3	D-W	-
II	Cleaning	55/131	10	FD-W	<ul style="list-style-type: none"><li>■ Concentrate, alkaline:<ul style="list-style-type: none"><li>- pH = 13</li><li>- &lt;5 % anionic surfactant</li></ul></li><li>■ working solution 0.5%<ul style="list-style-type: none"><li>- pH = 11*</li></ul></li></ul>
III	Intermediate rinse	>10/50	1	FD-W	-
IV	Thermal disinfecting	90/194	5	FD-W	-
V	Drying	-	-	-	According to the program for cleaning and disinfection device

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

\*Recommended: BBraun Helimatic Cleaner alkaline

- ▶ Check visible surfaces for residues after mechanical cleaning/disinfecting.

### 3.11 Inspection

- ▶ Allow the product to cool down to room temperature.
- ▶ Dry the product if it is wet or damp.

#### 3.11.1 Visual inspection

- ▶ Ensure that all soiling has been removed. In particular, pay attention to mating surfaces, hinges, shafts, recessed areas, drill grooves and the sides of the teeth on rasps.
- ▶ If the product is dirty: repeat the cleaning and disinfection process.
- ▶ Check the product for damage, e.g. insulation or corroded, loose, bent, broken, cracked, worn or severely scratched and fractured components.
- ▶ Check the product for missing or faded labels.
- ▶ Check the products with long, slim shapes (in particular rotating instruments) for deformities.
- ▶ Check the product for damage to the spiral element.
- ▶ Check the cutting edges for continuity, sharpness, nicks and other damage.
- ▶ Check the surfaces for rough spots.
- ▶ Check the product for burrs that could damage tissue or surgical gloves.
- ▶ Check the product for loose or missing parts.
- ▶ Immediately put aside damaged or inoperative products and send them to Aesculap Technical Service, see Technical service.

#### 3.11.2 Functional test

##### ⚠ CAUTION

**Damage (metal cold welding/friction corrosion) to the product caused by insufficient lubrication!**

- ▶ Prior to function checks, lubricate moving parts (e.g. joints, pusher components and threaded rods) with maintenance oil suitable for the respective sterilization process (e.g. for steam sterilization: STERILIT® I oil spray JG600 or STERILIT® I drip lubricator JG598).
- ▶ Assemble disassembled products, see Disassembly.
- ▶ Check that the product functions correctly.
- ▶ Check that all moving parts are working properly (e.g. hinges, locks/latches, sliding parts etc.).
- ▶ Check rotating products (e.g. reusable drills and cutters) for bends and deformities. To do this, roll the product on an even surface.
- ▶ Check for compatibility with associated products.
- ▶ Immediately put aside inoperative products and send them to Aesculap Technical Service, see Technical service.

### 3.12 Packaging

- ▶ Appropriately protect products with fine working tips.
- ▶ Place the product in its holder or on a suitable tray. Ensure that sharp edges are covered.
- ▶ Package trays appropriately for the sterilization process (e.g. in Aesculap sterile containers).
- ▶ Ensure that the packaging provides sufficient protection against contamination of the product during storage.

### 3.13 Steam sterilization

##### Note

The product can be sterilized either in disassembled or in assembled condition.

- ▶ Check to ensure that the sterilizing agent will come into contact with all external and internal surfaces (e.g., by opening any valves and faucets).
- ▶ Validated sterilization process
  - Steam sterilization using fractional vacuum process
  - Steam sterilizer according to DIN EN 285 and validated according to DIN EN ISO 17665
  - Sterilization using fractionated vacuum process at 134 °C/holding time 5 min
- ▶ If several devices are sterilized at the same time in the same steam sterilizer: Ensure that the maximum permitted load according to the manufacturers' specifications is not exceeded.

### 3.14 Storage

- ▶ Store sterile products in germ-proof packaging, protected from dust, in a dry, dark, temperature-controlled area.

## 4. Technical service

##### ⚠ CAUTION

Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

- ▶ Do not modify the product.
- ▶ For service and repairs, please contact your national B. Braun/Aesculap agency.

##### Service addresses

Aesculap Technischer Service

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78532 Tuttlingen / Germany

Phone: +49 7461 95-1601

Fax: +49 7461 16-2887

E-Mail: ats@aesculap.de

Other service addresses can be obtained from the address indicated above.

## 5. Disposal

##### ⚠ WARNING

Risk of infection due to contaminated products!

- ▶ Adhere to national regulations when disposing of or recycling the product, its components and its packaging.

##### ⚠ WARNING

Risk of injury due to sharp-edged and/or pointed products!

- ▶ When disposing of or recycling the product, ensure that the packaging prevents injury by the product.

##### Note

The user institution is obliged to reprocess the product before its disposal, see Validated reprocessing procedure.

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