en	Instructions for use/Technical description 3D Cages: CeSPACE® 3D, PROSPACE® 3D, PROSPACE® 3D Oblique, TSPACE® 3D
USA	Note for U.S. users This Instructions for Use is NOT intended for United States users. Please discard. The Instructions for Use for United States users can be obtained by visiting our website at www.aesculapimplantsystems.com and clicking the "Products" menu. If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-866-229-3002. A paper copy will be provided to you upon request at no additional cost.
de	Gebrauchsanweisung/Technische Beschreibung 3D Cages: CeSPACE® 3D, PROSPACE® 3D, PROSPACE® 3D Oblique, TSPACE® 3D
fr	Mode d'emploi/Description technique 3D Cages: CeSPACE® 3D, PROSPACE® 3D, PROSPACE® 3D Oblique, TSPACE® 3D
es	Instrucciones de manejo/Descripción técnica 3D Cages: CeSPACE® 3D, PROSPACE® 3D, PROSPACE® 3D Oblique, TSPACE® 3D
it	Istruzioni per l'uso/Descrizione tecnica 3D Cages: CeSPACE® 3D, PROSPACE® 3D, PROSPACE® 3D Oblique, TSPACE® 3D
pt	Instruções de utilização/Descrição técnica 3D Cages: CeSPACE® 3D, PROSPACE® 3D, PROSPACE® 3D Oblique, TSPACE® 3D
nl	Gebruiksaanwijzing/Technische beschrijving 3D Cages: CeSPACE® 3D, PROSPACE® 3D, PROSPACE® 3D Oblique, TSPACE® 3D
sv	Bruksanvisning/Teknisk beskrivning 3D Cages: CeSPACE® 3D, PROSPACE® 3D, PROSPACE® 3D Oblique, TSPACE® 3D
ru	Инструкция по применению/Техническое описание 3D Cages: CeSPACE® 3D, PROSPACE® 3D, PROSPACE® 3D Oblique, TSPACE® 3D
cs	Návod k použití/Technický popis 3D Cages: CeSPACE® 3D, PROSPACE® 3D, PROSPACE® 3D Oblique, TSPACE® 3D
pl	Instrukcja użytkowania/Opis techniczny 3D Cages: CeSPACE® 3D, PROSPACE® 3D, PROSPACE® 3D Oblique, TSPACE® 3D
sk	Návod na použitie/Technický opis 3D Cages: CeSPACE® 3D, PROSPACE® 3D, PROSPACE® 3D Oblique, TSPACE® 3D
tr	Kullanım Kılavuzu/Teknik açıklama 3D Cages: CeSPACE® 3D, PROSPACE® 3D, PROSPACE® 3D Oblique, TSPACE® 3D

Aesculap®

3D Cages: CeSPACE® 3D, PROSPACE® 3D, PROSPACE® 3D Oblique, TSPACE® 3D

1. Intended use

3D Cages are used as follows:

- CeSPACE® 3D: stabilization of the cervical spine C2-T1 through anterior approach, monosegmental and multisegmental.
- PROSPACE® 3D: stabilization of the lumbar and thoracic spine through posterior approach, monosegmental and multisegmental.
- PROSPACE® 3D Oblique: stabilization of the lumbar and thoracic spine through transforaminal approach, monosegmental and multisegmental.
- TSPACE® 3D: stabilization of the lumbar and thoracic spine through transforaminal approach, monosegmental and multisegmental.

2. Indications

- Degenerative instability
- Spondylolisthesis
- Post-discectomy syndrome
- Post-traumatic instabilities

Surgically installed implants are designed to support the normal healing process. They should neither replace normal structures of the body nor permanently bear the loads occurring in the case of incomplete healing.

3. Absolute contraindications

Do not use in the presence of:

- Severe damage to the bone structures of the spine that could prevent the stable implantation of the implant components; for example, osteopenia, severe osteoporosis, Paget's disease, bone tumors etc.
- Metabolic or degenerative metabolic bone diseases that could compromise the stable anchoring of the implant system
- Suspected allergy or sensitivity to the implant materials
- Acute or chronic vertebral infections of a local or systemic nature
- Cases not listed under indications

4. Relative contraindications

In the following circumstances, use of the implant system could represent an increased clinical risk and therefore requires precise, individual assessment by the surgeon:

- Medical or surgical conditions that could negatively impact the success of the implantation, including wound healing disorders
- Conditions that could subject the spine and implants to excessive pressure; for example, pregnancy, obesity, neuromuscular diseases or disorders
- Generally poor condition of the patient; for example, drug or alcohol addiction
- Poor patient compliance or limited ability to follow medical instructions, particularly in the post-op phase, including with regard to the restrictions on range of movement in terms of physical exercise and occupational activity

5. Side effects and interactions

In addition to surgery-related risks, potential complications in connection with intervertebral procedures can include, but are not limited to:

- Malpositioning, fracture, loosening, migration/dislocation of the implant
- Spondylolisthesis, pseudarthrosis, inadequate integration of the implant
- Loss of intervertebral disk height due to removal of healthy bone material
- Changes in bone density, degenerative changes in the region of the adjacent vertebral bodies
- Foreign body reactions, allergy
- Infection
- Neurological complications caused by overdistraction or trauma of the nerve roots or dura
- Persistent pain

6. Materials

The materials used in the implant are listed on the packaging:

- Titanium alloy Ti6Al4V ELI according to ASTM F3001 and in accordance with ASTM F136

7. Implantation 3D Cages

⚠ WARNING

The success of the implantation may be jeopardized if the implant bed is not adequately prepared beforehand or if the implant is aligned incorrectly!

- ▶ Make certain that the endplates of the neighboring vertebral bodies are not weakened, in order to minimize the risk of migration.
- ▶ Make certain that the implant bed is properly prepared to avoid damage to the implant when it is driven in.
- ▶ When inserting the implant into the intervertebral space, observe the markings to correctly align the implant.
 - CeSPACE® 3D: Marking points in the cranial direction
 - PROSPACE® 3D Oblique: Marking designates medial

⚠ WARNING

Surrounding structures may be injured due to the selection of an incorrect implant size or implant location!

- ▶ Always check the correct size and location every time with X-ray controlling by using trial implants.
- ▶ CeSPACE® 3D:
 - Use CeSPACE® 3D trial implant and insertion instrument with depth stop.
 - When inserting the CeSPACE® 3D trial implant, observe the markings to correctly align the trial implant.
- ▶ PROSPACE® 3D Oblique:
 - When inserting the PROSPACE® 3D Oblique trial implants, observe the markings to correctly align the trial implant.

⚠ WARNING

Damage to the implant due to excessive application of force!

- ▶ Observe the correct alignment of the implant when mounting on the instrument.
- ▶ Mount the implant on the insertion instrument hand-tight as far as it will go.
- ▶ When inserting the implant into the intervertebral space, avoid canting and levering, and take care to maintain an alignment parallel to the endplates.
- ▶ Do not use force during filling, mounting and implantation.

7.1 CeSPACE® 3D

Note

A cervical plate may be required for additional stabilization.

7.2 PROSPACE® 3D

- ▶ Always implant two implants per layer (PLIF technique).
- ▶ Always use PROSPACE® 3D in conjunction with an internal fixator.

7.3 PROSPACE® 3D Oblique

- ▶ Always use PROSPACE® 3D Oblique in conjunction with an internal fixator.

Note

PROSPACE® 3D Oblique can be implanted through an open or minimally invasive transforaminal access (oblique TULF technique).

7.4 TSPACE® 3D

- ▶ Always use TSPACE® 3D in conjunction with an internal fixator.

Note

TSPACE® 3D can be implanted through an open or minimally invasive transforaminal access.

8. General information

The operating surgeon shall devise an operation plan that specifies and accurately documents the following:

- Selection of the implant components and their dimensions
- Positioning of the implant components in the bone
- Location of intraoperative landmarks

The following conditions must be fulfilled prior to application:

- All requisite implant components are ready to hand.
- All requisite implantation instruments must be available and in working order, including specialized Aesculap implantation systems.
- The operating surgeon and operating room team are thoroughly conversant with the operating technique and with the available range of implants and instruments; information materials on these subjects must be complete and ready to hand.
- The manufacturer has been consulted if the preoperative situation was unclear and if implants were found in the area operated on.
- Aseptic operating conditions
- The operating surgeon is fully conversant with the rules governing medical practice, the current state of scientific knowledge, and the contents of relevant scientific articles by medical authors.

The patient has been informed of the procedure, taking into account the information provided in the sections on indications, contraindications, side effects and interactions, and has documented his/her understanding and consent regarding the following points:

- In the case of delayed or incomplete fusion, the implants can break and loosen due to high loads.
- The life-span of the implant depends on the patient's body weight.
- The implant components must not be overloaded by extreme strains, hard physical labor or sports.
- Corrective surgery may be necessitated by implant loosening, fracture or loss of correction.
- Smokers present an increased risk of bone fusion failure.
- The patient must undergo regular medical follow-up examinations of the implant components.

9. Safety notes

⚠ CAUTION

The use of high frequency surgical devices may damage the implant!

- ▶ Avoid using high frequency surgical devices in proximity to the implant.
- ▶ If the implant is damaged: Remove the implant.

- It is the operating surgeon's responsibility to ensure that the surgical procedure is performed properly.
- General risk factors associated with surgical procedures are not described in this documentation.
- The operating surgeon must have a thorough command of both the hands-on and conceptual aspects of the established operating techniques.
- The operating surgeon must be fully conversant with bone anatomy, including the pathways of nerves, blood vessels, muscles, and tendons.
- Aesculap is not responsible for any complications arising from incorrect diagnosis, choice of incorrect implant, incorrectly combined implant components and/or operating techniques, the limitations of treatment methods, or lack of asepsis.
- The user instructions for individual Aesculap implant components must be followed.
- The implant components were tested and approved in combination with Aesculap components. If other combinations are used, the responsibility for such action lies with the operating surgeon.
- Do not, under any circumstances, combine implant components from different manufacturers.
- Do not, under any circumstances, use damaged or surgically excised components.
- Implants that have been used before must not be reused.
- The attending physician shall make any decision with regard to the removal of implant components that have been used,

⚠ WARNING

The implant has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration or image artifact in the MR environment. Scanning a patient who has this implant may result in patient injury.

- The implant components applied, along with their article numbers, the name of the implant, as well as the batch number and serial number (if available) must be documented in all patient records.
- Postoperatively, individual patient information, as well as mobility and muscle training, is of particular importance.

10. Sterility

- The implant components come individually packed in protective packaging that is labeled according to its contents.
- The implant components are gamma-sterilized.
- ▶ Store implant components in their original packaging. Remove them from their original protective packaging only just prior to application.
- ▶ Prior to use, check the product expiry date and verify the integrity of the sterile packaging.
- ▶ Do not use implant components that are past their expiration date or whose packaging is damaged.

⚠ WARNING

Damage to implants caused by processing and resterilization!

- ▶ Do not reprocess or resterilize the implants.

Further information on Aesculap implant systems is always available from B. Braun/Aesculap or the relevant B. Braun/Aesculap office.

11. Disposal

- ▶ Adhere to national regulations when disposing of or recycling the product, its components and its packaging!

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