# Aesculap® Ennovate® Cervical

## Aesculap Spine

**en** Instructions for use/Technical description

Basic instruments

**USA** Note for U.S. users

This Instructions for Use is NOT intended for United States users. Please discard. The Instructions for Use for United States users can be obtained by visiting our website at [www.aesculapImplantsystems.com](http://www.aesculapImplantsystems.com). If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-866-229-3002. A paper copy will be provided to you upon request at no additional cost.

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





41

42

Ennovate Cervical Set Complementary

1 18 43 44 45 46 31 47 48 49

50 51 52 53 54 55 56

57 58 59 60 61 62 63 64 65 66

# Aesculap® Ennovate® Cervical

## Basic instruments

### Legend

#### Ennovate Cervical Set Basic

- 1 Ennovate Cervical ratchet handle straight SZ001R
- 2 Ennovate Cervical drill handle tear drop style SZ002R
- 3 Ennovate Cervical drill handle palm style SZ003R
- 4 Ennovate Cervical torque limiting T-handle SZ005R
- 5 Ennovate Cervical torque limiting handle straight SZ006R
- 6 Ennovate Cervical counter torque handle SZ008R
- 7 Ennovate Cervical drill  $\varnothing$  2.3 mm solid for universal drill guide SZ016SU/SZ016R
- 8 Ennovate Cervical drill  $\varnothing$  2.8 mm solid for universal drill guide SZ017SU/SZ017R
- 9 Ennovate Cervical drill  $\varnothing$  2.3 mm for variable/fixed drill guide SZ019SU/SZ019R
- 10 Ennovate Cervical drill  $\varnothing$  2.8 mm for variable/fixed drill guide SZ020SU/SZ020R
- 11 Ennovate Cervical center punch SZ050R
- 12 Ennovate Cervical fixed drill guide, 12 mm SZ048R
- 13 Ennovate Cervical fixed drill guide, 14 mm SZ049R
- 14 Ennovate Cervical variable drill guide for 3.5 mm / 3.6 mm and 4.0 mm screws SZ054R
- 15 Ennovate Cervical universal drill guide OPEN/NAV (incl. 3.5 mm / 3.6 mm and 4.0 mm tube) SZ055R
- 16 Ennovate Cervical MIS universal drill guide (incl. 3.5 mm / 3.6 mm and 4.0 mm tube) SZ143R
- 17 Ennovate Cervical tube for universal drill guide  $\varnothing$  4.5 mm screws SZ058R
- 18 Ennovate Cervical screwdriver for polyaxial screws SZ063R
- 19 Ennovate Cervical ball end screwdriver for polyaxial screws SZ064R
- 20 Ennovate Cervical screw body manipulator SZ065R
- 21 Ennovate Cervical sounder SZ066R
- 22 Ennovate Cervical set screw starter double-ended SZ068R
- 23 Ennovate Cervical set screw starter single-ended SZ069R
- 24 Ennovate Cervical set screw torque/removal driver shaft SZ070R
- 25 Ennovate Cervical / S4 Cervical rod holding forceps FW076R
- 26 Ennovate Cervical rod bending template 150 mm SZ073SU
- 27 Ennovate Cervical rod cutter SZ075R
- 28 Ennovate Cervical rod cutter CoCr SZ077R
- 29 Ennovate Cervical rod bender SZ076R
- 30 Spinal rod bending forceps FW024R
- 31 Ennovate Cervical persuader in-line SZ080R
- 32 Ennovate Cervical K-wire  $\varnothing$  1.0 mm for 3.6 mm screws SZ010SU/SZ010R
- 33 Ennovate Cervical K-wire  $\varnothing$  1.5 mm for 4.0 mm and 4.5 mm screws SZ011SU/SZ011R
- 34 Ennovate Cervical / S4 Cervical pedicle depth gauge FW042R
- 35 Ennovate Cervical bone probe blunt straight SZ022R
- 36 Ennovate Cervical bone probe blunt curved SZ023R
- 37 Ennovate Cervical bone probe straight thoracic cannulated SZ024R
- 38 Ennovate Cervical obturator for can. bone probe SZ025R
- 39 Ennovate Cervical bone probe thoracic straight SZ026R

- 40 Ennovate Cervical bone probe thoracic curved SZ027R
- 41 Ennovate Cervical / S4 Cervical bone probe straight FW674R
- 42 Ennovate Cervical / S4 Cervical bone probe curved FW675R

#### Ennovate Cervical Set Complementary

- 1 Ennovate Cervical ratchet handle straight SZ001R
- 18 Ennovate Cervical screwdriver for polyaxial screws SZ063R
- 43 Ennovate Cervical drill  $\varnothing$  3.3 mm solid for universal drill guide SZ018SU/SZ018R
- 44 Ennovate Cervical tap  $\varnothing$  3.5 mm /  $\varnothing$  3.6 mm SZ051R
- 45 Ennovate Cervical tap  $\varnothing$  4.0 mm SZ052R
- 46 Ennovate Cervical tap  $\varnothing$  4.5 mm SZ053R
- 31 Ennovate Cervical persuader in-line SZ080R
- 47 Ennovate Cervical rod persuader Basic SZ081R
- 48 Ennovate Cervical rod persuader pistol style SZ082R
- 49 Ennovate Cervical caliper for cross connectors SZ029R
- 50 Ennovate Cervical template for SQ020TS domino lateral 3.5 mm / 4.0 mm to 3.5 mm / 4.0 mm SZ030R
- 51 Ennovate Cervical template for SQ021TS domino lateral 3.5 mm / 4.0 mm to 5.5 mm SZ031R
- 52 Ennovate Cervical template for SQ022TS domino axial 3.5 mm / 4.0 mm to 3.5 mm / 4.0 mm SZ032R
- 53 Ennovate Cervical template for SQ023TS domino axial 3.5 mm / 4.0 mm to 5.5 mm SZ033R
- 54 Ennovate Cervical template for SQ024TS domino lateral 3.5 mm / 4.0 mm to 5.5 mm SZ034R
- 55 Ennovate Cervical template for SQ025TS domino lateral 3.5 mm / 4.0 mm to 3.5 mm / 4.0 mm SZ035R
- 56 Ennovate Cervical template for SQ010TS offset-connector SZ036R
- 57 Ennovate Cervical handle for connector templates SZ039R
- 58 Ennovate Cervical locking driver for head-to-head connector SZ089R
- 59 Ennovate Cervical rod bending template 60 mm SZ072SU
- 60 Ennovate Cervical rod bending template 290 mm SZ074SU
- 61 Ennovate Cervical in-situ rod bender SZ078R
- 62 Ennovate Cervical / S4 Cervical lamina preparator FW071R
- 63 Ennovate Cervical / S4 Cervical hook holding forceps curved FW528R
- 64 Ennovate Cervical / S4 Cervical compression forceps FW427R
- 65 Ennovate Cervical / S4 Cervical distraction forceps FW428R
- 66 Ennovate Cervical / S4 Cervical rod bending plates FW036R

The relative size of the instruments does not necessarily reflect their size in reality.

The depictions in this document are schematic representations only and may deviate from the actual product.

## Contents

1.	About this document	4	3.22	Insert laminoplasty plate	26
1.1	Scope	4	3.22.1	Prepare laminoplasty plate and hole	26
1.2	Safety messages	4	3.22.2	Fix laminoplasty plate	26
2.	Clinical use	4	3.22.3	Insert rod	27
2.1	Areas of use and limitations of use	4	4.	Validated reprocessing procedure	27
2.1.1	Intended use	4	4.1	General safety instructions	27
2.1.2	Indications	4	4.2	General information	27
2.1.3	Contraindications	4	4.3	Single-use products	28
2.2	Safety information	4	4.4	Reusable products	28
2.2.1	Clinical user	4	4.5	Preparations at the place of use	28
2.2.2	Product	5	4.6	Preparing for cleaning	28
2.2.3	Sterility	5	4.7	Disassembly	28
3.	Application	7	4.7.1	Depth gauge FW042R	28
3.1	Safe application	7	4.7.2	Taps SZ051R to SZ053R	28
3.2	Prepare screw hole with center punch 11 (SZ050R)	7	4.7.3	Variable drill guide SZ054R	29
3.3	Drill holes for the Ennovate Cervical screws	7	4.7.4	Universal drill guide SZ055R	29
3.3.1	Use variable drill guide 14 (SZ054R)	8	4.7.5	Rod benders, all instruments with springs and ratchets	29
3.3.2	Use universal drill guide 15 (SZ055R)	9	4.7.6	Bone probe SZ024R	29
3.4	Use bone probes	10	4.7.7	Rod bender SZ076R	29
3.5	Confirm depth	11	4.7.8	Linear rod persuader SZ080R	30
3.6	Use K-wires as guidance for screw placement (optional)	11	4.7.9	Persuador Pistol Style SZ082R	30
3.7	Cut thread with tap (optional)	11	4.8	Cleaning/Disinfection	30
3.8	Insert screw	12	4.8.1	Product-specific safety information on the reprocessing method	30
3.9	Align screw head	14	4.8.2	Validated cleaning and disinfection procedure	31
3.10	Select rod	14	4.9	Manual cleaning/disinfection	32
3.11	Cut rod	14	4.9.1	Manual cleaning with immersion disinfection	32
3.11.1	Cut rod with rod cutter 27 (SZ075R)	14	4.9.2	Manual cleaning with ultrasound and immersion disinfection	33
3.11.2	Cut rod with rod cutter 28 (SZ077R)	14	4.10	Mechanical cleaning/disinfection	34
3.12	Bend rod	15	4.10.1	Mechanical alkaline cleaning and thermal disinfecting	34
3.12.1	Bend rod with spinal rod bending forceps 30 (FW024R)	15	4.11	Mechanical cleaning/disinfection with manual pre-cleaning	35
3.12.2	Bend rod with Ennovate Cervical rod bender 29 (SZ076R)	15	4.11.1	Manual pre-cleaning with a brush	35
3.12.3	Bend rod with rod bending plates 66 (FW036R) or in situ rod benders 61 (SZ078R)	16	4.11.2	Manual pre-cleaning with ultrasound and brush	36
3.13	Insert rod with rod holding forceps 25 (FW076R)	16	4.11.3	Mechanical alkaline cleaning and thermal disinfecting	36
3.14	Insert set screw	17	4.12	Inspection	37
3.15	Persuade rod	17	4.12.1	Visual inspection	37
3.15.1	Use rod persuader in-line 31 (SZ080R)	17	4.12.2	Functional test	37
3.15.2	Use rod persuader pistol style 48 (SZ082R)	18	4.13	Assembly	37
3.15.3	Use rod persuader Basic 47 (SZ081R)	19	4.13.1	Depth gauge FW042R	37
3.16	Tighten set screw	19	4.13.2	Taps SZ051R to SZ053R	37
3.17	Remove Ennovate Cervical implants	21	4.13.3	Variable drill guide SZ054R	38
3.18	Place laminar hook (optional technique)	21	4.13.4	Universal drill guide SZ055R	38
3.18.1	Place laminar hook	21	4.13.5	Rod benders, all instruments with springs and ratchets	38
3.18.2	Insert set screw	21	4.13.6	Bone probe SZ024R	38
3.18.3	Tighten set screw	21	4.13.7	Rod bender SZ076R	38
3.19	Perform compression and distraction (optional technique)	21	4.13.8	Linear rod persuader SZ080R	38
3.20	Use connectors	22	4.13.9	Persuader Pistol Style SZ082R	39
3.20.1	Cross connectors: Rod-to-rod connectors	23	4.14	Packaging	39
3.20.2	Cross connectors: Head-to-head connectors	23	4.15	Steam sterilization	39
3.20.3	Domino connectors	24	4.16	Storage	39
3.20.4	Cable connectors	25	5.	Technical service	39
3.20.5	Lateral offset connectors (straight and L-shaped)	25	6.	Disposal	39
3.21	Use dual diameter rod SQ009TS	26			

## 1. About this document

### Note

General risk factors associated with surgical procedures are not described in these instructions for use.

### 1.1 Scope

These Instruction for Use apply to the Ennovate Cervical Basic technique.

### Note

The applicable CE mark for the product can be seen on the label or packaging of the product.

- ▶ For article-specific instructions for use as well as information on material compatibility and lifetime see B. Braun eFU at [eifu.bbraun.com](http://eifu.bbraun.com)

### 1.2 Safety messages

Safety messages make clear the dangers to patient, user and/or product that could arise during the use of the product. Safety messages are labeled as follows:

#### DANGER

Indicates a possible threat of danger. If not avoided, death or serious injury may result.

#### WARNING

Indicates a possible threat of danger. If not avoided, minor or moderate injury may result.

#### CAUTION

Indicates a possible threat of material damage. If not avoided, the product may be damaged.

## 2. Clinical use

### 2.1 Areas of use and limitations of use

#### 2.1.1 Intended use

The Ennovate Cervical instruments are used for the implantation of Ennovate Cervical system implants. The indications and contraindications are described in the Instruction for Use of the Ennovate Cervical implants TA015777.

#### 2.1.2 Indications

##### Note

The manufacturer is not responsible for any use of the product against the specified indications and/or the described applications.

For indications, see Chapter 2.1.1.

#### 2.1.3 Contraindications

For contraindications, see Chapter 2.1.1.

## 2.2 Safety information

### 2.2.1 Clinical user

#### General safety information

To prevent damage caused by improper setup or operation, and to not compromise the manufacturer warranty and liability:

- ▶ Use the product only according to these instructions for use.
- ▶ Follow the safety and maintenance instructions.
- ▶ Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge and experience.
- ▶ Store any new or unused products in a dry, clean, and safe place.
- ▶ Prior to use, check that the product is in good working order.
- ▶ Keep the instructions for use accessible for the user.

##### Note

The user is obligated to report all severe events in connection with the product to the manufacturer and the responsible authorities of the state in which the user is located.

#### Notes on surgical procedures

It is the user's responsibility to ensure that the surgical procedure is performed correctly.

Appropriate clinical training as well as a theoretical and practical proficiency of all the required operating techniques, including the use of this product, are prerequisites for the successful use of this product.

The user is required to obtain information from the manufacturer if there is an unclear preoperative situation regarding the use of the product.

## 2.2.2 Product

### Product-specific safety information

These Instruction for Use apply to the Ennovate Cervical Basic technique. For the operating manual of the Basic technique and for all other surgical techniques, please refer to following documents:

Art. no.	Designation
TA015777	IFU ENNOVATE CERVICAL IMPLANT SYSTEM FOR POSTERIOR, OCCIPITAL, CERVICAL AND THORACIC STABILIZATION
TA015762	IFU ENNOVATE CERVICAL NAVIGATED INSTRUMENTS
TA015760	IFU ENNOVATE CERVICAL INSTRUMENTS FOR OCCIPITAL FIXATION
TA015761	IFU ENNOVATE CERVICAL C1/C2 INSTRUMENTS
TA015766	IFU ENNOVATE CERVICAL MIS INSTRUMENTS
TA015755	IFU ENNOVATE CERVICAL TRAYS
TA013082	IFU SECURESPAN
TA013131	IFU SECURESPAN INSTRUMENTS
003102	ENNOVATE CERVICAL - SURGICAL TECHNIQUE
003202	ENNOVATE CERVICAL NAVIGATION - SURGICAL TECHNIQUE
003302	ENNOVATE CERVICAL MIS - SURGICAL TECHNIQUE
047602	SECURESPAN - SURGICAL TECHNIQUE

- ▶ Only combine Aesculap products with each other.
- ▶ Only use Aesculap Ennovate instruments and accessories.
- ▶ Always adhere to applicable standards.
- ▶ Observe the graphic/packing template when storing the instruments.
- ▶ Use the product only in applications according to its intended use, see Chapter 2.1.1.

Instruments can fall in situ or on the floor!

Handles can loosen during use!

- ▶ Check the seating of the coupling handles on the instruments for tightness.

Injury of spinal cord, nerve roots and blood vessels due to piercing of the cortex in case of incorrect application!

- ▶ Only carry out alignment and insertion of the instruments and polyaxial screws under radiographic control or with the aid of a navigation system.

## 2.2.3 Sterility

### Non-sterile, reusable products

Art. no.	Designation
<b>Ennovate Cervical Set Basic</b>	
SZ001R	Ennovate Cervical ratchet handle straight
SZ002R	Ennovate Cervical drill handle tear drop style
SZ003R	Ennovate Cervical drill handle palm style
SZ005R	Ennovate Cervical torque limiting T-Handle
SZ006R	Ennovate Cervical torque limiting handle straight
SZ008R	Ennovate Cervical counter torque handle
SZ050R	Ennovate Cervical center punch
SZ048R	Ennovate Cervical fixed drill guide, 12 mm
SZ049R	Ennovate Cervical fixed drill guide, 14 mm
SZ054R	Ennovate Cervical variable drill guide for 3.5 mm / 3.6 mm and 4.0 mm screws
SZ055R	Ennovate Cervical universal drill guide OPEN/NAV (incl. 3.5 mm and 4.0 mm tube)
SZ143R	Ennovate Cervical MIS universal drill guide (incl. 3.5 mm / 3.6 mm and 4.0 mm tube)
SZ058R	Ennovate Cervical tube for universal drill guide Ø 4.5 mm screws
SZ063R	Ennovate Cervical screwdriver for polyaxial screws
SZ064R	Ennovate Cervical ball end screwdriver for polyaxial screws
SZ065R	Ennovate Cervical screw body manipulator
SZ066R	Ennovate Cervical sounder
SZ068R	Ennovate Cervical set screw starter double-ended
SZ069R	Ennovate Cervical set screw starter single-ended
SZ070R	Ennovate Cervical set screw torque/removal driver shaft
FW076R	Ennovate Cervical / S4 Cervical rod holding forceps
SZ075R	Ennovate Cervical rod cutter
SZ077R	Ennovate Cervical rod cutter CoCr
SZ076R	Ennovate Cervical rod bender
FW024R	Spinal rod bending forceps
SZ080R	Ennovate Cervical persuader in-line
FW042R	Ennovate Cervical / S4 Cervical pedicle depth gauge
SZ022R	Ennovate Cervical bone probe blunt straight
SZ023R	Ennovate Cervical bone probe blunt curved
SZ024R	Ennovate Cervical bone probe straight thoracic cannulated
SZ025R	Ennovate Cervical obturator for can. bone probe
SZ026R	Ennovate Cervical bone probe thoracic straight
SZ027R	Ennovate Cervical bone probe thoracic curved
FW674R	Ennovate Cervical / S4 Cervical bone probe straight
FW675R	Ennovate Cervical / S4 Cervical bone probe curved

Art. no.	Designation
<b>Ennovate Cervical Set Complementary</b>	
SZ001R	Ennovate Cervical ratchet handle straight
SZ063R	Ennovate Cervical screwdriver for polyaxial screws
SZ080R	Ennovate Cervical persuader in-line
SZ081R	Ennovate Cervical rod persuader Basic
SZ082R	Ennovate Cervical rod persuader pistol style
SZ051R	Ennovate Cervical tap Ø 3.5 mm / Ø 3.6 mm
SZ052R	Ennovate Cervical tap Ø 4.0 mm
SZ029R	Ennovate Cervical caliper for cross connectors
SZ030R	Ennovate Cervical template for SQ020TS domino lateral 3.5 mm / 4.0 mm to 3.5 mm / 4.0 mm
SZ031R	Ennovate Cervical template for SQ021TS domino lateral 3.5 mm / 4.0 mm to 5.5 mm
SZ032R	Ennovate Cervical template for SQ022TS domino axial 3.5 mm / 4.0 mm to 3.5 mm / 4.0 mm
SZ033R	Ennovate Cervical template for SQ023TS domino axial 3.5 mm / 4.0 mm to 5.5 mm
SZ034R	Ennovate Cervical template for SQ024TS domino lateral 3.5 mm / 4.0 mm to 5.5 mm
SZ035R	Ennovate Cervical template for SQ025TS domino lateral 3.5 mm / 4.0 mm to 3.5 mm / 4.0 mm
SZ036R	Ennovate Cervical template for SQ010TS offset-connector
SZ039R	Ennovate Cervical handle for connector templates
SZ089R	Ennovate Cervical locking driver for head-to-head connector
SZ053R	Ennovate Cervical tap Ø 4.5 mm
SZ078R	Ennovate Cervical in-situ rod bender
FW071R	Ennovate Cervical / S4 Cervical lamina preparator
FW528R	Ennovate Cervical / S4 Cervical hook holding forceps curved
FW427R	Ennovate Cervical / S4 Cervical compression forceps
FW428R	Ennovate Cervical / S4 Cervical distraction forceps
FW036R	Ennovate Cervical / S4 Cervical rod bending plates

- ▶ Inspect the new product after removing its transport packaging and prior to first use to ensure it is in good working order.

The product is delivered in an unsterile condition.

- ▶ Clean the new product after removing its transport packaging and prior to its initial sterilization.

### Non-sterile single-use products

Art. no.	Designation
SZ016R	Ennovate Cervical drill Ø 2.3 mm solid for universal drill guide
SZ017R	Ennovate Cervical drill Ø 2.8 mm solid for universal drill guide
SZ018R	Ennovate Cervical drill Ø 3.3 mm solid for universal drill guide
SZ019R	Ennovate Cervical drill Ø 2.3 mm for variable/fixed drill guide
SZ020R	Ennovate Cervical drill Ø 2.8 mm for variable/fixed drill guide
SZ010R	Ennovate Cervical K-wire Ø 1.0 mm for 3.6 mm screws
SZ011R	Ennovate Cervical K-wire Ø 1.5 mm for 4.0 mm and 4.5 mm screws

The product is delivered in a non-sterile condition.

- ▶ Clean the new product after removing its transport packaging and prior to its initial sterilization.
- ▶ Do not reuse the product.

### Sterile, single-use products

Art. no.	Designation
SZ016SU	Ennovate Cervical drill Ø 2.3 mm solid for universal drill guide
SZ017SU	Ennovate Cervical drill Ø 2.8 mm solid for universal drill guide
SZ018SU	Ennovate Cervical drill Ø 3.3 mm solid for universal drill guide
SZ019SU	Ennovate Cervical drill Ø 2.3 mm for variable/fixed drill guide
SZ020SU	Ennovate Cervical drill Ø 2.8 mm for variable/fixed drill guide
SZ072SU	Ennovate Cervical rod bending template 60 mm
SZ073SU	Ennovate Cervical rod bending template 150 mm
SZ074SU	Ennovate Cervical rod bending template 290 mm
SZ010SU	Ennovate Cervical K-wire Ø 1.0 mm for 3.6 mm screws
SZ011SU	Ennovate Cervical K-wire Ø 1.5 mm for 4.0 mm and 4.5 mm screws

The product has been sterilized by irradiation and is supplied in sterile packaging.

- ▶ Do not use products from open or damaged sterile packaging.
- ▶ Do not use the product after its use-by date.
- ▶ Do not reuse the product.

The reprocessing of the product affects its functionality. Risk of injury, illness or death due to soiling and/or impaired functionality of the product.

- ▶ Do not reprocess the product.

### 3. Application

#### Note

All handles and handle surfaces are grey. All actuating elements that are actuated during the operation are gold colored. All torque wrench handles are purple.

#### 3.1 Safe application

When using K-wires, observe the following two warning notices:

##### ⚠ DANGER

Injury of the spinal cord, nerve roots and blood vessels due to piercing of the cortex in case of incorrect application!

- ▶ Only carry out alignment and insertion of the instruments, K-wires and polyaxial screws under radiographic control or with the aid of a navigation system.
- ▶ Carefully guide the instruments over the K-wires. When doing this, make sure that the K-wire remains in its position. It is recommended to hold the K-wire and to check the striped marking on the K-wires.
- ▶ Do not bend K-wires. Immediately replace bent K-wires.
- ▶ Inspect the cannulation of the instruments for blockage before every application.

##### ⚠ DANGER

Injury of spinal cord, nerve roots and blood vessels due to bent K-wire or K-wire not holding in the bone!

- ▶ Do not reuse K-wires.
- ▶ Carefully drive in K-wires.
- ▶ Do not bend K-wires. Immediately replace bent K-wire.
- ▶ Insert K-wire in the correct direction, with the roughened end forward.
- ▶ Only use K-wires SZ010SU/SZ010R, SZ011SU/SZ011R of the Ennovate Cervical system.

##### ⚠ WARNING

Risk of injury and/or malfunction!

- ▶ Prior to each use, inspect the product for loose, bent, broken, cracked, worn, or missing components.
- ▶ Always carry out a function test prior to each use of the product.
- ▶ Do not use the product if it is damaged or defective. Set aside the product if it is damaged.

##### ⚠ WARNING

Risk of injury to patients due to damaged instruments!

Ennovate Cervical instruments are precise and can be easily damaged!

- ▶ Handle Ennovate Cervical instruments with the greatest of care.
- ▶ Do not repair or modify the instruments (e.g. sharpening of the tips). In case of malfunction send them to Aesculap Technical Service .

#### 3.2 Prepare screw hole with center punch 11 (SZ050R)

##### ⚠ DANGER

Risk of damage to the spinal cord, nerve roots, adjacent intervertebral space or soft tissue in case of a too deep penetration of the center punch!

- ▶ Never insert the center punch beyond the raised edge.
- ▶ Do not use the center punch in case of poor bone quality.
- ▶ Determine the entry point and alignment of the polyaxial screw.


Fig. 1

- ▶ Use the center punch 11 (SZ050R) to open the cortex. A raised edge **a** is provided on the center punch, to indicate when the ideal depth has been reached, see Fig. 1.
- ▶ Alternatively, make a 1 mm to 2 mm drill hole using a small decortication burr.

#### 3.3 Drill holes for the Ennovate Cervical screws

The drill can be driven manually with a drill handle palm style (SZ003R) / drill handle tear drop style (SZ002R) or with a motor system and an AO small drill attachment (e.g. GB623R).

The drill must be applied with a drill guide. The following guides can be used:

- Variable drill guide 14 (SZ054R) for  $\varnothing$  3.5 mm and  $\varnothing$  4.0 mm screws
- Universal drill guide 15 (SZ055R) for navigated procedures and  $\varnothing$  4.5 mm screws

##### ⚠ DANGER

Risk of damage to the spinal cord, nerve roots, adjacent intervertebral space or soft tissue when drilling the holes!

- ▶ Always use the correct Ennovate Cervical drill guides when drilling the holes.
- ▶ Do not combine drills and guides incorrectly.

##### ⚠ DANGER

Injury to spinal cord and nerve roots caused by application of a drill that is too long!

- ▶ Use an appropriate imaging technique to select an appropriate drill length prior to the operation.
- ▶ Align and insert the drills only under radiographic monitoring and/or aided by a navigation system.
- ▶ Select the drill depth on the guide according to the intended hole depth.

**⚠ DANGER**

Risk of injury and damage to the drill if the rotational speed is too high!

- ▶ Make sure to set a rotational speed appropriate for the intended application.

**⚠ WARNING**

Risk of revision due to instable fixation or damage to the bone!

Damage to the implant or instrument!

Prepared hole is too shallow!

- ▶ Do not sharpen the drill, as this makes the reading on the depth gauge inaccurate and falsifies the reading.
- ▶ Replace blunt drills with new ones.

**Note**

When using the MIS universal drill guide SZ143R, refer to the instructions for use for Ennovate Cervical MIS instruments TA015766.

**3.3.1 Use variable drill guide 14 (SZ054R)**

For better identification the variable drill guide and corresponding drills are marked with a blue ring.

The positive stop of the drills are color-coded to match the color of the screw heads:

- The  $\varnothing$  2.3 mm drill 9 (SZ019SU/SZ019R) is required for the  $\varnothing$  3.5 mm /  $\varnothing$  3.6 mm polyaxial screws (silver).
- The  $\varnothing$  2.8 mm drill 10 (SZ020SU/SZ020R) is required for the  $\varnothing$  4.0 mm polyaxial screws (gold).

**⚠ DANGER**

Risk of damage to the spinal cord, nerve roots, adjacent intervertebral space or soft tissue when drilling a hole of an inaccurate size!

- ▶ Make sure to use the correct drill. Check the color coding and laser marking.

**⚠ DANGER**

Risk of damage to the spinal cord, nerve roots, adjacent intervertebral space or soft tissue when drilling a hole that is too deep or too shallow!

- ▶ Check the laser marking on the guide.
- ▶ Always check the drill lengths with the guide sleeve that is being used with a caliper prior to drilling.
- ▶ Select the appropriate drill.


Fig. 2

- ▶ Mount the drill on the chosen handle. To insert the drill 9/10 (SZ019SU/SZ019R, SZ020SU/SZ020R) in the handle 2/3 (SZ002R/SZ003R), see Fig. 2:
  - Pull back the locking sleeve a.
  - Insert the drill 9/10 up to the stop in the tool holder b of the attachment.
  - Release the sleeve a.
 The drill is attached.
- ▶ Pull at the drill 9/10 to check that it is securely coupled.


Fig. 3

- ▶ Set the desired drilling depth on the variable drill guide 14 (SZ054R) by rotating the inner sleeve a counter-clockwise, see Fig. 3.

**Note**

The variable drill guide offers a range of lengths up to 40 mm. Fixed 12 mm and 14 mm drill guides are available on special request: SZ048R (12 mm) and SZ049R (14 mm).


Fig. 4

- ▶ Check the depth: Insert the required drill into the guide up to the positive stop and verify the exposed length of the drill with a caliper (e.g. AA845R), see Fig. 4.


Fig. 5

- ▶ Position the guide at the desired entry site and advance the drill until the stop is reached, see Fig. 5.
- ▶ To uncouple the drill 9/10 from the handle 2/3:
  - Pull back the locking sleeve a, see Fig. 2.
  - Remove the drill 9/10.

**3.3.2 Use universal drill guide 15 (SZ055R)**

For a quick identification, the tubes are marked with a color-coded ring according to the screw head color:

- Silver ring for Ø 3.5 mm / Ø 3.6 mm tube
- Golden ring for Ø 4.0 mm tube
- Green ring for Ø 4.5 mm tube 17 (SZ058R)

**⚠ WARNING**

**Delay of surgery due to inaccurate hole position!**

- ▶ Make sure to select the right tube. Check the color coding.
- ▶ Select the appropriate drill guide tube.


Fig. 6

- ▶ Attach the tube a to drill guide body b, see Fig. 6. The positive stop of the drills is color-coded to match the color of the screw heads:
  - Ø 2.3 mm drill 7 (SZ016SU/SZ016R) for Ø 3.5 mm / Ø 3.6 mm polyaxial screws (silver)
  - Ø 2.8 mm drill 8 (SZ017SU/SZ017R) for Ø 4.0 mm polyaxial screws (gold)
  - Ø 3.3 mm drill 43 (SZ018SU/SZ018R) for Ø 4.5 mm polyaxial screws (green)

**⚠ DANGER**

**Risk of damage to the spinal cord, nerve roots, adjacent intervertebral space or soft tissue when drilling a hole of an inaccurate size!**

- ▶ Make sure to use the correct drill. Check the color coding and laser marking.
- ▶ Select the appropriate drill.
- ▶ Attach the drill to the handle, see Chapter 3.3.1, see Fig. 2.


Fig. 7

- ▶ Set the desired drilling depth by rotating the inner sleeve c of the drill guide: turn it counter-clockwise to increase the depth and clock-wise to decrease the depth (as indicated by the arrow on the inner sleeve c), see Fig. 7.

**⚠ DANGER**

**Risk of damage to the spinal cord, nerve roots, adjacent intervertebral space or soft tissue when drilling a hole that is too deep or too shallow!**

- ▶ Check the laser marking on the guide.


Fig. 8

- ▶ Insert the required drill into the guide up to the positive stop and verify the exposed length of the drill with a caliper (e.g. AA845R), see Fig. 8.

**⚠ DANGER**

**Risk of damage to the spinal cord, nerve roots, adjacent intervertebral space or soft tissue when drilling the holes!**

- ▶ Always check the drill lengths with the guide sleeve being used with a caliper before starting to drill.


Fig. 9

- Position the guide **15** at the desired entry site and advance the drill **7/8/43** until the stop is reached, see Fig. 9.

### 3.4 Use bone probes

Alternatively to drilling holes for the Ennovate Cervical screws, the pedicle preparation can be performed with bone probes. The probes are available in different designs: blunt, sharp, straight, curved. All probes are equipped with a scale to indicate the penetration depth in the bone:

- Cervical probes **35, 36, 41, 42** (SZ022R, SZ023R, FW674R, FW675R) with diameter 2.3 mm are dedicated for 3.5 mm / 3.6 mm screws.
- Thoracic probes **39, 40** (SZ026R, SZ027R) with a diameter of 2.8 mm are dedicated for 4.0 mm screws (marked with golden color coding ring). They can be used to prepare holes for 4.5 mm thoracic screws.


Fig. 10

#### Note

To determine the depth in the channel, the probes have ring markings **a** at 10 mm intervals, see Fig. 10.

#### ⚠ DANGER

Risk of damage to the spinal cord, nerve roots, adjacent intervertebral space or soft tissue!

Risk of revision in case of screw loosening or pedicle break when preparing a hole of an inaccurate size!

- Make sure to use the right probe. Check the color coding and laser marking.

#### ⚠ DANGER

Injury of spinal cord, nerve roots and blood vessels due to piercing the cortex in the case of incorrect application!

- Insert the bone probe carefully. Avoid levering and bending during insertion.
- Do not insert the bone probe too deep. Check the depth with the laser-marked scale. Verify suitable depth via radiographic control.

#### ⚠ WARNING

Delay of surgery if the curved bone probe is inserted in the wrong direction.

- Check the groove on the shaft.
- Prepare the channel with the chosen bone probe under radiographic control to the desired depth.


Fig. 11

- Make sure to insert the curved bone probes **36/40/42** (SZ023R/SZ027R/FW675R) in the right direction, use groove **b** for orientation.

#### Insert 4.0 mm and 4.5 mm screws with K-wire (optional)


Fig. 12

Ennovate Cervical offers a thoracic cannulated bone probe **37** (SZ024R) that can be assembled with an obturator **38** (SZ025R), as can be seen on the left side of the image, see Fig. 12.

- Insert the assembly in the same manner as the solid bone probe.
- Once the desired depth is reached, remove the obturator **38** by loosening the golden knob **a** counter-clockwise.
- Insert the K-wire **33** through the cannulation of the bone probe **37**.
- Remove the bone probe **37** while leaving the K-wire **33** in the prepared hole.

### 3.5 Confirm depth

#### Note

This step cannot be performed when using K-wires.


Fig. 13

- Use the sounder 21 (SZ066R) to confirm the depth of the prepared hole and the integrity of the wall, see Fig. 13.


Fig. 14

- Use the pedicle depth gauge 34 (FW042R) to measure the desired depth of the hole by using the retractable sleeve, see Fig. 14.  
The pedicle depth gauge 34 (FW042R) is marked in 2 mm intervals.

#### Note

The displayed depth reflects the actual screw thread length to be used as well as the depth of the hole, e.g. a 24 mm depth gauge reading represents a 24 mm drill depth as well as that a polyaxial screw with the length of 24 mm has to be selected.

### 3.6 Use K-wires as guidance for screw placement (optional)

K-wires can be used as guidance for screw placement.

- Insert the  $\varnothing$  4.0 mm and  $\varnothing$  4.5 mm screws over a 1.5 mm guidewire (SZ011SU/SZ011R, golden color coding).
- Insert the  $\varnothing$  3.6 mm screws over a 1.0 mm K-wire (SZ010SU/SZ010R).
- When tapping and inserting the screws over the K-wire, make sure that the K-wire does not advance unintendedly.
- For information on the safe use of K-wires, see Chapter 3.1.

### 3.7 Cut thread with tap (optional)

Ennovate Cervical screws are equipped with a fully threaded, tapered tip to reduce the need to tap. The use of a tap is recommended in case of hard bone. This step is optional.

The taps are equipped with a self-retracting sleeve which helps to avoid the risk of damaging surrounding tissue during tapping.

The taps, like the drills, are color-coded in accordance to the screw size:

- Silver: tap 44 (SZ051R, tap  $\varnothing$  3.5 mm /  $\varnothing$  3.6 mm)
- Gold: tap 45 (SZ052R, tap  $\varnothing$  4.0 mm)
- Green: tap 46 (SZ053R, tap  $\varnothing$  4.5 mm)


Fig. 15

The ratchet straight handle 1 (SZ001R) can be moved to the following positions by rotating the collar of the handle, as shown on the right side of the image, see Fig. 15:

- Forward (IN): for tapping
- Locked (●): blocks the ratchet
- Reverse (OUT): for removing tap from the bone

#### Note

The Ennovate Cervical taps are cannulated to allow guided tapping over a K-wire.

#### ⚠ DANGER

Risk of tissue injury or stripping bone threads when inserting the tap!

- Always check if the sleeve retracts correctly prior to applying the Ennovate Cervical tap.

**⚠ DANGER**

Risk of damage to the spinal cord, nerve roots, adjacent intervertebral space or soft tissue if the tap or K-wire goes too deep.

- ▶ Check the depth on the laser-marked scale.
- ▶ Hold the K-wire and check the striped marking on the wire to make sure that the K-wire does not advance unintendedly.

**⚠ DANGER**

Revision required in case of screw loosening or bone breakage!

Damage or breakage of instruments!

- ▶ Hold screw tap axially to the prepared hole and to the K-wire when turning in and out. Avoid levering, bending and strong applications of force when doing so.
- ▶ Make sure to use the right tap size. Check the color coding and laser marking.
- ▶ Select the correct tap. Check the color coding and laser marking.
- ▶ To attach the tap 44/45/46 to the ratchet straight handle 1 (SZ001R), push the ratchet straight handle 1 onto the hex connection of the corresponding tap 44/45/46 up to the stop, see Fig. 15.
- ▶ Align the rotation setting 'IN' with the laser marked line (or keep it locked to block the ratchet).
- ▶ Position the tap 44/45/46 on previously prepared entry point and cut thread to the required depth by rotating clockwise while maintaining the appropriate trajectory.


Fig. 16

- ▶ Read the depth from the scale a as the sleeve b retracts during tapping, see Fig. 16.
- ▶ Align the rotation setting 'OUT' with the laser marked line (or keep it locked to block the ratchet).
- ▶ Unscrew the tap counter-clockwise.
- ▶ Verify the integrity of the pre-drilled hole and cortex of the vertebral body with sounder 21 (SZ066R).
- ▶ Loosen the tap 44/45/46 from ratchet straight handle 1 by pulling on the black sleeve c.

### 3.8 Insert screw

**⚠ DANGER**

Incorrect application can result in trauma to the spinal cord and nerve roots!

Severe or fatal hemorrhage as a result of vascular injury!

- ▶ Check the insertion of screws by radiographic monitoring and/or with the aid of a navigation system.
- ▶ Take care not to breach the anterior and lateral cortices when screwing in the screw.
- ▶ Ensure that the vessels and nerves are not eroded, punctured during implantation or damaged as a result of implant breakage or migration after implantation.
- ▶ Take care to avoid pulsatile erosion of the vessels and nerves because of close apposition of the implants.

**⚠ WARNING**

Risk of injury to the patient if the polyaxial screw is damaged during removal from the packaging!

- ▶ When connecting the polyaxial screw with the screwdriver SZ063R never brace it against the packaging.

**⚠ WARNING**

Operation delay!

Damage to the polyaxial screw if the screwdriver is not correctly engaged in the screw!

Torque cannot be transferred to the screw!

- ▶ Connect the screwdriver completely to the screw.
- ▶ Check that the polyaxial screw is firmly attached to the screwdriver.
- ▶ Verify the axial alignment of the polyaxial screw to the screwdriver.
  - Only allow the axial load on the screw during screw insertion.
  - Avoid the induction of lateral forces.
- ▶ Hold the screwdriver axially to the prepared hole when turning in and out. Avoid levering, bending and strong applications of force when doing so.

**⚠ CAUTION**

Delay of surgery due to wrong screw selection!

- ▶ Check the laser marking and color-coding.


Fig. 17

**Note**

The screw length L is measured from the tip of the screw to the sphere of the screw shaft, see Fig. 17.

To ensure maximum polyaxiality with the screw, do not tighten the screw completely down to the bone. Leave a small gap below the head to allow rotation and angulation of the screw head.

- ▶ Select the appropriate screw (diameter, length, and type).


Fig. 18

- ▶ Set the rotation direction setting on the handle to 'IN' (or keep it locked to block the ratchet).
- ▶ Push the straight handle 1 (SZ001R) onto the hex connection a of the screwdriver 18 (SZ063R) up to the stop, see Fig. 18.


Fig. 19

- ▶ Place the tip b of the screwdriver into the head e of the screw until the tip b of the screwdriver sits in the socket of the polyaxial screw and thus axially aligns to the screw shaft, see Fig. 19.
- ▶ Screw in the screwdriver 18 into the thread of the polyaxial screw e up to the stop by rotating the golden knob d clockwise until the screwdriver thread engages the screw head e and tighten the golden knob d. Fastening is complete when the polyaxial alignment of the polyaxial screw is completely blocked.


Fig. 20

- ▶ Make sure that the screw is firmly seated on the screwdriver and the polyaxiality is blocked for insertion, see Fig. 20.
- ▶ Insert the screw while maintaining the appropriate trajectory. Use the outer rotational sleeve f to better hold the screwdriver while screwing.
- ▶ After insertion of the screw, rotate the golden knob d counter-clockwise to release the screwdriver.
- ▶ Pull the black sleeve c of the handle 1 while pulling the screwdriver 18 out of the handle.

**⚠ DANGER**

**Injury to the spinal cord, nerve root and blood vessels due to piercing of the cortex and faulty screw placement due to screw loosening when used improperly!**

- ▶ Ensure the axial alignment of the polyaxial screw to the K-wire when screwing in. When doing this, make sure that the K-wire remains in its position.
- ▶ When inserting the screw over the K-wire make sure that the K-wire does not advance forward. It is recommended to hold the K-wire and to check the striped marking on the wire.
- ▶ Remove the K-wire before the screw is completely screwed in to ensure that the K-wire does not advance forward.
- ▶ Do not remove the K-wire too early, otherwise the polyaxial screw guidance will fail.

**Note**

*The Ennovate Cervical screwdrivers are cannulated to allow guided screw insertion over a K-wire.*

- ▶ When using K-wires, observe the warnings described for using K-wires, see Chapter 3.1.

### 3.9 Align screw head

Once the screw is inserted, the position of the polyaxial head can be optimized for rod insertion. The polyaxial screw body can be rotated by 360° and angled up to 90°. The  $\varnothing$  4.0 mm favored angle screws provide additional angulation in the cephalad and caudal directions for a total of 110° with +/- 65° cone angulation.

Fig. 21

- ▶ To optimize the position of the head, use the screw body manipulator **20** (SZ065R), see Fig. 21.
  - Place the working end **a** of the screw body manipulator in the head of the polyaxial screw **b**.
  - Adjust the head angle as desired.

#### Note

If the screw is inserted too far, polyaxial movement of the screw body will be impeded due to bone contact.

If the screw is inserted too deep, turn the screw counter-clockwise using the ball end screwdriver **19** (SZ064R) or with the screwdriver **18** (SZ063R) attached to the ratchet handle straight **1** (SZ001R) (rotation setting set to 'OUT' or locked), until full polyaxial motion is achieved.

### 3.10 Select rod

The Ennovate Cervical rods are supplied in a straight design (excluding pre-bent rods for the occipital plate and curved MIS rods). Their curvature can be adjusted to the required spine lordosis using the bending instruments.

Exact contouring of the rods is essential to the successful assembly of the system.

- ▶ Determine the correct rod length and curvature that enable the rod to smoothly fit inside the heads of the polyaxial screws (or hooks).
- ▶ Use a rod bending template to aid in rod contouring or trimming to the required length (optional).

Rod templates exist in lengths of 60 mm, 150 mm and 290 mm (SZ072SU, SZ073SU, SZ074SU).

For differentiation the rods are color-coded:

- Silver: 3.5 mm rods
- Gold with blue orientation line: 4.0 mm titanium rods
- Silver with white orientation line: 4.0 mm cobalt chrome rods
- ▶ Select the appropriate rod (according to the diameter, length and material).

### 3.11 Cut rod

#### ⚠ DANGER

**Risk of damage to patient, personnel or equipment in the operating room due to flying fragments from cut rod!**

- ▶ Always hold the piece of rod to be cut at both ends.
- ▶ Use rod cutter **27** (SZ075R) to cut titanium rods.
- ▶ Use rod cutter **28** (SZ077R) to cut titanium and cobalt chrome rods and for a smoother rod cut result.
- ▶ If desired, file sharp edges left after cutting the rod with a suitable instrument, such as the Aesculap durogrip file OL429R.

#### 3.11.1 Cut rod with rod cutter 27 (SZ075R)

- ▶ To avoid projection, use the rod holding forceps **25** (FW076R) to grasp the small piece to be cut.


Fig. 22

- ▶ Cut rod with scissors **a** or center hole **b** of the rod cutter **27**, see Fig. 22.

#### 3.11.2 Cut rod with rod cutter 28 (SZ077R)

#### Note

It is not possible to cut the rods at the location where it is excessively bent, as it won't be possible to insert the rod into the rod cutter. Use the rod cutter SZ077R before bending the rods.


Fig. 23

- ▶ Make sure that the handles **a** of the rod cutter **28** (SZ077R) are in a free, relaxed and open position.
- ▶ Mark cutting location on the rod.
- ▶ Rotate the large knob **b** clockwise until both arrows **f** are aligned.
- ▶ Slide the rod **d** into the center hole **c** of the rod cutter.

- ▶ Repeatedly squeeze the handles until the rod is cut.  
The rod will be cut at the side of the rod cutter marked with "cutting line" e, as shown on the right side of the image, see Fig. 23.  
The cutting line is located approx. 8 mm from the upper surface of the instrument.

### 3.12 Bend rod

To bend cobalt chrome rods use the spinal rod bending forceps **30** (FW024R).

To bend titanium rods use the spinal rod bending forceps **30** or the Ennovate Cervical rod bender **29** (SZ076R).

Alternatively to the rod bender, Ennovate Cervical offers rod bending plates **66** (FW036R) as well as in situ rod benders **61** (SZ078R) which can be used to bend titanium and cobalt chrome rods.

Titanium is highly notch-sensitive and therefore care must be taken during rod contouring to ensure that surface damage to the rod is minimal. This is necessary to avoid potential fatigue failure of the implant.

When using cobalt chrome rods please be aware that due to the material characteristics sharp bends (e.g. not smaller than a radius of 50 mm or 130°) should be avoided. If sharper bends are needed bend the rod in several steps while changing positions along the rod.

#### Note

For long Ennovate Cervical rods, bending should be carried out in several steps in order to avoid excessive or insufficient lordosis.

#### ⚠ DANGER

Damage to and breakage of the Ennovate Cervical rods caused by excessive strain on the material!

Danger of injury from damage (breakage) of the rod if bent incorrectly!

- ▶ Always bend the rods in one direction only. Do not bend back the rods.
- ▶ Always use Ennovate Cervical bending instruments for bending the rods.
- ▶ Avoid small bending radii, notching or scratching of the rods.
- ▶ Always use Ennovate Cervical cutting instruments for cutting the rods.

#### ⚠ WARNING

Delay of surgery due to incorrect bending of the rod!

Rod cannot be inserted into the polyaxial screw head!

Rod cannot be fastened!

- ▶ Avoid severe rod bending in the region of the polyaxial screw.
- ▶ Adjust the curvature of the rod to the anatomical situation.

#### ⚠ WARNING

Delay of surgery due to wrong sequence in the surgery when using rod cutter SZ077R!

- ▶ Make sure to cut the rod before bending the extremities, as bent extremities will not fit into the rod cutter.

#### ⚠ CAUTION

Operation delay due to inability to bend cobalt chrome rod!

- ▶ Do not use rod bender **29** (SZ076R) to bend cobalt chrome rods.
- ▶ Use spinal rod bending forceps **30** (FW024R), or alternatively **61/66** (SZ078R/FW036R), to bend cobalt chrome rods.

#### Note

Rods with diameter of 4 mm have an orientation line that facilitates the alignment of the rod during the bending process, see Fig. 26.

#### 3.12.1 Bend rod with spinal rod bending forceps 30 (FW024R)


Fig. 24

- ▶ Select the radius (small, medium, large) by pulling and turning the knob **a**, see Fig. 24.
- ▶ Place the rod between the rollers **b** and the knob **a**.
- ▶ Repetitively squeeze the handles **c**.

#### 3.12.2 Bend rod with Ennovate Cervical rod bender 29 (SZ076R)


Fig. 25

- ▶ Make sure that the rod bender is in its starting position, with the central roller **a** moved completely back, see Chapter 4.13.7.
- ▶ Insert the rod between the central roller **a** and the two rollers **b**, see Fig. 25.
- ▶ Bend the rod to the desired radius by actuating the handle **c** of the rod bender.
- ▶ Repeatedly squeeze the handle until the desired contouring is achieved.
- ▶ To release the rod from the bender, press the release lever **d** and move the ratchet back to the starting position.


Fig. 26

As an optional bending guide, the rod bending template **26** (SZ073SU) can be inserted in the mounting on the distal roller for continuous visualization, as can be seen on the left side of the image, see Fig. 26.

- ▶ In case short rods need to be bent or the rod end needs to be contoured, insert the rod in one of the grooves **e** in front of the rollers.

### 3.12.3 Bend rod with rod bending plates **66** (FW036R) or in situ rod benders **61** (SZ078R)


Fig. 27

- ▶ To use the rod bending plates **66**, insert the rod in the distal hole of the bending plates. Contour the rod, see Fig. 27.
- ▶ To use the in situ rod benders **61**:
  - Position the in situ rod bender on the rod.
  - Contour the rod.

### 3.13 Insert rod with rod holding forceps **25** (FW076R)

#### ⚠ WARNING

Loss of the rod in situ due to incomplete coupling of the rod to the rod holding forceps or early release of the interlock and incorrect alignment of the rod slots of the polyaxial screw!

- ▶ Insert rod completely in jaw part of the rod holding forceps and clamp securely.
- ▶ Do not release clamping until the rod is correctly placed in the screw heads.
- ▶ Align rod slots/rod grooves of the polyaxial screw before inserting the rod.


Fig. 28

- ▶ Insert the rod using the rod holding forceps **25** (FW076R), see Fig. 28.
- ▶ Close the handles **a** to engage the locking mechanism.
- ▶ Insert the rod in the rod slots of the polyaxial screws.
- ▶ Relieve the lock **b** by squeezing the handles **a** together.

#### ⚠ WARNING

Inadequate fixation due to incorrect positioning of the rod in the screw head!

- ▶ Check that the screw head is at a right angle to the rod.
- ▶ Check that the rod is securely positioned in the groove base.


Fig. 29

- ▶ Check that the screw head is at a right angle to the rod.
- ▶ Check that the rod is securely positioned in the groove base.

### 3.14 Insert set screw

This step can only be performed where the rod is completely seated on the floor of the groove. If it is not, continue with the rod persuasion, see Chapter 3.15,

#### ⚠ WARNING

Revision required!

Inadequate fixation from screwing in at an angle/tilting of set screw (cross-threading)!

- ▶ Do not use the set screw to push the rod down.
- ▶ Make sure that the rod is correctly positioned on the floor of the groove before inserting the set screw.

#### ⚠ CAUTION

Intraoperative revision required!

Damage to the implant and instruments!

Operation delay!

- ▶ Load the set screw correctly and completely with the set screw starter.
- ▶ Always place the set screw starter axially to the head of the polyaxial screw and do not tilt it.


Fig. 30

- ▶ Load the set screw on one of the set screw starters (22 (SZ068R) double-ended or 23 (SZ069R) single-ended), see Fig. 30. Insert the working end b into the socket of the set screw a and push with light pressure.
- ▶ Verify solid seating of the set screw a on the set screw starter 22/23.


Fig. 31

- ▶ Place set screw a with set screw starter 22/23 aligned axially to the head of the polyaxial screw c and screw the set screw a into the head of the polyaxial screw c from the top until the set screw a contacts rod d, see Fig. 31.

### 3.15 Persuade rod

#### ⚠ DANGER

Risk of severe injury if too much force is applied when persuading the rod!

Polyaxial screw can tear out of the spinal column!

Notching can occur on the rod due to too high forces!

Revision required if rod breaks!

- ▶ Ensure that the forces applied when persuading the rod are not excessively high.
- ▶ Ensure that no notching occurs on the rod due to forces being too high.

#### ⚠ CAUTION

Operation delay!

Instrument slips from screw head!

Set screw cannot be screwed into screw head!

Rod cannot be fixated!

- ▶ Insert rod persuader correctly into the interface of the screw head. Check correct positioning.
- ▶ Push the rod down completely up to the stop in the screw head. Make certain that the rod is positioned securely in the screw head.

#### ⚠ CAUTION

Damage to and breakage of the rod persuader!

- ▶ Never use the rod persuader to bend the cervical rod.
- ▶ Never use the rod persuader to act as a counter torque for the Ennovate Cervical set screws.

#### 3.15.1 Use rod persuader in-line 31 (SZ080R)

#### ⚠ WARNING

Delay of surgery due to the persuader not being in the starting position!

- ▶ Make sure that the golden knob is pulled up so that the 3 lines are visible.

#### Note

The text and arrows on the instrument indicate the directions in which the knob can be moved to reach different positions, e.g. the release position is reached by pulling up the knob. The arrows do not point to the lines, see Fig. 32.

The starting position is reached when:

- The knob can freely move between the second and third line.
- The knob is pulled up, so that the 3 lines are visible.


Fig. 32

- ▶ Check that the instrument is in the starting position, see Chapter 4.13.
  - Pull the golden knob **a** up to the third line as shown on the left side of the image, see Fig. 32.
  - Place the persuader over the head of the screw.
- ▶ Attach the persuader to the screw head **c** by releasing the golden knob **d** down to the second line as shown at the center of the image, see Fig. 32.
 

The two arms on the distal end of the persuader **b** collapse and engage on the screw head **c**.
- ▶ Visually inspect the connection. In case of a faulty assembly, the next step cannot be performed.
- ▶ Persuade the rod by turning the golden knob **a** in the clockwise direction until the first line is reached as shown on the right side of the image, see Fig. 32.
 

The rod is fully seated in the screw head **c**.


Fig. 33

- ▶ Insert the set screw through the rod persuader in-line **31** (SZ080R) using the single-ended set screw starter **23** (SZ069R), see Fig. 33.


Fig. 34

- ▶ To disengage the persuader from the screw head, turn the golden knob **a** counter-clockwise as far as possible (until the second line can be seen) and pull it up to the third line, see Fig. 34.
 

Arms **b** spread and disengage from the screw head **c**.

### 3.15.2 Use rod persuader pistol style 48 (SZ082R)

#### ⚠ CAUTION

Delay of surgery if the persuader is not in the start position!

- ▶ Make sure that the handle is pressed up against the forward stop.


Fig. 35

- ▶ Release the bar **a** so that the handle **b** can move to the starting position, pressed up against the forward stop **c**.
 

To release the ratchet, press down the bar at location **e**, see Fig. 35.
- ▶ Squeeze the handle **b** to attach the rod persuader **48** to the screw head **f**.
- ▶ Persuade the rod by further squeezing the handle **b** until the rod is fully seated in the screw head.
- ▶ Insert the set screw **g** through the rod persuader **48** using the single-ended set screw starter **23** (SZ069R).
- ▶ To disengage the rod persuader **48**, release the bar **a** so that the handle can move back to the starting position.

#### Note

When applying the rod persuader pistol style, the set screw can only be inserted using the single-ended set screw starter **23**.

### 3.15.3 Use rod persuader Basic 47 (SZ081R)


Fig. 36

- ▶ Place the rod persuader **47** over the polyaxial screw head, see Fig. 36.
  - Make sure to engage the tip **b** of the persuader under the screw head **a**, as shown on the left side of the image, see Fig. 36.
  - If this is not possible (i.e. when the screw head is angulated or in contact with the bone), use rod persuader in-line **31** (SZ080R) or rod persuader pistol style **48** (SZ082R).
- ▶ Squeeze the handle **c** of the persuader to seat the rod **d** into the screw head **a**.


Fig. 37

- ▶ While the rod persuader is applied, insert the set screw through the working end with the set screw inserter **22** (SZ068R) or **23** (SZ069R), see Fig. 37.

### 3.16 Tighten set screw

#### ⚠ DANGER

Injury to the patient from initiation of high torques on the spinal column!

- ▶ Always use the counter torque handle.
- ▶ Always completely and correctly insert the counter torque and check that it is fully seated on the rod and screw head.
- ▶ Hold counter torque in position and do not twist or lever against the inserted rod or screw head.

#### ⚠ WARNING

Risk of correction loss/insufficient stabilization of the spinal column region to be treated due to under-tightening of the set screw!

- ▶ Set the set screw in place correctly.
- ▶ Make certain that the rods are correctly positioned on the floor of the groove.
- ▶ Fully tighten the set screw, using the torque wrench.

#### ⚠ WARNING

Risk of damage to the implant or instruments due to overtightening of the set screw!

Revision required due to damaged set screw.

- ▶ For final tightening of the set screws, use only 2.8 N-m torque handles **4/5** (SZ005R/SZ006R).
- ▶ Check the re-calibration/maintenance date on the torque wrench handle maintenance label prior to each use.
- ▶ When tightening or loosening the set screws, make certain that the working end of the torque wrench handle shaft is positioned entirely within the socket of the set screws.
- ▶ Apply torque slowly and without sudden movement until the torque wrench handle release is audibly triggered.

#### ⚠ WARNING

Inadequate fixation or damage to implants due to introduction of lateral forces!

- ▶ When tightening or loosening, only subject the set screws to axial pressure.
- ▶ Avoid exerting lateral forces when tightening or loosening the set screws.

#### ⚠ WARNING

Delay of surgery or insufficient stabilization of the spinal column region to be treated

- ▶ Do not loosen and retighten the set screw more than 3 times.

#### ⚠ CAUTION

Damage to the countering instrument due to improper usage!

- ▶ Always slide the countering instrument over the screw head (or over the rod next to the screw, when it does not fit on the screw head when connectors are used) onto the implant as far as it will go. When doing so, make sure that the implant entirely engages into the slot at the working tip. Make sure that the rod protrudes on both sides of the countering instrument.

The torque wrench must be used together with counter torque handle 6 (SZ008R).

**Note**

It is imperative to use the counter torque instrument to prevent applying torque directly to the patient's spine. Using the counter torque also ensures a perpendicular placement of the screwdriver.


Fig. 38

- ▶ Assemble the torque wrench by attaching a torque limiting handle 4/5 (SZ005R/SZ006R) to the set screw torque shaft 24 (SZ070R), see Fig. 38:
  - Pull back the sleeve b of handle 4/5 to insert the shaft 24 with compatible square coupling a up to the stop.
  - Check the instrument for firm seating.

Fig. 39

- ▶ Push counter torque handle 6 from the top on the head of the polyaxial screw c, completely and in correct position. The counter torque handle 6 can be placed in four different positions using the four cavities at the bottom of the instrument.


Fig. 40

- ▶ Insert the set screw torque shaft 24 into the cannulation of the counter torque handle 6 and carefully push down with light pressure into the socket of the set screw up to the stop, as shown on the left side of the image, see Fig. 40.
- ▶ Screw in the set screw by turning the 2.8 N-m torque handle 4/5 clockwise. The specified torque only applies for a clockwise rotation (direction of arrow).
- ▶ Hold the counter torque handle 6 securely in position until the acoustic signal of the triggering torque handle sounds. The acoustic signal indicates that 2.8 N-m has been reached.


Fig. 41

- ▶ Visually double check that the set screw is fully seated. When using a  $\varnothing$  3.5 mm rod, the set screw is located below the top of the polyaxial screw tulip after final tightening, as can be seen on the left side of the image, see Fig. 41. When using a  $\varnothing$  4 mm rod, the set screw is located plane with the top of the polyaxial screw tulip after final tightening, as can be seen on the right side of the image, see Fig. 41.
- ▶ Remove set screw torque shaft 24 together with counter torque handle 6.
- ▶ Repeat procedure for additional screws.

### 3.17 Remove Ennovate Cervical implants

#### ⚠ DANGER

Injury to patient from initiation of high torques on the spinal column when removing the set screw!

- ▶ Always use the counter torque handle.
- ▶ Always completely and correctly insert the counter torque handle and check that it is fully seated.

#### ⚠ WARNING

Delay of surgery or risk of insufficient stabilization due to damaged screw!

- ▶ Do not loosen and retighten the set screw more than 3 times.
- ▶ Remove the set screw using the set screw torque / removal driver shaft 24 (SZ070R) with a torque limiting handle 4/5 (SZ005R/SZ006R) in combination with counter torque handle 6 (SZ008R), see Fig. 40.
- ▶ Remove polyaxial screws with the ball end screwdriver 19 (SZ064R) or with the screwdriver 18 (SZ063R) attached to the ratchet handle straight 1 (SZ001R) (rotation setting set to 'OUT' or locked).

Implant removal can result in complications due to stuck implants, ingrown bone tissue, etc. Such complications can lead to damage to implants and/or instruments.

In case of complications when removing implant parts, use special instruments, such as the Cervical revision set, which can be ordered from Aesculap.

### 3.18 Place laminar hook (optional technique)

Ennovate Cervical offers four hook options. The choice of the hook used depends on the thickness of the lamina. The thin throat size is 4.5 mm and the thick throat is 6.0 mm. Hooks are color-coded: blue for left, green for right.

#### 3.18.1 Place laminar hook

#### ⚠ WARNING

Risk of instability due to an oversized hook.

- ▶ Check that the hook is firmly seated on the lamina.

#### ⚠ CAUTION

Delay of surgery due to an incorrectly sized hook.

- ▶ Check the color coding and laser marking.


Fig. 42

- ▶ Before using laminar hooks, prepare the lamina accordingly and dissect the ligamentum flavum.

Use lamina preparator 62 (FW071R) or the hook in the hook holding forceps 63 (FW528R).

- ▶ Select the correct lamina hook and position it on the lamina using the hook holding forceps 63, see Fig. 42.
- ▶ Check that the laminar hook is stable.

#### 3.18.2 Insert set screw

- ▶ Once the rod is in place load the set screw on one of the set screw starters 22/23 (SZ068R/SZ069R).
- ▶ Provisionally tighten the set screw into the lamina hook until it contacts the rod in the same manner as in the polyaxial screw, see Chapter 3.14.

#### 3.18.3 Tighten set screw

#### Note

The counter torque does not fit over the lamina hook.


Fig. 43

- ▶ Attach the counter torque 6 (SZ008R) next to the laminar hook a.
- ▶ Perform final tightening (2.8 N·m of torque), see Chapter 3.16.
- ▶ For risks related to the final tightening of the set screw, see Chapter 3.16.

### 3.19 Perform compression and distraction (optional technique)

If the surgeon sees a need for compression or distraction, straight and curved distraction and compression forceps 64/65 (FW427R/FW428R) are available.

#### ⚠ WARNING

Insufficient treatment or stabilization of the spinal column region to be treated!

Damage to the implant and/or instruments!

- ▶ For final tightening of the set screws, use only 2.8 N·m torque handles 4/5 (SZ005R/SZ006R).
- ▶ Before every corrective maneuver, verify correct seating of the attachments of the forceps on the rod and the polyaxial screws.

#### ⚠ WARNING

Injury to the patient due to initiation of high torques on the spinal column!

- ▶ When tightening the set screw always use the counter torque handle.
- ▶ Avoid excessive distraction or compression.

**⚠ WARNING**

Damage to the implant caused by applying the reduction instruments too high at the implant!

- ▶ Always apply the distraction and compression forceps on the rod at the implant.

**Note**

This technique may be useful especially for placing hooks in the cervical spine.

**Note**

Apply forceps **64/65** (FW427R/FW428R) before tightening the set screws. One set screw should be tightened before starting the process in order to create a fixed point for compression/distraction.


Fig. 44

- ▶ Select the desired starting point of distraction/compression.
- ▶ Make certain that the rod is completely inserted into the head of the polyaxial screw.
- ▶ Insert and finally tighten the first set screw, in order to be able to use it as an anchor point for the distraction or compression, see Chapter 3.14, see Chapter 3.16.
- ▶ Loosen or leave loose adjacent set screw while squeezing the handles **a** to start the distraction/compression process, see Fig. 44.
- ▶ Once the desired compression/distraction is achieved, fully tighten the remaining set screws, see Chapter 3.16.

### 3.20 Use connectors

**⚠ WARNING**

Inadequate fixation due to incorrect rod or connector position!

- ▶ Always position the connectors in such a way that the rod is completely inserted into the connectors.
- ▶ In case MIS rods are used, always position MIS rods in such a way that the teardrop-shaped tip and the conical tip are located outside the clamping area of the connector.
- ▶ Before tightening the set screws, ensure that the connector is positioned correctly.

**⚠ WARNING**

Risk of damage to the implant due to overtightening of the locking screw!

Insufficient fixation due to under-tightening or over-tightening of the set screw!

- ▶ Position the connectors correctly.
- ▶ Make certain that the rods are completely inserted into the connectors.
- ▶ Make certain that there is no soft tissue caught in the connectors.

**⚠ WARNING**

Injury to the patient due to initiation of high torques/forces on the spinal column!

Damage to the instruments!

- ▶ Hold counter torque handle in position and do not twist or lever against the connector.

**Note**

For risks related to the insertion and tightening of the set screw, see Chapter 3.14, see Chapter 3.16.

- ▶ Choose a connector based on the measured length, the individual anatomy of the patient and the rod diameter.
- ▶ Before installation of the domino and cross connector, ensure the following in order to prevent blocking during attachment on the rods:
  - Connector can be freely placed.
  - Integrated set screws of the connector are completely loosened.
- ▶ First correctly place connectors and tighten the integrated set screw.
  - Always use the 2.8 N·m torque handles **4/5** (SZ005R/SZ006R) and the counter torque handle **6** (SZ008R) to tighten and loosen set screws.
  - It is imperative to use the counter torque instrument to prevent applying torque directly to the patient's spine. For the straight lateral offset connector, this also ensures the perpendicular placement of the screwdriver thus simplifying correct tightening of the set screw.
  - Except for the straight lateral offset connector, the counter torque handle does not fit on the connectors. Place counter torque handle **6** next to the connector.

When using cross connectors (head-to-head connectors and rod-to-rod connectors), observe the following warnings.

**⚠ CAUTION**

Delay of surgery due to incorrectly sized connector!

- ▶ Make sure to select the correct connector size using the caliper **49** SZ029R. Check the color coding and laser marking.

**⚠ CAUTION**

Delay of surgery due to bent connector, connector cannot be adjusted!

- ▶ Make sure not to bend the connectors.

### 3.20.1 Cross connectors: Rod-to-rod connectors

Cross connectors are recommended in cases where additional torsional stability may be required.

#### ⚠ CAUTION

Connector is unusable if set screws are removed or not present!

- ▶ Do not unscrew set screws from the connectors under any circumstances.
- ▶ Do not swap set screws with other components of the connectors.
- ▶ Only use connectors if all the set screws are present.


Fig. 45

- ▶ Measure the distance between rods using the caliper 49 (SZ029R) to choose the connector size.
- ▶ Choose cross connector based on the measured length, the individual anatomy of the patient and the rod diameter (3.5 mm or 4.0 mm).

#### Note

For better differentiation, the cross connectors for the  $\varnothing 4.0$  mm rods are color coded in gold and blue whereas the cross connectors for the  $\varnothing 3.5$  mm rod have no color coding.

- ▶ Make sure to use the appropriate sized connector in rod diameter and length by checking the color coding and laser marking.
- ▶ Place the connector with the pre-assembled set screws on the two longitudinal rods at the predetermined position using forceps.


Fig. 46

- ▶ Once correctly placed, tighten the set screws, see Fig. 46, see Chapter 3.16.
- ▶ For risks related to the final tightening of the set screw, see Chapter 3.16.

### 3.20.2 Cross connectors: Head-to-head connectors

As an alternative to rod cross connectors and in case the screw heads are in close proximity to another head-to-head cross connector can be applied.

#### ⚠ WARNING

Risk of loosening the construct when not following the correct sequence!

- ▶ Always tighten the central screw of the head-to-head cross connector before tightening the lateral locking nuts.

#### ⚠ CAUTION

The connector is unusable if extended set screws are not present!

- ▶ Do not swap set screws with other components of the connectors.
- ▶ Only use connectors if all the set screws are present.


Fig. 47

#### Note

The head-to-head cross connectors allow angulation in the following planes: 45° lateral, 5° medial, 15° caudal and 15° cranial, as shown at the top of the image, see Fig. 47.

- ▶ Measure the distance e between the centers of the rods using the caliper 49 (SZ029R) to choose the connector size, as seen on the left side of the image, see Fig. 47.
- ▶ Place the extended set screw a and finally tightened it in the polyaxial screw head, see Fig. 47, see Chapter 3.16.
- ▶ Attach the cross connector b over the extended set screw a. If needed, loosen the center d to adjust the connector length.


Fig. 48

► Attach the lateral locking nuts **c** to locking driver **58** (SZ089R) and place them over the cross connector. Provisionally hand-tighten the components.

► Finally tighten the center screw **d** of the cross connector with the torque limiting handle **4/5** (SZ005R/SZ006R) and counter torque handle **6** (SZ008R).

The counter torque handle does not fit on the connector. Attach counter torque handle **6** (SZ008R) next to the cross connector, as shown at the center of the image, see Fig. 48.

► Attach the locking driver for head-to-head connector **58** (SZ089R) to a torque limiting handle **4/5** (SZ005R/SZ006R).

– Insert the locking driver for head-to-head connector **58** while pulling on the sleeve of the handle **4/5**.

– Finally tighten the locking nuts **c** using the locking driver for head-to-head cross connectors SZ089R and a torque limiting handle **4/5** (SZ005R/SZ006R) together with the counter torque handle **6** (SZ008R), as shown on the right side of the image, see Fig. 48.

– The counter torque handle does not fit on the connector. Attach counter torque handle **6** (SZ008R) next to the cross connector.

► For risks related to the final tightening of the set screw, see Chapter 3.16.

### 3.20.3 Domino connectors

Domino connectors allow to connect an existing spinal rod construct to a new spinal rod construct. Axial domino connectors connect 2 rods in the same alignment, lateral domino connectors connect 2 parallel laterally offset rods, see Fig. 49, see Fig. 50.

Cervical domino connectors allow to connect two  $\varnothing$  3.5 mm or  $\varnothing$  4.0 mm rods and are color coded in gold.

Thoracic domino connectors allow changing from a cervical  $\varnothing$  3.5 mm or  $\varnothing$  4.0 mm rod to a thoracic  $\varnothing$  5.5 mm rod. All thoracic rod options have a green color coding.

#### ⚠ WARNING

##### Risk of instability!

**The set screw is not properly seated as rods are not fully inserted!**

► **Make sure to insert the rods up to the hard stop in the axial domino connector.**

#### ⚠ CAUTION

**The connector is unusable if set screws are removed or not present!**

► Do not unscrew the set screws from the connectors under any circumstances.

► Do not swap set screws with other components of the connectors.

► Only use connectors if all the set screws are present.

#### ⚠ CAUTION

**Delay of surgery due to the rod being inserted into an incorrectly sized hole or due to wrong sized connector chosen!**

► Check the laser marking and color coding.

#### ⚠ CAUTION

**Delay of surgery due to incorrectly sized rod inserted into the connector!**

► Check the laser marking and color coding.


Fig. 49


Fig. 50

► Use connector templates **50** to **55** (SZ030R to SZ035R) to help in choosing the domino connector that fits the patient's anatomy, as shown on the left side of the image, see Fig. 50.

► Use handle **57** (SZ039R) to help insertion of the templates, see Fig. 50.

► Choose the correct connector. Check the laser coding and color-coding.

► Insert the connector into the implanted rod.

– For axial domino connectors, make sure to insert the rod until you reach the hard stop.

– For lateral domino connectors, make sure to insert the rod completely through the connector.

- ▶ Temporarily tighten the integrated set screw: place the set screw with the set screw starter **22/23** (SZ068R/SZ069R) aligned axially to the connector and screw in until the connector is fixed to the rod.
- ▶ Insert the second rod completely into the connector. Make sure to insert the rod completely.
- ▶ Temporarily tighten the integrated set screw, as described two steps prior.


Fig. 51

- ▶ Once the rod connector is positioned correctly, the set screws can be finally tightened, see Fig. 51, see Chapter 3.16.
- ▶ For risks related to the final tightening of the set screw, see Chapter 3.16.

**3.20.4 Cable connectors**

Cable connectors allow to connect a cerclage wire to the rod construct.


Fig. 52

- ▶ Place the connector **a** on rod **b**, insert and tighten the set screw **c**, see Fig. 52.
- ▶ For risks related to the final tightening of the set screw, see Chapter 3.16.

**3.20.5 Lateral offset connectors (straight and L-shaped)**

Lateral offset connectors can be used if a polyaxial screw or hook must be placed laterally to the longitudinal axis of the rod. Lateral offset connectors L-shaped are available on special request to offer variable placement of hooks or polyaxial screws.


Fig. 53

- ▶ Use connector template **56** (SZ036R) to check that the offset connector fits the patient's anatomy (optional).
- ▶ Place the offset connector **a** onto the rod **b**.
- ▶ Place the connector in the polyaxial screw head **c**, insert and provisionally tighten the set screws.


Fig. 54

- ▶ Once the connector is in correct position, finally tighten the set screws to 2.8 N·m, see Chapter 3.16.
- ▶ For risks related to the final tightening of the set screw, see Chapter 3.16.

### 3.21 Use dual diameter rod SQ009TS

The dual diameter rod can be used to form a construct that passes over the cervico-thoracic junction from an Ennovate Cervical construct to an Ennovate Thoracolumbar construct using only one rod.

#### ⚠ WARNING

Risk of breakage of the rod.

- ▶ Place the rod in such a way that the Ø 5.5 mm rod goes up as close as possible to the Ennovate Cervical screw. The Ø 5.5 mm rod is more stable than the Ø 4.0 mm rod.

#### Note

For all risks related to preparing and inserting the rod, see Chapter 3.10, see Chapter 3.12 and see Chapter 3.13.

### 3.22 Insert laminoplasty plate

#### ⚠ WARNING

Trauma to the spinal cord and nerve roots due to incorrect application!

- ▶ Select the appropriate size Ennovate Cervical laminoplasty plate length, polyaxial lateral mass screw length and diameter and SecureSpan screw length according to the individual indication, preoperative planning and bone situation found intraoperatively.

#### ⚠ CAUTION

Damage to or breakage of the plate caused by excessive strain on the material!

- ▶ Always bend the plate in one direction only.
- ▶ Do not bend back the plate.

Ennovate Cervical offers a special laminoplasty plate that allows combination of laminoplasty with unilateral lateral-mass (LM) screw fixation using the unilateral subaxial approach.

For the description of the technique and the risks related to it, refer to the following documents:

- SecureSpan system documents (O47602, TA013131, TA013082)
- Ennovate Cervical operating manual (OO3102)

In this section, only the steps specific to the Ennovate Cervical Laminoplasty plate are described.

#### 3.22.1 Prepare laminoplasty plate and hole

- ▶ Prepare the hinge in the lamina as described in operating manual O47602.
- ▶ Choose the appropriately sized Ennovate Cervical laminoplasty plate and screw length.
- ▶ Bend the plate if needed due to the patient's anatomy. For additional information on plate bending and plate holding, see O47602 and TA013131.
- ▶ Prepare the holes for the polyaxial lateral mass screws as described in section "Prepare screw hole with center punch 11 (SZ050R)" (see Chapter 3.2) to section "Cut thread with tap (optional)" (see Chapter 3.7).
- ▶ For the fixation of the stabilization block, see O47602 and TA013131.

#### 3.22.2 Fix laminoplasty plate

#### ⚠ WARNING

Risk of injury due to screw breakage!

- ▶ Keep screwdriver on the screw axis while picking up and inserting screws. Do not bend.

#### ⚠ WARNING

Injury to patient due to instable fixation of the plate.

- ▶ Check that the plate is well fixed after removal of the screwdriver.
- ▶ Make sure that all screws (Ennovate Cervical lateral mass screw and SecureSpan screws) are fully seated in the plate: Ensure that the Ennovate Cervical lateral mass screw head is in contact with the caudal border of the laminoplasty plate hole.

#### ⚠ WARNING

Risk of damage to soft tissue/neural structures during implantation.

- ▶ To avoid tearing out the screws from the bone, or damaging the screws or plates, never over-tighten.

#### ⚠ WARNING

Risk of loosening of the laminoplasty plate when the Ennovate Cervical lateral mass screw is not inserted first or if screws are not centered.

- ▶ Always fix the lateral mass screw before inserting the SecureSpan screws.
- ▶ Make sure that all SecureSpan screws are inserted in the center of the plate holes.


Fig. 55


Fig. 56

- ▶ Fix the laminoplasty plate with the polyaxial lateral mass screw, see Fig. 55.
  - For the insertion of the screw, see Chapter 3.8.)
  - While inserting the screw, make sure that the screw head **b** (not the tulip **c**) of the lateral mass screw directly contacts the caudal border **a** of the laminoplasty plate hole. This firmly fixes the laminoplasty plate to the lamina, see Fig. 56.
- ▶ To fix the plate on the Lamina, use the Aesculap SecureSpan screws
  - Follow the hole preparation and screw insertions described in O47602 and TA013131.
  - Do not bend the screws.
  - Do not over-tighten.
- ▶ Once the SecureSpan screws are fixed, re-tighten the polyaxial lateral mass screw, if necessary, as the contact of the polyaxial screw to the plate might slightly loosen when tightening the SecureSpan screws.
- ▶ Re-align the screw tulips **c** as needed, as can be seen on the right side of the image, see Fig. 56.

### 3.22.3 Insert rod


Fig. 57

- ▶ Insert and fix the rods, see Chapter 3.13, see Chapter 3.14, see Chapter 3.15.
- ▶ For risks related to the insertion and the final tightening of the set screw, see Chapter 3.16, see Chapter 3.14.

## 4. Validated reprocessing procedure

### 4.1 General safety instructions

#### Note

Adhere to national statutory regulations, national and international standards and directives, and local, clinical hygiene instructions for sterile processing.

#### Note

For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of products.

#### Note

Mechanical reprocessing should be favored over manual cleaning as it gives better and more reliable results.

#### Note

Successful processing of this medical device can only be ensured if the processing method is first validated. The operator/sterile processing technician is responsible for this.

#### Note

If there is no final sterilization, then a virucidal disinfectant must be used.

#### Note

For up-to-date information about reprocessing and material compatibility, see B. Braun eIFU at [eifu.bbraun.com](http://eifu.bbraun.com)

The validated steam sterilization procedure was carried out in the Aesculap sterile container system.

### 4.2 General information

Dried or affixed surgical residues can make cleaning more difficult or ineffective and lead to corrosion. Therefore the time interval between application and processing should not exceed 6 h; also, neither fixating pre-cleaning temperatures >45 °C nor fixating disinfecting agents (active ingredient: aldehydes/alcohols) should be used.

Excessive measures of neutralizing agents or basic cleaners may result in a chemical attack and/or to fading and the laser marking becoming unreadable visually or by machine for stainless steel.

Residues containing chlorine or chlorides e.g. in surgical residues, medicines, saline solutions and in the service water used for cleaning, disinfection and sterilization will cause corrosion damage (pitting, stress corrosion) and result in the destruction of stainless steel products. These must be removed by rinsing thoroughly with demineralized water and then drying.

Additional drying, if necessary.

Only process chemicals that have been tested and approved (e.g. VAH or FDA approval or CE mark) and which are compatible with the product's materials according to the chemical manufacturers' recommendations may be used for processing the product. All the chemical manufacturer's application specifications must be strictly observed. Failure to do so can result in the following problems:

- Optical changes of materials, e.g. fading or discoloration of titanium or aluminum. For aluminum, the application/process solution only needs to be of pH >8 to cause visible surface changes.
- Material damage such as corrosion, cracks, fracturing, premature aging or swelling.
- ▶ Do not use metal cleaning brushes or other abrasives that would damage the product surfaces and could cause corrosion.
- ▶ Further detailed advice on hygienically safe and material-/value-preserving reprocessing can be found at [www.a-k-i.org](http://www.a-k-i.org), link to "AKI-Brochures", "Red brochure".

### 4.3 Single-use products

For the list of single use products, see Chapter 2.2.3.

- ▶ Do not reuse the product.

The reprocessing of the product affects its functionality. Risk of injury, illness or death due to soiling and/or impaired functionality of the product.

- ▶ Do not reprocess sterile, single-use products.

Intraoperative contamination with blood, secretions, and other fluids may render the affected component unsuitable for resterilization!

- ▶ Handle the drills and K-wires with new gloves only.
- ▶ Keep the drill and K-wire storage trays covered or closed.
- ▶ Process drill and K-wire storage devices/trays separately from instrument trays.
- ▶ Clean drills and K-wires must not be cleaned together with contaminated drills and K-wires.
- ▶ Reprocess the drills and K-wires individually and separately if no implant system storage devices are available, When doing this, make certain that the implant components are not damaged.
- ▶ Do not reuse surgically contaminated drills and K-wires.

### 4.4 Reusable products

For the list of instruments that can be reused and must be sterilized prior to each use, see Chapter 2.2.3.

There is no set maximum number of uses and reprocessing cycles for the product.

The life of the product is limited by damage, normal wear and tear, type and duration of use, as well as handling, storage and transport of the product.

A careful visual and functional inspection before the next use is the best opportunity to recognize a product that is no longer functional.

### 4.5 Preparations at the place of use

- ▶ If applicable, rinse non-visible surfaces preferably with deionized water, with a disposable syringe for example.
- ▶ Remove any visible surgical residues to the extent possible with a damp, lint-free cloth.
- ▶ Transport the dry product in a sealed waste container for cleaning and disinfection within 6 hours.

### 4.6 Preparing for cleaning

- ▶ Do not disassemble fixation screws which permanently connect the various components.
- ▶ Disassemble the product prior to cleaning, see Chapter 4.7.
- ▶ Open up product with hinges.

## 4.7 Disassembly

### 4.7.1 Depth gauge FW042R


Fig. 58

- ▶ Remove the inner shaft **a**, see Fig. 58.

### 4.7.2 Taps SZ051R to SZ053R


Fig. 59

- ▶ Loosen and unscrew the nut **c** and remove it from the sleeve **a** in the direction of the arrow, see Fig. 59.
- ▶ Remove sleeve **a** together with the spring **b** in the direction of arrow.

4.7.3 Variable drill guide SZ054R

Fig. 60

- ▶ Turn the depth stop **a** clockwise until it is removed from the guide, see Fig. 60.

4.7.4 Universal drill guide SZ055R


Fig. 61

- ▶ Remove the tubes: turn the golden trigger **d** (see Fig. 61) and pull out the tube.
- ▶ Thread out the golden depth stop **a** by turning clockwise.
- ▶ Remove the inner shaft **b** of the handle **c**: rotate the golden knob **e** counter-clockwise.

4.7.5 Rod benders, all instruments with springs and ratchets

- ▶ Open springs and disengage ratchets.

4.7.6 Bone probe SZ024R


Fig. 62

- ▶ If used, remove the additional handle **a** from the grey handle **b**, see Fig. 62.

4.7.7 Rod bender SZ076R

*Note*

The central roller **a** is not intended to be completely loose.


Fig. 63

- ▶ Push the central roller **a** all the way forward until it is free to flip up. It is not intended to completely be loose, see Fig. 63.

#### 4.7.8 Linear rod persuader SZ080R

##### ⚠ CAUTION

Risk of damage to the instrument when detaching the golden knob a from the shaft b: spring arms can bend.

- ▶ Do not detach golden knob a from shaft b for cleaning.


Fig. 64

- ▶ Turn the golden knob a counter-clockwise until the persuasion shaft b is disengaged from the main body c. Do not remove the golden knob a from the persuasion shaft b, see Fig. 64.

#### 4.7.9 Persuador Pistol Style SZ082R


Fig. 65

- ▶ Release the ratchet c completely by pressing the ratchet c down while holding the trigger d, see Fig. 65.
- ▶ Remove the inner shaft e from the persuader.
- ▶ Force the springs b to disengage by pulling the trigger.

## 4.8 Cleaning/Disinfection

### 4.8.1 Product-specific safety information on the reprocessing method

Danger to the patient! SZ075R/SZ077R must not be manually reprocessed.

- ▶ Reprocess SZ075R/SZ077R only by manually pre-cleaning followed by cleaning in a mechanical washer/disinfector.

Damage to or destruction of the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!

- ▶ Following the manufacturer's instructions, use cleaning and disinfecting agents
  - that are approved for plastics and high-grade steel,
  - that do not attack softeners (e.g., in silicone).
- ▶ Observe specifications regarding concentration, temperature and exposure time.
- ▶ Do not exceed the maximum allowable disinfection temperature of 95 °C.

- ▶ Use suitable cleaning/disinfecting agents if the product is disposed of in a wet condition. To prevent foaming and degradation of the efficacy of the process chemicals: prior to mechanical cleaning and disinfection, rinse the product thoroughly with running water
- ▶ Manually pre-clean the product (with a cleansing brush) if there are residues of bone, tissue or ancillary materials (such as plaster or bone cement).
- ▶ If the products can be securely fixed in machines or storage devices in such a way that they will be cleaned thoroughly, clean and disinfect them mechanically.

#### 4.8.2 Validated cleaning and disinfection procedure

Validated procedure	Specific requirements	Reference
Manual cleaning with immersion disinfection FW024R, FW036R, FW042R, FW071R, FW076R, FW427R, FW528R, FW674R, FW675R, SZ008R, SZ010R, SZ011R, SZ016R to SZ020R, SZ022R to SZ027R, SZ029R to SZ036R, SZ039R, SZ048R to SZ054R, SZ058R, SZ064R to SZ066R, SZ068R to SZ070R, SZ076R, SZ078R	<ul style="list-style-type: none"> <li>■ Suitable cleaning brush</li> <li>■ 20 ml disposable syringe</li> <li>■ Drying phase: Use a lint-free cloth or medical compressed air</li> </ul>	<p>Chapter Manual cleaning/disinfection and subsection:</p> <ul style="list-style-type: none"> <li>■ Chapter Manual cleaning with immersion disinfection</li> </ul>
Manual cleaning with ultrasound and immersion disinfection FW428R, SZ001R to SZ003R, SZ005R, SZ006R, SZ055R, SZ063R, SZ080R to SZ082R, SZ089R, SZ143R	<ul style="list-style-type: none"> <li>■ Suitable cleaning brush</li> <li>■ 20 ml disposable syringe</li> <li>■ Drying phase: Use a lint-free cloth or medical compressed air</li> </ul>	<p>Chapter Manual cleaning/disinfection and subsection:</p> <ul style="list-style-type: none"> <li>■ Chapter Manual cleaning with ultrasound and immersion disinfection</li> </ul>
Mechanical alkaline cleaning and thermal disinfection FW024R, FW071R, FW076R, FW528R, FW674R, FW675R, SZ010R, SZ011R, SZ016R to SZ020R, SZ022R, SZ023R, SZ025R to SZ027R, SZ029R, SZ039R, SZ050R, SZ064R to SZ066R, SZ068R to SZ070R, SZ078R	<ul style="list-style-type: none"> <li>■ Place the product in a tray that is suitable for cleaning (avoiding rinsing blind spots).</li> <li>■ Place products in the tray with their hinges open.</li> </ul>	<p>Chapter Mechanical cleaning/disinfection and subsection:</p> <ul style="list-style-type: none"> <li>■ Chapter Mechanical alkaline cleaning and thermal disinfection</li> </ul>
Manual pre-cleaning with brush and subsequent mechanical alkaline cleaning and thermal disinfection FW036R, FW042R, FW427R, FW428R, SZ008R, SZ024R, SZ030R to SZ036R, SZ048R, SZ049R, SZ051R to SZ055R, SZ058R, SZ075R to SZ076R, SZ080R to SZ082R, SZ143R	<ul style="list-style-type: none"> <li>■ Suitable cleaning brush</li> <li>■ 20 ml disposable syringe</li> <li>■ Place the product in a tray that is suitable for cleaning (avoiding rinsing blind spots).</li> <li>■ Place products in the tray with their hinges open.</li> <li>■ FW042R, SZ008R, SZ024R, SZ048R, SZ049R, SZ051R to SZ055R, SZ058R, SZ080R, SZ082R, SZ143R: Connect components with lumens and channels directly to the rinsing port of the injector carriage.</li> </ul>	<p>Chapter Mechanical cleaning/disinfection with manual pre-cleaning and subsection:</p> <ul style="list-style-type: none"> <li>■ Chapter Manual pre-cleaning with a brush</li> <li>■ Chapter Mechanical alkaline cleaning and thermal disinfection</li> </ul>
Manual pre-cleaning with ultrasound and brush, and subsequent mechanical alkaline cleaning and thermal disinfection SZ001R to SZ003R, SZ005R, SZ006R, SZ063R, SZ089R	<ul style="list-style-type: none"> <li>■ 20 ml disposable syringe</li> <li>■ Place the product in a tray that is suitable for cleaning (avoiding rinsing blind spots).</li> <li>■ Place products in the tray with their hinges open.</li> <li>■ SZ001R, SZ002R, SZ003R, SZ005R, SZ006R: Rinse lumens using spray nozzle (3.8 bar water pressure) at least 3 times for 20 seconds each.</li> <li>■ SZ001R to SZ003R, SZ005R, SZ006R, SZ063R: Connect components with lumens and channels directly to the rinsing port of the injector carriage.</li> </ul>	<p>Chapter Mechanical cleaning/disinfection with manual pre-cleaning and subsection:</p> <ul style="list-style-type: none"> <li>■ Chapter Manual pre-cleaning with ultrasound and brush</li> <li>■ Chapter Mechanical alkaline cleaning and thermal disinfection</li> </ul>

## 4.9 Manual cleaning/disinfection

- ▶ Prior to manual disinfecting, allow water to drip off for a sufficient length of time to prevent dilution of the disinfecting solution.
- ▶ After manual cleaning/disinfection, check visible surfaces visually for residues.
- ▶ Repeat the cleaning/disinfection process if necessary.

### 4.9.1 Manual cleaning with immersion disinfection

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Disinfecting cleaning	RT (cold)	>15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
II	Intermediate rinse	RT (cold)	1	-	D-W	-
III	Disinfection	RT (cold)	5	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
IV	Final rinse	RT (cold)	1	-	FD-W	-
V	Drying	RT	-	-	-	-

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

RT: Room temperature

\*Recommended: BBraun Stabimed fresh

- ▶ Note the information on appropriate cleaning brushes and disposable syringes, see Chapter 4.8.2.

#### Phase I

- ▶ Fully immerse the product in the cleaning/disinfectant for at least 15 min. Ensure that all accessible surfaces are moistened.
- ▶ Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- ▶ If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
- ▶ Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- ▶ Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

#### Phase II

- ▶ Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- ▶ Drain any remaining water fully.

#### Phase III

- ▶ Fully immerse the product in the disinfectant solution.
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- ▶ Rinse lumens at least 5 times at the beginning of the exposure time using an appropriate disposable syringe. Ensure that all accessible surfaces are moistened.

#### Phase IV

- ▶ Rinse/flush the product thoroughly (all accessible surfaces).
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during final rinse.
- ▶ Rinse lumens with an appropriate disposable syringe at least five times.
- ▶ Drain any remaining water fully.

#### Phase V

- ▶ Dry the product in the drying phase with suitable equipment (e.g. cloth, compressed air), see Chapter 4.8.2.

#### 4.9.2 Manual cleaning with ultrasound and immersion disinfection

Phase	Step	D [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	<b>Ultrasonic cleaning</b>	RT (cold)	>15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
II	<b>Intermediate rinse</b>	RT (cold)	1	-	D-W	-
III	<b>Disinfection</b>	RT (cold)	5	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
IV	<b>Final rinse</b>	RT (cold)	1	-	FD-W	-
V	<b>Drying</b>	RT	-	-	-	-

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

RT: Room temperature

\*Recommended: BBraun Stabimed fresh

- ▶ Note the information on appropriate cleaning brushes and disposable syringes, see Chapter 4.8.2.

##### Phase I

- ▶ Clean the product in an ultrasonic cleaning bath (frequency 35 kHz) for at least 15 min. Ensure that all accessible surfaces are immersed and acoustic shadows are avoided.
- ▶ Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- ▶ If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
- ▶ Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- ▶ Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

##### Phase II

- ▶ Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- ▶ Drain any remaining water fully.

##### Phase III

- ▶ Fully immerse the product in the disinfectant solution.
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- ▶ Rinse lumens at least five times at the beginning of the exposure time with an appropriate disposable syringe. Ensure that all accessible surfaces are moistened.

##### Phase IV

- ▶ Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during final rinse.
- ▶ Rinse lumens with an appropriate disposable syringe at least five times.
- ▶ Drain any remaining water fully.

##### Phase V

- ▶ Dry the product in the drying phase with suitable equipment (e.g. cloth, compressed air), see Chapter 4.8.2.

## 4.10 Mechanical cleaning/disinfection

### Note

The cleaning and disinfection device must be of tested and approved effectiveness (e.g. FDA approval or CE mark according to DIN EN ISO 15883).

### Note

The cleaning and disinfection device used for processing must be serviced and checked at regular intervals.

### 4.10.1 Mechanical alkaline cleaning and thermal disinfecting

Machine type: single-chamber cleaning/disinfection device without ultrasound

Phase	Step	T [°C/°F]	t [min]	Water quality	Chemistry/Note
I	Pre-rinse	<25/77	3	D-W	-
II	Cleaning	55/131	10	FD-W	<ul style="list-style-type: none"> <li>■ Concentrate, alkaline:               <ul style="list-style-type: none"> <li>- pH = 13</li> <li>- &lt;5 % anionic surfactant</li> </ul> </li> <li>■ 0.5 % working solution               <ul style="list-style-type: none"> <li>- pH ~ 11*</li> </ul> </li> </ul>
III	Intermediate rinse	>10/50	1	FD-W	-
IV	Thermal disinfection	90/194	5	FD-W	-
V	Drying	-	-	-	According to the program for cleaning and disinfection device

DW: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

\*Recommended: BBraun Helimatic Cleaner alkaline

- ▶ Check visible surfaces for residues after mechanical cleaning/disinfecting.

## 4.11 Mechanical cleaning/disinfection with manual pre-cleaning

### Note

The cleaning and disinfection device must be of tested and approved effectiveness (e.g. FDA approval or CE mark according to DIN EN ISO 15883).

### Note

The cleaning and disinfection device used for processing must be serviced and checked at regular intervals.

### 4.11.1 Manual pre-cleaning with a brush

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Disinfecting cleaning	RT (cold)	>15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
II	Rinsing	RT (cold)	1	-	D-W	-

D-W: Drinking water

RT: Room temperature

\*Recommended: BBraun Stabimed fresh

- ▶ Note the information on appropriate cleaning brushes and disposable syringes, see Chapter 4.8.2.

#### Phase I

- ▶ Fully immerse the product in the cleaning/disinfectant for at least 15 min. Ensure that all accessible surfaces are moistened.
- ▶ Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- ▶ If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
- ▶ Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- ▶ Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

#### Phase II

- ▶ Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.

#### 4.11.2 Manual pre-cleaning with ultrasound and brush

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Ultrasonic cleaning	RT (cold)	>15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
II	Rinsing	RT (cold)	1	-	D-W	-

D-W: Drinking water

RT: Room temperature

\*Recommended: BBraun Stabimed fresh

- Note the information on appropriate cleaning brushes and disposable syringes, see Chapter 4.8.2.

##### Phase I

- Clean the product in an ultrasonic cleaning bath (frequency 35 kHz) for at least 15 min. Ensure that all accessible surfaces are immersed and acoustic shadows are avoided.
- Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
- Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

##### Phase II

- Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.

#### 4.11.3 Mechanical alkaline cleaning and thermal disinfecting

Machine type: single-chamber cleaning/disinfecting machine without ultrasound

Phase	Step	T [°C/°F]	t [min]	Water quality	Chemical
I	Prerinse	<25/77	3	D-W	-
II	Cleaning	55/131	10	FD-W	<ul style="list-style-type: none"> <li>■ Concentrate, alkaline:               <ul style="list-style-type: none"> <li>- pH ~ 13</li> <li>- &lt;5 % anionic surfactant</li> </ul> </li> <li>■ 0.5 % working solution               <ul style="list-style-type: none"> <li>- pH = 11*</li> </ul> </li> </ul>
III	Intermediate rinse	>10/50	1	FD-W	-
IV	Thermal disinfection	90/194	5	FD-W	-
V	Drying	-	-	-	According to the program for cleaning and disinfection device

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

\*Recommended: BBraun Helimatic Cleaner alkaline

- Check visible surfaces for residues after mechanical cleaning/disinfecting.

## 4.12 Inspection

- ▶ Allow the product to cool down to room temperature.
- ▶ Dry the product if it is wet or damp.

### 4.12.1 Visual inspection

- ▶ Make sure all dirt has been removed. Pay particular attention to e.g. fitting surfaces, hinges, shanks, recesses, drilled grooves.
- ▶ For contaminated products: Repeat the cleaning and disinfection procedure.
- ▶ Check instruments having interfaces to the implant for damage, burrs, or deformation at the corresponding functional geometries as these may damage the implants. This applies in particular to:
  - Star bit geometries (for example of screwdrivers)
  - Surfaces in contact with screw heads (for example for screw drivers, screw body manipulator)
  - Surfaces in contact with nuts (locking driver 58)
  - Thread geometries for connection to the female thread of the screw head (for example screw drivers)
  - Surfaces in contact with the rod (for example rod benders, rod holding forceps, rod persuaders, counter torque handle, caliper for cross connectors, templates for rod connectors)
- ▶ Check the product for damage, e.g. insulation, corroded, loose, bent, broken, cracked, worn, heavily scratched or broken parts.
- ▶ Check the product for missing or faded labels and/or paint marks (e.g. color of handle, gold plating).
- ▶ Check products having long, narrow geometries (in particular rotating instruments) for deformations.
- ▶ Check cutting instruments for:
  - Sharpness (example tips of punches and bone probes)
  - Continuous cutting edge (may not have any notches, dents, or other damage)
- ▶ Check surfaces for variations in roughness.
- ▶ Check the product for burrs that may damage tissue or surgical gloves.
- ▶ Check the product for loose or missing parts.
- ▶ Check the cannulation of the instruments for continuity.
- ▶ Immediately put aside damaged or inoperative products and send them to Aesculap Technical Service, see Chapter 5..

### 4.12.2 Functional test

#### ⚠ CAUTION

**Damage (metal seizure/friction corrosion) to the product caused by insufficient lubrication!**

- ▶ **Oil any moving parts (e. g. joints, pusher components and threaded rods) prior to the functional test using oil suitable for the applied sterilization process (e. g. STERILIT® I oil spray JG600 for steam sterilization or STERILIT® I drip-feed lubricator JG598).**
- ▶ Assemble disassembled products, see Chapter 4.13.
- ▶ Check instruments having interfaces to other instruments for ease of movement and compatibility.
- ▶ Make sure that instruments having a square, hexagonal or AO connection are self-locking in Ennovate handles.
- ▶ Check rod persuader in line 31 for the following criteria:
  - Ease of operation of the thread
  - Can be assembled and disassembled
- ▶ Single out rod persuader if it is difficult to operate.
- ▶ Check that the rod persuader pistol style 48 can be assembled and disassembled.
- ▶ Check that the depth stop can be smoothly screwed in the variable and universal drill guides 14/15. Check that the depth stop can be moved in and out of the MIS universal drill guide 16.

- ▶ Check that the tubes can be attached and removed from the universal drill guide 15/16: check that once the tube is attached it does not disengage from the drill guide.
- ▶ Check self-retaining features of instruments (as for example for 22/23/24, 58).
- ▶ Check that the product functions correctly.
- ▶ Check all moving parts (e.g. hinges, locks/interlocks, sliding parts, etc.) for full and smooth movement.
- ▶ Check locks for secure latching function.
- ▶ Check for compatibility with associated products.
- ▶ Immediately put aside inoperative products and send them to Aesculap Technical Service, see Chapter 5..

## 4.13 Assembly

### 4.13.1 Depth gauge FW042R


Fig. 66

- ▶ Insert the inner shaft a: turn the shaft until the orientation of the groove enables you to insert the shaft, see Fig. 66.

### 4.13.2 Taps SZ051R to SZ053R


Fig. 67

- ▶ Slide the sleeve a with spring b over the tap c in the direction of the arrow and fasten it by turning the nut d clockwise, see Fig. 67.

4.13.3 Variable drill guide SZ054R

Fig. 68

- ▶ Insert the depth stop **a** into the guide and turn it counter-clockwise, see Fig. 68.

4.13.4 Universal drill guide SZ055R


Fig. 69

- ▶ Insert the depth stop **a** and turn it counter-clockwise, see Fig. 69.
- ▶ Insert the inner shaft **b** into the handle: turn the golden knob **e** clockwise and tighten.

4.13.5 Rod benders, all instruments with springs and ratchets

- ▶ Engage the ratchets.

4.13.6 Bone probe SZ024R


Fig. 70

- ▶ For the radiographic control additional handle **a** can be mounted to the bone probe. Turn the additional handle **a** clockwise into the grey handle **b** of the bone probe, see Fig. 70.

4.13.7 Rod bender SZ076R


Fig. 71

- ▶ Push the central roller **a** back so that it lays flat on the surface **b**, as shown at the top of the image, see Fig. 71.
- ▶ Push the central roller **a** back in the direction of the arrow, as shown at the bottom of the image, see Fig. 71.
- ▶ Press on the golden trigger **c** and move the central roller **a** completely back, as shown on the right side of the image, see Fig. 71.

4.13.8 Linear rod persuader SZ080R


Fig. 72

- ▶ Check that the golden knob **a** is well inserted into the persuasion shaft **b** and that the spring arms are not deformed.
- ▶ Hold the main body **c** vertically so that the arms **d** are flush with the main body **c**. Insert the persuasion shaft **b** into the main body **c**, see Fig. 72.

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