Aesculap® Ennovate® Cervical

Aesculap Spine

en	Instructions for use/Technical description MIS instruments	lv	Lietošanas instrukcijas/tehniskais apraksts MIS instrumentu
USA	Note for U.S. users This Instructions for Use is NOT intended for United States users. Please discard. The Instructions for Use for United States users can be obtained by visiting our website at www.aesculapImplantsystems.com . If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-866-229-3002. A paper copy will be provided to you upon request at no additional cost.	lt	Naudojimo instrukcija/techninis aprašas MIS instrumentai
de	Gebrauchsanweisung/Technische Beschreibung MIS Instrumente	ru	Инструкция по применению/Техническое описание MIS инструменты
fr	Mode d'emploi/Description technique Instruments MIS	cs	Návod k použití/Technický popis Nástroje MIS
es	Instrucciones de manejo/Descripción técnica Juego de MIS	pl	Instrukcja użytkowania/Opis techniczny Instrumenty MIS
it	Istruzioni per l'uso/Descrizione tecnica Strumenti MIS	sk	Návod na použitie/Technický opis Nástroje MIS
pt	Instruções de utilização/Descrição técnica MIS Instrumentos	hu	Használati útmutató/Műszaki leírás MIS műszerek
nl	Gebruiksaanwijzing/Technische beschrijving MIS instrumenten	sl	Navodila za uporabo/Tehnični opis MIS instrumenti
da	Brugsanvisning/Teknisk beskrivelse MIS instrumenter	hr	Upute za uporabu/Tehnički opis MIS instrumenti
sv	Bruksanvisning/Teknisk beskrivning MIS instrument	ro	Manual de utilizare/Descriere tehnică Instrumente MIS
fi	Käyttöohje/Tekninen kuvaus MIS-instrumentit	bg	Упътване за употреба/Техническо описание Инструменти MIS
		tr	Kullanım Kılavuzu/Teknik açıklama MIS alet
		el	Οδηγίες χρήσης/Τεχνική περιγραφή Εργαλεία MIS

1

2

3

4

5

6

7

8

9

10

11

12

13

14


15

16c

16b

16a

16

Aesculap® Ennovate® Cervical

MIS instruments

Legend

Ennovate Cervical set MIS

- 1 Ennovate Cervical MIS dilator handle (SZ145R)
- 2 Ennovate Cervical drill cannulated Ø 2.3 mm (SZ013SU/SZ013R)
- 3 Ennovate Cervical drill cannulated Ø 2.8 mm (SZ014SU/SZ014R)
- 4 Ennovate Cervical drill cannulated Ø 3.3 mm (SZ015SU/SZ015R)
- 5 Ennovate Cervical MIS obturator for universal drill guide (SZ144R)
- 6 Ennovate Cervical MIS and C1/C2 K-wire long Ø 1.5 mm (SZ137SU/SZ137R)
- 7 Ennovate Cervical MIS K-wire long Ø 1.0 mm (SZ158SU/SZ158R)
- 8 Ennovate Cervical MIS dilator small (SZ146R)
- 9 Ennovate Cervical MIS dilator large (SZ147R)
- 10 Ennovate Cervical MIS counter torque (SZ148R)
- 11 Ennovate Cervical MIS extended tab sleeve for MIS screw (SZ149R)
- 12 Ennovate Cervical MIS screwdriver (SZ150R)
- 13 Ennovate Cervical MIS rod inserter (SZ151R)
- 14 Ennovate Cervical MIS tab breaker (SZ153R)
- 15 Ennovate Cervical MIS drill removal tool (SZ157R)
- 16 Ennovate Cervical MIS universal drill guide (inclusive 3.5 mm / 3.6 mm and 4.0 mm tube) (SZ143R) consisting of:
 - 16a Drill guide body
 - 16b Drill guide 3.5 mm / 3.6 mm tube
 - 16c Drill guide 4.0 mm tube

The relative size of the instruments does not necessarily reflect their size in reality.

The depictions in this document are schematic representations only and may deviate from the actual product.

Contents

1.	About this document.	3
1.1	Scope	3
1.2	Safety messages.	3
2.	Clinical use.	3
2.1	Areas of use and limitations of use	3
2.1.1	Intended use.	3
2.1.2	Indications	3
2.1.3	Contraindications.	3
2.2	Safety information.	3
2.2.1	Clinical user	3
2.2.2	Product	4
2.2.3	Sterility.	4
3.	Application	5
3.1	Safe application.	5
3.2	MIS streamlined workflow	5
3.2.1	Prepare MIS universal drill guide 16 (SZ143R)	5
3.2.2	Insert MIS universal drill guide.	6
3.2.3	Open cortex: create pilot hole	6
3.2.4	Drill hole.	7

3.2.5	Remove drill guide body.	7
3.2.6	Insert dilators	7
3.2.7	Remove drill.	8
3.2.8	Insert K-Wire.	8
3.2.9	Remove small dilator and drill tube	8
3.2.10	Cut thread with tap (optional).	8
3.2.11	Prepare screw	8
3.2.12	Perform navigated screw placement.	9
3.2.13	Prepare rod	9
3.2.14	Insert rod	9
3.2.15	Insert set screw	10
3.2.16	Tighten set screw	10
3.2.17	Remove long tabs	10
3.3	MIS extended workflow	11
3.3.1	Insert K-wire into the drill removal tool	11
3.3.2	Insert drill removal tool	11
3.3.3	Align K-wire to the drill tip	11
3.3.4	Remove drill.	12
4.	Validated reprocessing procedure	12
4.1	General safety instructions	12
4.2	General information	12
4.3	Single-use products	13
4.4	Reusable products.	13
4.5	Preparations at the place of use	13
4.6	Preparing for cleaning	13
4.7	Disassembly	13
4.7.1	MIS universal drill guide (SZ143R)	13
4.7.2	MIS rod inserter (SZ151R)	13
4.7.3	Drill removal tool (SZ157R)	13
4.8	Cleaning/Disinfection	14
4.8.1	Product-specific safety information on the reprocessing method.	14
4.8.2	Validated cleaning and disinfection procedure	14
4.9	Manual cleaning/disinfection	15
4.9.1	Manual cleaning with immersion disinfection.	15
4.9.2	Manual cleaning with ultrasound and immersion disinfection	16
4.10	Mechanical cleaning/disinfection	17
4.10.1	Mechanical alkaline cleaning and thermal disinfecting	17
4.11	Mechanical cleaning/disinfection with manual pre-cleaning.	18
4.11.1	Manual pre-cleaning with a brush	18
4.11.2	Manual pre-cleaning with ultrasound and brush	19
4.11.3	Mechanical alkaline cleaning and thermal disinfecting	19
4.12	Inspection	20
4.12.1	Visual inspection.	20
4.12.2	Functional test	20
4.13	Assembly	20
4.13.1	MIS universal drill guide (SZ143R)	20
4.13.2	MIS rod inserter (SZ151R)	20
4.13.3	Drill removal tool (SZ157R)	20

4.14	Packaging	21
4.15	Steam sterilization.....	21
4.16	Storage	21
5.	Technical service	21
6.	Disposal	21
7.	Symbols on product and packaging	21

1. About this document

Note

General risk factors associated with surgical procedures are not described in these instructions for use.

1.1 Scope

These instructions for use are limited to the specificities of the MIS technique and apply only to the instruments listed on the foldout page.

Note

The applicable CE mark for the product can be seen on the label or packaging of the product.

- ▶ For article-specific instructions for use as well as information on material compatibility and lifetime see B. Braun eFU at eifu.bbraun.com

1.2 Safety messages

Safety messages make clear the dangers to patient, user and/or product that could arise during the use of the product. Safety messages are labeled as follows:

DANGER

Indicates a possible threat of danger. If not avoided, death or serious injury may result.

WARNING

Indicates a possible threat of danger. If not avoided, minor or moderate injury may result.

CAUTION

Indicates a possible threat of material damage. If not avoided, the product may be damaged.

2. Clinical use

2.1 Areas of use and limitations of use

2.1.1 Intended use

The Ennovate Cervical MIS instruments are used for the minimally invasive implantation of Ennovate Cervical system implants. The indications and contraindications are described in the instruction for use of the Ennovate Cervical implants (TA015777).

2.1.2 Indications

Note

The manufacturer is not responsible for any use of the product against the specified indications and/or the described applications.

For indications, see Chapter 2.1.1.

2.1.3 Contraindications

For contraindications, see Chapter 2.1.1.

2.2 Safety information

2.2.1 Clinical user

General safety information

To prevent damage caused by improper setup or operation, and to not compromise the manufacturer warranty and liability:

- ▶ Use the product only according to these instructions for use.
- ▶ Follow the safety and maintenance instructions.
- ▶ Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge and experience.
- ▶ Store any new or unused products in a dry, clean, and safe place.
- ▶ Prior to use, check that the product is in good working order.
- ▶ Keep the instructions for use accessible for the user.

Note

The user is obligated to report all severe events in connection with the product to the manufacturer and the responsible authorities of the state in which the user is located.

Notes on surgical procedures

It is the user's responsibility to ensure that the surgical procedure is performed correctly.

Appropriate clinical training as well as a theoretical and practical proficiency of all the required operating techniques, including the use of this product, are prerequisites for the successful use of this product.

The user is required to obtain information from the manufacturer if there is an unclear preoperative situation regarding the use of the product.

2.2.2 Product

Product-specific safety information

These instructions for use are limited to the specificities of Ennovate Cervical MIS instruments. Several steps are described in the instructions for use for the basic technique. The MIS technique is always performed with the Brainlab navigation system. Refer to the operating manual and instructions for use dedicated to navigation. Refer to the following documents:

Art. no.	Designation
TA015777	IFU ENNOVATE CERVICAL IMPLANT SYSTEM FOR POSTERIOR, OCCIPITAL, CERVICAL AND THORACIC STABILIZATION
TA015763	IFU ENNOVATE CERVICAL BASIC INSTRUMENTS
TA015762	IFU ENNOVATE CERVICAL NAVIGATED INSTRUMENTS
TA015755	IFU ENNOVATE CERVICAL TRAYS
O03102	ENNOVATE CERVICAL - SURGICAL TECHNIQUE
O03202	ENNOVATE CERVICAL NAVIGATION - SURGICAL TECHNIQUE
O03302	ENNOVATE CERVICAL MIS - SURGICAL TECHNIQUE

- ▶ Only combine Aesculap products with each other.
- ▶ Only use Aesculap Ennovate instruments and accessories.
- ▶ **When disposing of or recycling the product, ensure that the packaging prevents injury by the product.**
- ▶ Observe the graphic/packing template when storing the instruments.
- ▶ Use the product only in applications according to its intended use, see Chapter 2.1.1.

Instruments can fall in situ or on the floor!

Handles can loosen during use!

- ▶ Check the seating of the coupling handles on the instruments for tightness.

Injury of spinal cord, nerve roots and blood vessels due to piercing of the cortex in case of incorrect application!

- ▶ Only carry out alignment and insertion of the instruments and polyaxial screws with the aid of a navigation system.

2.2.3 Sterility

Non-sterile, reusable products

Art. no.	Designation
SZ145R	Ennovate Cervical MIS dilator handle
SZ144R	Ennovate Cervical MIS obturator for universal drill guide
SZ146R	Ennovate Cervical MIS dilator small
SZ147R	Ennovate Cervical MIS dilator large
SZ148R	Ennovate Cervical MIS counter torque
SZ149R	Ennovate Cervical MIS extended tab sleeve for MIS screw
SZ150R	Ennovate Cervical MIS screwdriver
SZ151R	Ennovate Cervical MIS rod inserter
SZ153R	Ennovate Cervical MIS tab breaker
SZ157R	Ennovate Cervical MIS drill removal tool
SZ143R	Ennovate Cervical MIS universal drill guide

- ▶ Inspect the new product after removing its transport packaging and prior to first use to ensure it is in good working order.

The product is delivered in a non-sterile condition.

- ▶ Clean the new product after removing its transport packaging and prior to its initial sterilization.

Non-sterile, single-use products

Art. no.	Designation
SZ013R	Ennovate Cervical drill cannulated Ø 2.3 mm
SZ014R	Ennovate Cervical drill cannulated Ø 2.8 mm
SZ015R	Ennovate Cervical drill cannulated Ø 3.3 mm
SZ137R	Ennovate Cervical MIS and C1/C2 K-wire long Ø 1.5 mm
SZ158R	Ennovate Cervical MIS K-wire long Ø 1.0 mm

The product is delivered in a non-sterile condition.

- ▶ Clean the new product after removing its transport packaging and prior to its initial sterilization.
- ▶ Do not reuse the product.

Sterile, single-use products

Art. no.	Designation
SZ013SU	Ennovate Cervical drill cannulated Ø 2.3 mm
SZ014SU	Ennovate Cervical drill cannulated Ø 2.8 mm
SZ015SU	Ennovate Cervical drill cannulated Ø 3.3 mm
SZ137SU	Ennovate Cervical MIS and C1/C2 K-wire long Ø 1.5 mm
SZ158SU	Ennovate Cervical MIS K-wire long Ø 1.0 mm

The product has been sterilized by irradiation and is supplied in sterile packaging.

- ▶ Do not use products from open or damaged sterile packaging.
- ▶ Do not use the product after its use-by date.
- ▶ Do not reuse the product.

The reprocessing of the product affects its functionality. Risk of injury, illness or death due to soiling and/or impaired functionality of the product.

- ▶ Do not reprocess the product.

3. Application

Note

All handles and handle surfaces are grey. All actuating elements that are actuated during the operation are gold colored.

3.1 Safe application

When using K-wires, observe the following two warning notices:

DANGER

Injury of spinal cord, nerve roots and blood vessels due to piercing of the cortex in the case of incorrect application!

- ▶ Only carry out the alignment and insertion of the instruments, K-wires and polyaxial screws under radiographic control or with the aid of a navigation system.
- ▶ Carefully guide the instruments over the K-wire. When doing this, make sure that the K-wire remains in its position. It is recommended to hold the K-wire and to check the striped marking on the K-wire.
- ▶ Do not bend K-wires. Immediately replace bent K-wires.
- ▶ Inspect the cannulation of the instruments for blockage before every application.

DANGER

Injury of spinal cord, nerve roots and blood vessels due to bent K-wires or K-wires not holding in the bone!

- ▶ Do not reuse K-wires.
- ▶ Carefully drive in K-wires.
- ▶ Do not bend K-wires. Immediately replace bent K-wires.
- ▶ Insert the K-wires in the correct direction, with the roughened end pointing forward.
- ▶ Only use K-wires SZ137SU/SZ137R, SZ158SU/SZ158R of the Ennovate Cervical system.

WARNING

Risk of injury and/or malfunction!

- ▶ Prior to each use, inspect the product for loose, bent, broken, cracked, worn, or missing components.
- ▶ Always carry out a function test prior to each use of the product.

WARNING

Risk of injury to the patients due to damaged instruments!

Ennovate Cervical instruments are precise and can be easily damaged.

- ▶ Handle Ennovate Cervical instruments with the greatest of care.
- ▶ Do not repair or modify the instruments (e.g. sharpening of the tips). In case of malfunction send them to Aesculap Technical Service.

3.2 MIS streamlined workflow

For patient registration, navigation set-up and instrument calibration refer to the IFU for Ennovate Cervical navigated instruments (TA015762).

3.2.1 Prepare MIS universal drill guide 16 (SZ143R)

The tubes and the positive stop of the drills are color-coded to match the color of the screw heads:

- \varnothing 3.5 mm / 3.6 mm tube **16b** and \varnothing 2.3 mm drill **2** (SZ013SU/SZ013R) for \varnothing 3.5 mm / 3.6 mm polyaxial screws (silver).
- \varnothing 4.0 mm tube **16c** and \varnothing 2.8 mm drill **3** (SZ014SU/SZ014R) for \varnothing 4.0 mm polyaxial screws (gold).
- \varnothing 4.5 mm tube (SZ058R) and \varnothing 3.3 mm drill **4** (SZ015SU/SZ015R) for \varnothing 4.5 mm polyaxial screws (green)

DANGER

Risk of damage to the spinal cord, nerve roots, adjacent intervertebral space or soft tissue when drilling the holes!

- ▶ Always use the correct Ennovate Cervical drill guides when drilling holes.
- ▶ Do not combine drills and guides incorrectly.

DANGER

Risk of damage to the spinal cord, nerve roots, adjacent intervertebral space or soft tissue when drilling a hole that is too deep or too shallow.

- ▶ Use an appropriate imaging technique to select an appropriate drill length prior to the operation.
- ▶ Align and insert the drills only under radiographic monitoring and/or aided by a navigation system.
- ▶ Select a drill of a length equivalent to the intended drill hole depth.
- ▶ Always check the drill length with a caliper before starting to drill.
- ▶ Check the laser marking on the drill guide.

DANGER

Risk of damage to the spinal cord, nerve roots, adjacent intervertebral space or soft tissue when drilling a hole of an inaccurate size!

- ▶ Make sure to attach the appropriate drill tube. Check the color coding and laser marking.
- ▶ Make sure to use the right drill. Check the color coding and laser marking.

DANGER

Risk of injury and damage to the drill if the rotational speed is too high!

- ▶ Make sure to set a rotational speed appropriate for the intended application.

WARNING

Risk of revision due to instable fixation or damage to the bone!

Damage to the implant or instrument!

Prepared hole is too shallow!

- ▶ Do not sharpen the drill, as this makes the reading on the depth gauge inaccurate and falsifies the reading.
- ▶ Replace blunt drills with new ones.


Fig. 1

- ▶ Select the appropriate drill guide tube **16b/16c** or \varnothing 4.5 mm tube SZ058R.
- ▶ Attach the tube to the drill guide body **16a**, see Fig. 1.
- ▶ Select the appropriate drill **2/3/4**.
- ▶ Set the desired depth on the MIS universal drill guide **16** (SZ143R): hold down the golden square button **d** and adjust the depth stop **e** to the desired position. Release the golden button **d** to lock the depth stop **e** in place.
- ▶ Insert the required drill into the MIS universal drill guide **16** up to the positive stop and verify the exposed length with a caliper (e.g. AA845R). The universal drill guide offers a range of up to 40 mm.
- ▶ Attach the Aesculap navigation adapter with 4 marker spheres (SZ091R) to the universal drill guide, see TA015762.
- ▶ Perform the calibration, see TA015762.

3.2.2 Insert MIS universal drill guide

⚠ DANGER

The patient can suffer from severe complications if the MIS universal drill guide slips out of its position or/and is positioned incorrectly!

The position of the MIS universal drill guide affects all subsequent operative steps!

- ▶ When removing instruments (for example the obturator) during the procedure, make sure that the MIS universal drill guide remains firmly in place.
- ▶ When removing the obturator, make sure that the spring of the obturator does not collide with the depth stop of the drill guide.
- ▶ Always check the navigation to maintain the desired trajectory and positioning of the guide.

⚠ WARNING

Delay of surgery if the obturator is not inserted correctly!

- ▶ Check that the tip of the obturator is protruding out of the guide.
- ▶ Make a small stab incision through the skin, subcutaneous tissue and dorsal fascia. The incision can be dilated with a haemostatic forceps.


Fig. 2

- ▶ Introduce the MIS obturator **5** (SZ144R) into the MIS universal drill guide **16** (SZ143R). Check that the tip of the obturator **5** protrudes out of the guide, as shown on the left side of the image, see Fig. 2.
- ▶ Insert the assembly through the lateral stab incision into the surgical site. The flexibility of the soft tissue enables you to adjust its position in order to place it at the precise planned entry site for the pedicle screw.
- ▶ Remove the obturator **5** while holding the MIS universal drill guide **16**, thereby maintaining its trajectory. Before removing the obturator **5** align the flattened section **a** of the obturator handle with the opening **b** of the depth stop **c** on the drill guide **16**, as shown on the right side of the image, see Fig. 2. This ensures that the spring **d** on the shaft of the obturator **5** does not collide with the depth stop **c** of the drill guide **16**.

3.2.3 Open cortex: create pilot hole

⚠ WARNING

Delay of surgery due to the introduction of a too small cortical punch in the right tube!

- ▶ Make sure to select the right cortical punch. Check the color coding and laser marking

Cortical punches (SZ059R/SZ060R/SZ061R) are included in the Ennovate Cervical navigation set (see TA015762). The cortical punches are marked with a color coded ring according to the screw head color.

- ▶ Choose the appropriate cortical punch.
- ▶ Introduce the cortical punch into the MIS universal drill guide **16** (SZ143R).
- ▶ Position the cortical punch together with the MIS universal drill guide **16** on the desired entry point and advance the cortical punch into the bone until the stop is reached.

Note

The cortical punch has a positive stop in the drill guide **16** and protrudes 6 mm (independently of the depth stop).

- ▶ Remove the cortical punch.

3.2.4 Drill hole

⚠ WARNING

Operation delay if hole position is lost!

- ▶ Do not remove the drill after drilling the hole. Only remove the drill handle or the motor system. The drill remains in the bone for the next step (dilation).

The drill can be driven manually with a drill handle palm style (SZ003R) / drill handle teardrop style (SZ002R) or with a motor system and an AO small drill attachment (e.g. GB623R).


Fig. 3

- ▶ Attach the appropriate drill 2/3/4 (SZ013SU/SZ013R, SZ014SU/SZ014R, SZ015SU/SZ015R) to the desired drill handle a (SZ002R or SZ003R) or to the motor system as described in the IFU for Ennovate Cervical Basic instruments (TA015763).
- ▶ Verify that the depth stop b of the MIS universal drill guide 16 (SZ143R) is still set at the desired depth.
- ▶ Insert the drill into the MIS universal drill guide 16 and advance the drill to the predefined depth under trajectory control with the Brainlab navigation system until the stop is reached.
- ▶ Leave the drill in place and only remove the drill handle a (SZ002R or SZ003R) or the motor system.

3.2.5 Remove drill guide body

⚠ WARNING

Operation delay if hole position is lost!

- ▶ Do not remove drill from the bone until the dilators have been placed. Only remove the drill guide body 16a.
- ▶ When removing the drill guide body 16a, make sure that the drill does not get caught on the depth stop and gets removed from the bone.


Fig. 4

- ▶ Remove the drill guide body 16a while the drill tube 16b/16c or Ø 4.5 mm tube SZ058R remains in situ to protect the soft tissue:
 - Slide the golden lever arm d at the distal end of the drill guide body in a clockwise direction, as shown on the left side of the image, see Fig. 4
 - Hold the golden lever d in this position and unlock the drill guide body 16a from the drill tube 16b/16c or Ø 4.5 mm tube SZ058R by pulling the drill guide body upwards. Make sure that the drill does not get caught on the depth stop e.
 - Remove the drill guide body 16a by sliding it over the drill.

3.2.6 Insert dilators


Fig. 5

The fascia and muscles need to be dilated to allow for screw placement.

- ▶ Slide the dilators over the drill 2/3/4 (SZ013SU/SZ013R, SZ014SU/SZ014R, SZ015SU/SZ015R) and drill tube 16b/16c or Ø 4.5 mm tube SZ058R starting with the small dilator 8 (SZ146R) and continuing with the large dilator 9 (SZ147R). Push and twist to dilate tissue.
- ▶ If preferred, a dilator handle 1 (SZ145R) is available: slide the pins of the handle into the slots a of the dilators. Turn the dilator handle 1 clockwise.

The dilators should be docked on the bony anatomy to minimize the risk of soft tissue irritations.

3.2.7 Remove drill

⚠ CAUTION

Operation delay (change to open technique) if the hole is lost (and cannot be relocated) when removing the drill!

- ▶ Hold the drill tube firmly in place while removing the drill.
- ▶ Attach a desired handle (SZ002R/SZ003R) and remove the drill 2/3/4 (SZ013SU/SZ013R, SZ014SU/SZ014R, SZ015SU/SZ015R) by turning counter-clockwise.
- ▶ Once the drill is removed keep the drill tube firmly in place and immediately continue with the next surgical step inserting the K-wire.

3.2.8 Insert K-Wire

Two K-wires are available:

- Ø 1.0 mm K-wire long 7 (SZ158SU/SZ158R) for Ø 3.6 mm MIS screws
- Ø 1.5 mm K-wire long 6 (SZ137SU/SZ137R) for Ø 4.0 mm and 4.5 mm MIS screws

⚠ DANGER

Risk of injury if the K-wire is not properly fixed in the bone and thereby loses its position!

- ▶ Make sure to properly insert the K-wire in the pilot hole.
- ▶ Always check the position of the K-wire.

⚠ CAUTION

Operation delay (change to open technique) if the hole is lost (and cannot be relocated) before inserting the K-wire!

- ▶ Hold the drill tube firmly in place once the drill is removed until the K-wire is fully inserted.


Fig. 6

- ▶ Select the appropriate K-wire 6/7.
- ▶ Insert the K-wire 6/7 through the drill tube 16b/16c or Ø 4.5 mm tube SZ058R into the pilot hole. Make sure that the threaded tip of the K-wire 6/7 is located at the distal end.
- ▶ Insert the K-wire 6/7 until the thread is fully inserted in the bone. The end of thread indicates the end position of the screw.

3.2.9 Remove small dilator and drill tube


Fig. 7

- ▶ Remove the small dilator 8 (SZ146R) and the drill tube 16b/16c or Ø 4.5 mm tube SZ058R by hand. Constantly hold the K-wire 6/7 in place to make sure it does not move. Leave the large dilator 9 (SZ147R) in place as it serves as a working port during the subsequent surgical steps.

3.2.10 Cut thread with tap (optional)

Ennovate Cervical MIS screws are equipped with a fully threaded, tapered tip reducing the need to tap. The use of a tap is recommended in case of a hard bone. This step is optional.

- ▶ Perform navigated tapping, see TA015762.

3.2.11 Prepare screw

Cannulated long tab screws are available and dedicated to the MIS procedure. Refer to the surgical technique 003302 for more details.

⚠ WARNING

Operation delay if the long tabs of the MIS screw prematurely break or if implants get damaged due to loosening of the screw from the extended sleeve!

- ▶ Proceed with care when inserting the MIS long tabs into the extended tab sleeve.
- ▶ Visually inspect correct seating of the extended sleeve pin within one of the MIS sleeve holes.


Fig. 8

- ▶ Select the correct screw.
- ▶ Slide the extended tab sleeve 11 (SZ149R) over the MIS long tabs a in the direction of the arrow. Slightly press the long tabs a together to ease the sliding of the extended sleeve.
A click is heard when the small pin c at the distal end of the extended tab sleeve clicks into one of the small holes b of the MIS long tabs.

3.2.12 Perform navigated screw placement


Fig. 9

Note

The screw length L is measured from the tip of the screw to the sphere of the screw shaft, see Fig. 9.

To ensure maximum polyaxiality with the screw, do not tighten the screw completely down to the bone. Leave a small gap below the head to allow rotation and angulation of the screw head.

Prepare the MIS screwdriver 12 (SZ150R), see TA015762. In TA015762 the steps are described with the screwdriver SZ063R. Proceed exactly in the same way with SZ150R.

- ▶ Perform the navigated screw insertion, see TA015762.
- ▶ To check if the final screw position has been reached:
 - Check the screw depth information of the navigation system.
 - Check the relative position of the MIS screwdriver 12 in relation to the dilator 9 as described in the next step.


Fig. 10

- ▶ Check the relative position of the two lines a, b on the MIS screwdriver in relation to the upper edge of the dilator 9:
 - Make sure that the dilator is fully seated on the bone.
 - If the first line a on the screwdriver lines up with the upper edge of the dilator, the threaded part of the screw is seated in the bone.
 - If the second line b on the screwdriver lines up with the upper edge of the dilator, the full screw length is seated in the bone.
 - To ensure maximum polyaxiality of the screw, do not tighten the screw completely down to the second line b. Leave a small gap below the head to allow rotation and angulation of the screw head.
- ▶ Once the screw is placed, remove the large dilator 9.

3.2.13 Prepare rod

The MIS rods have a teardrop shaped geometry to engage with the rod inserter 13 (SZ151R). Rods are available in 10 mm increments. Do not cut the rods in order to maintain the blunt tip for screw insertion.

- ▶ Define the appropriate rod length with the navigation system or using rod templates.
- ▶ MIS rods have a predefined lordotic shape. If further bending is needed, bend the rod as described in TA015763.

Note

Be aware that an extremely bent rod might be difficult to insert through the MIS incision.

3.2.14 Insert rod

⚠ WARNING

Loss of the rod in situ may require the expansion of the skin incision!

- ▶ Do not decouple the rod inserter from the rod until the final position has been reached and the rod is completely fastened by at least one set screw.

⚠ WARNING

Inadequate fixation due to incorrect rod position!

- ▶ Check that the screw head is oriented orthogonally to the rod.
- ▶ Always position the MIS rods in such a way that the teardrop shaped tip and the conical tip are outside of the screw head.
- ▶ Select suitable rod lengths. When doing this, consider any possible corrective maneuver to be performed (e.g. distraction).
- ▶ Check the rod positioning before fastening with the set screw.

⚠ WARNING

Operation delay due to incomplete coupling of the rod on the rod inserter or incorrect alignment of the rod with the grooves of the extended sleeve.

- ▶ Insert the rod in the correct direction: check that the laser marked arrow on the rod inserter points to the rod tip.
- ▶ Verify a solid fixation of the rod in the rod inserter.
- ▶ Only tighten/loosen the golden rotary knob on the rod inserter by hand. Do not use any additional tools.
- ▶ Align the grooves of the extended sleeve so that the rod can be inserted.


Fig. 11

- ▶ Unscrew the golden knob **a** of the rod inserter **13** (SZ151R) and slide the teardrop shape **b** of the rod into the distal opening of the rod inserter so that the arrow **c** on the rod inserter points to the rod tip **b**.
- ▶ Turn the golden knob **a** clockwise to firmly tighten the rod in place.
- ▶ Turn the polyaxial screw so that the long slot of the extended tab sleeve **11** (SZ149R) is turned towards the rod. This can be done by hand or using the rod inserter **13**.
- ▶ Guide the rod down through the long slots of the extended tab sleeve **11**.

3.2.15 Insert set screw


Fig. 12

- ▶ Use the MIS counter torque **10** (SZ148R) to persuade the rod in the polyaxial screw head.
- ▶ Load the set screw on one of the set screw starters **a** (SZ068R or SZ069R), see TA015763.
- ▶ Insert the set screw through the long tabs of the MIS screw into the polyaxial body. Provisionally tighten the rod with the set screw, see TA015763.

3.2.16 Tighten set screw

The final tightening of the set screw is performed as described in TA015763. The only difference is that the MIS counter torque handle **10** (SZ148R) is used instead of the counter torque handle SZ008R.

- ▶ Refer to TA015763 for all risks related to this step.
- ▶ Once the final tightening is achieved, detach the rod inserter **13** (SZ151R) from the rod by turning the golden knob **a** counter-clockwise, see Fig. 11. Remove the connection by pulling the rod inserter **13** out of the rod, e.g. by levering the rod inserter towards the construct.
- ▶ Once all set screws are fully tightened, visually double-check that the set screws are fully seated.

The set screws are located plane with the top of the polyaxial screw tulips after final tightening with the \varnothing 4.0 mm MIS rod.

3.2.17 Remove long tabs

⚠ WARNING

Risk of correction loss / insufficient stabilization of the spinal column due to under-tightening of the set screw!

- ▶ **Break off the long tabs only once the set screws have been fully tightened.**


Fig. 13

- ▶ Remove the extended tab sleeve **11** (SZ149R): slightly press or pull the flute lateral to the MIS screw long tabs **a** to disengage the connection and slide off the extended tab sleeve.
- ▶ Remove the long tabs **a** with the MIS tab breaker **14** (SZ153R).

3.3 MIS extended workflow

The extended workflow offers an enhanced control of the screw entry point. The surgical workflow is equal to the streamlined workflow until the pilot hole has been drilled. In the extended workflow the cannulated drill stays in the pilot hole while the K-wire is inserted. A dedicated drill removal tool **15** (SZ157R) allows to remove the drill while the K-wire safely stays in place.

The only difference to the streamlined workflow concerns the removal of the drill (see Chapter 3.2.7) and the insertion of the K-wire (see Chapter 3.2.8). All other steps are performed in the exact same manner and are therefore not described in this section.

- ▶ Perform all steps of the streamlined workflow until you finish the insertion of the dilators (see Chapter 3.2.6).

Two K-wires are available:

- Ø 1.0 mm K-wire long **7** (SZ158SU/SZ158R) for Ø 3.6 mm MIS screws
- Ø 1.5 mm K-wire long **6** (SZ137SU/SZ137R) for Ø 4.0 mm and 4.5 mm MIS screws

3.3.1 Insert K-wire into the drill removal tool

- ▶ Select the appropriate K-wire **6/7**.


Fig. 14

- ▶ Insert the K-wire into the drill removal tool **15** (SZ157R):
 - Turn the golden K-wire lock **a** counter-clockwise to allow K-wire insertion.
 - Hold the K-wire with the threaded tip pointing down. After the insertion the threaded tip should be located at the distal end of the drill removal tool to be able to engage with the bone.
- ▶ Slide the K-wire into the drill removal tool **15** until the central black marking **b** is partly covered by the K-wire lock **a**.
- ▶ Re-tighten the K-wire by turning the K-wire lock **a** clockwise.
- ▶ Press the golden ratchet lever **c** and move the slider to the distal starting position.

3.3.2 Insert drill removal tool


Fig. 15

- ▶ Slide the assembly through the cannulated drill **2/3/4** (SZ013SU/SZ013R, SZ014SU/SZ014R, SZ015SU/SZ015R).
- ▶ Swing the drill removal tool **15** (SZ157R) slightly lateral so that the slider **a** passes the drill stop.
- ▶ Slide the assembly down until its distal part **b** sits on the large dilator **9** (SZ147R). This allows the slider to engage with the drill stop in the next steps.

3.3.3 Align K-wire to the drill tip


Fig. 16

- ▶ Turn the golden K-wire lock **a** counter-clockwise to release the K-wire.
- ▶ Advance the K-wire manually until the central black marking **b** on the K-wire is fully inserted into the cannulated drill **2/3/4** (SZ013SU/SZ013R, SZ014SU/SZ014R, SZ015SU/SZ015R). This aligns the K-wire tip with the drill tip.
- ▶ Re-tighten the K-wire by turning the K-wire lock **a** clockwise.

Note

To avoid a protruding K-wire do not advance it beyond the laser marking.

3.3.4 Remove drill

⚠ WARNING

Operation delay due to difficulty to remove the drill with drill removal tool (bone material stuck in the cannulated drill)!

- ▶ Change to the streamlined workflow if the drill cannot be removed with the drill removal tool.


Fig. 17

- ▶ Turn the golden hand wheel **c** clockwise. This extracts the drill **2/3/4** (SZ013SU/SZ013R, SZ014SU/SZ014R, SZ015SU/SZ015R) while holding the K-wire **6/7** (SZ137SU/SZ137R or SZ158SU/SZ158R) in place, as can be seen on the left side of the image, see Fig. 17.
- ▶ Once the drill is fully removed from the bone, release the K-wire lock **a** by turning it counter-clockwise, as can be seen on the right side of the image, see Fig. 17.
- ▶ Remove the drill removal tool **15** (SZ157R) and the drill **2/3/4**. Ensure that the K-wire is not pulled out of the pilot hole: hold the K-wire in place.

Note

Due to bone material in the tip of the cannulated drill it may happen that the black marking on the K-wire cannot be fully inserted into the cannulated drill and that the drill can not be removed with the drill removal tool. Please change to the streamlined workflow, see see Chapter 3.2.

- ▶ Perform all subsequent steps starting from the removal of the small dilator (see Chapter 3.2.9) to the end.

4. Validated reprocessing procedure

4.1 General safety instructions

Note

Adhere to national statutory regulations, national and international standards and directives, and local, clinical hygiene instructions for sterile processing.

Note

For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of products.

Note

Mechanical reprocessing should be favored over manual cleaning as it gives better and more reliable results.

Note

Successful processing of this medical device can only be ensured if the processing method is first validated. The operator/sterile processing technician is responsible for this.

Note

If there is no final sterilization, then a virucidal disinfectant must be used.

Note

For up-to-date information about reprocessing and material compatibility, see B. Braun eIFU at eifu.bbraun.com

The validated steam sterilization procedure was carried out in the Aesculap sterile container system.

4.2 General information

Dried or affixed surgical residues can make cleaning more difficult or ineffective and lead to corrosion. Therefore the time interval between application and processing should not exceed 6 h; also, neither fixating pre-cleaning temperatures >45 °C nor fixating disinfecting agents (active ingredient: aldehydes/alcohols) should be used.

Excessive measures of neutralizing agents or basic cleaners may result in a chemical attack and/or to fading and the laser marking becoming unreadable visually or by machine for stainless steel.

Residues containing chlorine or chlorides e.g. in surgical residues, medicines, saline solutions and in the service water used for cleaning, disinfection and sterilization will cause corrosion damage (pitting, stress corrosion) and result in the destruction of stainless steel products. These must be removed by rinsing thoroughly with demineralized water and then drying.

Additional drying, if necessary.

Only process chemicals that have been tested and approved (e.g. VAH or FDA approval or CE mark) and which are compatible with the product's materials according to the chemical manufacturers' recommendations may be used for processing the product. All the chemical manufacturer's application specifications must be strictly observed. Failure to do so can result in the following problems:

- Optical changes of materials, e.g. fading or discoloration of titanium or aluminum. For aluminum, the application/process solution only needs to be of pH >8 to cause visible surface changes.
- Material damage such as corrosion, cracks, fracturing, premature aging or swelling.
- ▶ Do not use metal cleaning brushes or other abrasives that would damage the product surfaces and could cause corrosion.
- ▶ Further detailed advice on hygienically safe and material-/value-preserving reprocessing can be found at www.a-k-i.org, link to "AKI-Brochures", "Red brochure".

4.3 Single-use products

For the list of single use products, see Chapter 2.2.3.

- ▶ Do not reuse the product.

The reprocessing of the product affects its functionality. Risk of injury, illness or death due to soiling and/or impaired functionality of the product.

- ▶ Do not reprocess sterile, single-use products.

Intraoperative contamination with blood, secretions, and other fluids may render the affected component unsuitable for resterilization!

- ▶ Handle the drills and K-wires with new gloves only.
- ▶ Keep the drill and K-wire storage trays covered or closed.
- ▶ Process drill and K-wire storage devices/trays separately from instrument trays.
- ▶ Clean drills and K-wires must not be cleaned together with contaminated drills and K-wires.
- ▶ Reprocess the drills and K-wires individually and separately if no implant system storage devices are available, When doing this, make certain that the implant components are not damaged.
- ▶ Do not reuse surgically contaminated drills and K-wires.

4.4 Reusable products

For the list of instruments that can be reused and must be sterilized prior to each use, see Chapter 2.2.3.

There is no set maximum number of uses and reprocessing cycles for the product.

The life of the product is limited by damage, normal wear and tear, type and duration of use, as well as handling, storage and transport of the product.

A careful visual and functional inspection before the next use is the best opportunity to recognize a product that is no longer functional.

4.5 Preparations at the place of use

- ▶ If applicable, rinse non-visible surfaces preferably with deionized water, with a disposable syringe for example.
- ▶ Remove any visible surgical residues to the extent possible with a damp, lint-free cloth.
- ▶ Transport the dry product in a sealed waste container for cleaning and disinfection within 6 hours.

4.6 Preparing for cleaning

- ▶ Do not disassemble fixation screws which permanently connect the various components.
- ▶ Disassemble the product prior to cleaning, see Chapter 4.7.
- ▶ Open up product with hinges.

4.7 Disassembly

4.7.1 MIS universal drill guide (SZ143R)


Fig. 18

- ▶ Remove the tubes 16b/16c.
- ▶ Pull the depth stop a all the way up so that it is loose and the inner part of the drill guide can be rinsed.
- ▶ Unscrew and remove the inner shaft b from the handle c.

4.7.2 MIS rod inserter (SZ151R)


Fig. 19

- ▶ Unscrew and remove the golden locking knob a.

4.7.3 Drill removal tool (SZ157R)

Fig. 20

- ▶ Turn the K-wire locking knob a counter-clockwise until it is loose..

4.8 Cleaning/Disinfection

4.8.1 Product-specific safety information on the reprocessing method

Damage to or destruction of the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!

- ▶ Following the manufacturer's instructions, use cleaning and disinfecting agents
 - that are approved for plastics and high-grade steel,
 - that do not attack softeners (e.g., in silicone).
- ▶ Observe specifications regarding concentration, temperature and exposure time.
- ▶ Do not exceed the maximum allowable disinfection temperature of 95 °C.

- ▶ Use suitable cleaning/disinfecting agents if the product is disposed of in a wet condition. To prevent foaming and degradation of the efficacy of the process chemicals: prior to mechanical cleaning and disinfection, rinse the product thoroughly with running water
- ▶ Manually pre-clean the product (with a cleansing brush) if there are residues of bone, tissue or ancillary materials (such as plaster or bone cement).
- ▶ If the products can be securely fixed in machines or storage devices in such a way that they will be cleaned thoroughly, clean and disinfect them mechanically.

4.8.2 Validated cleaning and disinfection procedure

Validated procedure	Specific requirements	Reference
Manual cleaning with immersion disinfection ■ SZ013R to SZ015R ■ SZ137R ■ SZ144R to SZ149R ■ SZ151R ■ SZ153R ■ SZ158R	■ Suitable cleaning brush ■ 20 ml disposable syringe ■ Drying phase: Use a lint-free cloth or medical compressed air	Chapter Manual cleaning/disinfection and subsection: ■ Chapter Manual cleaning with immersion disinfection
Manual cleaning with ultrasound and immersion disinfection ■ SZ143R ■ SZ150R ■ SZ157R	■ Suitable cleaning brush ■ 20 ml disposable syringe ■ Drying phase: Use a lint-free cloth or medical compressed air	Chapter Manual cleaning/disinfection and subsection: ■ Chapter Manual cleaning with ultrasound and immersion disinfection
Mechanical alkaline cleaning and thermal disinfection ■ SZ013R to SZ015R ■ SZ137R ■ SZ144R to SZ147R ■ SZ158R	■ Place the product in a tray that is suitable for cleaning (avoiding rinsing blind spots). ■ Place products in the tray with their hinges open. ■ SZ145R to SZ147R: Rinse the product under running water. Move the moving parts while doing so. ■ SZ145R to SZ147R: Connect components with lumens and channels directly to the rinsing port of the injector carriage.	Chapter Mechanical cleaning/disinfection and subsection: ■ Chapter Mechanical alkaline cleaning and thermal disinfecting
Manual pre-cleaning with brush and subsequent mechanical alkaline cleaning and thermal disinfection ■ SZ143R ■ SZ148R to SZ149R ■ SZ151R ■ SZ153R ■ SZ157R	■ Suitable cleaning brush ■ 20 ml disposable syringe ■ Place the product in a tray that is suitable for cleaning (avoiding rinsing blind spots). ■ Place products in the tray with their hinges open. ■ SZ143R, SZ148R and SZ151R: Connect components with lumens and channels directly to the rinsing port of the injector carriage.	Chapter Mechanical cleaning/disinfection with manual pre-cleaning and subsection: ■ Chapter Manual pre-cleaning with a brush ■ Chapter Mechanical alkaline cleaning and thermal disinfecting
Manual pre-cleaning with ultrasound and brush, and subsequent mechanical alkaline cleaning and thermal disinfection ■ SZ150R	■ 20 ml disposable syringe ■ Place the product in a tray that is suitable for cleaning (avoiding rinsing blind spots). ■ Place products in the tray with their hinges open. ■ Connect components with lumens and channels directly to the rinsing port of the injector carriage.	Chapter Mechanical cleaning/disinfection with manual pre-cleaning and subsection: ■ Chapter Manual pre-cleaning with ultrasound and brush ■ Chapter Mechanical alkaline cleaning and thermal disinfecting

4.9 Manual cleaning/disinfection

- ▶ Prior to manual disinfecting, allow water to drip off for a sufficient length of time to prevent dilution of the disinfecting solution.
- ▶ After manual cleaning/disinfection, check visible surfaces visually for residues.
- ▶ Repeat the cleaning/disinfection process if necessary.

4.9.1 Manual cleaning with immersion disinfection

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Disinfecting cleaning	RT (cold)	>15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
II	Intermediate rinse	RT (cold)	1	-	D-W	-
III	Disinfection	RT (cold)	5	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
IV	Final rinse	RT (cold)	1	-	FD-W	-
V	Drying	RT	-	-	-	-

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

RT: Room temperature

*Recommended: BBraun Stabimed fresh

- ▶ Note the information on appropriate cleaning brushes and disposable syringes, see Chapter 4.8.2.

Phase I

- ▶ Fully immerse the product in the cleaning/disinfectant for at least 15 min. Ensure that all accessible surfaces are moistened.
- ▶ Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- ▶ If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
- ▶ Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- ▶ Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

Phase II

- ▶ Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- ▶ Drain any remaining water fully.

Phase III

- ▶ Fully immerse the product in the disinfectant solution.
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- ▶ Rinse lumens at least 5 times at the beginning of the exposure time using an appropriate disposable syringe. Ensure that all accessible surfaces are moistened.

Phase IV

- ▶ Rinse/flush the product thoroughly (all accessible surfaces).
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during final rinse.
- ▶ Rinse lumens with an appropriate disposable syringe at least five times.
- ▶ Drain any remaining water fully.

Phase V

- ▶ Dry the product in the drying phase with suitable equipment (e.g. cloth, compressed air), see Chapter 4.8.2.

4.9.2 Manual cleaning with ultrasound and immersion disinfection

Phase	Step	D [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Ultrasonic cleaning	RT (cold)	>15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
II	Intermediate rinse	RT (cold)	1	-	D-W	-
III	Disinfection	RT (cold)	5	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
IV	Final rinse	RT (cold)	1	-	FD-W	-
V	Drying	RT	-	-	-	-

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

RT: Room temperature

*Recommended: BBraun Stabimed fresh

- ▶ Note the information on appropriate cleaning brushes and disposable syringes, see Chapter 4.8.2.

Phase I

- ▶ Clean the product in an ultrasonic cleaning bath (frequency 35 kHz) for at least 15 min. Ensure that all accessible surfaces are immersed and acoustic shadows are avoided.
- ▶ Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- ▶ If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
- ▶ Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- ▶ Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

Phase II

- ▶ Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- ▶ Drain any remaining water fully.

Phase III

- ▶ Fully immerse the product in the disinfectant solution.
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- ▶ Rinse lumens at least five times at the beginning of the exposure time with an appropriate disposable syringe. Ensure that all accessible surfaces are moistened.

Phase IV

- ▶ Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during final rinse.
- ▶ Rinse lumens with an appropriate disposable syringe at least five times.
- ▶ Drain any remaining water fully.

Phase V

- ▶ Dry the product in the drying phase with suitable equipment (e.g. cloth, compressed air), see Chapter 4.8.2.

4.10 Mechanical cleaning/disinfection

Note

The cleaning and disinfection device must be of tested and approved effectiveness (e.g. FDA approval or CE mark according to DIN EN ISO 15883).

Note

The cleaning and disinfection device used for processing must be serviced and checked at regular intervals.

4.10.1 Mechanical alkaline cleaning and thermal disinfecting

Machine type: single-chamber cleaning/disinfection device without ultrasound

Phase	Step	T [°C/°F]	t [min]	Water quality	Chemistry/Note
I	Pre-rinse	<25/77	3	D-W	-
II	Cleaning	55/131	10	FD-W	<ul style="list-style-type: none"> ■ Concentrate, alkaline: <ul style="list-style-type: none"> - pH = 13 - <5 % anionic surfactant ■ 0.5 % working solution <ul style="list-style-type: none"> - pH ~ 11*
III	Intermediate rinse	>10/50	1	FD-W	-
IV	Thermal disinfection	90/194	5	FD-W	-
V	Drying	-	-	-	According to the program for cleaning and disinfection device

DW: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

*Recommended: BBraun Helimatic Cleaner alkaline

- ▶ Check visible surfaces for residues after mechanical cleaning/disinfecting.

4.11 Mechanical cleaning/disinfection with manual pre-cleaning

Note

The cleaning and disinfection device must be of tested and approved effectiveness (e.g. FDA approval or CE mark according to DIN EN ISO 15883).

Note

The cleaning and disinfection device used for processing must be serviced and checked at regular intervals.

4.11.1 Manual pre-cleaning with a brush

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Disinfecting cleaning	RT (cold)	>15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
II	Rinsing	RT (cold)	1	-	D-W	-

D-W: Drinking water

RT: Room temperature

*Recommended: BBraun Stabimed fresh

- ▶ Note the information on appropriate cleaning brushes and disposable syringes, see Chapter 4.8.2.

Phase I

- ▶ Fully immerse the product in the cleaning/disinfectant for at least 15 min. Ensure that all accessible surfaces are moistened.
- ▶ Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- ▶ If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
- ▶ Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- ▶ Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

Phase II

- ▶ Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.

4.11.2 Manual pre-cleaning with ultrasound and brush

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Ultrasonic cleaning	RT (cold)	>15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
II	Rinsing	RT (cold)	1	-	D-W	-

D-W: Drinking water

RT: Room temperature

*Recommended: BBraun Stabimed fresh

- Note the information on appropriate cleaning brushes and disposable syringes, see Chapter 4.8.2.

Phase I

- Clean the product in an ultrasonic cleaning bath (frequency 35 kHz) for at least 15 min. Ensure that all accessible surfaces are immersed and acoustic shadows are avoided.
- Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
- Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

Phase II

- Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.

4.11.3 Mechanical alkaline cleaning and thermal disinfecting

Machine type: single-chamber cleaning/disinfecting machine without ultrasound

Phase	Step	T [°C/°F]	t [min]	Water quality	Chemical
I	Prerinse	<25/77	3	D-W	-
II	Cleaning	55/131	10	FD-W	<ul style="list-style-type: none"> ■ Concentrate, alkaline: <ul style="list-style-type: none"> - pH ~ 13 - <5 % anionic surfactant ■ 0.5 % working solution <ul style="list-style-type: none"> - pH = 11*
III	Intermediate rinse	>10/50	1	FD-W	-
IV	Thermal disinfection	90/194	5	FD-W	-
V	Drying	-	-	-	According to the program for cleaning and disinfection device

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

*Recommended: BBraun Helimatic Cleaner alkaline

- Check visible surfaces for residues after mechanical cleaning/disinfecting.

4.12 Inspection

- ▶ Allow the product to cool down to room temperature.
- ▶ Dry the product if it is wet or damp.

4.12.1 Visual inspection

- ▶ Make sure all dirt has been removed. Pay particular attention to e.g. fitting surfaces, hinges, shanks, recesses, drilled grooves.
- ▶ For contaminated products: Repeat the cleaning and disinfection procedure.
- ▶ Check instruments having interfaces to the implant for damage, burrs, or deformation at the corresponding functional geometries as these may damage the implants before use. This applies in particular to:
 - Star bit geometries (for example of screwdriver 12)
 - Surfaces in contact with screw heads (for example for screwdriver 12)
 - Surfaces in contact with screw tabs (tab breaker 14, tab sleeve 11)
 - Thread geometries for connection to the female thread of the screw head (for example screw driver 12)
 - Surfaces in contact with the rod (for example counter torque handle 10, rod inserter 13)
- ▶ Check product for damage, e.g. insulation, corroded, loose, bent, broken, cracked, worn, heavily scratched or broken parts.
- ▶ Check product for missing or faded labels and/or paint marks (e.g. color of handle, gold plating).
- ▶ Check products having long, narrow geometries (in particular rotating instruments) for deformations.
- ▶ Check surfaces for variations in roughness.
- ▶ Check product for burrs that may damage tissue or surgical gloves.
- ▶ Check product for loose or missing parts.
- ▶ Check the cannulation of the instruments for continuity.
- ▶ Immediately put aside damaged or inoperative products and send them to Aesculap Technical Service, see Chapter 5..

4.12.2 Functional test

⚠ CAUTION

Damage (metal seizure / friction corrosion) to the product caused by insufficient lubrication!

- ▶ **Oil any moving parts (e. g. joints, pusher components and threaded rods) prior to the functional test using oil suitable for the applied sterilization process (e. g. STERILIT® I oil spray JG600 for steam sterilization or STERILIT® I drip-feed lubricator JG598). Do not forget to oil the Ennovate Cervical MIS rod inserter SZ151R at the marked location.**
- ▶ Assemble disassembled products, see Chapter 4.13.
- ▶ Check instruments having interfaces to other instruments for ease of movement and compatibility.
- ▶ Make sure that instruments having a square, hexagonal or an AO connection are self-locking in the Ennovate handles.
- ▶ Check that the dilators 8/9 can be mounted and dismounted from the dilator handle 1.
- ▶ Check that the depth stop can be smoothly screwed into the universal drill guide 16.
- ▶ Check that tubes 16b/16c can be attached and removed from the universal drill guide body 16a: Check that once the tube is attached it does not disengage from the drill guide.

- ▶ Check that the obturator 5 can be inserted through the universal drill guide 16. Check that once it is inserted it does not disengage from the drill guide.
- ▶ Check that the product functions correctly.
- ▶ Check all moving parts (e.g. hinges, locks/interlocks, sliding parts, etc.) for full and smooth movement.
- ▶ Check for compatibility with associated products.
- ▶ Immediately put aside inoperative products and send them to Aesculap Technical Service, see Chapter 5..

4.13 Assembly

4.13.1 MIS universal drill guide (SZ143R)


Fig. 21

- ▶ Insert the depth stop a.
- ▶ Insert the shaft b into the handle c.

4.13.2 MIS rod inserter (SZ151R)


Fig. 22

- ▶ Insert and screw in the golden locking knob a.

4.13.3 Drill removal tool (SZ157R)

Fig. 23

- ▶ Engage the K-wire locking knob a by turning it clockwise.

4.14 Packaging

- ▶ Appropriately protect products with fine working tips.
- ▶ Store products with ratchet locks fully opened or locked no further than in the first notch.
- ▶ Place the product in its holder or on a suitable tray. Ensure that sharp edges are covered.
- ▶ Package trays appropriately for the sterilization process (e.g. in Aesculap sterile containers).
- ▶ Ensure that the packaging provides sufficient protection against contamination of the product during storage.

4.15 Steam sterilization

- ▶ Check to ensure that the sterilizing agent will come into contact with all external and internal surfaces (e.g., by opening any valves and faucets).
- ▶ Validated sterilization process
 - Steam sterilization using fractional vacuum process
 - Steam sterilizer according to DIN EN 285 and validated according to DIN EN ISO 17665
 - Sterilization using fractionated vacuum process at 134 °C / holding time 5 min
- ▶ If several devices are sterilized at the same time in the same steam sterilizer: Ensure that the maximum permitted load according to the manufacturers' specifications is not exceeded.

4.16 Storage

- ▶ Store sterile products in germ-proof packaging, protected from dust, in a dry, dark, temperature-controlled area.
- ▶ Store sterile single-use products in germ-proof packaging in a dust-protected, dry, dark and temperature-controlled room.

5. Technical service

⚠ CAUTION

Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

- ▶ Do not modify the product.
- ▶ For service and repairs, please contact your national B. Braun/Aesculap agency.

Service addresses

Aesculap Technischer Service

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78532 Tuttlingen / Germany

Phone: +49 7461 95-1601

Fax: +49 7461 16-2887

E-Mail: ats@aesculap.de

Other service addresses can be obtained from the address indicated above.

6. Disposal

⚠ WARNING

Risk of infection due to contaminated products!

- ▶ Adhere to national regulations when disposing of or recycling the product, its components and its packaging.

⚠ WARNING

Risk of injury due to sharp-edged and/or pointed products!

- ▶ When disposing of or recycling the product, ensure that the packaging prevents injury by the product.

Note

The user institution is obliged to reprocess the product before its disposal, see Chapter 4..

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