# Aesculap® Ennovate

## Aesculap Spine

<b>en</b>	<b>Instructions for use/Technical description</b> Instruments for MIS access	<b>lv</b>	<b>Lietošanas instrukcijas/tehniskais apraksts</b> MIS ķirurģijas instrumenti
<b>USA</b>	<b>Note for U.S. users</b> This Instructions for Use is NOT intended for United States users. Please discard. The Instructions for Use for United States users can be obtained by visiting our website at <a href="http://www.aesculapImplantsystems.com">www.aesculapImplantsystems.com</a> . If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-866-229-3002. A paper copy will be provided to you upon request at no additional cost.	<b>it</b>	<b>Naudojimo instrukcija/techninis aprašas</b> MIS (minimalios invazinės chirurgijos) instrumentai
<b>de</b>	<b>Gebrauchsanweisung/Technische Beschreibung</b> Instrumente für MIS-Zugang	<b>ru</b>	<b>Инструкция по применению/Техническое описание</b> Инструменты для MIS-доступа Ennovate
<b>fr</b>	<b>Mode d'emploi/Description technique</b> Instruments pour accès MIS	<b>cs</b>	<b>Návod k použití/Technický popis</b> Nástroje pro přístup MIS
<b>es</b>	<b>Instrucciones de manejo/Descripción técnica</b> Instrumental para abordaje mínimamente invasivo	<b>pl</b>	<b>Instrukcja użytkowania/Opis techniczny</b> Instrumenty do dostępu MIS
<b>it</b>	<b>Istruzioni per l'uso/Descrizione tecnica</b> Strumentario MIS	<b>sk</b>	<b>Návod na použitie/Technický opis</b> Nástroje pre prístup MIS
<b>pt</b>	<b>Instruções de utilização/Descrição técnica</b> Instrumentos para acesso MIS	<b>hu</b>	<b>Használati útmutató/Műszaki leírás</b> Eszközök minimálisan invazív beavatkozásokhoz (MIS) történő hozzáféréshoz
<b>nl</b>	<b>Gebruiksaanwijzing/Technische beschrijving</b> Instrumenten voor MIS-toegang	<b>sl</b>	<b>Navodila za uporabo/Tehnični opis</b> Instrumenti za minimalno invazivno operacijo
<b>da</b>	<b>Brugsanvisning/Teknisk beskrivelse</b> Instrumenter til MIS-adgang	<b>hr</b>	<b>Upute za uporabu/Tehnički opis</b> Instrumenti za pristup za minimalno-invazivnu kirurgiju (MIS)
<b>sv</b>	<b>Bruksanvisning/Teknisk beskrivning</b> Instrument för MIS-åtkomst	<b>ro</b>	<b>Manual de utilizare/Descriere tehnică</b> Instrumente pentru acces în endoprotetică de șold minim invazivă (MIS)
<b>fi</b>	<b>Käyttöohje/Tekninen kuvaus</b> MIS-sisäänvienti-instrumentit	<b>bg</b>	<b>Упътване за употреба/Техническо описание</b> Инструменти за достъп до МИС
		<b>tr</b>	<b>Kullanım Kılavuzu/Teknik açıklama</b> MIS erişimi için aletler
		<b>el</b>	<b>Οδηγίες χρήσης/Τεχνική περιγραφή</b> Εργαλεία ελάχιστα επεμβατικής χειρουργικής (MIS)

1

2

3

4

5

6

7

8

9

10

11

12

13

14

15

16

17

18

19

20

21

22

23

24

25

26

27

28

29

30

31

32

33

34

35

36

37

38

39

# Aesculap® Ennovate

## Instruments for MIS access

### Legend

- 1 Ennovate parallel distractor forceps SZ231R
- 2 Ennovate parallel compressor forceps SZ232R
- 3 Ennovate screw tap for Ø 4.5 mm to 10.5 mm SZ254R to SZ260R
- 4 Ennovate MIS bone access needle, short SZ364SU
- 5 Ennovate french rod bender, adjustable SZ270R
- 6 Ennovate MIS K-wire trocar, short SZ365R
- 7 Ennovate MIS K-wire sleeve, short SZ366R
- 8 Ennovate MIS K-wire forceps SZ367R
- 9 Ennovate MIS slotted hammer SZ368R
- 10 Ennovate MIS K-wire steel SZ369S, Ennovate MIS K-wire steel, sterile SZ369SU, Ennovate MIS K-wire nitinol SZ370, Ennovate MIS K-wire nitinol, sterile SZ370SU
- 11 Ennovate MIS skin incision guide SZ371R
- 12 Ennovate MIS dilator, S SZ372T
- 13 Ennovate MIS dilator, M SZ373T
- 14 Ennovate MIS dilator, L SZ374T
- 15 Ennovate MIS dilator handle SZ375T
- 16 Ennovate MIS pedicle probe lumbar, straight SZ376R
- 17 Ennovate MIS screw length measuring device SZ377R
- 18 Ennovate MIS downtube, short SZ378R/Ennovate FRI downtube SZ405R
- 19 Ennovate MIS locking wrench for downtube SZ379R
- 20 Ennovate MIS removal key, short SZ380R
- 21 Ennovate MIS screwdriver, short SZ381R
- 22 Ennovate MIS rod length caliper, short SZ382R
- 23 Ennovate MIS rod inserter, short SZ384R
- 24 Ennovate MIS rod inserter, long SZ385R
- 25 Ennovate MIS rod indicator, short SZ387R
- 26 Ennovate MIS parallel sleeve SZ388R
- 27 Ennovate MIS rod pusher, short SZ389R
- 28 Ennovate MIS rod pusher handle SZ390R
- 29 Ennovate set screw driver, short SZ391R
- 30 Ennovate MIS counter torque handle SZ392R
- 31 Ennovate PolyBlocker SZ393R
- 32 Ennovate MIS compression/distraction device SZ394T
- 33 Ennovate MIS C-Ring for compression/distraction device SZ395R
- 34 Ennovate MIS cage sleeve SZ396R
- 35 Ennovate MIS quick rod pusher, short SZ397R
- 36 Ennovate MIS alignment working port, short SZ398P
- 37 Ennovate MIS alignment tool, short SZ399R
- 38 Ennovate injection cannula MIS SR138SU
- 39 Ennovate cannula sleeve MIS SR139R

The depictions in this document are schematic representations only and may deviate from the actual product.

### Contents

1.	About this document .....	3
1.1	Scope .....	3
1.2	Warnings .....	3
2.	Clinical application .....	4
2.1	Areas of application and restrictions on use .....	4
2.1.1	Intended purpose .....	4
2.1.2	Indications .....	4
2.1.3	Absolute contraindications .....	4
2.1.4	Relative contraindications .....	4
2.2	Safety instructions .....	4
2.2.1	Clinical user .....	4
2.2.2	Product .....	4
2.2.3	Sterility .....	4
3.	Safe operation .....	5
3.1	Place K-wire .....	5
3.1.1	Place K-wire with bone access needle SZ364SU .....	5
3.1.2	Place K-wire with K-wire targeting device (SZ365R and SZ366R) .....	6
3.2	Perform skin incision with skin incision guide SZ371R .....	6
3.3	Perform dilation with dilator SZ372T to SZ374T and dilator handle SZ375T .....	6
3.4	Enlarge pedicle canal with pedicle probe SZ376R .....	7
3.5	Cut thread with screw tap SZ254R to SZ260R (optional) .....	7
3.6	Determine screw length with screw length measuring device SZ377R .....	8
3.7	Couple downtube SZ378R/SZ405R to pedicle screw .....	8
3.8	Screw in pedicle screw .....	8
3.8.1	Screw in pedicle screw with screwdriver SZ381R .....	8
3.8.2	Screw in fenestrated pedicle screw with injection cannula SR138SU and cannula sleeve SR139R and fix with bone cement .....	9
3.9	Determine rod length with rod length caliper SZ382R .....	11
3.10	Adjust rod length with rod cutter FW206R .....	11
3.11	Bend rod with french rod bender SZ270R .....	11
3.12	Insert rod with rod inserter SZ384R/SZ385R .....	12
3.13	Verify rod with rod indicator SZ387R .....	13
3.14	Push rod into screw head .....	13
3.14.1	Push rod into screw head with quick rod pusher SZ397R (option 1) .....	13
3.14.2	Push rod into screw head with rod pusher SZ389R (option 2) .....	13
3.15	Insert and tighten set screw with the insertion instrument for the set screw SZ391R .....	14
3.15.1	Take up set screw .....	14
3.15.2	Insert and tighten set screw with MIS quick rod pusher SZ397R (option 1) .....	15
3.15.3	Insert and tighten set screw with rod pusher SZ389R (option 2) .....	15
3.16	Remove downtube SZ378R/SZ405R with removal key SZ380R .....	15

3.17	Carry out corrective measure (optional) .....	16	4.13	Assembling .....	32
3.17.1	Prepare distraction/compression .....	16	4.13.1	Downtube SZ378R/SZ405R .....	32
3.17.2	Perform parallel distraction/compression with parallel distractor forceps SZ231R and parallel compressor forceps SZ232R (option 1) .....	16	4.13.2	Screwdriver SZ381R .....	32
3.17.3	Perform distraction with compression/distraction device SZ394T (option 2) .....	17	4.13.3	Rod length measuring instrument for short construction SZ382R .....	32
3.18	Interbody fusion management: distraction/compression with alignment tool SZ399R and cage sleeve SZ396R (optional) ..	18	4.13.4	Rod inserter for short construction SZ384R .....	32
3.18.1	Remove downtube SZ378R/SZ405R .....	18	4.13.5	Rod inserter for long construction SZ385R .....	33
3.18.2	Perform distraction/compression .....	19	4.13.6	Alignment working port SZ398P .....	33
3.19	Recouple downtube after cage implantation .....	19	4.13.7	Alignment tool SZ399R .....	33
3.19.1	Insert alignment working port SZ398P .....	19	4.14	Packaging .....	33
3.19.2	Recouple downtube SZ378R/SZ405R with alignment tool SZ380R .....	20	4.15	Steam sterilization .....	33
3.20	Recouple downtube SZ378R/SZ405R after accidental decoupling .....	21	4.16	Storage .....	33
3.20.1	Couple alignment tool SZ399R .....	21	5.	Maintenance and Service .....	33
3.20.2	Recouple downtube SZ378R/SZ405R .....	21	5.1	Technical Service .....	33
4.	Validated reprocessing procedure .....	22	5.2	Accessories/Spare parts .....	33
4.1	General safety instructions .....	22	6.	Disposal .....	33
4.2	General information .....	22	7.	Symbols on product and packaging .....	33
4.3	Single-use products .....	22			
4.4	Reusable products .....	22			
4.5	Preparations at the place of use .....	22			
4.6	Preparation before cleaning .....	22			
4.7	Disassembling .....	23			
4.7.1	Parallel distractor forceps SZ231R and parallel compressor forceps SZ232R .....	23			
4.7.2	K-wire trocar SZ365R and K-wire sleeve SZ366R .....	23			
4.7.3	Skin incision guide SZ371R .....	23			
4.7.4	Downtube SZ378R/SZ405R .....	23			
4.7.5	Screwdriver SZ381R .....	23			
4.7.6	Rod length measuring instrument for short construction SZ382R .....	24			
4.7.7	Rod inserter for short construction SZ384R .....	24			
4.7.8	Rod inserter for long construction SZ385R .....	24			
4.7.9	Alignment working port SZ398P .....	24			
4.7.10	Alignment tool SZ399R .....	24			
4.7.11	Injection cannula SR138SU and cannula sleeve SR139R .....	24			
4.8	Cleaning/disinfection .....	25			
4.8.1	Product-specific safety notes on the reprocessing procedure ..	25			
4.8.2	Validated cleaning and disinfection procedure .....	25			
4.9	Manual cleaning/disinfection .....	27			
4.9.1	Manual cleaning with immersion disinfection .....	27			
4.9.2	Manual cleaning with ultrasound and immersion disinfection ..	28			
4.10	Mechanical cleaning/disinfecting .....	29			
4.10.1	Mechanical alkaline cleaning and thermal disinfecting .....	29			
4.11	Mechanical cleaning/disinfection with manual pre-cleaning ..	30			
4.11.1	Manual pre-cleaning with a brush .....	30			
4.11.2	Manual pre-cleaning with ultrasound and brush .....	31			
4.11.3	Mechanical alkaline cleaning and thermal disinfecting .....	31			
4.12	Inspection .....	32			
4.12.1	Visual inspection .....	32			
4.12.2	Function test .....	32			

## 1. About this document

### Note

General risks of surgical interventions are not described in these instructions for use.

### 1.1 Scope

These instructions for use apply to all products listed in the Key.

### Note

The applicable CE mark for the product can be found on the label or packaging of the product.

► For article-specific instructions for use and information on material compatibility, see Aesculap Extranet at <https://extranet.bbraun.com>

### 1.2 Warnings

Warnings indicate hazards to patients, users and/or products that may emerge during the use of the product. Warnings are marked as follows:

#### DANGER

Indicates a potential danger. If it is not prevented, death or extremely severe injuries may result.

#### WARNING

Indicates a potential danger. If it is not prevented, minor or moderate injuries may result.

#### CAUTION

Indicates a potential damage to property. If these are not prevented, the product may be damaged.

## 2. Clinical application

### 2.1 Areas of application and restrictions on use

#### 2.1.1 Intended purpose

Ennovate instruments are used for the implantation of Ennovate Spinal System implants.

#### 2.1.2 Indications

##### Note

*Use of the product contrary to the indications mentioned above and/or described applications is outside the manufacturer's liability.*

Indications are described in the instructions for use of Ennovate implants TA014887.

#### 2.1.3 Absolute contraindications

Absolute contraindications are described in the instructions for use of Ennovate implants TA014887.

#### 2.1.4 Relative contraindications

Relative contraindications are described in the instructions for use of Ennovate implants TA014887.

## 2.2 Safety instructions

### 2.2.1 Clinical user

#### General safety information

To prevent damage caused by improper setup or operation, and to not compromise the manufacturer warranty and liability:

- ▶ Use the product only according to these instructions for use.
- ▶ Follow the safety and maintenance instructions.
- ▶ Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge and experience.
- ▶ Store any new or unused products in a dry, clean, and safe place.
- ▶ Prior to use, check that the product is in good working order.
- ▶ Keep the instructions for use accessible for the user.

##### Note

*The user is obligated to report all severe events in connection with the product to the manufacturer and the responsible authorities of the state in which the user is located.*

#### Notes on surgical procedures

It is the user's responsibility to ensure that the surgical procedure is performed correctly.

Appropriate clinical training as well as a theoretical and practical proficiency of all the required operating techniques, including the use of this product, are prerequisites for the successful use of this product.

The user is required to obtain information from the manufacturer if there is an unclear preoperative situation regarding the use of the product.

### 2.2.2 Product

#### Product specific safety information

- ▶ See usage instructions for the Ennovate handles TA015007 and Ennovate Implants TA014887 and operating manual O00702.
- ▶ Use only Aesculap Ennovate instruments and accessories.
- ▶ Hints on inserting the instruments into the supports: Observe the graphic/packing template.

### 2.2.3 Sterility

#### Non-sterile packaged products

The product will be delivered non-sterile.

- ▶ Clean brand-new product after removal of the transport packaging and before the first sterilization.

#### Sterile packaged products

Art. no.	Designation
BB540	Carbon steel scalpel blades
SR138SU	Ennovate injection cannula MIS
SZ364SU	Ennovate MIS bone access needle
SZ369SU	Ennovate MIS K-wire SS, sterile
SZ370SU	Ennovate MIS K-wire nitinol, sterile

The product has been sterilized using irradiation or ethylene oxide and is supplied in sterile packaging.

- ▶ Do not use products from open or damaged sterile packaging.
- ▶ Do not use product after expiration date.

Processing the product will impair its functionality. Contamination and/or impaired function of the products may result in injury or illness and death.

- ▶ Do not process product.

#### Single use products

Art. no.	Designation
BB540	Carbon steel scalpel blades
SR138SU	Ennovate injection cannula MIS
SZ364SU	Ennovate MIS bone access needle
SZ369SU	Ennovate MIS K-wire SS, sterile
SZ369S	Ennovate MIS K-wire SS, unsterile
SZ370SU	Ennovate MIS K-wire nitinol, sterile
SZ370	Ennovate MIS K-wire nitinol, unsterile

- ▶ Do not reuse product.

### 3. Safe operation

#### ⚠ DANGER

Injury of spinal cord, nerve roots and blood vessels due to piercing the cortex in the case of incorrect application!

- ▶ Alignment and insertion of the instruments, k-wires and pedicle screws may only be carried out under radiographic control or aided by a navigation system.
- ▶ Carefully guide instruments over the K-wires. While doing so, ensure that the K-wire remains in its position.
- ▶ Do not bend K-wires. Immediately replace bent K-wire.
- ▶ Inspect the cannulization of the instruments for passage before every application.

#### ⚠ DANGER

Injury of spinal cord, nerve roots and blood vessels due to bent K-wire or K-wire not holding in the bone!

- ▶ Do not reuse K-wires.
- ▶ Use K-wires only with intended K-wire forceps.
- ▶ Carefully drive in K-wires.
- ▶ Do not bend K-wires. Immediately replace bent K-wire.
- ▶ Insert K-wire in the correct direction, with the roughened end forward.
- ▶ Use only K-wires SZ369S, SZ369SU, SZ370 or SZ370SU of the Ennovate system.

#### ⚠ WARNING

Risk of injury and/or malfunction!

- ▶ Check product for loose, bent, broken, cracked, worn or missing parts before every application.
- ▶ Always carry out a function check prior to using the product.

#### ⚠ CAUTION

Instrument can fall in situ or on the floor!

Loosening of the handle during application!

- ▶ Inspect the seating of the coupling handles on the instruments.

#### Note

All actuation elements that are actuated during the operation are gold colored.

#### PolyLock®

PolyLock® can be used to temporarily block the polyaxiality of a polyaxial pedicle screw using the specialized instruments for the Ennovate system, and without the use of a rod or a locking screw. Then the polyaxial screw temporarily behaves like a monoaxial screw.

The procedure for activating and deactivating PolyLock® is described in the relevant chapters of these instructions for use.

PolyLock® can be activated and deactivated multiple times during instrumentation.

#### Note

PolyLock® does not necessarily have to be loosened before the final interlocking of the polyaxial screw using the rod and locking screw.

### 3.1 Place K-wire

#### 3.1.1 Place K-wire with bone access needle SZ364SU


Fig. 1

- ▶ Read and observe the enclosed instructions for use of the aspiration needle SZ364SU.
- ▶ Determine entry point and alignment of the pedicle screws.
- ▶ Perforate cortex with bone access needle 4.
- ▶ Insert the bone access needle 4 under radiographic control to the desired depth.
- ▶ Rotate top part a of the bone access needle 4 clockwise by 90° and remove from the sleeve b.


Fig. 2

- ▶ Insert K-wire 10 in sleeve b of the bone access needle 4 and insert to the desired depth with K-wire forceps 8 and slotted hammer 9 under radiographic control.
- ▶ After reaching the desired depth: Remove the sleeve b for the aspiration needle 4. Use slotted hammer 9 if necessary. In doing so, ensure that the position of the K-wire 10 does not change.

### 3.1.2 Place K-wire with K-wire targeting device (SZ365R and SZ366R)

- ▶ Determine entry point and alignment of the pedicle screws.


Fig. 3

- ▶ Push K-wire trocar **6** in K-wire sleeve **7** from above up to the stop and screw completely together clockwise using the gold colored sleeve **a**. Ensure that the trocar and sleeve are securely connected with the gold colored sleeve **a**. Together, the K-wire trocar **6** and K-wire sleeve **7** from the K-wire targeting device.

#### **⚠ DANGER**

**Injury of spinal cord, nerve roots and blood vessels due to piercing the cortex in the case of incorrect application!**

- ▶ Insert K-wire targeting device max. 30 mm deep.
- ▶ Perforate cortex with K-wire targeting device.
- ▶ Insert the K-wire targeting device under radiographic control to the desired depth.

#### Note

*In order that the hand does not enter the beam path when holding the targeting device, the PEEK rod **b** can be mounted on the trocar and the trocar held at a certain distance from the beam path.*

- ▶ Completely loosen the gold colored sleeve **a** on the K-wire trocar **6**. While doing so, brace on the handle of the K-wire trocar **6**.
- ▶ Remove K-wire trocar **6** from trocar sleeve **7**. Use slotted hammer **9** if necessary.


Fig. 4

- ▶ Insert K-wire **10** through K-wire sleeve **7** and insert to the desired depth with K-wire forceps **8** and slotted hammer **9** under radiographic control. When doing so, ensure that the K-wire is clamped in the recesses **c** on the mouth part of the K-wire forceps **8**.
- ▶ After reaching the desired depth: remove K-wire sleeve **7** using K-wire pliers **8** and slotted hammer **9**.
  - In doing so, ensure that the position of the K-wire **10** does not change.
  - To do so, position K-wire pliers **8** using the provided ledge **d** behind the articulation on the provided square **e** on the K-wire sleeve **7** beneath the fastening thread and clamp by actuating the handle parts.
- ▶ Place slotted hammer **9** on the clamped K-wire forceps **8**, knock out K-wire sleeve **7** and remove over the K-wire **10**.

### 3.2 Perform skin incision with skin incision guide SZ371R


Fig. 5

#### **⚠ WARNING**

**Risk of injury due to use of unsuitable scalpel blades!**

- ▶ Mount only scalpel blades BB540 on the skin incision guide.
- ▶ Mount two scalpel blades BB540 **a** on the skin incision guide **11**.
- ▶ Make skin incision matched to the diameter of the downtube: guide skin incision guide **11** over K-wire **10** and make skin incision. When doing so, ensure the correct depth of the skin incision.
- ▶ Remove skin incision guide **11** over K-wire **10**.

### 3.3 Perform dilation with dilator SZ372T to SZ374T and dilator handle SZ375T

Fasciae and muscles are bluntly dilated to expand access for inserting the screws.


Fig. 6

- ▶ Push dilators **12** to **14** successively over the K-wire **10**. While doing so, ensure that the K-wire is always held in its position.
- ▶ Place dilators directly on the bone in order to prevent soft tissue irritations.

#### Note

The dilator handle **15** can be optionally used to facilitate the insertion of the dilators. To do so, the dilators **12** to **14** are inserted into the corresponding opening **a** to **c** of the dilator handle **15**.

### 3.4 Enlarge pedicle canal with pedicle probe SZ376R

If necessary, the pedicle canal can be opened further with the pedicle probe.

#### ⚠ WARNING

Revision required in case of screw loosening or pedicle break!

- ▶ Use pedicle probe only for screws with  $\varnothing$  5.5 mm and larger.
- ▶ Do not use pedicle probe deeper than 80 mm. Coordinate use to determined screw length.
- ▶ Carefully insert the pedicle preparator. Avoid levering and bending during insertion.
- ▶ Carefully insert the pedicle preparator. Avoid levering and bending during insertion

#### Note

The pedicle probe is either guided through the medium dilator **13** or the large dilator **14** or used without dilators.

- ▶ Remove dilators **12** (and **13** and **14** as applicable) over the K-wire, see Fig. 6.


Fig. 7

- ▶ Guide pedicle probe **16** over the K-wire **10** and prepare to required depth. When using dilators, ensure that pedicle probe into the dilator **13/14**.
- ▶ Remove pedicle probe **16** from the site.

### 3.5 Cut thread with screw tap SZ254R to SZ260R (optional)

The Ennovate pedicle screws are self-cutting, with a continuous threaded tip. No thread cutting is usually required as a result. screw taps are optionally available in various diameters, however.

The screw taps have a smaller diameter than the corresponding screw. The color of the screw tap corresponds to the color of the corresponding screw.

#### ⚠ WARNING

Revision required in case of screw loosening or pedicle break!

Damage or breakage of instruments!

- ▶ Hold screw tap axially to the pre-prepared pedicle canal and to the K-wire when turning in and out. Avoid levering, bending and strong applications of force when doing so.
- ▶ Coordinate selection of screw diameter to the pedicle diameter. Use only screw tap in accordance with the screw diameter.

#### Note

In the case of large screw diameters, the thread must be cut to the nominal diameter by gradually raising the screw tap.

#### Note

The handles have a rotation direction setting in which the three positions of screw in (IN), locked (:), and screw out (OUT) can be set.


Fig. 8

- ▶ Push preferred handle **a** (SZ222R to SZ225R) onto the hex connection **c** of the corresponding screw tap **3** up to the stop.
- ▶ Set rotation direction setting on the handle to "IN".
- ▶ Guide screw tap **3** over the K-wire **10** and cut threads to required depth by rotating clockwise.
- ▶ Set rotation direction setting on the handle to "OUT".
- ▶ Unscrew screw tap **3** counter-clockwise and remove over the K-wire **10**.
- ▶ Loosen screw tap **3** from the handle **a**. To do so, push back the black sleeve **b** on the handle and remove screw tap **3**.

### 3.6 Determine screw length with screw length measuring device SZ377R

#### Note

The displayed screw length is approx. 2 mm shorter than the K-wire in the screw length measuring device. This increases the security that the screw cannot be screwed in further than the K-wire, injuring underlying structures. The readout is only an approximate value, which is dependent on the insertion depth of the K-wire. The final selection of screw length must be made by the operator.

- ▶ Make sure that the position of the K-wire **10** approximately corresponds with the final screw depth.


Fig. 9

- ▶ Guide screw length measuring device **17** over the K-wire **10** until it lies on the vertebral body. When doing so, ensure the correct orientation of the screw length measuring device.
- ▶ Read screw length between the two widest markings **a** on the K-wire **10**.

### 3.7 Couple downtube SZ378R/SZ405R to pedicle screw

#### ⚠ CAUTION

Operation delay!

Damage to the pedicle screw, the set screw and the downtube from loosening of the screw interlock from the downtube!

- ▶ Check that the pedicle screw is solidly fastened.
- ▶ Visually inspect correct seating of the support arms in the screw head.
- ▶ Tighten gold colored sleeve on downtube only by hand or with locking wrench. Do not use additional any additional tools.


Fig. 10

- ▶ Pull gold color sleeve **a** of the downtube **18** back against the spring mechanism up to the stop so that the retaining lugs **b** are visible.
- ▶ Push the downtube **18** with the retaining lugs **b** into the interface **c** on the screw head of the pedicle screw **d** in this position until an acoustic signal can be heard. This signals the latching of the downtube into the screw head.
- ▶ Visually inspect the connection between the downtube **18** and the pedicle screw **d**.
- ▶ Release the gold colored sleeve **a**. The downtube **18** is moved forward over the pedicle screw **d** via the spring mechanism.
- ▶ Fasten pedicle screw **d** to the downtube: tighten gold colored sleeve **a** counter-clockwise (left-hand thread) to the stop and simultaneously hold the downtube **18** securely against it. Fastening is complete when the marking (line and lock) is visible.

#### Note

If necessary, the locking wrench for downtube **19** can be pushed over the downtube in order to additionally tighten the connection and secure it against loosening.

### 3.8 Screw in pedicle screw

#### ⚠ DANGER

Injury to spinal cord, nerve root and blood vessels due to piercing of the cortex and faulty screw placement or inspection due to screw loosening when used improperly!

- ▶ Ensure axial alignment of the pedicle screw to the K-wire when screwing in. While doing so, ensure that the K-wire remains in its position.
- ▶ Remove K-wire before the screw is completely screwed in.
- ▶ Do not remove K-wire too early, as otherwise the pedicle screw guide will fail.

#### 3.8.1 Screw in pedicle screw with screwdriver SZ381R

To screw in the Ennovate pedicle screw, it is connected with the screwdriver. To do so, the pedicle screw must already be fastened completely in the downtube.

#### ⚠ CAUTION

Operation delay!

Damage to the pedicle screw, the set screw and by loosening the screwdriver from the screw!

Torque cannot be transferred!

- ▶ Completely connect the screwdriver with the screw so that the marking line on the screwdriver terminates with the end of the downtube.
- ▶ Check that the pedicle screw is securely fastened to the screwdriver.
- ▶ Verify the axial alignment of the pedicle screw to the screwdriver.
- ▶ Do not hold the gold colored sleeve of the screwdriver when screwing in the pedicle screw. Hold the downtube for guidance.


Fig. 11

- ▶ Insert screwdriver **21** into the downtube **18** until the tip **a** of the screwdriver is located in the inner profile of the screw head **b** and thus axially aligned to the screw shaft.
- ▶ Turn the screwdriver **21** into the thread of the screw body by rotating the gold colored sleeve **c** of the screwdriver clockwise. Securely hold the gold colored sleeve **d** of the downtube **18** while doing so. Fastening is complete when the marking line **e** on the screwdriver **21** terminates with the end of the downtube **18** and the polyaxial alignment of the screw (when using polyaxial screws) is completely blocked.
- ▶ Push preferred handle **f** (SZ222R to SZ225R) onto the hex connection of the screwdriver **21** up to the stop.
- ▶ Set rotation direction setting on the handle to "IN".
- ▶ Carefully guide the connected instruments over the K-wire **10** and screw in pedicle screw by rotating clockwise several rotations under radiographic control. Do not hold the screwdriver **21** by the gold colored sleeve **c** while doing so, as otherwise the screwdriver will release from the screw.
- ▶ Once the pedicle screw has sufficient hold in the bone, remove K-wire **10** to prevent further advancing of the K-wire. The K-wire cutter **8** can be used together with the slotted hammer **9** for this purpose.
- ▶ Insert the pedicle screw under radiographic control up to the desired depth.
  - When using polyaxial screws, ensure that the polyaxial alignment of the pedicle screw is ensured and the screw body does not strike against any anatomical structures.
  - When using monoaxial screws, ensure that the downtubes are aligned to each other so that the other instrumentation steps can be carried out.

#### ⚠ CAUTION

##### Operation delay!

Damage to the pedicle screw, the set screw and the downtube and from loosening of the screw interlock from the downtube!

- ▶ When loosening the screwdriver, do not hold the gold colored sleeve of the downtube, or do not release the gold colored sleeve of the downtube.
- ▶ After removing the screwdriver, verify solid fastening of the pedicle screw with the downtube and tighten as applicable.
- ▶ Release screwdriver **21** from the downtube **18**: release screwdriver **21** from the pedicle screw by turning the gold colored sleeve **c** counter-clockwise and remove from the downtube **18**.

- ▶ Remove handle **f** by actuating the black sleeve on the handle of the screwdriver **21**.
- ▶ Verify position and size of the pedicle screw before placing the rod with intraoperative imaging.

#### 3.8.2 Screw in fenestrated pedicle screw with injection cannula SR138SU and cannula sleeve SR139R and fix with bone cement

For the following steps, the fenestrated Ennovate pedicle screw must already be completely fixed in the downtube.

The injection cannula and the cannula sleeve are used to screw in the fenestrated pedicle screw, instead of the screwdriver. These are first connected to each other.

To increase anchoring stability, the fenestrated Ennovate pedicle screw can be fixed with bone cement.

#### ⚠ WARNING

Risk of injury due to material incompatibility!

- ▶ Only use application systems that have been approved for use with bone cement.

#### Screw in fenestrated pedicle screw

#### ⚠ WARNING

Risk of injury from cement leakage!

- ▶ Always apply imaging techniques for control during the preparation and positioning of the implants.
- ▶ Make certain that the slots of the fenestrated Ennovate pedicle screw are positioned completely in the bone.
- ▶ Take care not to perforate the pedicle cortex.
- ▶ Do not penetrate the opposite cortex and ensure that the opposite cortex is undamaged.
- ▶ Do not insert or cement the fenestrated Ennovate pedicle screw in fractured or partially resected vertebral bodies. Carry out preoperative thin-section CT control.
- ▶ Insert and screw the injection cannula and cannula sleeve into the screw head without exerting appreciable force. Make certain that the axes of the pedicle screw and the injection cannula are aligned with each other.
- ▶ Cement application must always be accompanied by high-quality, image-producing X-ray screening. Use no more than 2 ml of cement for each pedicle screw. Take the dead space volume of the injection cannula into consideration.

#### ⚠ CAUTION

Damage to the pedicle screw, the injection cannula, and delayed surgery due to loosening of the injection cannula from the pedicle screw and resulting replacement of the injection cannula!

Torque cannot be transferred!

- ▶ Completely connect the injection cannula and cannula sleeve to the pedicle screw so that the marking line on the cannula sleeve terminates with the end of the downtube.
- ▶ Check that the pedicle screw is solidly fastened to the injection cannula.
- ▶ Verify the axial alignment of the pedicle screw to the injection cannula.
- ▶ Do not hold the gold colored sleeve for the injection cannula when screwing in the pedicle screw. Hold the downtube for guidance.

### ⚠ CAUTION

**Implant damage caused by the application of excessive shear forces with the injection cannula!**

- ▶ Do not subject the injection cannula to transverse loads.
- ▶ Do not bend the injection cannula.


Fig. 12

- ▶ Connect injection cannula **38** to cannula sleeve **39**: insert injection cannula **38** at a slight incline through the tip of the cannula sleeve **39** and pivot into the sleeve at the recessed location.
- ▶ Insert connected instrument into the downtube **18** until the tip of the injection cannula **38** sits in the inner profile of the screw head and is thus axially aligned to the screw shaft.
- ▶ Screw in cannula sleeve **39** into the thread of the screw body by rotating the gold colored sleeve **a** clockwise. Securely hold the gold colored sleeve **c** of the downtube **18** while doing so.

Fastening is complete when the marking line **b** on the cannula sleeve **39** terminates with the end of the downtube **18** and the polyaxial alignment of the screw (when using polyaxial screws) is completely blocked.


Fig. 13

- ▶ Push preferred handle **d** (SZ222R to SZ225R) onto the hex connection of the injection cannula **38** up to the stop.
- ▶ Set rotation direction setting on the handle to "IN".

- ▶ Carefully guide the connected instruments over the K-wire **18/38/39** and screw in pedicle screw by rotating clockwise several rotations under radiographic control. When doing so, do not hold the sleeve for the injection cannula **39**, as otherwise the sleeve and injection cannula **38** will release from the pedicle screw.
- ▶ Once the pedicle screw has sufficient hold in the bone, remove K-wire **10** to prevent further advancing of the K-wire. The K-wire cutter **8** can be used together with the slotted hammer **9** for this purpose.
- ▶ Insert the pedicle screw under radiographic control up to the desired depth.
  - When using polyaxial screws, ensure that the polyaxial alignment of the pedicle screw is ensured and the screw body does not strike against any anatomical structures.
  - When using monoaxial screws, ensure that the downtubes are aligned to each other so that the other instrumentation steps can be carried out.
- ▶ Remove the handle **d** by retracting the black sleeve **e** of the injection cannula **38**.

### ⚠ CAUTION

**Damage to the pedicle screw, the set screw and the downtube and operation delay from loosening of the screw interlock from the downtube!**

- ▶ When loosening the sleeve for the injection cannula, do not hold the gold colored sleeve of the downtube, or do not release the gold colored sleeve of the downtube.
- ▶ After removing the injection cannula, verify solid fastening of the pedicle screw with the downtube and tighten as applicable.

### Cement fenestrated pedicle screw

#### ⚠ WARNING

**Revision and loss of anchoring stability by alignment (axial twisting) of the downtube after cementing!**

- ▶ Align rod slots/rod grooves of the downtube before cementing the pedicle screw.

#### ⚠ WARNING

**Wound contamination due to insufficient curing of cement!**

**Risk of set screw thread adhesion, resulting in system instability!**

- ▶ Observe the correct curing time.
- ▶ For monitoring the cement texture, prepare a reference volume of cement together with the cement to be applied.
- ▶ Make certain the implant system is not exposed to any load before the cement has reached its final strength.

#### ⚠ CAUTION

**Operation delay!**

**Uncontrolled escape of cement on Luer connection!**

- ▶ Ensure that the Luer lock adapter screw is fully screwed in on the injection cannula.

#### ⚠ CAUTION

**Inadequate fixation during removal of the application system due to incorrect viscosity of the cement!**

**Contamination of the interlock thread!**

- ▶ Make certain the cement has the correct viscosity.

#### Note

*Cementing takes place after all the screws have been aligned and the connection rods have been adapted (see Chapter 3.9), but before inserting the rods (see Chapter 3.12) and before locking the fenestrated pedicle screw.*

- ▶ Screw the Luer lock adapter **f** completely on the injection cannula **38**.
- ▶ Check that the pedicle screw is solidly fastened to the injection cannula **38**.
- ▶ Verify the position of the pedicle screw by a radiographic process.
- ▶ Prepare the cement. Follow the usage instructions for the cement and applicator used.
- ▶ Cement pedicle screw only with a cement suitable for screw augmentation.

#### Note

The dead space volume of the injection cannula equals 0.8 ml.

- ▶ Screw applicator into the Luer connection of the injection cannula **38** from above. Ensure that the connection is secure.
- ▶ Regularly verify with radiographic process during application of the bone cement.
- ▶ Carefully actuate applicator and apply cement under low pressure. Ensure the axial alignment of the applicator to the injection cannula **38**. Avoid levering and inclined placement.
- ▶ Apply cement until it escapes from the lateral slots of the pedicle screw. Avoid uncontrolled cement escape by continually observing under radiographic control.


Fig. 14

- ▶ Continue to apply cement until a cement plume is visible. If necessary, press the remaining cement in the cannula with the plunger **g**. To do so, push plunger **g** into the injection cannula **38** from the top and press down.
- ▶ Cement additional pedicle screws with the same procedure.
- ▶ Before removing the injection cannula **38**: allow the cement to cure completely.
- ▶ After the cement cures completely: remove injection cannula **38** and cannula sleeve **39** from the downtube **18**. To do so, release cannula sleeve **39** from the pedicle screw by turning the gold colored sleeve **a** counter-clockwise and remove from the downtube **18**, see Fig. 12.

### 3.9 Determine rod length with rod length caliper SZ382R

The necessary rod length is determined with the rod length caliper SZ382R.

#### Note

The determination of rod length is only an estimate. The final selection of rod length must be made by the operator.

#### ⚠ WARNING

Screw in or unscrew polyaxial screw by alignment (axial twisting) of the clamp sleeve after activation of PolyLock®!

- ▶ Do not axially twist clamp sleeve if the PolyLock® is activated.
- ▶ Align rod slots/rod grooves for the clamp sleeve before activating the PolyLock®.

#### ⚠ WARNING

Screwing in or unscrewing of monoaxial screw by alignment (axial twisting) of the downtube!

- ▶ Align rod slots/rod grooves on the downtube under radiographic control.

#### ⚠ CAUTION

Rod length cannot be determined!

- ▶ Align downtubes parallel to each other for measurement of the rod length with the rod length caliper. When using monoaxial screws, the downtubes cannot be subsequently aligned parallel.
- ▶ Take any subsequent distraction into account when determining the rod length.


Fig. 15

- ▶ Guide rod length caliper **22** with both hollow cylinders **a** over both outer downtubes **18** (cranial and caudal downtube). Ensure that both downtubes are aligned parallel to each other.
- ▶ Read the rod length through the window **b** in the slide unit.

#### Note

The length indicated on the scale indicates the actual rod length. If the measured rod length is between two existing lengths, the two closest rod lengths are indicated by small arrows in the window **b**. The operator has to select the rod length.

Insertion tip and rod hex are already considered in the rod length specification.

- ▶ Remove rod length caliper **22** from the downtubes **18**.
- ▶ Align downtubes **18** with the grooves in the rod direction. To do so, set rod length caliper **22** inclined by 90° in recesses **c** of the downtubes **18**.

### 3.10 Adjust rod length with rod cutter FW206R

- ▶ If necessary, cut the rod to desired length with rod cutter FW206R to the proper length.

### 3.11 Bend rod with french rod bender SZ270R

If necessary, the rod can be bent with the french rod bender SZ270R. Three different bending radii can be set on the french rod bender.

#### Note

Alternatively, the french rod bender FW024R can also be used with titanium rods.

**⚠ WARNING**

Danger of injury from damage (breakage) of the rod if bent incorrectly!

Rod cannot be inserted into the downtube!

Rod cannot be fastened!

- ▶ Avoid severe rod bending in the region of the pedicle screw.
- ▶ Only bend rod with provided french rod bender.
- ▶ Do not bend rod multiple times at the same location.
- ▶ Do not bend rod back


Fig. 16

- ▶ Adjust desired bending radius on french rod bender 5: pull rotary knob a and set the desired radius by rotating.
- ▶ Place rod b between rotary knob a and the two rollers c.

**Note**

The rod b has an orientation line that facilitates the alignment of the rod during the bending process

- ▶ Bend rod b to the desired radius by actuating the handle d of the french rod bender 5.

**3.12 Insert rod with rod inserter SZ384R/SZ385R**

**⚠ WARNING**

Loss of the rod in situ or expansion of the skin incision required!

- ▶ Do not decouple rod inserter from the rod until the rod is completely fastened by at least one set screw.

**⚠ WARNING**

Revision required due to faulty fastening of the rod!

- ▶ Always place hex and point of the rod outside of the screw head.
- ▶ Select suitable rod lengths in consideration of any possible corrective maneuver to be performed (e.g. distraction).
- ▶ Verify rod positioning under radiographic control before fastening with the set screw.

**⚠ CAUTION**

Operation delay due to incomplete coupling of the rod on the rod inserter or incorrect alignment of the rod slot/rod grooves of the downtube!

- ▶ Insert rod up to the stop in the rod inserter and completely tighten fixation.
- ▶ Verify solid fixation of the rod in the rod inserter.
- ▶ Only tighten/loosen gold colored rotary knob on the rod inserter by hand. Do not use additional any additional tools.
- ▶ Align rod slots/rod grooves of the downtube before insertion of the rod.

**Note**

The MIS rods have a conical tip on one side for easier implantation through the soft tissue. On the opposite side of the rod, there is a hex that forms the solid interface to the rod inserter.

**Note**

When using the long rod inserter SZ385R, an additional distal incision can be required.


Fig. 17


Fig. 18

- ▶ Open the mouth part of the rod inserter **23/24** by turning the gold colored rotary knob **a** counter-clockwise (up to the line mark with SZ384R).
- ▶ Insert rod **b** with the hex into the mouth part **c** of the rod inserter **23/24** up to the stop.

#### Note

The rod can be aligned to the instrument with the orientation lines on the rod and on the rod inserter.

- ▶ Rotate the gold colored rotary knob **a** clockwise and tighten. Verify solid fixation of the rod in the rod inserter and tighten if necessary.
- ▶ Guide the clamped rod **b** with the rod inserter **23/24** through the longitudinal slot of the downtube **18**.
- ▶ Do not decouple rod inserter **23/24** from the rod until the rod **b** is completely fastened by at least one set screw, see Chapter 3.15.

### 3.13 Verify rod with rod indicator SZ387R

The rod indicator SZ387R is used in order to ensure that the rod was placed correctly by the downtubes.

#### ⚠ CAUTION

Faulty display of the rod indicator due to tissue or instrument within the downtube!

- ▶ Before inserting the rod indicator, ensure that no tissue is present in the downtube and that no other instruments are located in the downtube.


Fig. 19

- ▶ Insert rod indicator **25** from the top into the downtube **18**. Ensure that the rod indicator **25** lies flush with the large diameter on the head of the downtube.
- ▶ Read rod position on the window **a** of the rod indicator **25**.
  - If the "rod" label is visible, the rod was placed correctly. If the "no rod" label is visible, the rod is not in the downtube.
- ▶ If the rod is not located in the downtube: correct position of the rod. To do so, remove rod indicator **25** from the downtube and insert rod again, see Chapter 3.12.

### 3.14 Push rod into screw head

#### 3.14.1 Push rod into screw head with quick rod pusher SZ397R (option 1)

Once the rod has been placed through all downtubes, the rod can be pressed into the screw head with the quick rod pusher.


Fig. 20

- ▶ Insert quick rod pusher **35** from the top into the downtube **18**.

#### Note

The quick rod pusher can be pushed into the downtube in two positions, so that the handle points in the direction of the rod axis to the right or left of the downtube.

#### ⚠ WARNING

Operation delay!

Set screw cannot be screwed into screw head!

Rod cannot be fixated!

- ▶ Push quick rod pusher down completely until the line marking terminates flush with the head of the downtube.
- ▶ If the pressure force by hand is insufficient: use rod pusher SZ389R.
- ▶ Push rod **a** into the screw head: push quick rod pusher **35** down on the handle until the line marking **b** terminates flush with the head of the downtube **18**.
  - Rod **a** is located completely in the screw head.
- ▶ Continue, see Chapter 3.15.

#### 3.14.2 Push rod into screw head with rod pusher SZ389R (option 2)

Once the rod has been positioned through all clamp sleeves, the rod can be pressed into the screw head using the rod pusher.

#### ⚠ DANGER

Risk of severe injury!

Pedicle screw can tear out of the spinal column if too much force is applied when pushing the rod!

- ▶ Make sure that no excessive force is applied when turning the gold-colored sleeve for the rod pusher to the marking "0". When the marking "0" is reached, the forces are no longer transmitted to the screw, but rather converted into the locking mechanism of the polyaxiality (PolyLock®).

**CAUTION**

Operation delay!

Set screw cannot be screwed into screw head!

Rod cannot be fixated!

- ▶ Turn rod pusher all the way down until the "0" marking on the gold-colored sleeve of the rod persuader instrument terminates flush with the head of the clamp sleeve.

**Note**

The rod pusher handle 28 can be used for easier rotation.


Fig. 21

- ▶ Guide rod pusher 27 from the top into the downtube 18 and turn golden sleeve a of the rod pusher by hand or with rod pusher handle 28 into the downtube 18 clockwise.

Rod b is pushed down. The scale c on the rod pusher indicates the remaining reduction path.

**Note**

The counter torque handle can be placed on the downtube in eight different positions.

- ▶ If the forces are high when screwing in, use counter torque handle 30: push counter torque handle 30 from the top on the head of the downtube 18 up to the stop. Hold the bracing handle 30 securely in position while screwing in.

The rod b is completely in the screw head when the scale on the rod pusher reaches the "0" marking.

- ▶ Repeat procedure for additional pedicle screws. Rod b is located completely in the screw head.
- ▶ Continue, see Chapter 3.15.

### 3.15 Insert and tighten set screw with the insertion instrument for the set screw SZ391R

**WARNING**

Insufficient stabilization of the spinal column region to be treated due to undertightening the set screw!

Risk of damage to the implant or instruments due to overtightening of the set screw!

- ▶ For the final tightening of the locking screws and when activating PolyLock®, use only 10 Nm torque handle SZ228R.
- ▶ Check the recalibration/maintenance date on the torque wrench handle maintenance label prior to each use.
- ▶ When tightening or loosening the set screws, make certain that the working end of the torque wrench handle shaft is positioned entirely within the socket of the set screws.
- ▶ Apply the torque slowly and without sudden movement until the torque wrench handle release is audibly triggered.

**WARNING**

Injury to patient from initiation of high torques on the spinal column!

- ▶ Always use the counter holder for tightening the locking screw and when activating PolyLock®.
- ▶ Hold brace in position and do not twist against the downtube or the inserted rod.

**CAUTION**

Operation delay!

Damage to the implant and/or instruments!

- ▶ Do not allow set screw driver to fall into the downtube.
- ▶ Receive set screw correctly and completely up to the line marking with the insertion instrument for the set screw.
- ▶ Insert set screw driver axially in the rod pusher. Ensure that the insertion instrument for the set screw does not tilt.

#### 3.15.1 Take up set screw


Fig. 22

- ▶ Push torque wrench handle SZ228R a onto the hex connection b of the insertion instrument for the set screw 29 up to the stop.
- ▶ Completely take up set screw SY001T/SY001TS c with the insertion instrument for the set screw 29. To do so, insert the working end of the insertion instrument for the set screw 29 into the socket of the set screw c and push with light pressure until the line marking d terminates flush with the surface of the set screw c.
- ▶ Verify secure seat of the locking screw c on the insertion instrument for the locking screw 29.

### 3.15.2 Insert and tighten set screw with MIS quick rod pusher SZ397R (option 1)


Fig. 23

- ▶ Insert set screw driver 29 with mounted set screw into the cannulization of the quick rod pusher 35 and carefully guide down onto the screw head.
- ▶ Screw in the set screw by turning the torque wrench handle SZ228R a clockwise.  
Hold the handle of the quick rod pusher 35 securely in position while screwing in until the acoustic signal of the triggering torque wrench handle sounds. The acoustic signal indicates that 10 Nm has been reached.
- ▶ Remove set screw driver 29 from the quick rod pusher 35.
- ▶ Repeat procedure for additional screws.
- ▶ Decouple rod inserter 23/24 from the rod d by turning the golden rotary knob c counter-clockwise and remove.

### 3.15.3 Insert and tighten set screw with rod pusher SZ389R (option 2)

#### Note

The counter torque handle can be placed on the downtube in eight different positions.


Fig. 24

- ▶ Push counter torque handle 30 from the top on the head of the downtube 18 up to the stop.
- ▶ Insert set screw driver 29 with mounted set screw into the cannulization of the rod pusher 27 and carefully guide down onto the screw head.

- ▶ Screw in the set screw by turning the torque wrench handle SZ228R a clockwise.  
Hold the counter torque handle 30 securely in position while screwing in until the acoustic signal of the triggering counter torque handle 30 sounds. The acoustic signal indicates that 10 Nm has been reached.
- ▶ Remove set screw driver 29 from the rod pusher 27.
- ▶ Repeat procedure for additional screws.
- ▶ Remove rod pusher 27 from the downtube 18 by turning the gold colored sleeve b counter-clockwise.
- ▶ Decouple rod inserter 23/24 from the rod d by turning the golden rotary knob c counter-clockwise and remove.

### 3.16 Remove downtube SZ378R/SZ405R with removal key SZ380R

#### ⚠ WARNING

Operation delay!

Downtube cannot be removed!

- ▶ Loosen gold colored sleeve of the downtube one to two clockwise rotations before inserting the removal key.
- ▶ Ensure that all instruments have been removed from the downtube 18.


Fig. 25

- ▶ Loosen gold colored sleeve a on downtube 18 one to two clockwise revolutions until the line marking is no longer visible.

#### Note

If necessary, the locking wrench for downtube 19 can be used to make loosening easier.

- ▶ Insert removal key 20 from the top into the downtube 18 up to the stop so that the line marking b terminates flush with the head of the downtube 18.
- ▶ Turn removal key 20 clockwise by 90° towards the clamp sleeve and remove from site in conjunction with the clamp sleeve 18 by pulling on the removal key 20.
- ▶ Repeat procedure for additional downtubes.

### 3.17 Carry out corrective measure (optional)

#### ⚠ WARNING

Insufficient treatment or stabilization of the spinal column region to be treated!

Damage to the implant and/or instruments!

- ▶ For final tightening of the locking screws and for activating PolyLock® use only 10 Nm torque handle SZ228R.
- ▶ Before every corrective maneuver: verify correct seating of the attachments on the parallel distractor forceps or on the parallel compressor forceps and to the downtube.

#### ⚠ WARNING

Injury to patient from initiation of high torques on the spinal column!

- ▶ When tightening the locking screw and when activating the PolyLock®, counteract any movement by tightly gripping the parallel distractor handle or the parallel compressor handle.
- ▶ Hold the handles of the parallel distractor forceps or parallel compressor forceps in position and do not twist against the downtube or the inserted rod.

#### ⚠ CAUTION

Screw in or unscrew polyaxial screw by alignment (axial twisting) of the clamp sleeve after activation of PolyLock®!

- ▶ Do not axially twist clamp sleeve if the PolyLock® is activated.
- ▶ Always align rod slots/rod grooves of the clamp sleeve before PolyLock® is activated.

#### 3.17.1 Prepare distraction/compression

If necessary, the distraction or compression can be performed with the parallel distractor forceps 1 or with the parallel compressor forceps 2 and the parallel sleeve 26 or C-ring attachments for the compression/distraction device 33.

Alternatively, the french rod bender compression/distraction device 32 can also be used for distraction.

- ▶ Perform operation steps up to pushing the rod, see Chapter 3.14.
- ▶ Select starting point of distraction/compression.


Fig. 26

- ▶ Ensure that the rod pusher 27 is completely inserted into the downtube.
- ▶ Completely tighten set screw, see Chapter 3.15.3.  
The set screw is used as an anchor point for the distraction/compression.
- ▶ Loosen or leave loose adjacent set screw in order to be able to start the distraction/compression process. Ensure that the downtubes 18 are aligned parallel to each other.

#### 3.17.2 Perform parallel distraction/compression with parallel distractor forceps SZ231R and parallel compressor forceps SZ232R (option 1)

#### Note

Only straight rods can be used for the parallel distraction/compression.


Fig. 27

- ▶ Guide parallel sleeves 26 or C-rings for compression/distraction device 33 over heads of the downtubes 18 that are to be distracted/compressed. Ensure correct alignment with parallel sleeves 26 when doing so ("UP" text points up).
- ▶ Couple parallel distractor forceps 1 or parallel compressor forceps 2 to parallel sleeves 26 or C-rings for compression/distraction device 33.


Fig. 28

- ▶ For polyaxial screws: Activate PolyLock®:
  - Slide torque handle SZ228R a all the way onto the hex connection of the PolyBlocker 31.
  - Set torque handle ST228R a including mounted PolyBlocker 31 from above into the gold colored sleeve of the rod pusher 27 and turn the torque handle SZ228R a clockwise until the acoustic signal of the triggering torque handle is heard. Counteract any movement by tightly gripping the parallel distractor handle 1 or the parallel compressor handle 2. The acoustic signal indicates that the PolyLock® was activated.
  - Remove torque handle SZ228R a including mounted PolyBlocker 31 from the rod pusher 27.
- ▶ Perform distraction/compression: push handle parts of the parallel distractor forceps 1 or parallel compressor forceps 2 together.
- ▶ Secure distraction/compression with gold colored nut b on the parallel distractor forceps 1 or parallel compressor forceps 2.
- ▶ Push torque wrench handle SZ228R a onto the hex connection of the insertion instrument for the set screw 29 up to the stop.

- ▶ Guide insertion instrument for the set screw **29** through the cannulization of the rod pusher **27** under light pressure completely into the socket of the previously loosened set screw until the line marking terminates flush with the surface.
- ▶ Check for secure seating of the insertion instrument for the set screw **29**, see Chapter 3.15.1.
- ▶ Turn torque wrench handle **SZ228R a** clockwise until the acoustic signal of the triggering torque wrench handle sounds. Counter-hold parallel distractor forceps **1** or parallel compressor forceps **2** securely on the handles.

The acoustic signal indicates that 10 Nm has been reached.

- ▶ Relieve parallel distractor forceps **1** or parallel compressor forceps **2** by loosening the gold colored nut **b** and opening up the locking device **c**.
- ▶ Retighten both set screws as described above, while bracing with torque wrench handle **SZ228R a**.
- ▶ For polyaxial screws: Deactivate PolyLock®:
  - Slide torque handle **SZ228R a** all the way onto the hex connection of the PolyBlocker **31**.
  - Insert torque handle **SZ228R a** with mounted PolyBlocker **31** from the top into the gold colored sleeve of the rod pusher **27** and turn the torque handle **SZ228R a** counter-clockwise until the rod persuader instrument **27** can be unscrewed without significant resistance. Counteract any movement by tightly gripping the parallel distractor handle **1** or the parallel compressor handle **2**.
  - Remove torque handle **SZ228R a** including the mounted PolyBlocker **31** from the rod pusher **27**.
- ▶ Loosen rod pusher **27** by turning the gold colored sleeve clockwise and pulling out of the downtube **18**.
- ▶ Remove parallel distractor forceps **1** or parallel compressor forceps **2** from the downtubes **18**.
- ▶ If needed, use the process for

### 3.17.3 Perform distraction with compression/distraction device **SZ394T (option 2)**


Fig. 29

- ▶ Guide C-rings for compression/distraction device **33** over heads of the downtubes **18** that are to be distracted/compressed.
- ▶ Couple parallel compressor forceps **2** to C-rings for compression/distraction device **33**.


Fig. 30

- ▶ For polyaxial screws: Activate PolyLock®:
  - Slide torque handle **SZ228R a** all the way onto the hex connection of the PolyBlocker **31**.
  - Insert torque handle **SZ228R a** including mounted PolyBlocker **31** from above into the gold colored sleeve of the rod pusher **27** and turn the torque handle **SZ228R a** clockwise until the acoustic signal of the triggering torque handle is heard. Counteract any movement by tightly gripping the handles of the parallel distractor **2** in doing so. The acoustic signal indicates that the PolyLock® was activated.
  - Remove torque handle **SZ228R a** including mounted PolyBlocker **31** from the rod pusher **27**.
- ▶ Place compression/distraction device **32** between the two downtubes **18** and rotate 90° to form the rotation point. When doing this, ensure that the compression/distraction device **32** is placed as close to the body as possible.
- ▶ Perform distraction: push handle parts of the parallel compressor forceps **2** together.
- ▶ Secure distraction/compression with gold colored nut **b** on the parallel compressor forceps **2**.
- ▶ Push torque wrench handle **SZ228R a** onto the hex connection of the insertion instrument for the set screw **29** up to the stop.
- ▶ Guide insertion instrument for the set screw **29** through the cannulization of the rod pusher **27** under light pressure completely into the socket of the previously loosened set screw until the line marking terminates flush with the surface.
- ▶ Check for secure seating of the insertion instrument for the set screw **29**, see Chapter 3.15.1.
- ▶ Turn torque wrench handle **SZ228R a** clockwise until the acoustic signal of the triggering torque wrench handle sounds. Counter-hold parallel distractor forceps **2** securely on the handles when doing so. The acoustic signal indicates that 10 Nm has been reached.
- ▶ Relieve parallel compressor forceps **2** by loosening the gold colored nut **b** and opening up the locking device **c**.
- ▶ Remove compression/distraction device **32** between the two downtubes **18**.
- ▶ Retighten both set screws as described above, while bracing with torque wrench handle **SZ228R a**.

- ▶ For polyaxial screws: Deactivate PolyLock®:
  - Slide torque handle SZ228R a all the way onto the hex connection of the PolyBlocker 31.
  - Insert torque handle SZ228R a with mounted PolyBlocker 31 from the top into the gold colored sleeve of the rod pusher 27 and turn the torque handle SZ228R a counter-clockwise until the rod persuader instrument 27 can be unscrewed without significant resistance. Counteract any movement by tightly gripping the handles of the parallel distractor 2 in doing so.
  - Remove torque handle SZ228R a including the mounted PolyBlocker 31 from the rod pusher 27.
- ▶ Loosen rod pusher 27 by turning the gold colored sleeve clockwise and pulling out of the downtube 18.
- ▶ Remove parallel compressor forceps 2 from downtubes 18.
- ▶ If needed, activate the process for

### 3.18 Interbody fusion management: distraction/compression with alignment tool SZ399R and cage sleeve SZ396R (optional)

In order to create more space in the case of a cage implantation, the downtubes can be temporarily removed from the pedicle screws and reattached later. Furthermore, the intervertebral disc space can be distracted for inserting the cage and compressed as required after insertion. The alignment tool 37 is used for this.

#### ⚠ WARNING

**Insufficient treatment or stabilization of the spinal column region to be treated!**

**Damage to the implant and/or instruments!**

- ▶ Never use alignment tool without inner rod.
- ▶ To activate the PolyLock® with the alignment instrument, completely connect the clamp sleeve to the pedicle screw and apply counter-action force.
- ▶ To activate PolyLock® with the alignment instrument, use only the 10 Nm torque handle SZ228R.

#### ⚠ WARNING

**Injury to patient from initiation of high torques on the spinal column!**

- ▶ To activate the PolyLock® with the alignment instrument, completely connect clamp sleeve to the pedicle screw and apply counter-action force.
- ▶ Hold brace in position and do not twist against the downtube or the inserted rod.

#### ⚠ WARNING

**Screw in or unscrew polyaxial screw by twisting the alignment instrument in an axial direction after activation of PolyLock®!**

- ▶ Do not axially twist alignment instrument with activated PolyLock® without bracing using the mounted clamp sleeve.

#### ⚠ WARNING

**Insufficient treatment or stabilization of the spinal column region to be treated!**

**Damage to the implant and/or instruments!**

- ▶ Before every corrective maneuver verify correct seating of the attachments on the parallel distractor forceps or on the parallel compressor forceps.

#### ⚠ CAUTION

**Operation delay!**

**Damage to the implant and downtube!**

- ▶ Loosen gold colored sleeve of downtube one to two rotations clockwise before the outer sleeve of the alignment tool is twisted by 90° in order to loosen the downtube from the pedicle screw.

#### ⚠ CAUTION

**Operation delay!**

**Damage to the implant and/or instruments!**

- ▶ Do not allow alignment tool to fall into the downtube.
- ▶ Do not allow outer sleeve to fall over the hex shaft onto the implant.

#### 3.18.1 Remove downtube SZ378R/SZ405R

- ▶ Perform operation steps up to screwing in the pedicle screw, see Chapter 3.8.


Fig. 31

- ▶ Insert alignment instrument 37, consisting of outer sleeve a, hexagonal shaft b and inner rod c, completely mounted and correctly positioned, from above through the clamp sleeve 18 and set into the screw head.

The directional arrow on the working end of the alignment tool 37 indicates the correct insertion direction.

- ▶ Screw in hexagonal end of the alignment tool 37 clockwise by hand into the screw head in order to anchor the instrument in the screw head.

The instrument is fully screwed in when the upper marking line d on the outer sleeve a of the alignment tool 37 terminates flush with the end of the downtube 18.

#### Note

*The counter torque handle can be placed in eight different positions.*

- ▶ Push counter torque handle 30 from the top on the head of the downtube 18 up to the stop.
- ▶ Push torque wrench handle SZ228R e onto the hex connection of the alignment tool 37 up to the stop.
- ▶ Align downtubes 18 parallel to each other.
- ▶ Tighten the alignment tool 37 by turning the torque wrench handle SZ228R e clockwise. Hold the counter torque handle 30 securely in position until the acoustic signal of the triggering torque wrench handle sounds.

The acoustic signal indicates that 10 Nm have been reached and PolyLock® is activated.

**Note**

PolyLock® provides for the parallel distraction of the vertebral bodies in the following.


Fig. 32

- ▶ Loosen gold colored sleeve **f** on downtube **18** one to two clockwise revolutions until the line marking is no longer visible.

**Note**

If necessary, the locking wrench for downtube **19** can be used to make loosening easier.

- ▶ Pull outer sleeve **a** of the alignment tool **37** out of the downtube until the circumferential line **g** on the outer sleeve terminates flush with the end of the downtube **18**.
- ▶ Twist outer sleeve **a** of the alignment tool **37** 90° to the downtube and remove the two instruments from the site together by pulling on the outer sleeve.
- ▶ Repeat the steps for adjacent downtube **18**.

The two hex shafts **b** of the alignment tool **37** are now parallel to each other in the two adjacent polyaxially locked pedicle screws.

**3.18.2 Perform distraction/compression**

The distraction/compression can now be performed via the hex shafts of the alignment tool.


Fig. 33

- ▶ Guide one cage sleeve **34** over the two hex shafts **b** of the alignment tool **37**.
- ▶ To perform a distraction: couple parallel distractor forceps **1** to cage sleeve **34**.

- ▶ To perform a compression: couple parallel compressor forceps **2** to cage sleeve **34**.
- ▶ Perform desired distraction/compression by pushing handle parts of parallel distractor forceps **1** or parallel compressor forceps **2** together.

**Note**

The distraction/compression can additionally be secured by the gold colored nut **f** on the parallel distractor forceps **1** or parallel compressor forceps **2**.

The cages can now be inserted in accordance with the corresponding cage surgery technology.

**3.19 Recouple downtube after cage implantation****⚠ CAUTION**

Operation delay!

Difficult coupling of the downtube!

Damage to the instruments!

- ▶ Remove alignment working port before the downtube is completely engaged on the screw.

**⚠ CAUTION**

Operation delay due to faulty recoupling of downtube to the pedicle screw!

Damage to the pedicle screw, the set screw and the downtube!

- ▶ Check that the pedicle screw is solidly fastened.
- ▶ Visually inspect correct seating of the support arms in the screw head.
- ▶ Axially pull or push downtube as applicable until correct seating of the downtube to the pedicle screw is attained.
- ▶ Tighten gold colored sleeve on downtube only by hand or with locking wrench. Do not use additional any additional tools.


**3.19.1 Insert alignment working port SZ398P**

Fig. 34

- ▶ Push completely mounted alignment working port **36** over the hex shaft **b** of the alignment tool **37** into the existing skin incision of the previously removed downtube **18**. When doing so, ensure that the tissue is dilated far enough that there is free access to the pedicle screw.
- ▶ Guide the outer sleeve **a** of the alignment tool **37** over the fixed hex shaft **b** through the alignment working port **36** up to the stop on the pedicle screw.

The correct position of the outer sleeve **a** has been reached when it sits on the pedicle screw, rotationally fixed, and the circumferential marking line **c** is visible on the hex shaft **b**.

### 3.19.2 Recouple downtube SZ378R/SZ405R with alignment tool SZ380R

#### ⚠ WARNING

#### Operation delay!

#### Damage or deformation of the K-wire!

- ▶ Axially align screw and downtube.
- ▶ Completely loosen gold colored sleeve **d** of the downtube **18** clockwise (left-hand thread) so that it can be pulled back against the spring mechanism up to the stop.
- ▶ Release the gold colored sleeve **d**.
- ▶ Place downtube **18** from the top on the outer sleeve **a** of the alignment tool **37** and push down onto the pedicle screw.


Fig. 35

- ▶ Open alignment working port **36** by loosening the gold colored sleeve **e**.
- ▶ Remove gold colored sleeve **e** over the downtube **18**.
- ▶ Individually remove remaining black half shells **f** of the alignment working port **36** from the site.


Fig. 36

- ▶ Push downtube **18** completely down on the pedicle screw and let it click into place.
  - The acoustic signal indicates that the downtube **18** has clicked into place in the screw head. The downtube is correctly latched in place when the upper marking line **g** on the outer sleeve of the alignment tool **37** terminates flush with the head of the downtube **18**.
- ▶ If this is not the case, push or pull downtube **18** in the corresponding direction until the marking line terminates flush with the head of the downtube **18**.
- ▶ Check support arms in screw head for correct seating.
- ▶ Fasten pedicle screw to the downtube: tighten gold colored sleeve **d** counter-clockwise (left-hand thread) to the stop. Securely hold the head of the downtube **18** while doing so.
  - Fastening is complete when the marking **h** (line and lock) is visible.

#### Note

If necessary, the locking wrench for downtube **19** can be pushed over the downtube in order to additionally tighten the connection and secure it against loosening.

#### Note

The counter torque handle can be placed on the downtube in eight different positions.

- ▶ Push counter torque handle **30** from the top on the head of the downtube **18** up to the stop.
- ▶ Push torque wrench handle SZ228R **i** onto the hex connection of the alignment tool **37** up to the stop.
- ▶ Loosen the alignment tool **37** by turning the torque wrench handle SZ228R **i** counter-clockwise. Hold the counter torque handle **30** securely in position until the alignment tool **37** can be unscrewed without significant resistance.


Fig. 37

- ▶ Remove all instruments from the downtube **18**.
- ▶ Repeat the steps for adjacent downtube.

### 3.20 Recouple downtube SZ378R/SZ405R after accidental decoupling

If the pedicle screw is already screwed in and the downtube SZ378R/SZ405R was accidentally removed, the downtube can be reattached to the pedicle screw with the alignment tool SZ399R.

#### 3.20.1 Couple alignment tool SZ399R

- ▶ Push completely mounted alignment working port **36** into the existing skin incision of the previously removed downtube **18**. When doing so, ensure that the tissue is dilated far enough that there is free access to the pedicle screw.

#### Note

The dilators **12** to **14** can be used for support when inserting.

#### Note

For clarity, the alignment working port **36** is not shown in the following figures. For figure for inserting and removing the alignment working port, see Chapter 3.19.


Fig. 38

- ▶ Remove outer sleeve **a** from the alignment tool **37**.  
The alignment tool **37** now consists of the hex shaft **b** and inner rod **c**.
- ▶ Carefully guide inner rod **c** together with hex shaft **b** through the alignment working port **36** onto the pedicle screw **d**.
- ▶ Palpate screw head **d** with inner rod **c** and position correctly in the rod slot. Carefully push hex shaft downward against the spring and screw into the screw head slightly by rotating clockwise.
- ▶ Guide the outer sleeve **a** of the alignment tool **37** over the fixed hex shaft **b** up to the stop on the pedicle screw.  
The correct position of the outer sleeve **a** has been reached when it sits on the pedicle screw, rotationally fixed.
- ▶ Carefully push hex shaft **b** downward against the spring and screw into the screw head completely by rotating clockwise. Brace outer sleeve **a** while doing so.

#### 3.20.2 Recouple downtube SZ378R/SZ405R

#### ⚠ WARNING

#### Operation delay!

#### Damage or deformation of the K-wire!

- ▶ Axially align screw and downtube.
- ▶ Completely loosen gold colored sleeve **e** of the downtube **18** clockwise (left-hand thread) so that it can be pulled back against the spring mechanism up to the stop.
- ▶ Release the gold colored sleeve **c**.


Fig. 39

- ▶ Place downtube **18** from the top on the outer sleeve **a** of the alignment tool **37** and push down onto the pedicle screw.
- ▶ Open alignment working port **36** by loosening the gold colored sleeve.
- ▶ Remove gold colored sleeve of the alignment working port **36** over the downtube **18**.
- ▶ Individually remove remaining black half shells of the alignment working port **36** from the site.
- ▶ Push downtube **18** completely down on the pedicle screw and let it click into place.  
The acoustic signal indicates that the downtube **18** has clicked into place in the screw head. The downtube is correctly latched in place when the upper marking line **f** on the outer sleeve of the alignment tool **37** terminates flush with the head of the downtube **18**.
- ▶ If this is not the case, push or pull downtube **18** in the corresponding direction until the marking line terminates flush with the head of the downtube **18**.
- ▶ Check support arms in screw head for correct seating.
- ▶ Fasten pedicle screw to the downtube: tighten gold colored sleeve **e** counter-clockwise (left-hand thread) to the stop. Securely hold the head of the downtube **18** while doing so.  
Fastening is complete when the marking **g** (line and lock) is visible.

#### Note

If necessary, the locking wrench for downtube **19** can be pushed over the downtube in order to additionally tighten the connection and secure it against loosening.

- ▶ Loosen alignment tool **37** by rotating the hex shaft **b** counter-clockwise by hand or with chosen handle from SZ222R to SZ225R.
- ▶ Remove alignment tool **37** from the downtube **18**.

## 4. Validated reprocessing procedure

### 4.1 General safety instructions

#### Note

Adhere to national statutory regulations, national and international standards and directives, and local, clinical hygiene instructions for sterile processing.

#### Note

For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of products.

#### Note

Mechanical reprocessing should be favored over manual cleaning as it gives better and more reliable results.

#### Note

Successful processing of this medical device can only be ensured if the processing method is first validated. The operator/sterile processing technician is responsible for this.

The recommended chemistry was used for validation.

#### Note

If there is no final sterilization, then a virucidal disinfectant must be used.

#### Note

For up-to-date information about reprocessing and material compatibility, see also the Aesculap Extranet at <https://extranet.bbraun.com>

The validated steam sterilization procedure was carried out in the Aesculap sterile container system.

### 4.2 General information

Dried or affixed surgical residues can make cleaning more difficult or ineffective and lead to corrosion. Therefore the time interval between application and processing should not exceed 6 h; also, neither fixating pre-cleaning temperatures >45 °C nor fixating disinfecting agents (active ingredient: aldehydes/alcohols) should be used.

Excessive measures of neutralizing agents or basic cleaners may result in a chemical attack and/or to fading and the laser marking becoming unreadable visually or by machine for stainless steel.

Residues containing chlorine or chlorides e.g. in surgical residues, medicines, saline solutions and in the service water used for cleaning, disinfection and sterilization will cause corrosion damage (pitting, stress corrosion) and result in the destruction of stainless steel products. These must be removed by rinsing thoroughly with demineralized water and then drying.

Additional drying, if necessary.

Only process chemicals that have been tested and approved (e.g. VAH or FDA approval or CE mark) and which are compatible with the product's materials according to the chemical manufacturers' recommendations may be used for processing the product. All the chemical manufacturer's application specifications must be strictly observed. Failure to do so can result in the following problems:

- Optical changes of materials, e.g. fading or discoloration of titanium or aluminum. For aluminum, the application/process solution only needs to be of pH >8 to cause visible surface changes.
- Material damage such as corrosion, cracks, fracturing, premature aging or swelling.
- ▶ Do not use metal cleaning brushes or other abrasives that would damage the product surfaces and could cause corrosion.
- ▶ Further detailed advice on hygienically safe and material-/value-preserving reprocessing can be found at [www.a-k-i.org](http://www.a-k-i.org), link to Publications, Red Brochure – Proper maintenance of instruments.

### 4.3 Single-use products

Art. no.	Designation
BB540	Carbon steel scalpel blades
SR138SU	Ennovate injection cannula MIS
SZ364SU	Ennovate MIS bone access needle
SZ369SU	Ennovate MIS K-wire SS, sterile
SZ370SU	Ennovate MIS K-wire nitinol, sterile

Processing the product will impair its functionality. Contamination and/or impaired function of the products may result in injury or illness and death.

- ▶ Do not process product.

Art. no.	Designation
BB540	Carbon steel scalpel blades
SR138SU	Ennovate injection cannula MIS
SZ364SU	Ennovate MIS bone access needle
SZ369SU	Ennovate MIS K-wire SS, sterile
SZ369S	Ennovate MIS K-wire SS, unsterile
SZ370SU	Ennovate MIS K-wire nitinol, sterile
SZ370	Ennovate MIS K-wire nitinol, unsterile

- ▶ Do not reuse product.

### 4.4 Reusable products

There are no known effects of processing resulting in damage to the product.

The manufacturer provided proof of biocompatibility and processability for the product after 200 processing cycles.

Careful visual and functional testing prior to next use is the best way to identify a malfunctioning product see Chapter 4.12.

### 4.5 Preparations at the place of use

- ▶ If applicable, rinse non-visible surfaces preferably with deionized water, with a disposable syringe for example.
- ▶ Remove any visible surgical residues to the extent possible with a damp, lint-free cloth.
- ▶ Transport the dry product in a sealed waste container for cleaning and disinfection within 6 hours.

### 4.6 Preparation before cleaning

- ▶ Dismantle the product prior to cleaning, see Chapter 4.7
- ▶ Open up instruments with hinges.

## 4.7 Disassembling

### 4.7.1 Parallel distractor forceps SZ231R and parallel compressor forceps SZ232R


Fig. 40

- ▶ Remove attachments **26**, **33**, **34**, or SZ411R from the parallel distractor forceps **1** or from the parallel compressor forceps **2**. To do so, push and hold gold colored lever **a** and remove attachment in the direction of the arrow.


Fig. 41

- ▶ To clean the parallel distractor forceps **1** and the parallel compressor forceps **2**, move the locking device **b** to the center position.

### 4.7.2 K-wire trocar SZ365R and K-wire sleeve SZ366R


Fig. 42

- ▶ Remove K-wire sleeve **7** from the K-wire trocar **6**. To do so, loosen gold colored sleeve **b** on the K-wire trocar **6** by turning in the direction of the arrow and remove K-wire sleeve **7**.
- ▶ Unscrew beige handle **c** on the K-wire trocar **6** by unscrewing in the direction of the arrow.

### 4.7.3 Skin incision guide SZ371R

#### Note

The scalpel blades BB540 **a** are single-use products and must be disposed of after use.

- ▶ Remove scalpel blades BB540 **a** from skin incision guide **11** and dispose.


Fig. 43

### 4.7.4 Downtube SZ378R/SZ405R

#### Note

The downtube **18** cannot be completely dismantled.


Fig. 44

- ▶ Unscrew gold colored rotation sleeve **a** of downtube **18** from the head **b** by turning in the direction of the arrow and push in the direction of the arrow up to the stop.

### 4.7.5 Screwdriver SZ381R


Fig. 45

- ▶ Press and hold push button **a** on screwdriver **21**.
- ▶ Pull out screwdriver shaft **c** from screwdriver sleeve **b**.

**4.7.6 Rod length measuring instrument for short construction SZ382R**


Fig. 46

- ▶ Push movable arm **a** in the direction of the arrow and remove from guide arm **b** of the rod length measuring instrument **22**.

**4.7.7 Rod inserter for short construction SZ384R**


Fig. 47

- ▶ Unscrew gold colored rotary knob **a** out of the rod inserter **23** by turning in the direction of the arrow and pull out.

**4.7.8 Rod inserter for long construction SZ385R**


Fig. 48

- ▶ Turn gold colored rotary knob **a** in the direction of the arrow and pull rod **b** out of shaft **c** of the rod inserter **24**.

**4.7.9 Alignment working port SZ398P**


Fig. 49

- ▶ Rotate gold colored ring **a** of the alignment working port **36** in the direction of the arrow and remove.
- ▶ Push half shells **b** apart.

**4.7.10 Alignment tool SZ399R**


Fig. 50

- ▶ Pull inner shaft **c** out of threaded shaft **b** of the alignment tool **37**, see Fig. 50.
- ▶ Pull threaded shaft **b** out of outer sleeve **a**.

**4.7.11 Injection cannula SR138SU and cannula sleeve SR139R**

*Note*

*The injection cannula **38** is a single-use product and must be disposed of after use.*


Fig. 51

- ▶ Pivot injection cannula **38** out of the cannula sleeve **39** and remove by pulling in the direction of the arrow.

## 4.8 Cleaning/disinfection

### 4.8.1 Product-specific safety notes on the reprocessing procedure

Damage or destruction of the product due to inappropriate cleaning agents/disinfectants and/or excessive temperatures!

- ▶ Use cleaning agents and disinfectants according to the manufacturer's instructions,
  - which are approved for plastics and high-grade steel,
  - that do not attack softeners (e.g., in silicone).
- ▶ Observe information concerning concentration, temperature and exposure time.
- ▶ Do not exceed a disinfection temperature of 95 °C.

- ▶ Manually pre-clean the product (with a cleansing brush) if there are residues of bone, tissue or ancillary materials (such as plaster or bone cement).
- ▶ Use suitable cleaning/disinfecting agents if the product is put away in a wet condition. To prevent foam formation and reduced effectiveness of the process chemicals: Prior to mechanical cleaning and disinfection, rinse the product thoroughly with running water.
- ▶ Clean and disinfect microsurgical products mechanically if they can be placed securely in the machine or on the positioning aids.

### 4.8.2 Validated cleaning and disinfection procedure

Validated procedure	Specific requirements	Reference
Manual cleaning with immersion disinfection SR139R, SZ254R-SZ260R, SZ270R, SZ365R-SZ368R, SZ369S, SZ370, SZ371R, SZ372T-SZ375T, SZ376R-SZ377R, SZ379R-SZ380R, SZ382R, SZ387R-SZ388R, SZ390R-SZ393R, SZ394T, SZ395R-SZ397R, SZ398P, SZ399R	<ul style="list-style-type: none"> <li>■ Suitable cleaning brush</li> <li>■ 20 ml disposable syringe</li> <li>■ Keep working ends open for cleaning.</li> <li>■ When cleaning products with movable hinges, ensure that these are in an open position and, if applicable, move the joint while cleaning.</li> <li>■ Drying phase: Use a lint-free cloth or medical compressed air</li> <li>■ SZ369S/SZ370: the product may not be reprocessed following contamination with blood (directly or indirectly).</li> </ul>	Chapter Manual cleaning/disinfection and sub-chapter: <ul style="list-style-type: none"> <li>■ Chapter Manual cleaning with immersion disinfection</li> </ul>
Manual cleaning with ultrasound and immersion disinfection SZ231R-SZ232R, SZ378R, SZ381R, SZ384R-SZ385R, SZ389R, SZ405R	<ul style="list-style-type: none"> <li>■ Suitable cleaning brush</li> <li>■ 20 ml disposable syringe</li> <li>■ Keep working ends open for cleaning.</li> <li>■ When cleaning products with movable hinges, ensure that these are in an open position and, if applicable, move the joint while cleaning.</li> <li>■ Drying phase: Use a lint-free cloth or medical compressed air</li> <li>■ SZ378R/SZ389R/SZ405R:               <ul style="list-style-type: none"> <li>– <b>Before phase I:</b> Clean the product with a suitable cleaning brush under running water until all discernible residues have been removed from the surface. Mobilize non-rigid components, such as set screws, links, etc. during cleaning.</li> </ul> </li> <li>■ SZ381R/SZ378R/SZ389R/SZ405R:               <ul style="list-style-type: none"> <li>– Rinse lumen with spray nozzle (3.8 bar water pressure) at least 3 times for 20 seconds each.</li> </ul> </li> </ul>	Chapter Manual cleaning/disinfection and sub-chapter: <ul style="list-style-type: none"> <li>■ Chapter Manual cleaning with ultrasound and immersion disinfection</li> </ul>
Mechanical alkaline cleaning and thermal disinfection SZ270R, SZ368R, SZ369S, SZ370, SZ375T, SZ379R, SZ382R, SZ387R, SZ390R-SZ392R, SZ394T, SZ398P	<ul style="list-style-type: none"> <li>■ Place the product in a tray that is suitable for cleaning (avoiding rinsing blind spots).</li> <li>■ Keep working ends open for cleaning.</li> <li>■ Place products in the tray with their hinges open.</li> <li>■ SZ369S/SZ370: the product may not be reprocessed following contamination with blood (directly or indirectly).</li> <li>■ SZ387R/SZ391R/SZ392R/SZ394T: Manually pre-clean the product using a cleansing brush if there are residues of bone/tissue or medical materials and supplies (e.g. casting compound or bone cement).</li> <li>■ SZ379R/SZ390R:               <ul style="list-style-type: none"> <li>– Rinse product under running water before cleaning.</li> <li>– Connect components with lumens and channels directly to the rinsing port of the injector carriage.</li> </ul> </li> </ul>	Chapter Mechanical cleaning/disinfecting and sub-chapter: <ul style="list-style-type: none"> <li>■ Chapter Mechanical alkaline cleaning and thermal disinfecting</li> </ul>

Validated procedure	Specific requirements	Reference
<p>Manual pre-cleaning with brush and subsequent mechanical alkaline cleaning and thermal disinfection</p> <p>SR139R, SZ254R-SZ260R, SZ365R-SZ367R, SZ371R, SZ372T-SZ374T, SZ376R-SZ377R, SZ380R, SZ384R-SZ385R, SZ388R, SZ393R, SZ395R-SZ396R, SZ397R, SZ399R</p>	<ul style="list-style-type: none"> <li>■ Suitable cleaning brush</li> <li>■ 20 ml disposable syringe</li> <li>■ Place the product in a tray that is suitable for cleaning (avoiding rinsing blind spots).</li> <li>■ Keep working ends open for cleaning.</li> <li>■ Place products in the tray with their hinges open.</li> <li>■ SZ384R/SZ385R: Manual pre-cleaning               <ul style="list-style-type: none"> <li>– <b>Phase II:</b> Rinse lumen with spray nozzle (3.8 bar water pressure) at least 3 times for 20 seconds each.</li> </ul> </li> <li>■ SZ254R-SZ260R/SZ366R/SZ371R/SZ372T-SZ374T/SZ376R/SZ377R/SZ384R/SZ385R/SZ388R/SZ393R/SZ396R/SZ397R/SZ399R:               <ul style="list-style-type: none"> <li>– Connect components with lumens and channels directly to the rinsing port of the injector carriage.</li> </ul> </li> </ul>	<p>Chapter Mechanical cleaning/disinfection with manual pre-cleaning and sub-chapter:</p> <ul style="list-style-type: none"> <li>■ Chapter Manual pre-cleaning with a brush</li> <li>■ Chapter Mechanical alkaline cleaning and thermal disinfecting</li> </ul>
<p>Manual pre-cleaning with ultrasound and brush, and subsequent mechanical alkaline cleaning and thermal disinfection</p> <p>SZ231R-SZ232R, SZ378R, SZ381R, SZ389R, SZ405R</p>	<ul style="list-style-type: none"> <li>■ Suitable cleaning brush</li> <li>■ 20 ml disposable syringe</li> <li>■ Place the product in a tray that is suitable for cleaning (avoiding rinsing blind spots).</li> <li>■ Keep working ends open for cleaning.</li> <li>■ Place products in the tray with their hinges open.</li> <li>■ SZ231R/SZ232R: Deviating manual pre-cleaning               <ul style="list-style-type: none"> <li>– <b>Phase I:</b> Place product in cold water for at least 5 min. and move movable parts twice.</li> <li>– <b>Phase II:</b> Clean the product with a suitable cleaning brush under running water until all discernible residues have been removed from the surface. Mobilize non-rigid components, such as set screws, links, etc. during cleaning.</li> <li>– <b>Phase III:</b> Clean the product in an ultrasound cleaning bath (frequency of 35 kHz) for at least 5 minutes. Ensure that all accessible surfaces are immersed and acoustic shadows are avoided.</li> </ul> </li> <li>■ SZ378R/SZ389R/SZ405R: Deviating manual pre-cleaning               <ul style="list-style-type: none"> <li>– <b>Phase I:</b> Clean the product using a suitable cleaning brush under running water until no residue is visible on the surface. If applicable, use a suitable cleaning brush to brush any surfaces that are not visible for at least 1 min. Move non-rigid components, such as adjusting screws, joints, etc., during cleaning.</li> <li>– <b>Phase II:</b> Clean the product in an ultrasound cleaning bath (frequency of 35 kHz) for at least 15 minutes. Ensure that all accessible surfaces are immersed and acoustic shadows are avoided.</li> <li>– <b>Phase III:</b> See Phase I</li> <li>– <b>Phase IV:</b> Rinse lumen with spray nozzle (3.8 bar water pressure) at least 3 times for 20 seconds each.</li> </ul> </li> <li>■ SZ381R: Rinse lumens using spray nozzle (3.8 bar water pressure) at least 3 times for 20 seconds each.</li> <li>■ SZ378R/SZ381R/SZ389R/SZ405R: Mechanical cleaning/disinfection               <ul style="list-style-type: none"> <li>– Connect the lumens and channels of individual parts directly to the injector unit's special flushing connector.</li> </ul> </li> </ul>	<p>Chapter Mechanical cleaning/disinfection with manual pre-cleaning and sub-chapter:</p> <ul style="list-style-type: none"> <li>■ Chapter Manual pre-cleaning with ultrasound and brush</li> <li>■ Chapter Mechanical alkaline cleaning and thermal disinfecting</li> </ul>

## 4.9 Manual cleaning/disinfection

- ▶ Prior to manual disinfecting, allow water to drip off for a sufficient length of time to prevent dilution of the disinfecting solution.
- ▶ After manual cleaning/disinfection, check visible surfaces visually for residues.
- ▶ Repeat the cleaning /disinfection process if necessary.

### 4.9.1 Manual cleaning with immersion disinfection

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	<b>Disinfecting cleaning</b>	RT (cold)	>15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
II	<b>Intermediate rinse</b>	RT (cold)	1	-	D-W	-
III	<b>Disinfection</b>	RT (cold)	15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
IV	<b>Final rinse</b>	RT (cold)	1	-	FD-W	-
V	<b>Drying</b>	RT	-	-	-	-

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

RT: Room temperature

\*Recommended: BBraun Stabimed

- ▶ Note the information on appropriate cleaning brushes and disposable syringes, see Chapter 4.8.2.

#### Phase I

- ▶ Fully immerse the product in the cleaning/disinfectant for at least 15 min. Ensure that all accessible surfaces are moistened.
- ▶ Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- ▶ If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
- ▶ Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- ▶ Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

#### Phase II

- ▶ Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- ▶ Drain any remaining water fully.

#### Phase III

- ▶ Fully immerse the product in the disinfectant solution.
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- ▶ Rinse lumens at least 5 times at the beginning of the exposure time using an appropriate disposable syringe. Ensure that all accessible surfaces are moistened.

#### Phase IV

- ▶ Rinse/flush the product thoroughly (all accessible surfaces).
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during final rinse.
- ▶ Rinse lumens with an appropriate disposable syringe at least five times.
- ▶ Drain any remaining water fully.

#### Phase V

- ▶ Dry the product in the drying phase with suitable equipment (e.g. cloth, compressed air), see Chapter 4.8.2.

#### 4.9.2 Manual cleaning with ultrasound and immersion disinfection

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	<b>Ultrasonic cleaning</b>	RT (cold)	>15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
II	<b>Intermediate rinse</b>	RT (cold)	1	-	D-W	-
III	<b>Disinfection</b>	RT (cold)	15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
IV	<b>Final rinse</b>	RT (cold)	1	-	FD-W	-
V	<b>Drying</b>	RT	-	-	-	-

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

RT: Room temperature

\*Recommended: BBraun Stabimed

- ▶ Note the information on appropriate cleaning brushes and disposable syringes, see Chapter 4.8.2.

##### Phase I

- ▶ Clean the product in an ultrasonic cleaning bath (frequency 35 kHz) for at least 15 min. Ensure that all accessible surfaces are immersed and acoustic shadows are avoided.
- ▶ Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- ▶ If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
- ▶ Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- ▶ Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

##### Phase II

- ▶ Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- ▶ Drain any remaining water fully.

##### Phase III

- ▶ Fully immerse the product in the disinfectant solution.
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- ▶ Rinse lumens at least five times at the beginning of the exposure time with an appropriate disposable syringe. Ensure that all accessible surfaces are moistened.

##### Phase IV

- ▶ Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during final rinse.
- ▶ Rinse lumens with an appropriate disposable syringe at least five times.
- ▶ Drain any remaining water fully.

##### Phase V

- ▶ Dry the product in the drying phase with suitable equipment (e.g. cloth, compressed air), see Chapter 4.8.2.

## 4.10 Mechanical cleaning/disinfecting

### Note

The cleaning and disinfection device must be of tested and approved effectiveness (e.g. FDA approval or CE mark according to DIN EN ISO 15883).

### Note

The cleaning and disinfection device used for processing must be serviced and checked at regular intervals.

### 4.10.1 Mechanical alkaline cleaning and thermal disinfecting

Machine type: single-chamber cleaning/disinfection device without ultrasound

Phase	Step	T [°C/°F]	t [min]	Water quality	Chemical/Note
I	Prerinse	<25/77	3	D-W	-
II	Cleaning	55/131	10	FD-W	<ul style="list-style-type: none"> <li>■ Concentrate, alkaline:               <ul style="list-style-type: none"> <li>- pH = 13</li> <li>- &lt;5 % anionic surfactant</li> </ul> </li> <li>■ 0.5 % working solution               <ul style="list-style-type: none"> <li>- pH = 11*</li> </ul> </li> </ul>
III	Intermediate rinse	>10/50	1	FD-W	-
IV	Thermal disinfecting	90/194	5	FD-W	-
V	Drying	-	-	-	According to the program for cleaning and disinfection device

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

\*Recommended: BBraun Helimatic Cleaner alkaline

- ▶ Check visible surfaces for residues after mechanical cleaning/disinfecting.

## 4.11 Mechanical cleaning/disinfection with manual pre-cleaning

### Note

The cleaning and disinfection device must be of tested and approved effectiveness (e.g. FDA approval or CE mark according to DIN EN ISO 15883).

### Note

The cleaning and disinfection device used for processing must be serviced and checked at regular intervals.

### 4.11.1 Manual pre-cleaning with a brush

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Disinfectant cleaning	RT (cold)	>15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
II	Rinsing	RT (cold)	1	-	D-W	-

D-W: Drinking water

RT: Room temperature

\*Recommended: BBraun Stabimed

- ▶ Note the information on appropriate cleaning brushes and disposable syringes, see Chapter 4.8.2.

#### Phase I

- ▶ Fully immerse the product in the cleaning/disinfectant for at least 15 min. Ensure that all accessible surfaces are moistened.
- ▶ Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- ▶ If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
- ▶ Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- ▶ Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

#### Phase II

- ▶ Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.

#### 4.11.2 Manual pre-cleaning with ultrasound and brush

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	<b>Ultrasonic cleaning</b>	RT (cold)	>15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
II	<b>Rinsing</b>	RT (cold)	1	-	D-W	-

D-W: Drinking water

RT: Room temperature

\*Recommended: BBraun Stabimed

- Note the information on appropriate cleaning brushes and disposable syringes, see Chapter 4.8.2.

##### Phase I

- Clean the product in an ultrasonic cleaning bath (frequency 35 kHz) for at least 15 min. Ensure that all accessible surfaces are immersed and acoustic shadows are avoided.
- Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
- Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

##### Phase II

- Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.

#### 4.11.3 Mechanical alkaline cleaning and thermal disinfecting

Machine type: single-chamber cleaning/disinfection device without ultrasound

Phase	Step	T [°C/°F]	t [min]	Water quality	Chemical
I	<b>Prerinse</b>	<25/77	3	D-W	-
II	<b>Cleaning</b>	55/131	10	FD-W	<ul style="list-style-type: none"> <li>■ Concentrate, alkaline:               <ul style="list-style-type: none"> <li>- pH = 13</li> <li>- &lt;5 % anionic surfactant</li> </ul> </li> <li>■ 0.5 % working solution               <ul style="list-style-type: none"> <li>- pH = 11*</li> </ul> </li> </ul>
III	<b>Intermediate rinse</b>	>10/50	1	FD-W	-
IV	<b>Thermal disinfecting</b>	90/194	5	FD-W	-
V	<b>Drying</b>	-	-	-	According to the program for cleaning and disinfection device

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

\*Recommended: BBraun Helimatic Cleaner alkaline

- Check visible surfaces for residues after mechanical cleaning/disinfecting.

## 4.12 Inspection

- ▶ Allow the product to cool down to room temperature.
- ▶ Dry the product if it is wet or damp.

### 4.12.1 Visual inspection

- ▶ Make sure all contamination has been removed. Pay particular attention to fitting surfaces, hinges, shanks, recesses, drilled grooves and the sides of teeth on rasps.
- ▶ For contaminated products: Repeat cleaning and disinfection procedure.
- ▶ Check instruments having interfaces to the implant for damage, burrs, and deformation at the corresponding functional geometries as these may damage the implants before use. This applies in particular to:
  - Torx geometry for instruments **21** and **29**
  - Contact geometries of the rod for instruments **5, 18, 23, 24, 25, 27** and **35**
  - Contact geometries for screw heads on instruments **18, 21, 27, 35, 37** and **39**
  - Thread geometries for connection to the female thread of the screw head (threads of the locking screw) for instruments **21, 37** and **39**
- ▶ Check product for damage, e.g. insulation, corroded, loose, bent, broken, cracked, worn, heavily scratched or missing parts.
- ▶ Check product for missing or faded labels and/or paint marks (e.g. color of handle, gold plating).
- ▶ Check products having long, narrow geometries (in particular rotating instruments) for deformations.
- ▶ Check cutting edges:
  - Sharpness
  - Continuous cutting edge (may not have any notches, dents and other damage)
- ▶ Check surfaces for variations in roughness.
- ▶ Check product for burrs that may damage tissue or surgical gloves.
- ▶ Check product for loose or missing parts.
- ▶ Check the cannulation of the instruments for continuity.
- ▶ Immediately separate out any damaged or worn products and return them to the Aesculap Technical Service, see Chapter 5.1.

### 4.12.2 Function test

#### CAUTION

**Damage (metal seizure/friction corrosion) to the product caused by insufficient lubrication!**

- ▶ Oil any moving parts (e. g. joints, spool parts and threaded rods) prior to the functional test using oil suitable for the applied sterilization process (e. g. STERILIT® I oil spray JG600 for steam sterilization or STERILIT® I drip-feed lubricator JG598).
- ▶ Assemble dismountable products, see Chapter 4.13

- ▶ Check any instruments having interfaces to other instruments for ease of movement and compatibility.
  - Check the secure coupling of the K-wire trocar **6** to the K-wire sleeve **7**, see Fig. 42. The thread must turn easily.
  - Check the connection of the K-wire cutters **8** to the K-wire sleeve **7**.
  - The rod pusher **27** can be fully and smoothly inserted into the clamp sleeve, see Fig. 21.
  - The quick rod pusher **35** can be fully inserted into the clamping sleeve, see Fig. 20.
  - Insertion instrument for locking screw **29** can be fully pushed into place by the quick rod pusher **35** and the rod pusher **27**.
  - Check the coupling of the counter holder **30** on the clamp sleeve **18**, see Fig. 21.
  - Check the connection of the PolyBlocker **31** to the rod pusher **27**, see Fig. 28.
  - Secure coupling incl. catch function of the sleeves. Check **26, 33** and **34** using the distraction/compression handles **1** and **2**, see Fig. 40.
  - Check connection of sleeves **26** and **33** to the clamp sleeve **18**, see Fig. 27/see Fig. 29.
  - Check the connection of the sleeve **34** to the alignment instrument **37**, see Fig. 31/see Fig. 33.
- ▶ Make sure that instruments having square or hexagonal connections are self-locking in Ennovate handles, see Fig. 8.
- ▶ Check clamp sleeve **18** for the criteria listed below, see Fig. 44:
  - Ease of operation of the thread
  - Can be assembled and disassembled
  - Moveability on the support arms
- ▶ Separate out clamp sleeve **18** if it is difficult to operate.
- ▶ Turning knobs on the rod insert instruments **23** and **24** can be easily twisted into and out of the respective instruments, see Fig. 47/see Fig. 48.
- ▶ Check that the product functions correctly.
- ▶ Check that all moving parts are working properly (e.g. hinges, locks/latches, sliding parts etc.).
- ▶ Check locks for secure latching function.
- ▶ Check for compatibility with associated products.
- ▶ Immediately single out the defective product and send it to the Aesculap Technical Service, see Chapter 5.1.

## 4.13 Assembling

### 4.13.1 Downtube SZ378R/SZ405R

- ▶ Push gold colored rotation sleeve **a** of downtube **18** opposite the direction of the arrow onto the head **b** and screw on by turning opposite the direction of the arrow until it can move freely in the direction of the head **b**, see Fig. 44.

### 4.13.2 Screwdriver SZ381R

- ▶ Press and hold push button **a** on screwdriver **21**.
- ▶ Push screwdriver shaft **c** opposite the direction of the arrow up to the stop in screwdriver sleeve **b** and latch in place, see Fig. 45.

### 4.13.3 Rod length measuring instrument for short construction SZ382R

- ▶ Push movable arm **a** opposite the direction of the arrow on guide arm **b** of the rod length measuring instrument **22**, see Fig. 46.

### 4.13.4 Rod inserter for short construction SZ384R

- ▶ Push gold colored rotary knob **a** opposite the direction of the arrow up to the stop in the shaft of the rod inserter **23** and screw in two to four rotations by turning opposite the direction of the arrow, see Fig. 47.

#### 4.13.5 Rod inserter for long construction SZ385R

- ▶ Push rod **b** opposite the direction of the arrow into the shaft **c** of the rod inserter **24** and turn back by turning the gold colored rotary knob **a** opposite the direction of the arrow up to the stop, see Fig. 48.

#### 4.13.6 Alignment working port SZ398P

- ▶ Push half shells **b** together, see Fig. 49.
- ▶ Rotate gold colored ring **a** of the alignment working port **36** opposite the direction of the arrow on half shells **b**.

#### 4.13.7 Alignment tool SZ399R

- ▶ Push inner shaft **c** opposite the direction of the arrow up to the stop in the threaded shaft **b** of the alignment tool **37**, see Fig. 50.
- ▶ Push threaded shaft **b** with inserted inner shaft **b** opposite the direction of the arrow up to the stop in the outer sleeve **a**.

### 4.14 Packaging

- ▶ Place the product in its holder or on a suitable tray. Ensure that all cutting edges are protected.
- ▶ Pack trays appropriately for the intended sterilization process (e.g. in Aesculap sterile containers).
- ▶ Ensure that the packaging provides sufficient protection against recontamination of the product during storage.

### 4.15 Steam sterilization

#### Note

SZ384R/SZ385R: The products can only be sterilized in the disassembled state.

- ▶ Check to ensure that the sterilizing agent will come into contact with all external and internal surfaces (e.g. by opening any valves and faucets).
- ▶ Validated sterilization process
  - Disassemble the product if necessary
  - Steam sterilization using fractional vacuum process
  - Steam sterilizer according to DIN EN 285 and validated according to DIN EN ISO 17665
  - Sterilization using fractionated vacuum process at 134 °C/holding time 5 min
- ▶ When sterilizing several products at the same time in a steam sterilizer, ensure that the maximum load capacity of the steam sterilizer specified by the manufacturer is not exceeded.

### 4.16 Storage

- ▶ Store sterile products in germ-proof packaging, protected from dust, in a dry, dark, temperature-controlled area.
- ▶ Store sterile packed single-use products dust-protected in a dry, dark and temperature-controlled room.

## 5. Maintenance and Service

### 5.1 Technical Service

#### ⚠ CAUTION

Modifications to medical equipment may result in loss of guarantee/warranty claims and any approvals.

- ▶ Do not modify the product.
- ▶ Contact national B. Braun/Aesculap representative for service and repair.

#### Service addresses

Aesculap Technischer Service

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78532 Tuttlingen / Germany

Phone: +49 7461 95-1601

Fax: +49 7461 16-2887

E-Mail: ats@aesculap.de

Other service addresses can be obtained from the address indicated above.

### 5.2 Accessories/Spare parts

Art. no.	Designation
BB540	Carbon steel scalpel blades

## 6. Disposal

#### ⚠ WARNING

Risk of infection from contaminated products!

- ▶ Observe national regulations when disposing or recycling the product, its components and their packaging.

#### ⚠ WARNING

Risk of injury from sharp and/or pointed products!

- ▶ When disposing of or recycling the product, make sure the packaging prevents injuries from the product.

#### Note

The operator has to process the product before disposal, see Chapter 4.

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