AESCULAP® Ennovate® Thoracolumbar

en Instructions for use/Technical description

Navigated Instruments

USA Note for U.S. users

This Instructions for Use is NOT intended for United States users. Please discard. The Instructions for Use for United States users can be obtained by visiting our website at www.aesculapimplantsystemsifus.com. If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-866-229-3002. A paper copy will be provided to you upon request at no additional cost.

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
Contents

1.	About this document	2
1.1	Scope	3
1.2	Safety messages	3
2.	Clinical use	3
2.1	Components required for use	3
2.2	Areas of use and limitations of use	3
2.2.1	Intended use	3
2.2.2	Indications	3
2.2.3	Contraindications	3
2.3	Safety information	4
2.3.1	Clinical user	4
2.3.2	Product	4
2.3.3	Sterility	4
3.	Application	4
3.1	Safe Application	4
3.2	Set-up Navigation	5
3.2.1	Mounting the navigation adapter on the instruments	5
3.2.2	Manual calibration of Ennovate Thoracolumbar instruments	5
3.3	Preparing the pedicle	5
3.3.1	For the OPEN technique	5
3.3.2	For the MIS technique	5
3.4	Using the drill guide	5
3.5	Cut thread with taps using navigation (optional)	6
3.6	Perform navigated screw insertion	6

4.	Validated reprocessing procedure	6
4.1	General safety instructions	6
4.2	General information	6
4.3	Single-use products	7
4.4	Reusable products	7
4.5	Preparations at the place of use	7
4.6	Preparing for cleaning	7
4.7	Disassembly	7
4.8	Cleaning/Disinfection	7
4.8.1	Product-specific safety information on the reprocessing method	7
4.8.2	Validated cleaning and disinfection procedure	8
4.9	Manual cleaning/disinfection	9
4.9.1	Manual cleaning with immersion disinfection	9
4.9.2	Manual cleaning with ultrasound and immersion disinfection	10
4.10	Mechanical cleaning/disinfection	11
4.10.1	Mechanical alkaline cleaning and thermal disinfecting	11
4.11	Mechanical cleaning/disinfection with manual pre-cleaning	12
4.11.1	Manual pre-cleaning with a brush	12
4.11.2	Manual pre-cleaning with ultrasound and brush	13
4.11.3	Mechanical alkaline cleaning and thermal disinfecting	13
4.12	Inspection	14
4.12.1	Visual inspection	14
4.12.2	Functional test	14
4.13	Assembly	14
4.14	Packaging	14
4.15	Steam sterilization	14
4.16	Storage	14
5.	Technical service	14
6.	Disposal	15

1. About this document

Note

General risk factors associated with surgical procedures are not described in these instructions for use.

1.1 Scope

This instruction for use is applied to the Ennovate Thoracolumbar Navigation Instruments listed in the table below and may only be used in combination with Brainlab Navigation software and hardware.

Art. no.	Designation	Sterility
SZ464R	Ennovate navigation pedicle awl	non-sterile
SZ465R	Ennovate navigation K-wire trocar	non-sterile
SZ466R	Ennovate navigation K-wire sleeve	non-sterile
SZ468R	Ennovate navigation lumbar pedicle probes straight	non-sterile
SZ469R	Ennovate navigation lumbar pedicle probes straight	non-sterile
SZ470R	Ennovate navigation cannulated pedicle probe	non-sterile
SZ471R	Ennovate navigation thoracic pedicle probe, straight	non-sterile
SZ472R	Ennovate navigation thoracic pedicle probe, straight	non-sterile
SZ474R	Ennovate navigation drill guide handle	non-sterile
SZ475R	Ennovate navigation drill guide tube, Ø 4.5 mm / 5.5 mm screws	non-sterile
SZ476R	Ennovate navigation drill Ø 4.5 mm / 5.5 mm screws	non-sterile
SZ476SU	Ennovate drill, Ø 4.5 mm / 5.5 mm screws, single-use	sterile
SZ477R	Ennovate navigation drill guide tube, Ø 6.5 mm / 7.5 mm screws	non-sterile
SZ478R	Ennovate drill Ø 6.5 mm / 7.5 mm screws, reusable	non-sterile
SZ478SU	Ennovate drill, Ø 6.5 mm / 7.5 mm screws, single-use	sterile
SZ479R	Ennovate navigation handle tightening key	non-sterile
SZ480R	Ennovate cortical reamers	non-sterile
SZ481R	Ennovate navigation screwdriver	non-sterile
SZ482R	Ennovate navigation ratchet handle, T-shape	non-sterile
SZ483R	Ennovate navigation ratchet handle, straight	non-sterile
SZ488R	Ennovate trocar for drill guide tube, Ø 4.5 mm / 5.5 mm screws	non-sterile
SZ489R	Ennovate trocar for drill guide tube, Ø 6.5 mm / 7.5 mm screws	non-sterile
SZ090R	Aesculap navigation adapter , 3 marker spheres	non-sterile
SZ092R	Aesculap navigation adapter , 3 marker spheres, size L	non-sterile

► For article-specific instructions for use as well as information on material compatibility and lifetime see B. Braun eIFU at eifu.bbraun.com

1.2 Safety messages

Safety messages make clear the dangers to patient, user and/or product that could arise during the use of the product. Safety messages are labeled as follows:

DANGER

Indicates a possible threat of danger. If not avoided, death or serious injury may result.

WARNING

Indicates a possible threat of danger. If not avoided, minor or moderate injury may result.

CAUTION

Indicates a possible threat of material damage. If not avoided, the product may be damaged.

2. Clinical use

2.1 Components required for use

All instruments listed in the legend are required for use. Additionally, instruments from the Ennovate Lumbar basic set, Ennovate Cervical navigated instruments as well as from the Brainlab instruments are required for use.

2.2 Areas of use and limitations of use

2.2.1 Intended use

Ennovate Thoracolumbar Navigation Instruments are used for navigated bone preparation workflows and the implantation of Ennovate PentaCore pedicle screws.

2.2.2 Indications

Note

The manufacturer is not responsible for any use of the product against the specified indications and/or the described applications.

The indications are described in the instruction for use of the Aesculap Ennovate Spinal System implants (TA014887).

2.2.3 Contraindications

The contraindications are described in the instruction for use of the Aesculap Ennovate Spinal System implants (TA014887).

2.3 Safety information

2.3.1 Clinical user

General safety information

To prevent damage caused by improper setup or operation, and to not compromise the manufacturer warranty and liability:

- ▶ Use the product only according to these instructions for use.
- ▶ Follow the safety and maintenance instructions.
- ▶ Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge and experience.
- ▶ Store any new or unused products in a dry, clean, and safe place.
- ▶ Prior to use, check that the product is in good working order.
- ▶ Keep the instructions for use accessible for the user.

Note

The user is obligated to report all severe events in connection with the product to the manufacturer and the responsible authorities of the state in which the user is located.

Notes on surgical procedures

It is the user's responsibility to ensure that the surgical procedure is performed correctly.

Appropriate clinical training as well as a theoretical and practical proficiency of all the required operating techniques, including the use of this product, are prerequisites for the successful use of this product.

The user is required to obtain information from the manufacturer if there is an unclear preoperative situation regarding the use of the product.

2.3.2 Product

Product-specific safety information

All the risks related to the handling of the instruments are described in detail in the respective instructions for use.

- ▶ Thoroughly read the instructions for use before performing navigated surgery.

These instructions for use apply to the Ennovate Thoracolumbar navigated instruments listed in the legend and may only be used with the Brainlab navigation system. The surgical technique and its related risks are described in detail in following documents:

Art. no.	Designation
TA014887	Aesculap Ennovate Spinal System implants
TA014986	IFU Ennovate instruments for open access
TA014987	IFU Ennovate instruments for MIS access
TA015762	IFU Ennovate Cervical navigated instruments
TA016113	IFU Ennovate navigation trays
003702	Surgical manual Ennovate Thoracolumbar navigation

For more information on the safe handling of Brainlab instruments as well as on the Navigation Software Spine & Trauma 3D, refer to the relevant Brainlab instrument and software user guides.

Visit <https://www.brainlab.com> to access online user guides or to contact the Brainlab customer support.

Aesculap and Brainlab accept absolutely no responsibility if instruments other than those named below are used.

Follow the combinations described in this document.

- ▶ Only combine Aesculap products (and Brainlab marker spheres) with each other.
- ▶ Only use Aesculap Ennovate instruments and accessories.
- ▶ Observe the graphic/packing template when storing the instruments.

- ▶ Use the product only in applications according to its intended use, see Intended use.

Instruments can fall in situ or on the floor!

Handles can loosen during use!

- ▶ Check the seating of the coupling handles on the instruments for tightness.

Injury of spinal cord, nerve roots and blood vessels due to piercing of the cortex in case of incorrect application!

- ▶ Only carry out alignment and insertion of the instruments and Ennovate PentaCore pedicle screws under radiographic control or with the aid of the Brainlab navigation system.

2.3.3 Sterility

Non-sterile reusable products

The product is delivered in an unsterile condition, see Scope.

- ▶ Inspect the new product after removing its transport packaging and prior to first use to ensure it is in good working order.
- ▶ Clean the new product after removing its transport packaging and prior to its initial sterilization.

Sterile packed products

The product is sterilized by irradiation and is supplied in sterile packaging, see Scope.

- ▶ Do not use products from open or damaged sterile packaging.
- ▶ Do not use the product after its use-by date.
- ▶ Do not reuse the product.
- ▶ Store sterile products in germ-proof packaging, protected from dust, in a dry, dark, temperature-controlled area.

3. Application

3.1 Safe Application

⚠ WARNING

Risk of injury and/or malfunction!

- ▶ Prior to each use, inspect the product for loose, bent, broken, cracked, worn, or fractured components. Specifically check the navigation adapters, including pins, are not bent.
- ▶ Always carry out a function test prior to each use of the product.
- ▶ Do not use the product if it is damaged or defective. Set aside the product if it is damaged.

⚠ WARNING

Risk of injury to the patient due to interruption of the navigation!

- ▶ Prior to the operation, plan the configuration of the operating room, the assembly of the instruments and the alignment of the navigation adapter.
- ▶ Make sure that the navigation camera has an unrestricted view of the reflective marker spheres of the instruments.

⚠ WARNING

Inaccurate navigation due to bent instruments!

- ▶ Once the instrument (e.g. probe, drill) is inserted in the bone, do not correct the axis anymore, as this may cause the instrument to bend and the navigation to be inaccurate. If a correction of the axis is required, remove the instrument and insert it again.

3.2 Set-up Navigation

3.2.1 Mounting the navigation adapter on the instruments

- ▶ To ensure a safe mounting of the reflective Brainlab marker spheres onto Aesculap adapters 13, please refer to the Brainlab user guides.
- ▶ To ensure a safe assembly and combination of the Ennovate Navigation Instruments and Adapters, please refer to the surgical manual O03702.

3.2.2 Manual calibration of Ennovate Thoracolumbar instruments

⚠ WARNING

Risk of injury to the patient due to incorrect navigation!

- ▶ Make sure that the navigation adapter is properly fixed on the instruments.
- ▶ Make sure to perform a correct calibration of the instruments before each use.
- ▶ Make sure that each instrument is correctly assembled.

To ensure a safe calibration and navigation workflow, please refer to the Brainlab User Guides

3.3 Preparing the pedicle

The preparation of the pedicle depends on the used operation technique.

3.3.1 For the OPEN technique

- ▶ Use the pedicle awl (SZ464R) with adapter SZ090R.
- ▶ Use the appropriate pedicle probes (SZ468R, SZ469R, SZ471R and SZ472R) with adapter SZ092R or SZ096R. For the size, pay attention on the indications on the pedicle probes.

3.3.2 For the MIS technique

- ▶ Use the pedicle K-wire aim instrument consisting of the trocar SZ465R inserted in the sleeve SZ466R (see surgical manual O03702 for assembly) with the adapter SZ092R.
- ▶ Use the cannulated pedicle probe (SZ470R) with the adapter SZ090R only for screw diameters starting from 5.5 mm.

Note

Alternatively, navigated pedicle preparation can be performed with the navigated drill guide for both approaches (Open or MIS).

3.4 Using the drill guide

⚠ CAUTION

Delay of surgery due to incorrect hole size!

- ▶ Check the color coding and laser marking to ensure the correct drill tube, trocar and drill are used.


Fig. 1

- ▶ Mount the appropriate drill tube 6 to the drill guide handle 5 by turning the knob until the stop. Use the drill guide with the adapter SZ092R.
- ▶ Choose the appropriate sized trocar (SZ488R or SZ489R).


Fig. 2

- ▶ Set the depth sleeve 5a at the intended drilling depth.
- ▶ Introduce the trocar into the drill guide until stop.
- ▶ Position the drill guide together with the trocar on the desired entry point and pull back the trocar while pushing the drill guide down.
- ▶ Check that the depth stop of the guide is set at the intended drilling depth, select the corresponding drill (colour coding) and attach it to an appropriate motor system for spine surgeries.
- ▶ Leave the drill guide 5+6 at the desired entry site and advance the drill 7 carefully to the predefined depth under trajectory control with the Brainlab navigation system until the stop is reached.

3.5 Cut thread with taps using navigation (optional)

⚠ CAUTION

Operation delay!

- ▶ Always tighten handle prior to calibration.
- ▶ Make sure that the assembly is fixed and correctly assembled.
- ▶ Do not retighten the handle to the shaft after having performed the calibration.
- ▶ Mount the Aesculap navigation adapter 13 (SZ092R) on the navigation handle 10/11 (SZ482R or SZ483R), see surgical manual 003702.


Fig. 3

3.6 Perform navigated screw insertion

⚠ DANGER

Risk of injury to the patient due to inaccurate navigation!

- ▶ Always repeat the calibration procedure when changing screws (even if using screws of the same dimensions).
- ▶ Always use the Aesculap navigation adapter with four marker spheres (SZ091R) for calibration of the Ennovate screwdriver SZ481R.

⚠ CAUTION

Operation delay!

- ▶ Always tighten handle prior to calibration.
- ▶ Make sure that the assembly is fixed and correctly assembled.
- ▶ Do not retighten the handle to the shaft after having performed the calibration.

Fig. 4

- ▶ Insert the screw according to the instructions of the navigation system. Refer to TA014986 (open access) and TA014987 (MIS access) for a detailed description of screw insertion and the risks related to it.

4. Validated reprocessing procedure

4.1 General safety instructions

Note

Adhere to national statutory regulations, national and international standards and directives, and local, clinical hygiene instructions for sterile processing.

Note

For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of products.

Note

Mechanical reprocessing should be favored over manual cleaning as it gives better and more reliable results.

Note

Successful processing of this medical device can only be ensured if the processing method is first validated. The operator/sterile processing technician is responsible for this.

Note

If there is no final sterilization, then a virucidal disinfectant must be used.

Note

For up-to-date information about reprocessing and material compatibility, see B. Braun eIFU at eifu.bbraun.com

The validated steam sterilization procedure was carried out in the Aesculap sterile container system.

4.2 General information

Dried or affixed surgical residues can make cleaning more difficult or ineffective and lead to corrosion. Therefore the time interval between application and processing should not exceed 6 h; also, neither fixating pre-cleaning temperatures >45 °C nor fixating disinfecting agents (active ingredient: aldehydes/alcohols) should be used.

Excessive measures of neutralizing agents or basic cleaners may result in a chemical attack and/or to fading and the laser marking becoming unreadable visually or by machine for stainless steel.

Residues containing chlorine or chlorides e.g. in surgical residues, medicines, saline solutions and in the service water used for cleaning, disinfection and sterilization will cause corrosion damage (pitting, stress corrosion) and result in the destruction of stainless steel products. These must be removed by rinsing thoroughly with demineralized water and then drying.

Additional drying, if necessary.

Only process chemicals that have been tested and approved (e.g. VAH or FDA approval or CE mark) and which are compatible with the product's materials according to the chemical manufacturers' recommendations may be used for processing the product. All the chemical manufacturer's application specifications must be strictly observed. Failure to do so can result in the following problems:

- Optical changes of materials, e.g. fading or discoloration of titanium or aluminum. For aluminum, the application/process solution only needs to be of pH >8 to cause visible surface changes.
- Material damage such as corrosion, cracks, fracturing, premature aging or swelling.
- ▶ Do not use metal cleaning brushes or other abrasives that would damage the product surfaces and could cause corrosion.
- ▶ Further detailed advice on hygienically safe and material-/value-preserving reprocessing can be found at www.a-k-i.org, link to "AKI-Brochures", "Red brochure".

4.3 Single-use products

- ▶ Do not reuse the product, see Scope.

The reprocessing of the product affects its functionality. Risk of injury, illness or death due to soiling and/or impaired functionality of the product.

- ▶ Do not reprocess the product.

4.4 Reusable products

Influences of the reprocessing which lead to damage to the product are not known.

A careful visual and functional inspection before the next use is the best opportunity to recognize a product that is no longer functional, see Inspection.

4.5 Preparations at the place of use

- ▶ If applicable, rinse non-visible surfaces preferably with deionized water, with a disposable syringe for example.
- ▶ Remove any visible surgical residues to the extent possible with a damp, lint-free cloth.
- ▶ Transport the dry product in a sealed waste container for cleaning and disinfection within 6 hours.

4.6 Preparing for cleaning

- ▶ Do not disassemble fixation screws which permanently connect the various components.
- ▶ Disassemble the product prior to cleaning, see Disassembly.

4.7 Disassembly

SZ482R and SZ483R

Fig. 5

- ▶ Take off PEEK sleeve a.
- ▶ Unscrew pusher with the help of tightening key 14.

SZ474R


Fig. 6

- ▶ Release depth sleeve 5a and handle lock 5c from the handle 5b.

4.8 Cleaning/Disinfection

4.8.1 Product-specific safety information on the reprocessing method

Damage to or destruction of the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!

- ▶ Following the manufacturer's instructions, use cleaning and disinfecting agents
 - that are approved for plastics and high-grade steel,
 - that do not attack softeners (e.g., in silicone).
- ▶ Observe specifications regarding concentration, temperature and exposure time.
- ▶ Do not exceed the maximum allowable disinfection temperature of 95 °C.

4.8.2 Validated cleaning and disinfection procedure

Validated procedure	Specific requirements	Reference
<p>Manual cleaning with immersion disinfection</p> <ul style="list-style-type: none"> ■ SZ464R to SZ466R, SZ468R to SZ472R, SZ475R to SZ480R, SZ488R, SZ489R 	<ul style="list-style-type: none"> ■ Suitable cleaning brush ■ Disposable syringe 20 ml ■ Drying phase: Use a lint-free cloth or medical compressed air <p>SZ466R:</p> <ul style="list-style-type: none"> ■ Before phase I: Clean the product with a suitable cleaning brush under running water until all discernible residues have been removed from the surface. Mobilize non-rigid components, such as set screws, links, etc. during cleaning. ■ Rinse lumen with spray nozzle (3.8 bar water pressure) at least 3 times for 20 seconds each. 	<p>Chapter Manual cleaning/disinfection and subsection:</p> <ul style="list-style-type: none"> ■ Chapter Manual cleaning with immersion disinfection
<p>Manual cleaning with ultrasound and immersion disinfection</p> <ul style="list-style-type: none"> ■ SZ474R ■ SZ481R ■ SZ090R ■ SZ092R 	<ul style="list-style-type: none"> ■ Suitable cleaning brush ■ Disposable syringe 20 ml ■ Drying phase: Use a lint-free cloth or medical compressed air <p>SZ474R, SZ481R:</p> <ul style="list-style-type: none"> ■ Before phase I: Clean the product with a suitable cleaning brush under running water until all discernible residues have been removed from the surface. Mobilize non-rigid components, such as set screws, links, etc. during cleaning. ■ Rinse lumen with spray nozzle (3.8 bar water pressure) at least 3 times for 20 seconds each. 	<p>Chapter Manual cleaning/disinfection and subsection:</p> <ul style="list-style-type: none"> ■ Chapter Manual cleaning with ultrasound and immersion disinfection
<p>Mechanical alkaline cleaning and thermal disinfection</p> <ul style="list-style-type: none"> ■ SZ464R, SZ465R, SZ468R, SZ469R, SZ471R, SZ472R, SZ476R to SZ479R, SZ488R, SZ489R 	<ul style="list-style-type: none"> ■ Place the product on a tray that is suitable for cleaning (avoid rinsing blind spots). ■ Place products in the tray with their hinges open. 	<p>Chapter Mechanical cleaning/disinfection and subsection:</p> <ul style="list-style-type: none"> ■ Chapter Mechanical alkaline cleaning and thermal disinfecting
<p>Manual pre-cleaning with brush and subsequent mechanical alkaline cleaning and thermal disinfection</p> <ul style="list-style-type: none"> ■ SZ466R, SZ470R, SZ474R, SZ475R, SZ477R, SZ480R, SZ482R, SZ483R, SZ090R, SZ092R 	<ul style="list-style-type: none"> ■ Suitable cleaning brush ■ Disposable syringe 20 ml ■ Place the product on a tray that is suitable for cleaning (avoid rinsing blind spots). ■ Place products in the tray with their hinges open. <p>SZ466R, SZ470R, SZ474R, SZ475R, SZ477R, SZ480R, SZ482R, SZ483R:</p> <ul style="list-style-type: none"> ■ Connect the lumens and channels of individual parts directly to the injector unit's special flushing connector. 	<p>Chapter Mechanical cleaning/disinfection with manual pre-cleaning and subsection:</p> <ul style="list-style-type: none"> ■ Chapter Manual pre-cleaning with a brush ■ Chapter Mechanical alkaline cleaning and thermal disinfecting
<p>Manual pre-cleaning with ultrasound and brush, and subsequent mechanical alkaline cleaning and thermal disinfection</p> <ul style="list-style-type: none"> ■ SZ481R 	<ul style="list-style-type: none"> ■ Suitable cleaning brush ■ Disposable syringe 20 ml ■ Rinse lumen with spray nozzle (3.8 bar water pressure) at least 3 times for 20 seconds each. ■ Place the product on a tray that is suitable for cleaning (avoid rinsing blind spots). ■ Place products in the tray with their hinges open. ■ Connect the lumens and channels of individual parts directly to the injector unit's special flushing connector. 	<p>Chapter Mechanical cleaning/disinfection with manual pre-cleaning and subsection:</p> <ul style="list-style-type: none"> ■ Chapter Manual pre-cleaning with ultrasound and brush ■ Chapter Mechanical alkaline cleaning and thermal disinfecting

4.9 Manual cleaning/disinfection

- ▶ Prior to manual disinfecting, allow water to drip off for a sufficient length of time to prevent dilution of the disinfecting solution.
- ▶ After manual cleaning/disinfection, check visible surfaces visually for residues.
- ▶ Repeat the cleaning/disinfection process if necessary.

4.9.1 Manual cleaning with immersion disinfection

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Disinfecting cleaning	RT (cold)	>15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
II	Intermediate rinse	RT (cold)	1	-	D-W	-
III	Disinfection	RT (cold)	5	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
IV	Final rinse	RT (cold)	1	-	FD-W	-
V	Drying	RT	-	-	-	-

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

RT: Room temperature

*Recommended: BBraun Stabimed fresh

- ▶ Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

Phase I

- ▶ Fully immerse the product in the cleaning/disinfectant for at least 15 min. Ensure that all accessible surfaces are moistened.
- ▶ Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- ▶ If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
- ▶ Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- ▶ Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

Phase II

- ▶ Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- ▶ Drain any remaining water fully.

Phase III

- ▶ Fully immerse the product in the disinfectant solution.
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- ▶ Rinse lumens at least 5 times at the beginning of the exposure time using an appropriate disposable syringe. Ensure that all accessible surfaces are moistened.

Phase IV

- ▶ Rinse/flush the product thoroughly (all accessible surfaces).
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during final rinse.
- ▶ Rinse lumens with an appropriate disposable syringe at least five times.
- ▶ Drain any remaining water fully.

Phase V

- ▶ Dry the product in the drying phase with suitable equipment (e.g. cloth, compressed air), see Validated cleaning and disinfection procedure.

4.9.2 Manual cleaning with ultrasound and immersion disinfection

Phase	Step	D [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Ultrasonic cleaning	RT (cold)	>15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
II	Intermediate rinse	RT (cold)	1	-	D-W	-
III	Disinfection	RT (cold)	5	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
IV	Final rinse	RT (cold)	1	-	FD-W	-
V	Drying	RT	-	-	-	-

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

RT: Room temperature

*Recommended: BBraun Stabimed fresh

- ▶ Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

Phase I

- ▶ Clean the product in an ultrasonic cleaning bath (frequency 35 kHz) for at least 15 min. Ensure that all accessible surfaces are immersed and acoustic shadows are avoided.
- ▶ Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- ▶ If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
- ▶ Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- ▶ Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

Phase II

- ▶ Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- ▶ Drain any remaining water fully.

Phase III

- ▶ Fully immerse the product in the disinfectant solution.
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- ▶ Rinse lumens at least five times at the beginning of the exposure time with an appropriate disposable syringe. Ensure that all accessible surfaces are moistened.

Phase IV

- ▶ Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during final rinse.
- ▶ Rinse lumens with an appropriate disposable syringe at least five times.
- ▶ Drain any remaining water fully.

Phase V

- ▶ Dry the product in the drying phase with suitable equipment (e.g. cloth, compressed air), see Validated cleaning and disinfection procedure.

4.10 Mechanical cleaning/disinfection

Note

The cleaning and disinfection device must be of tested and approved effectiveness (e.g. FDA approval or CE mark according to DIN EN ISO 15883).

Note

The cleaning and disinfection device used for processing must be serviced and checked at regular intervals.

4.10.1 Mechanical alkaline cleaning and thermal disinfecting

Machine type: single-chamber cleaning/disinfection device without ultrasound

Phase	Step	T [°C/°F]	t [min]	Water quality	Chemistry/Note
I	Pre-rinse	<25/77	3	D-W	-
II	Cleaning	55/131	10	FD-W	<ul style="list-style-type: none"> ■ Concentrate, alkaline: <ul style="list-style-type: none"> - pH = 13 - <5 % anionic surfactant ■ 0.5 % working solution <ul style="list-style-type: none"> - pH ~ 11*
III	Intermediate rinse	>10/50	1	FD-W	-
IV	Thermal disinfection	90/194	5	FD-W	-
V	Drying	-	-	-	According to the program for cleaning and disinfection device

DW: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

*Recommended: BBraun Helimatic Cleaner alkaline

- ▶ Check visible surfaces for residues after mechanical cleaning/disinfecting.

4.11 Mechanical cleaning/disinfection with manual pre-cleaning

Note

The cleaning and disinfection device must be of tested and approved effectiveness (e.g. FDA approval or CE mark according to DIN EN ISO 15883).

Note

The cleaning and disinfection device used for processing must be serviced and checked at regular intervals.

4.11.1 Manual pre-cleaning with a brush

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Disinfecting cleaning	RT (cold)	>15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
II	Rinsing	RT (cold)	1	-	D-W	-

D-W: Drinking water

RT: Room temperature

*Recommended: BBraun Stabimed fresh

- ▶ Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

Phase I

- ▶ Fully immerse the product in the cleaning/disinfectant for at least 15 min. Ensure that all accessible surfaces are moistened.
- ▶ Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- ▶ If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
- ▶ Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- ▶ Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

Phase II

- ▶ Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.

4.11.2 Manual pre-cleaning with ultrasound and brush

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Ultrasonic cleaning	RT (cold)	>15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
II	Rinsing	RT (cold)	1	-	D-W	-

D-W: Drinking water

RT: Room temperature

*Recommended: BBraun Stabimed fresh

- ▶ Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

Phase I

- ▶ Clean the product in an ultrasonic cleaning bath (frequency 35 kHz) for at least 15 min. Ensure that all accessible surfaces are immersed and acoustic shadows are avoided.
- ▶ Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- ▶ If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
- ▶ Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- ▶ Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

Phase II

- ▶ Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.

4.11.3 Mechanical alkaline cleaning and thermal disinfecting

Machine type: single-chamber cleaning/disinfecting machine without ultrasound

Phase	Step	T [°C/°F]	t [min]	Water quality	Chemical
I	Prerinse	<25/77	3	D-W	-
II	Cleaning	55/131	10	FD-W	<ul style="list-style-type: none"> ■ Concentrate, alkaline: <ul style="list-style-type: none"> - pH ~ 13 - <5 % anionic surfactant ■ 0.5 % working solution <ul style="list-style-type: none"> - pH = 11*
III	Intermediate rinse	>10/50	1	FD-W	-
IV	Thermal disinfection	90/194	5	FD-W	-
V	Drying	-	-	-	According to the program for cleaning and disinfection device

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

*Recommended: BBraun Helimatic Cleaner alkaline

- ▶ Check visible surfaces for residues after mechanical cleaning/disinfecting.

4.12 Inspection

- ▶ Allow the product to cool down to room temperature.
- ▶ Dry the product if it is wet or damp.

4.12.1 Visual inspection

- ▶ Ensure that all soiling has been removed. In particular, pay attention to mating surfaces, hinges, shafts, recessed areas, drill grooves and the sides of the teeth on rasps.
- ▶ If the product is dirty: repeat the cleaning and disinfection process.
- ▶ Check the product for damage, e.g. insulation or corroded, loose, bent, broken, cracked, worn or severely scratched and fractured components.
- ▶ Check the product for missing or faded labels.
- ▶ Check the cutting edges for continuity, sharpness, nicks and other damage.
- ▶ Check the surfaces for rough spots.
- ▶ Check the product for burrs that could damage tissue or surgical gloves.
- ▶ Check the product for loose or missing parts.
- ▶ Immediately put aside damaged or inoperative products and send them to Aesculap Technical Service, see Technical service.

4.12.2 Functional test

CAUTION

Damage (metal cold welding/friction corrosion) to the product caused by insufficient lubrication!

- ▶ **Prior to function checks, lubricate moving parts (e.g. joints, pusher components and threaded rods) with maintenance oil suitable for the respective sterilization process (e.g. for steam sterilization: STERILIT® I oil spray JG600 or STERILIT® I drip lubricator JG598).**
- ▶ Assemble disassembled products, see Assembly.
- ▶ Check that the product functions correctly.
- ▶ Check that all moving parts are working properly (e.g. hinges, locks/latches, sliding parts etc.).
- ▶ Check rotating products (e.g. reusable drills and cutters) for bends and deformities. To do this, roll the product on an even surface.
- ▶ Check for compatibility with associated products.
- ▶ Immediately put aside inoperative products and send them to Aesculap Technical Service, see Technical service.

4.13 Assembly

SZ482R and SZ483R

Slide the PEEK sleeve **a** over the ring with adapter until the click. Screw the pusher with the help of SZ390R, see Fig. 5.

SZ474R

Assemble the depth sleeve **5a** and the handle lock **5c** in the handle **5b**, see Fig. 6.

4.14 Packaging

- ▶ Place the product in its holder or on a suitable tray. Ensure that sharp edges are covered.
- ▶ Package trays appropriately for the sterilization process (e.g. in Aesculap sterile containers).
- ▶ Ensure that the packaging provides sufficient protection against contamination of the product during storage.

4.15 Steam sterilization

Note

The product may only be sterilized in disassembled condition.

- ▶ Check to ensure that the sterilizing agent will come into contact with all external and internal surfaces (e.g., by opening any valves and faucets).
- ▶ Validated sterilization process
 - Disassemble the product
 - Steam sterilization using fractional vacuum process
 - Steam sterilizer according to DIN EN 285 and validated according to DIN EN ISO 17665
 - Sterilization using fractionated vacuum process at 134 °C/holding time 5 min
- ▶ If several devices are sterilized at the same time in the same steam sterilizer: Ensure that the maximum permitted load according to the manufacturers' specifications is not exceeded.

4.16 Storage

- ▶ Store sterile single-use products in germ-proof packaging in a dust-protected, dry, dark and temperature-controlled room.

5. Technical service

CAUTION

Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

- ▶ **Do not modify the product.**
- ▶ **For service and repairs, please contact your national B. Braun/Aesculap agency.**

Service addresses

Aesculap Technischer Service

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78532 Tuttlingen / Germany

Phone: +49 7461 95-1601

Fax: +49 7461 16-2887

E-Mail: ats@aesculap.de

Other service addresses can be obtained from the address indicated above.

6. Disposal

WARNING

Risk of infection due to contaminated products!

- ▶ Adhere to national regulations when disposing of or recycling the product, its components and its packaging.

WARNING

Risk of injury due to sharp-edged and/or pointed products!

- ▶ When disposing of or recycling the product, ensure that the packaging prevents injury by the product.

Note

The user institution is obliged to reprocess the product before its disposal, see Validated reprocessing procedure.

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