**GB** Instructions for use/Technical description  
Modulift VBR – Size Small Implants

**USA** Instructions for use/Technical description  
Modulift VBR – Size Small Implants  
Note for U.S. users

This Instructions for Use is NOT intended for United States users. Please discard. The Instructions for Use for United States users can be obtained by visiting our website at [www.aesculapimplantsystems.com](http://www.aesculapimplantsystems.com) and clicking the "Products" menu. If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-866-229-3002. A paper copy will be provided to you upon request at no additional cost.

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## Aesculap® Modulift VBR – Size Small Implants

### Intended use

The Aesculap Modulift VBR System is used for partial or total replacement of a collapsed, damaged or unstable vertebral body in the cervical spine (C3 – C7). The Aesculap Modulift VBR System is used with supplemental spinal fixation systems. The Aesculap Modulift VBR System may be used with bone graft.

### Materials

The materials used in the implant are listed on the packaging.

The sequence of materials is:

- Ti6Al4V acc. to ASTM F136/ISO 5832-3
- CoCrMo acc. to ASTM F1537

The titanium implants are coated with a colored oxide. Slight changes in color may occur, but do not affect the implant quality.

### Indications

Surgically installed implants serve to support normal healing processes. They are not intended for use either as replacement for natural body parts or to bear loads over the long term in the event the healing process is incomplete or slow.

Use for the following indications:

- Fractures
- Anterior correction of spinal deformities
- Degenerative disc disease (e.g., spondylolisthesis, spondylolysis, spinal stenosis)
- Deformities
- Degenerative instability
- Post-trauma instability
- Spinal tumors

WARNING

#### Risk of implant breaking or improper healing!

- ▶ Support the implant with supplemental fixation such as a posterior rod or anterior plating system.

### Absolute contraindications

Do not use in the presence of:

- Severe damage to the bone structures of the spine that could prevent the stable implantation of the implant components; for example, osteopenia, severe osteoporosis, Paget's disease, bone tumors etc.
- Metabolic or degenerative metabolic bone diseases that could compromise the stable anchoring of the implant system
- Suspected allergy or sensitivity to the implant materials
- Acute or chronic vertebral infections of a local or systemic nature
- Poor patient compliance or limited ability to follow medical instructions, particularly in the post-op phase, including with regard to the restrictions on range of movement in terms of physical exercise and occupational activity
- Cases not listed under indications

### Relative contraindications

In the following circumstances, use of the implant system could represent an increased clinical risk and therefore requires precise, individual assessment by the surgeon:

- Medical or surgical conditions that could negatively impact the success of the implantation, including wound healing disorders
- Conditions that could subject the spine and implants to excessive pressure; for example, pregnancy, obesity, neuromuscular diseases or disorders
- Generally poor condition of the patient; for example, drug or alcohol addiction

### Side effects and interactions

The general risks of surgery are assumed known and are therefore not described in these instructions for use. Potential risks and clinical consequences associated with use of the implant system include, but are not limited to:

- Implant failure resulting from excessive load
  - Warping or bending
  - Loosening
  - Breakage
- Inadequate fixation
- Failed or delayed fusion
- Infection
- Fractured vertebral body or bodies
- Injuries to
  - nerve roots
  - spinal cord
  - blood vessels
  - organs
- Changes of the normal spine lordosis
- Impairment of the gastrointestinal, urological and/or reproductive systems
- Pain or indisposition
- Bursitis
- Decreased bone density due to load avoidance
- Bone atrophy/fracture above or below the operated spine section
- Limited physical capacity
- Persistence of symptoms that were to be treated by the implantation
- Death

### Safety notes

- It is the operating surgeon's responsibility to ensure that the surgical procedure is performed properly.
- General risk factors associated with surgical procedures are not described in this documentation.
- The operating surgeon must have thorough command of both the hands-on and conceptual aspects of the established operating techniques.
- The operating surgeon must be fully conversant with the bone anatomy and the pathways of nerves and blood vessels, muscles, and tendons.
- Aesculap is not responsible for any complications arising from erroneous indication, wrong choice of implant, incorrect combination of implant components and operating technique, the limitations of the treatment method, or inadequate asepsis.
- The user instructions for individual Aesculap implant components must be followed.
- The implant components were tested and approved in combination with Aesculap components. If other combinations are used, the responsibility for such action lies with the operating surgeon.
- Do not, under any circumstances, combine implant components from different manufacturers.
- Do not, under any circumstances, use damaged or surgically excised components.
- Implants that have been used before must not be reused.

WARNING

**The implant has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration or image artifact in the MR environment. Scanning a patient who has this implant may result in patient injury.**

- Delayed healing can cause implant breakage due to metal fatigue.
- The attending physician shall make any decision with regard to the removal of implants that have been locked down.
- The implant components applied, along with their article numbers, the name of the implant, as well as the batch number and serial number (if available) must be documented in all patient records.
- Postoperatively, individual patient information, as well as mobility and muscle training, is of particular importance.
- Damage to the load-bearing structures of the implant can lead to loosening of components, dislocation, migration and other severe complications.

### Sterility

- The implants come individually packed in protective packaging that is labeled according to its contents.
- The implants are sterilized by irradiation.
- ▶ Store implant components in their original packaging. Remove them from their original protective packaging only just prior to application.
- ▶ Prior to use, check the product expiration date and verify the integrity of the sterile packaging.
- ▶ Do not use implants that are past their expiration date or whose packaging is damaged.

WARNING

#### Damage to implants caused by processing and resterilization!

- ▶ Do not reprocess or resterilize the implants.
- ▶ Never resterilize used, explanted, or damaged implants.

### Application

The operating surgeon shall devise an operation plan that specifies and accurately documents the following:

- Selection of the implant components and their dimensions.
- Positioning of the implant components in the bone.
- Location of intraoperative landmarks.

The following conditions must be fulfilled prior to application:

- All requisite implant components are ready to hand.
- Operating conditions are highly aseptic.
- All requisite implantation instruments must be available and in working order, including specialized Aesculap implantation systems.
- The operating surgeon and operating room team are thoroughly conversant with the operating technique and with the available range of implants and instruments; information materials on these subjects must be complete and ready to hand.
- The operating surgeon is fully conversant with the rules governing medical practice, the current state of scientific knowledge, and the contents of relevant scientific articles by medical authors.
- The manufacturer has been consulted if the preoperative situation was unclear and if implants were found in the area operated on.

The surgical procedure and following information has been explained to the patient, and the patient's consent has been documented:

- The patient is aware of the risks involved in neurosurgery, general surgery, orthopedic surgery, and general anesthesia.
- The patient has been informed of the advantages and disadvantages of a Modulift VBR implant and has been made aware of possible alternative methods of treatment.
- The special instrument set for the installation of implants and for distraction, compression, and reductions of the cervical, lumbar, and thoracic spine have been prepared and is ready.
- In the case of delayed or incomplete fusion, the implants can break and loosen due to high loads.
- The life-span of the implant depends on the patient's body weight.
- The implant components must not be overloaded by extreme strains, hard physical labor, or exercise.
- Corrective surgery may be necessitated by implant loosening, fracture, or loss of correction.
- Smokers present an increased risk of bone fusion failure.
- The patient must undergo medical check-ups of the implants at regular intervals.

### Preparation of surgical site

- ▶ To avoid internal stress on, and weakening of, the implant: avoid scoring or scratching of the implant components.
- ▶ Do not bend the implants.
- ▶ Using standard instruments, remove disc material and all cartilage from the adjacent vertebral end plates.
- ▶ Carry out a partial or complete corpectomy of the damaged vertebral body.

### Verification of Implant size

- ▶ Select the appropriate Modulift VBR implants according to the individual indication, preoperative planning and bone situation found intraoperatively.
- ▶ When selecting Modulift VBR implants, make certain to select the correct widths, height and correction angle.
- ▶ To select the proper implant, use the footplate sizer for footplate size, the trial implants for body height, and the trial implant caps for correction angle.
- ▶ Always use the largest possible trial implant for selecting the implant size.
- ▶ Use radiographic control to ensure that the trial implants are positioned centrally.

#### Note

The height of the trial implants corresponds to the maximum length of the implant not including the spikes.

## Implant Inserter

- ▶ Prepare and assemble the implant inserter according to its Instructions for Use.
- ▶ Connect the desired implant to the implant inserter assembly using the turn knob and clamp. Proper alignment is guided by using the clamp jaw tip geometry (mating pins, and dove tail).
- ▶ The orientation of the handle can be adjusted by loosening the tightening knob at the end of the handle. Retighten the knob after the desired orientation has been chosen.

## Bone Graft

WARNING

Bone chips may fall out while transferring implants from bench block to surgical site!

- ▶ Make sure that bone chips are firmly packed within the central column of the implant.
- ▶ Do not swing, shake, or otherwise move the implant in a way that would encourage bone chips to leave the implant.

WARNING

Damage to the implant caused by overpacking with bone material!

- ▶ Do not overpack implants.

- ▶ Using the bench block, graft funnel, and bone tamp, it is possible to add bone graft or other media into the central column of the chosen implant.
- ▶ It is up to the operating surgeon whether or not to add bone graft to the implant.

## Insertion and Expansion

WARNING

Incorrect application can result in trauma to the spinal cord and nerve roots!

- ▶ Insert, align, and lock down the implant only under X-ray control.
- ▶ Ensure that the implant spikes fully engage the vertebral end plates.

WARNING

Potential damage to implant or delay of surgery!

- ▶ Use the guide marks on the implant inserter ("Open", "Clamp") to attach implant to implant inserter.
- ▶ Only use hand strength to tighten implant inserter.

WARNING

Potential damage to implant or delay of surgery!

- ▶ Make sure that the expansion driver is properly seated in the implant inserter (held in by ball pluners and able to spin freely).
- ▶ Only use hand strength to expand the implant.

- ▶ For positioning the implants, always use the implant insertion instrument intended for this purpose. An implant tamp is also provided for minor adjustments after withdrawing the implant insertion instrument.
- ▶ Expand the implant so that the implant spikes engage with the end plates of the adjacent vertebral bodies.
- ▶ Always expand the implant using the provided expansion driver.
- ▶ Use the implant insertion instrument as a guide for the expansion driver so that it retains consistent engagement with the implant.

## Final Implant Verification

WARNING

Risk of injury due to improper tightening of the set screw!

- ▶ Make certain that the set screw driver is fully seated within the set screw before tightening.
- ▶ Fully tighten the set screw using the provided torque wrench.
- ▶ Always use the implant inserter as a counter torque during tightening.

WARNING

Improper positioning of the implant may lead to improper healing!

- ▶ Never loosen the set screw after they have been tightened if the implant will stay in the patient.
- ▶ Tighten the set screw only after the implant is expanded to the correct height and positioned within the patient.

WARNING

Damage to the implant caused by overtightening of the set screw!

- ▶ Always tighten the clamping screw with the exact torque specified for this purpose.

WARNING

Damage to the implant caused by tightening of the set screw without using the countering instrument!

- ▶ Always use the countering instrument (implant insertion instrument) when tightening the set screw.

WARNING

Damage to the implant caused by expanding the implant after tightening the set screws!

- ▶ Always expand and position the implant before tightening the set screws.
- ▶ Do not expand the implant after tightening the set screws.

- ▶ After confirming the final position of the implant, remove the expansion driver and tighten the set screws. Always use the torque handle, set screw driver, and implant insertion instrument (which acts as a counter-torque) provided for this purpose.
- ▶ If the implant must be removed, loosen the set screws and contract the implant. Always use the provided set screw removal driver and implant insertion instrument as a countertorque.

Further information on Aesculap implant systems is always available from B. Braun/Aesculap or the appropriate B. Braun/Aesculap office.

## Disposal

- ▶ Adhere to national regulations when disposing of or recycling the product, its components and its packaging!

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