AESCULAP® S4® Spinal System

en	Instructions for use/Technical description Augmentation screw
USA	Note for U.S. users This Instructions for Use is NOT intended for United States users. Please discard. The Instructions for Use for United States users can be obtained by visiting our website at www.aesculapimplantsystem-sifus.com . If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-866-229-3002. A paper copy will be provided to you upon request at no additional cost.
de	Gebrauchsanweisung/Technische Beschreibung Augmentationschraube
fr	Mode d'emploi/Description technique Vis d'augmentation
es	Instrucciones de manejo/Descripción técnica Tornillo para sustitución
it	Istruzioni per l'uso/Descrizione tecnica Vite di augmentazione
pt	Instruções de utilização/Descrição técnica Parafuso de augmentação
nl	Gebruiksaanwijzing/Technische beschrijving Augmentatieschroef
da	Bruksanvisning/Teknisk beskrivelse Ekspansionskrue
nb	Bruksanvisning/Teknisk beskrivelse Forsterkningskrue
sv	Bruksanvisning/Teknisk beskrivning Förstärkningskruv
fi	Käyttöohje/Tekninen kuvaus Augmentaatioruuvi
et	Kasutusjuhend/Tehniline kirjeldus Tugikruvi
lv	Lietošanas instrukcijas/tehniskais apraksts Pastiprināšanas skrūve
lt	Naudojimo instrukcija/techninis aprašas Augmentacijos sraigtas
ru	Инструкция по применению/Техническое описание Усилительный винт
cs	Návod k použití/Technický popis Augmentační šroub
pl	Instrukcja użytkowania/Opis techniczny Śruba augmentacyjna
sk	Návod na použitie/Technický opis Augmentačná skrutka
hu	Használati útmutató/Műszaki leírás Augmentációs csavar
sl	Navodila za uporabo/Tehnični opis Avgmentacijski vijak
hr	Upute za uporabu/Tehnički opis Vijak za augmentaciju
ro	Manual de utilizare/Descriere tehnică Șurub de augmentare
bg	Упътване за употреба/Техническо описание Аугментационен винт
tr	Kullanım Kılavuzu/Teknik açıklama Büyüme vidası
el	Οδηγίες χρήσης/Τεχνική περιγραφή Βίδα ενίσχυσης
zh	使用说明 / 技术说明 S4® 胸腰椎后路固定系统 - 骨水泥增强型螺钉

B | BRAUN

Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany
Phone +49 (0) 7461 95-0 | Fax +49 (0) 7461 95-26 00 | www.bbraun.com

AESCULAP® – a B. Braun brand

TA012865 2022-06 Change No. AE0061814

AESCULAP® S⁴ Spinal System

Augmentation screw

1. About this document

Note

General risk factors associated with surgical procedures are not described in these instructions for use.

Note

Generally, the S⁴ Spinal System augmentation screw is covered by the instructions for use of the S⁴ Spinal System – Lumbar/Deformity TAO11187. This information on the augmentation screw is intended to complement the respective information contained in the instructions for use of the S⁴ Spinal System – Lumbar/Deformity.

1.1 Scope

These instructions for use apply to all S⁴ Spinal System augmentation screws.

- ▶ For article-specific instructions for use as well as information on material compatibility and lifetime see B. Braun eIFU at eifu.bbraun.com

1.2 Safety messages

Safety messages make clear the dangers to patient, user and/or product that could arise during the use of the product. Safety messages are labeled as follows:

⚠ WARNING

Indicates a possible threat of danger. If not avoided, minor or moderate injury may result.

⚠ CAUTION

Indicates a possible threat of material damage. If not avoided, the product may be damaged.

2. Clinical use

Note

The short summary of safety and clinical performance of the product is available in the European Database for Medical Products (EUDAMED).

2.1 Product description

Additional information on Aesculap implant systems is available from B. Braun/Aesculap or from your local B. Braun/Aesculap agency.

2.1.1 Materials

The materials used in the implant are listed on the packaging:

- ISOTAN®_F titanium forged alloy Ti6Al4V according to ISO 5832-3
- ISOTAN®_P pure titanium according to ISO 5832-2

The titanium implants are anodized with a colored oxide layer. Slight changes in coloration may occur, but do not affect the implant quality.

ISOTAN® is a registered trademarks of Aesculap AG, 78532 Tuttlingen / Germany.

2.2 Areas of use and limitations of use

2.2.1 Intended use

The S⁴ Spinal System implants are used for dorsal monosegmental and multisegmental stabilization of the lumbar and thoracic spine.

2.2.2 Indications

Note

The manufacturer is not responsible for any use of the product against the specified indications and/or the described applications.

Use for severe disorders of the spine that cannot be treated by conservative therapies:

- Fractures
- Pseudarthrosis or delayed healing
- Degenerative instability
- Post-discectomy syndrome
- Spondylolisthesis
- Hyper-/hypokyphosis
- Hyper-/hypolordosis
- Scoliosis
- Stenosis
- Post-trauma instability
- Tumors
- Osteoporosis (with the simultaneous use of bone cement and fenestrated screws)

2.2.3 Absolute contraindications

Do not use in the presence of:

- Severe damage to the bone structures of the spine that could prevent the stable implantation of the implant components; for example, Paget's disease, bone tumors etc.
- Metabolic or degenerative metabolic bone diseases that could compromise the stable anchoring of the implant system
- Suspected allergy or sensitivity to the implant materials
- Acute or chronic vertebral infections of a local or systemic nature
- Poor patient compliance or limited ability to follow medical instructions, particularly in the post-op phase, including with regard to the restrictions on range of movement in terms of physical exercise and occupational activity
- Cases not listed under indications

2.2.4 Relative contraindications

In the following circumstances, use of the implant system could represent an increased clinical risk and therefore requires precise, individual assessment by the surgeon:

- Medical or surgical conditions that could negatively impact the success of the implantation, including wound healing disorders
- Conditions that could subject the spine and implants to excessive pressure; for example, pregnancy, obesity, neuromuscular diseases or disorders
- Generally poor condition of the patient; for example, drug or alcohol addiction

2.3 Risks, adverse effects and interactions

As part of the legal duty to inform, the following typical risks, interactions and side effects associated with the use of surgical instruments are pointed out.

The general risks of surgery are assumed known and are therefore not described in these instructions for use. Potential risks and clinical consequences associated with use of the implant system include, but are not limited to:

- Implant failure resulting from excessive load
 - Warping or bending
 - Loosening
 - Breakage
- Poor stabilization, requiring subsequent revision surgery or surgery to remove the implant system
- Failed or delayed fusion

- Cement leakages (cement emerging from the vertebral body)
- Infection
- Fractured vertebral body or bodies
- Injuries to
 - Nerve roots
 - Spinal cord
 - Blood vessels
 - Organs
- Changes of the normal spine lordosis
- Neurological disorders, including impairment of the gastrointestinal, urological and/or reproductive system
- Pain or indisposition
- Bursitis
- Decreased bone density due to load avoidance
- Bone atrophy/fracture above or below the spine section provided for
- Limited physical ability
- Persistence of symptoms that were to be treated by the implantation
- Paresis
- Subsequent surgery to remove the implant system
- Metal/foreign body allergy or hypersensitivity
- Dural tears
- Meningitis

2.4 Safety information

2.4.1 Clinical user

General safety information

To prevent damage caused by improper setup or operation, and to not compromise the manufacturer warranty and liability:

- ▶ Use the product only according to these instructions for use.
- ▶ Always follow the safety advice and information given in the instructions for use.
- ▶ Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge and experience.
- ▶ Store any new or unused products in a dry, clean, and safe place.
- ▶ Keep the instructions for use accessible for the user.

Note

The user is obligated to report all severe events in connection with the product to the manufacturer and the responsible authorities of the state in which the user is located.

Notes on surgical procedures

It is the user's responsibility to ensure that the surgical procedure is performed correctly.

Appropriate clinical training as well as a theoretical and practical proficiency of all the required operating techniques, including the use of this product, are prerequisites for the successful use of this product.

Aesculap is not responsible for complications caused by:

- incorrect indication or implant selection
- incorrect surgical technique
- incorrect combination of implant components
- combination not approved by Aesculap with components from other manufacturers
- exceeding the limitations of the treatment method or non-observance of essential medical precautions

The user is required to obtain information from the manufacturer if there is an unclear preoperative situation regarding the use of the product.

Note

There are special S⁴ instruments provided for the implantation of these system components and for the augmentation, distraction, compression, and reduction of the lumbar and thoracic spine.

2.4.2 Product-specific safety information

⚠ WARNING

The implant has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration or image artifact in the MR environment. Scanning a patient who has this implant may result in patient injury.

2.4.3 Sterility

The product has been sterilized by irradiation and is supplied in sterile packaging.

The safe sterile provision of the product is only guaranteed if the sterile packaging is undamaged and unopened and the use-by date has not passed.

- ▶ Protect products from direct sunlight.
- ▶ Store products in their original packaging in a dust-protected, dry and temperature-controlled room.
- ▶ Remove products from their original protective packaging only just prior to application.
- ▶ Visually inspect the product packaging to ensure that the sterile barrier system is intact.
- ▶ Do not use products from open or damaged sterile packaging.
- ▶ Do not use the product after its use-by date.

The product is intended for single-use.

- ▶ Do not reuse the product.

The processing of the product affects its functionality. Risk of injury, illness or death due to soiling and/or impaired functionality of the product.

- ▶ Do not process the product.

2.4.4 Storage

- ▶ Store sterile packed implant components dust-protected in a dry, dark and temperature-controlled room.

2.5 Patient education

Within the framework of the patient education, the relevant circumstances needed for consent must be explained to the patient in accordance with their level of understanding, pre-existing knowledge and need for information. This includes:

- Diagnosis, procedure and risk clarification
- Operative procedure
- Advantages and disadvantages of the procedure
- All alternative procedures that can be considered

The patient must be properly informed about the procedure and in particular about the following information:

- Delayed healing or incomplete fusion can cause the implant to fracture or loosen as a result of the extreme load to which it is subjected.
- The life-span of the implant depends on the patient's body weight.
- The implant components must not be overloaded by extreme strains, hard physical labor or sports.
- Corrective surgery may be necessitated by implant loosening, fracture or loss of correction.
- Smokers present an increased risk of bone fusion failure.
- The patient must undergo medical check-ups of the implant components at regular intervals.
- Postoperatively, individual patient information, as well as mobility and muscle training, is of particular importance.
- The attending physician shall make any decision with regard to the removal of implant components that have been used.
- Damage to the load-bearing structures of the implant can lead to loosening of components, dislocation, migration, and other severe complications.

2.6 Application

2.6.1 Documentation

The user shall devise an operation plan that specifies and accurately documents the following:

- Selection of the implant components and their dimensions
- Positioning of the implant components in the bone
- Location of intraoperative landmarks

Important information concerning the implanted product and the operation can be noted down on the patient ID. The patient ID can be ordered separately from the manufacturer.

Each package contains additional labels showing the designation, article and lot number and – if applicable – the individual serial number of the product.

- ▶ Use these labels for documentation in the patient's file (for the hospital) and the patient ID (for the patient).

The surgical procedure has been explained to the patient, and the patient's consent has been documented.

2.6.2 Implantation

The S⁴ Spinal System augmentation screw can be fixed with bone cement to increase anchoring stability. In this case, the injection cannula is inserted in the S⁴ Spinal System augmentation screw for application of the bone cement.

- S⁴ Monoaxial/polyaxial screws (augmentation screw), supplied in sterile condition
- S⁴ Element monoaxial/polyaxial screws (augmentation screw), supplied in sterile condition
- Cement injection cannula (sterile), see TA013132
- for percutaneous application with S⁴ Element monoaxial/polyaxial screws (augmentation screws): S⁴ Element Augmentation instruments, see TA014315

⚠ WARNING

Risk of injury and/or malfunction!

- ▶ Prior to each use, inspect instruments for loose, bent, broken, cracked, worn, or fractured components.
- ▶ Always carry out a function test prior to each use of the instruments.

To implant the S⁴ Spinal System – augmentation screw, proceed as follows:

- ▶ Select the appropriate S⁴ Spinal System variant and anchoring implants according to the individual indication, preoperative planning, and bone situation found intraoperatively.
- ▶ When selecting S⁴ Spinal System screw types, make certain to select the correct diameter, length, and canal orientation.

⚠ WARNING

Incorrect application can result in trauma to the spinal cord and nerve roots!

- ▶ Align and screw in the instruments, the pedicle screw, and the polyaxial screw only under X-ray control or using a navigation system.
- ▶ To avoid internal stress on, and weakening of, the implant: avoid scoring or scratching of the implant components.
- ▶ Only allow axial load on the screws during screw insertion. Avoid induction of lateral forces during screw insertion.

⚠ WARNING

Danger of screw breakage due to hard bone or contact with lateral cortex!

- ▶ Prepare the vertebral body over the entire screw length with thread cutter.

⚠ WARNING

Risk of connecting rod breakage!

- ▶ Do not bend back or bend excessively the connecting rods and cross connectors.

⚠ CAUTION

Damage to the hexagon socket due to incorrect application of the screwdriver tip!

- ▶ Make certain the screwdriver tip is fully inserted in the hexagon socket of the screw head when turning the screw.
- ▶ If the screw head loosens during revision of the larger polyaxial screws (∅ 7.5 mm / ∅ 8.5 mm), remove the remaining screw component with the size 3.5 screwdriver.

Fixing the augmentation screw with bone cement

To increase anchoring stability, the S⁴ Spinal System – augmentation screw can be fixed with bone cement. The following must be observed when cementing:

Note

Cementing is carried out only after all screws have been aligned and the connecting rods have been adapted, and before locking the S⁴ Spinal System augmentation screw.

⚠ WARNING

Risk of breaking the injection cannula by inducing excessive shear forces!

- ▶ Do not subject the injection cannula to transverse loads.
- ▶ Do not bend the injection cannula.

Note

Once the S⁴ Spinal System augmentation screw has been cemented, it will be impossible to change the implant position, use the reduction lever, or make corrections!

Note

With the S⁴ Spinal System augmentation screw finally cemented, explanting the screw can become extremely difficult, especially in cases of osteopenia!

- ▶ Prior to applying the bone cement, fully screw in the S⁴ injection cannula at the thread in the bone screw head, without applying perceptible force.
- ▶ Check the bone cement for proper texture and viscosity before positioning the applicator on the injection cannula. Observe the manufacturer's instructions for the bone cement!
- ▶ Cement S⁴ Spinal System augmentation screw with a cement suitable for screw augmentation.

⚠ WARNING

Risk of injury from cement leakages!

- ▶ Always apply imaging techniques for control during the preparation and positioning of the implants.
- ▶ Make certain the slots of the S⁴ Spinal System augmentation screw are positioned fully inside the bone.
- ▶ Make sure that the pedicle cortex is not perforated.
- ▶ Do not penetrate the opposite cortex and ensure that the opposite cortex is undamaged.
- ▶ Do not insert or cement the S⁴ Spinal System augmentation screw in a fractured or partially resected vertebral body. Perform a thin-slice CT scan preoperatively.
- ▶ Screw the S⁴ injection cannula into the screw head without applying any perceptible force. Make certain the axes of the pedicle screw and the injection cannula are aligned with each other.
- ▶ Apply cement only under high-quality radiographic control (max. volume of cement per polyaxial screw: 2 ml).

⚠ WARNING

Risk of injury due to material incompatibility!

- ▶ Only use application systems that have been approved for use with bone cement.

⚠ WARNING

Contamination of the locking thread/insufficient fixation when removing the application system, due to incorrect viscosity of the cement!

- ▶ Make certain the cement has the correct viscosity.
- ▶ Do not allow the cement to cure completely before removing the injection cannula. Remove the cement application system incl. cannula as soon as the cement shows a doughy texture (allow for a certain amount of cement for reference purposes when preparing the cement).

⚠ WARNING

Wound contamination due to insufficient curing of cement! Risk of caulking of the locking screw thread, resulting in system instability!

- ▶ Observe the correct curing period after removing the injection cannula.
- ▶ For monitoring the cement texture, prepare a reference volume of cement together with the cement to be applied.
- ▶ Make certain the system is not exposed to any load before the cement has reached its final strength.

⚠ WARNING

Risk of screw breakage caused by pedicle screws used for spondylolisthesis!

- ▶ Support the stabilization measure with interbody fusion (ALIF, PLIF or TLIF).

3. Disposal

⚠ WARNING

Risk of infection due to contaminated products!

- ▶ Adhere to national regulations when disposing of or recycling the product, its components and its packaging.

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