
Aesculap Spine

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## Aesculap®

S<sup>4</sup> Injection cannula for augmentation screws SR146SU,

S<sup>4</sup> Injection cannula for augmentation screws SR148SU

### Legend

- 1 Injection cannula for open access SR146SU
- 2 Injection cannula for percutaneous access SR148SU
- 2a Cannula
- 2b Plunger

## 1. About this document

*Note*  
General risk factors associated with surgical procedures are not described in these instructions for use.

### 1.1 Scope

These instructions for use apply for the following products:

Art. no.	Designation
SR146SU	S <sup>4</sup> Injection cannula for augmentation screws
SR148SU	S <sup>4</sup> Injection cannula percutaneous for augmentation screws

► For article specific instructions for use and material compatibility and lifetime information, see B. Braun eFU at [eifu.bb.raun.com](http://eifu.bb.raun.com)

### 1.2 Safety messages

Safety messages make clear the dangers to patient, user and/or product that could arise during the use of the product. Safety messages are labeled as follows:

#### ⚠ WARNING

Indicates a possible threat of danger. If not avoided, minor or moderate injury may result.

#### ⚠ CAUTION

Indicates a possible threat of material damage. If not avoided, the product may be damaged.

## 2. Clinical use

### 2.1 Areas of use and limitations of use

#### 2.1.1 Intended use

The S<sup>4</sup> Spinal System injection cannula is inserted into the S<sup>4</sup> Spinal System augmentation screw to apply bone cement.

#### 2.1.2 Indications

*Note*  
The manufacturer is not responsible for any use of the product against the specified indications and/or the described applications.

For indications, see Intended use.

#### 2.1.3 Contraindications

No known contraindications.

## 2.2 Safety information

### 2.2.1 Clinical user

#### General safety information

To prevent damage caused by improper setup or operation, and to not compromise the manufacturer warranty and liability:

- Use the product only according to these instructions for use.
- Follow the safety and maintenance instructions.
- Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge and experience.
- Store any new or unused products in a dry, clean, and safe place.
- Prior to use, check that the product is in good working order.
- Keep the instructions for use accessible for the user.

#### *Note*

The user is obligated to report all severe events in connection with the product to the manufacturer and the responsible authorities of the state in which the user is located.

#### Notes on surgical procedures

It is the user's responsibility to ensure that the surgical procedure is performed correctly.

Appropriate clinical training as well as a theoretical and practical proficiency of all the required operating techniques, including the use of this product, are prerequisites for the successful use of this product.

The user is required to obtain information from the manufacturer if there is an unclear preoperative situation regarding the use of the product.

#### 2.2.2 Product specific safety information

For the S<sup>4</sup> Spinal System augmentation screw the instructions for use for the S<sup>4</sup> Spinal System lumbar/deformity TA011187 and S<sup>4</sup> Spinal System augmentation screw TA012865 generally apply.

This information on the S<sup>4</sup> Spinal System Injection cannula supplements the individual information in the instructions for use for the S<sup>4</sup> Spinal System lumbar/deformity and augmentation screw.

#### 2.2.3 Sterility

The product is sterilized by irradiation and sterile packaged.

- Do not use any product from open or damaged sterile packaging.
- Do not use product after expiration date.

The S<sup>4</sup> Spinal System Injection cannula is for single use only in combination with the S<sup>4</sup> Spinal System augmentation screw.

- Do not reuse product.

Processing the product will impair its functionality. Contamination and/or impaired function of the products may result in injury or illness and death.

- Do not process product.

## 2.3 Application

#### ⚠ WARNING

Risk of injury and/or malfunction!

- Prior to each use, inspect the product for loose, bent, broken, cracked, worn, or fractured components.
- Always carry out a function test prior to each use of the product.

#### ⚠ WARNING

Risk of injury from leakage of cement!

- Screw the injection cannula into the screw head without applying force.
- Make sure that the axes of the pedicle screw and the injection cannula are aligned with each other.
- Cement application only under high-quality radiographic control (max. volume of cement per pedicle screw: 2ml in the pedicle).

#### ⚠ WARNING

Contamination of the locking thread of the S<sup>4</sup> Spinal System augmentation screw and/or insufficient fixation when removing the application system, due to incorrect viscosity of the cement!

- Make sure the cement has the correct viscosity.

#### ⚠ WARNING

Risk of injury and/or malfunction!

- Screw the injection cannula hand-tight into the thread of the pedicle screw.
- Screw the injection cannula into the pedicle screw and leave in place until the cement has completely set.

#### ⚠ CAUTION

Risk of breakage of injection cannula if excessive shear forces are applied!

- Do not subject the cement injection cannula to excessive shear forces.
- Do not bend the injection cannula.

### 2.3.1 Injection cannula SR146SU for open access

#### ⚠ WARNING

Risk of injury due to application of the product if used outside the field of view!

- Only use the product under visual control.

- Screw in the pedicle screw in accordance with the instructions. Leave the K-wire in the pedicle screw.
- Remove the injection cannula SR146SU 1 from the sterile packaging and rinse using sterile saline solution.
- Screw the injection cannula SR146SU 1 using K-wire, hand-tight, in a line with the bone thread of the pedicle screw into the screw head.
- Remove K-wire.
- Check that the injection cannula SR146SU 1 is positioned properly.
- Connect the cement applicator to the Luer lock connector.
- Apply the cement.
- Remove the cement applicator from the Luer lock connector.
- Unscrew the injection cannula SR146SU 1 from the pedicle screw once the cement has completely set.
- Dispose of the injection cannula SR146SU 1 and do not reuse or process the latter.
- Continue the surgery and insert the rod.

### 2.3.2 Injection cannula SR148SU for percutaneous access

The SR148SU injection cannula can be used percutaneously only in conjunction with the S<sup>4</sup> clamp sleeve percutaneous and the S<sup>4</sup> augmentation screws.

#### ⚠ WARNING

Risk of injury and/or malfunction!

- Only insert the plunger 2b into the injection cannula SR148SU 2a.

#### ⚠ WARNING

Risk of injury and/or malfunction!

- Insert injection cannula SR148SU 2a using the K-wire or plunger 2b.
- Screw in the pedicle screw in accordance with the instructions. Leave the K-wire in the pedicle screw.
- Remove the injection cannula SR148SU 2a from the sterile packaging and rinse using sterile saline solution.
- Screw the injection cannula SR148SU 2a using the K-wire, hand-tight, in a line with the bone thread of the pedicle screw, into the screw head and remove the K-wire.
- or -
- Insert the plunger into the injection cannula. Insert the cannula into the screw head via the plunger.
- To check for proper positioning of the injection cannula SR148SU 2a, insert the plunger 2b into the injection cannula SR148SU 2a as far as it will go. When properly seated, the marking on the plunger 2b disappears completely inside the injection cannula SR148SU 2a.
- If the marking remains visible, check the alignment of the injection cannula SR148SU 2a. If necessary, screw in the injection cannula SR148SU 2a once again.
- Remove the plunger 2b from the injection cannula SR148SU 2a.
- Connect the cement applicator to the Luer lock connector.
- Apply the cement.
- Remove the cement applicator from the Luer lock connector. If necessary, use the plunger 2b to press cement into the screw.
- Unscrew the injection cannula SR148SU 2a from the pedicle screw once the cement has completely set.
- Dispose of the injection cannula SR148SU 2a and do not reuse or resterilize the latter.
- Continue the surgery and insert the rod.

## 3. Single use products

The S<sup>4</sup> Spinal System injection cannula is for single use only in combination with the S<sup>4</sup> Spinal System augmentation screw.

- Do not reuse product.

Processing the product will impair its functionality. Contamination and/or impaired function of the products may result in injury or illness and death.

- Do not process product.

## 4. Storage

- Store sterile packaged single-use products in a dry, dark and evenly tempered room.

## 5. Technical service

#### ⚠ CAUTION

Modifications of medical equipment may result in the voiding of any guarantee/warranty claims and any approvals.

- Do not modify the product.
- For service and repairs, please contact your national B. Braun/Aesculap agency.

#### Service addresses

Aesculap Technischer Service

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78532 Tuttlingen / Germany

Phone: +49 7461 95-1601

Fax: +49 7461 16-2887

E-Mail: [ats@aesculap.de](mailto:ats@aesculap.de)

Other service addresses can be obtained from the address indicated above.

## 6. Disposal

### WARNING

Risk of infection due to contaminated products!

- ▶ Adhere to national regulations when disposing of or recycling the product, its components and its packaging.

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