en **Instructions for use/Technical description**  
S<sup>4</sup>® percutaneous instruments

USA **Note for U.S. users**  
This Instructions for Use is NOT intended for United States users. Please discard. The Instructions for Use for United States users can be obtained by visiting our website at [www.aesculapImplantsystems.com](http://www.aesculapImplantsystems.com). If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-866-229-3002. A paper copy will be provided to you upon request at no additional cost.

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
## Legend

- 1 Cleaning instrument FW692R
- 2 Clamp sleeve FW693R
- 3 Percutaneous outer sleeve FW735R
- 4 Inner sleeve FW694R
- 5 Inner sleeve FW734R
- 6 Screwdriver FW695R (for polyaxial S<sup>4</sup> bolts)
- 7 Screwdriver FW696R (for monoaxial S<sup>4</sup> bolts)
- 8 Spacer axial FW143P
- 9 Insertion instrument for locking screws FW697R
- 10 Counterholder for percutaneous outer sleeve FW736R

## 1. About this document

### Note

General risk factors associated with surgical procedures are not described in these instructions for use.

### 1.1 Scope

These instructions for use apply for the following products:

Art. no.	Designation
FW692R	Cleaning instrument
FW693R	Clamp sleeve
FW735R	Percutaneous outer sleeve
FW694R	Inner sleeve
FW734R	Inner sleeve
FW695R	Screwdriver for polyaxial S <sup>4</sup> bolts
FW696R	Screwdriver for monoaxial S <sup>4</sup> bolts
FW143P	Spacer axial
FW697R	Insertion instrument for set screws
FW736R	Counterholder for percutaneous outer sleeve

### Note

The applicable CE mark for the product can be found on the label or packaging of the product.

- ▶ For article specific instructions for use and material compatibility and lifetime information, see B. Braun eIFU at [ifu.bbraun.com](http://ifu.bbraun.com)

### 1.2 Safety messages

Safety messages make clear the dangers to patient, user and/or product that could arise during the use of the product. Safety messages are labeled as follows:

#### ⚠ DANGER

Indicates a potential danger. If it is not prevented, death or severe injuries may result.

#### ⚠ WARNING

Indicates a possible threat of danger. If not avoided, minor or moderate injury may result.

#### ⚠ CAUTION

Indicates a possible threat of material damage. If not avoided, the product may be damaged.

## 2. Clinical use

### 2.1 Areas of use and limitations of use

#### 2.1.1 Intended use

S<sup>4</sup> percutaneous instruments are used to implant cannulated S<sup>4</sup> implants of the S<sup>4</sup> Spinal System using a percutaneous operating technique.

#### 2.1.2 Indications

The manufacturer is not responsible for any use of the product against the specified indications and/or the described applications.

For indications, see Intended use.

#### 2.1.3 Contraindications

No known contraindications.

## 2.2 Safety information

### 2.2.1 Clinical user

#### General safety information

To prevent damage caused by improper setup or operation, and to not compromise the manufacturer warranty and liability:

- ▶ Use the product only according to these instructions for use.
- ▶ Follow the safety and maintenance instructions.
- ▶ Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge and experience.
- ▶ Store any new or unused products in a dry, clean, and safe place.
- ▶ Prior to use, check that the product is in good working order.
- ▶ Keep the instructions for use accessible for the user.

#### Note

The user is obligated to report all severe events in connection with the product to the manufacturer and the responsible authorities of the state in which the user is located.

#### Notes on surgical procedures

It is the user's responsibility to ensure that the surgical procedure is performed correctly.

Appropriate clinical training as well as a theoretical and practical proficiency of all the required operating techniques, including the use of this product, are prerequisites for the successful use of this product.

The user is required to obtain information from the manufacturer if there is an unclear preoperative situation regarding the use of the product.

### 2.2.2 Product specific safety information

- ▶ Read, observe and safely store the OR Manuals S<sup>4</sup> Spinal System Instrumentation, S<sup>4</sup> Modular MIS Pedicle Screw System - Surgical technique 069002.

### 2.2.3 Sterility

The product is delivered in an unsterile condition.

- ▶ Clean the new product after removing its transport packaging and prior to its initial sterilization.

## 2.3 Application

### ⚠ DANGER

If the guidewire is pushed into the aorta this can be life-threatening for patients!

- ▶ When inserting S<sup>4</sup> implants or S<sup>4</sup> instruments using the guidewire, observe the markings on the guidewire.
- ▶ Use forceps to hold the guidewire.
- ▶ Remove the guidewire in time.

### ⚠ WARNING

Risk of injury and/or malfunction!

- ▶ Prior to each use, inspect the product for loose, bent, broken, cracked, worn, or fractured components.
- ▶ Always carry out a function test prior to each use of the product.

### 2.3.1 Cleaning instrument FW692R

The cleaning instrument is used to clean/clear the cannulation if a cannulated instrument is blocked intra-operatively by soft tissue or bone chips.

#### ⚠ WARNING

Risk of injury if the cleaning instrument punctures the cannulated instrument!

- ▶ Carefully remove the obstruction by applying a slight pressure.
- ▶ Always insert the cleaning instrument in the direction of the working end.
- ▶ Do not use your hand to close the outlet port.

#### ⚠ CAUTION

Incorrect use of the cleaning instrument may damage cannulated instruments!

- ▶ Always insert the cleaning instrument in the direction of the working end.
- ▶ Check the cleaning instrument for bends, kinks and other damage before use.
- ▶ Disassemble the instrument where possible on a separate sterile table, if applicable see Disassembly
- ▶ Carefully insert the cleaning instrument into the cannulation from the handle side towards the working end.
- ▶ Remove the obstruction by using the handle part to apply slight pressure. Be careful not to obstruct the cannulation.
- ▶ Install the instrument.

### 2.3.2 Clamp sleeve FW693R and screwdriver FW695R/FW696R

Screwdrivers facilitate the percutaneous application of pedicle screws of the S<sup>4</sup> Spinal System. Depending on the type of bolt, the corresponding screwdriver for monoaxial or polyaxial bolts must be chosen.

- ▶ Assemble the clamp sleeve, see Assembly.

#### Insert the pedicle screws in the clamp sleeve

- ▶ Turn clamping bolt C in the direction of the arrow until the thread in the direction of the working end becomes visible, see Fig. 1.
- ▶ Pull the outer sleeve A back as far as possible towards the clamping bolt C, see Fig. 2.
- ▶ Insert screwdriver F into clamp sleeve 2 from the rear end. Make sure the side lugs slide into the grooves all the way to the bottom, see Fig. 3.
- ▶ Insert the working ends of the pedicle screws into the tension sleeve B, see Fig. 4.
- ▶ For polyaxial pedicle screws: Ensure that the hexagon of the screwdriver is fully inserted in the hexagon of the polyaxial pedicle bolt.
- ▶ For monoaxial pedicle bolts: Ensure the dihedral is seated securely in the bolt head.
- ▶ Push the outer sleeve A towards the pedicle bolt and clamp it into the clamp sleeve 2 using the clamping bolt C, see Fig. 5.
- ▶ Check to ensure that the implant is properly seated in the clamp sleeve 2.
- ▶ For implanting the pedicle screw and inserting the rod, see operating technique or operating manual for the S<sup>4</sup> Spinal System, S<sup>4</sup> FRI and S<sup>4</sup> SRI.

#### Remove the clamp sleeve from the inserted pedicle screw

- ▶ Turn clamping bolt C in the direction of the arrow until the thread in the direction of the working end becomes visible, see Fig. 6.
- ▶ Pull back the outer sleeve A towards the clamping bolt C.
- ▶ Remove Complete clamp sleeve 2, see Fig. 7.

### 3. Validated reprocessing procedure

#### 3.1 General safety information

##### Note

Adhere to national statutory regulations, national and international standards and directives, and local, clinical hygiene instructions for sterile processing.

##### Note

For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of products.

##### Note

Mechanical reprocessing should be favored over manual cleaning as it gives better and more reliable results.

##### Note

Successful processing of this medical device can only be ensured if the processing method is first validated. The operator/sterile processing technician is responsible for this.

##### Note

If there is no final sterilization, then a virucidal disinfectant must be used.

##### Note

For up-to-date information about reprocessing and material compatibility, see B. Braun eFU at [eifu.bbraun.com](http://eifu.bbraun.com). The validated steam sterilization procedure was carried out in the Aesculap sterile container system.

#### 3.2 General information

Dried or affixed surgical residues can make cleaning more difficult or ineffective and lead to corrosion. Therefore the time interval between application and processing should not exceed 6 h; also, neither fixing pre-cleaning temperatures >45 °C nor fixing disinfecting agents (active ingredient: aldehydes/alcohols) should be used.

Excessive measures of neutralizing agents or basic cleaners may result in a chemical attack and/or to fading and the laser marking becoming unreadable visually or by machine for stainless steel.

Residues containing chlorine or chlorides e.g. in surgical residues, medicines, saline solutions and in the service water used for cleaning, disinfection and sterilization will cause corrosion damage (pitting, stress corrosion) and result in the destruction of stainless steel products. These must be removed by rinsing thoroughly with demineralized water and then drying.

Additional drying, if necessary.

Only process chemicals that have been tested and approved (e.g. VAH or FDA approval or CE mark) and which are compatible with the product's materials according to the chemical manufacturers' recommendations may be used for processing the product. All the chemical manufacturer's application specifications must be strictly observed. Failure to do so can result in the following problems:

- Optical changes of materials, e.g. fading or discoloration of titanium or aluminum. For aluminum, the application/process solution only needs to be of pH >8 to cause visible surface changes.
- Material damage such as corrosion, cracks, fracturing, premature aging or swelling.
- Do not use metal cleaning brushes or other abrasives that would damage the product surfaces and could cause corrosion.
- Further detailed advice on hygienically safe and material-/value-preserving reprocessing can be found at [www.a-k-i.org](http://www.a-k-i.org), link to "AKI-Brochures", "Red brochure".

#### 3.3 Reusable products

Influences of the reprocessing which lead to damage to the product are not known.

A careful visual and functional inspection before the next use is the best option to determine if a product is no longer functional, see Inspection.

#### 3.4 Preparations at the place of use

- If applicable, rinse non-visible surfaces preferably with deionized water, with a disposable syringe for example.
- Remove any visible surgical residues to the extent possible with a damp, lint-free cloth.
- Transport the dry product in a sealed waste container for cleaning and disinfection within 6 hours.

#### 3.5 Preparing for cleaning

- Do not disassemble anchoring screws/bolts that permanently connect the various components.
- Disassemble the product prior to cleaning, see Disassembly.

#### 3.6 Disassembly

##### 3.6.1 Clamp sleeve FW693R

- Pull the outer sleeve A off over the collet B towards the working end.
- Unscrew clamping bolt C towards the rear end and remove from collet B, see Fig. 8.

##### 3.6.2 Inner sleeve FW694R and FW734R

- Unscrew nut D from inner sleeve E, see Fig. 9.

#### 3.7 Cleaning/Disinfection

##### 3.7.1 Product-specific safety information on the reprocessing method

Hazardous to patients! The product may not be processed manually.

- Process the product exclusively by machine.
- Process the product only by mechanical cleaning preceded by manual pre-cleaning.

Mechanical cleaning/disinfection may cause damage to or destruction of the product!

- Only clean/disinfect the product manually.
- Do not sterilize the product under any circumstances.

Unsuitable cleaning/disinfecting agents may cause damage to or destruction of the product!

- Only use cleaning/disinfecting agents approved for surface cleaning in accordance with the manufacturer's instructions.

Damage to or destruction of the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!

- Use cleaning agents and disinfectants in accordance with the manufacturer's instructions.
- Observe specifications regarding concentration, temperature and exposure time.
- Do not exceed the maximum allowable disinfection temperature of 95 °C.

#### 3.7.2 Validated cleaning and disinfection procedure

Validated procedure	Specific requirements	Reference
Manual cleaning with immersion disinfection	<ul style="list-style-type: none"> <li>Cleaning brush: e.g. TA007747 and/or TE654202</li> <li>Disposable syringe 20ml</li> <li>Keep working ends open for cleaning purposes.</li> <li>Clean products having movable hinges in the open position or while moving the joints.</li> <li>Drying phase: Use a lint-free cloth or medical compressed air</li> </ul>	Chapter Manual cleaning/disinfection and subsection: <ul style="list-style-type: none"> <li>Chapter Manual cleaning with immersion disinfection</li> </ul>
Alkaline machine cleaning and thermal disinfection	<ul style="list-style-type: none"> <li>Spacer FW143P</li> <li>Cleaning instrument FW692R</li> <li>Clamp sleeve FW693R</li> <li>Insertion instrument for set screws FW697R</li> <li>Percutaneous outer sleeve FW735R</li> <li>Counterholder for percutaneous outer sleeve FW736R</li> </ul>	Chapter Machine cleaning/disinfection and subsections: <ul style="list-style-type: none"> <li>Chapter Mechanical alkaline cleaning and thermal disinfection</li> </ul>
Manual pre-cleaning with brush and subsequent mechanical alkaline cleaning and thermal disinfection	<ul style="list-style-type: none"> <li>Cleaning brush: e.g. TA007747 and/or TE654202</li> <li>Disposable syringe 20ml</li> <li>Place product on a sterile screen basket suitable for cleaning (make sure all areas will be reached by water jets).</li> <li>Connect the lumens and channels of individual parts directly to the injector unit's special flushing connector.</li> <li>Keep working ends open for cleaning purposes.</li> <li>Open the links and joints before placing the product on the screen basket.</li> </ul>	Chapter Mechanical cleaning/disinfection with manual pre-cleaning and subsection: <ul style="list-style-type: none"> <li>Chapter Manual pre-cleaning with a brush</li> <li>Chapter Mechanical alkaline cleaning and thermal disinfecting</li> </ul>

#### 3.8 Manual cleaning/disinfection

- Prior to manual disinfecting, allow water to drip off for a sufficient length of time to prevent dilution of the disinfecting solution.
- After manual cleaning/disinfection, check visible surfaces visually for residues.
- Repeat the cleaning/disinfection process if necessary.

##### 3.8.1 Manual cleaning with immersion disinfection

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Disinfecting cleaning	RT (cold)	>15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
II	Intermediate rinse	RT (cold)	1	-	D-W	-
III	Disinfection	RT (cold)	5	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
IV	Final rinse	RT (cold)	1	-	FD-W	-
V	Drying	RT	-	-	-	-

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

RT: Room temperature

\*Recommended: BBraun Stabimed fresh

- Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

##### Phase I

- Fully immerse the product in the cleaning/disinfectant for at least 15 min. Ensure that all accessible surfaces are moistened.
- Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
- Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

##### Phase II

- Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- Drain any remaining water fully.

##### Phase III

- Fully immerse the product in the disinfectant solution.
- Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- Rinse lumens at least 5 times at the beginning of the exposure time using an appropriate disposable syringe. Ensure that all accessible surfaces are moistened.

##### Phase IV

- Rinse/flush the product thoroughly (all accessible surfaces).
- Mobilize non-rigid components, such as set screws, joints, etc. during final rinse.
- Rinse lumens with an appropriate disposable syringe at least five times.
- Drain any remaining water fully.

## Phase V

- ▶ Dry the product in the drying phase with suitable equipment (e.g. cloth, compressed air), see Validated cleaning and disinfection procedure.

### 3.9 Machine cleaning/disinfection

#### Note

The cleaning and disinfection device must be of tested and approved effectiveness (e.g. FDA approval or CE mark according to DIN EN ISO 15883).

#### Note

The cleaning and disinfection device used for processing must be serviced and checked at regular intervals.

#### 3.9.1 Mechanical alkaline cleaning and thermal disinfection

Machine type: single-chamber cleaning/disinfection device without ultrasound

Phase	Step	T [°C/°F]	t [min]	Water quality	Chemical/Note
I	Pre-rinse	<25/77	3	D-W	-
II	Cleaning	55/131	10	FD-W	<ul style="list-style-type: none"><li>■ Concentrate, alkaline:<ul style="list-style-type: none"><li>- pH = 13</li><li>- &lt;5 % anionic surfactant</li></ul></li><li>■ 0.5 % working solution<ul style="list-style-type: none"><li>- pH = 11*</li></ul></li></ul>
III	Intermediate rinse	>10/50	1	FD-W	-
IV	Thermal disinfecting	90/194	5	FD-W	-
V	Drying	-	-	-	According to the program for cleaning and disinfection device

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

\*Recommended: BBraun Helimatic Cleaner alkaline

- ▶ Check visible surfaces for residues after mechanical cleaning/disinfecting.

### 3.10 Mechanical cleaning/disinfection with manual pre-cleaning

#### Note

The cleaning and disinfection device must be of tested and approved effectiveness (e.g. FDA approval or CE mark according to DIN EN ISO 15883).

#### Note

The cleaning and disinfection device used for processing must be serviced and checked at regular intervals.

#### 3.10.1 Manual pre-cleaning with a brush

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Disinfectant cleaning	RT (cold)	>15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
II	Rinsing	RT (cold)	1	-	D-W	-

D-W: Drinking water

RT: Room temperature

\*Recommended: BBraun Stabimed fresh

- ▶ Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

#### Phase I

- ▶ Fully immerse the product in the cleaning/disinfectant for at least 15 min. Ensure that all accessible surfaces are moistened.
- ▶ Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- ▶ If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
- ▶ Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- ▶ Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

#### Phase II

- ▶ Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.

### 3.10.2 Mechanical alkaline cleaning and thermal disinfecting

Machine type: single-chamber cleaning/disinfection device without ultrasound

Phase	Step	T [°C/°F]	t [min]	Water quality	Chemical
I	Prerinse	<25/77	3	D-W	-
II	Cleaning	55/131	10	FD-W	<ul style="list-style-type: none"><li>■ Concentrate, alkaline:<ul style="list-style-type: none"><li>- pH ~ 13</li><li>- &lt;5 % anionic surfactant</li></ul></li><li>■ working solution 0.5%<ul style="list-style-type: none"><li>- pH = 11*</li></ul></li></ul>
III	Intermediate rinse	>10/50	1	FD-W	-
IV	Thermal disinfecting	90/194	5	FD-W	-
V	Drying	-	-	-	In accordance with the program for the cleaning and disinfecting machine

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

\*Recommended: BBraun Helimatic Cleaner alkaline

- ▶ Check visible surfaces for residues after mechanical cleaning/disinfecting.

### 3.11 Inspection

- ▶ Allow the product to cool down to room temperature.
- ▶ Dry the product if it is wet or damp.

#### 3.11.1 Visual inspection

- ▶ Ensure that all soiling has been removed. In particular, pay attention to mating surfaces, hinges, shafts, recessed areas, drill grooves and the sides of the teeth on rasps.
- ▶ If the product is dirty: repeat the cleaning and disinfection process.
- ▶ Check the product for damage, e.g. insulation or corroded, loose, bent, broken, cracked, worn or severely scratched and fractured components.
- ▶ Check the product for missing or faded labels.
- ▶ Check products having long, narrow geometries (in particular rotating instruments) for deformations.
- ▶ Check the cutting edges for continuity, sharpness, nicks and other damage.
- ▶ Check the surfaces for rough spots.
- ▶ Check the product for burrs that could damage tissue or surgical gloves.
- ▶ Check the product for loose or missing parts.
- ▶ Immediately put aside damaged or inoperative products and send them to Aesculap Technical Service, see Technical service.

#### 3.11.2 Functional test

##### ⚠ CAUTION

The product may become damaged (metal pitting/fretting corrosion) if not sufficiently oiled!

- ▶ Oil any moving parts (e.g. joints, spool parts and threaded rods) prior to the functional test using oil suitable for the applied sterilization process (e.g. STERILIT® I oil spray JG600 for steam sterilization or STERILIT® I drip-feed lubricator JG598).
- ▶ Assemble disassembled products, see Assembly.
- ▶ Check that the product functions correctly.
- ▶ Check that all moving parts are working properly (e.g. hinges, locks/latches, sliding parts etc.).
- ▶ Check rotating products (e.g. reusable drills and cutters) for bends and deformities. To do this, roll product, for example, on a flat surface.
- ▶ Check for compatibility with associated products.
- ▶ Immediately put aside inoperative products and send them to Aesculap Technical Service, see Technical service.

### 3.12 Assembly

#### 3.12.1 Clamp sleeve FW693R

- ▶ Screw the clamping bolt C onto the collet B until the thread of the collet is just covered.
- ▶ Slide outer sleeve A over the collet B from the front, see Fig. 10. Make sure the slots are aligned with the T-nuts. The ribbed end of the clamp sleeve should point towards the clamping bolt C for ease of use.

#### 3.12.2 Inner sleeve FW694R and FW734R

- ▶ Bolt nut D onto inner sleeve E, see Fig. 11.

### 3.13 Packaging

- ▶ Place the product in its holder or on a suitable tray. Ensure that sharp edges are covered.
- ▶ Package trays appropriately for the sterilization process (e.g. in Aesculap sterile containers).
- ▶ Ensure that the packaging provides sufficient protection against contamination of the product during storage.

### 3.14 Steam sterilization

#### Note

The product may only be sterilized in a disassembled state.

- ▶ Check to ensure that the sterilizing agent will come into contact with all external and internal surfaces (e.g., by opening any valves and faucets).
- ▶ Validated sterilization process
  - Disassemble the product
  - Steam sterilization in a fractionated vacuum process
  - Steam sterilizer in accordance with EN285 and validated in accordance with EN ISO17665
  - Sterilization in the fractionated vacuum process at 134 °C, holding time 5 min
- ▶ If several devices are sterilized at the same time in the same steam sterilizer: Ensure that the maximum permitted load according to the manufacturers' specifications is not exceeded.

### 3.15 Storage

- ▶ Store sterile products in germ-proof packaging, protected from dust, in a dry, dark, temperature-controlled area.

## 4. Technical service

### CAUTION

Modifications to medical equipment may result in the voiding of any guarantee/warranty claims and any approvals.

- ▶ Do not modify the product.
- ▶ For service and repairs, please contact your national B. Braun/Aesculap agency.

### Service addresses

Aesculap Technischer Service

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78532 Tuttlingen / Germany

Phone: +49 7461 95-1601

Fax: +49 7461 16-2887

E-Mail: [ats@aesculap.de](mailto:ats@aesculap.de)

Other service addresses can be obtained from the address indicated above.

## 5. Disposal

### WARNING

Risk of infection due to contaminated products!

- ▶ Adhere to national regulations when disposing of or recycling the product, its components and its packaging.

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