
# Aesculap® SecureSpan

## Aesculap Spine

en	<b>Instructions for use/Technical description</b> Instruments	lv	<b>Lietošanas instrukcijas/tehniskais apraksts</b> Instrumenti
USA	<b>Note for U.S. users</b> This Instructions for Use is NOT intended for United States users. Please discard. The Instructions for Use for United States users can be obtained by visiting our website at <a href="http://www.aesculapimplantsystems.com">www.aesculapimplantsystems.com</a> . If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-866-229-3002. A paper copy will be provided to you upon request at no additional cost.	lt	<b>Naudojimo instrukcija/techninis aprašas</b> Instrumentai
de	<b>Gebrauchsanweisung/Technische Beschreibung</b> Instrumente	ru	<b>Инструкция по применению/Техническое описание</b> Инструменты
fr	<b>Mode d'emploi/Description technique</b> Instruments	cs	<b>Návod k použití/Technický popis</b> Nástroje
es	<b>Instrucciones de manejo/Descripción técnica</b> Instrumental	pl	<b>Instrukcja użytkowania/Opis techniczny</b> Instrumenty
it	<b>Istruzioni per l'uso/Descrizione tecnica</b> Strumenti	sk	<b>Návod na použitie/Technický opis</b> Nástroje
pt	<b>Instruções de utilização/Descrição técnica</b> Instrumentos	hu	<b>Használati útmutató/Műszaki leírás</b> Műszerek
nl	<b>Gebruiksaanwijzing/Technische beschrijving</b> Instrumenten	sl	<b>Navodila za uporabo/Tehnični opis</b> Instrumenti
da	<b>Brugsanvisning/Teknisk beskrivelse</b> Instrumenter	hr	<b>Upute za uporabu/Tehnički opis</b> Instrumenti
sv	<b>Bruksanvisning/Teknisk beskrivning</b> Instrument	ro	<b>Manual de utilizare/Descriere tehnică</b> Instrumente
fi	<b>Käyttöohje/Tekninen kuvaus</b> Instrumentit	bg	<b>Упътване за употреба/Техническо описание</b> Инструменти
		tr	<b>Kullanım Kılavuzu/Teknik açıklama</b> Aletler
		el	<b>Οδηγίες χρήσης/Τεχνική περιγραφή</b> Όργανα


**A**


**B**

**C**

**D**


**E**


**F**


**G**

**H**


I


J

K

L

M

# Aesculap® SecureSpan

## Instruments

### Legend

A	Trial implants FG836R-FG840R
B	Lamina lifter FG847R
C	Locking screwdriver FG848R
D	Implant holder FG849R
E	Plate holder FG851R
F	Drill handle FG861R
G	Depth gauge FG866R
H	Bone awl FG867R
I	Adjustable single drill guide FG868R
J	Fixed depth drills FG856SU to FG860SU & FG842SU
K	Drill FG869SU for use with drill guide
L	Plate bending pliers LX185R
M	Plate cutters LX159R

### Contents

1.	About this document.	2
1.1	Scope	2
1.2	Safety messages.	2
2.	Clinical use.	3
2.1	Areas of use and limitations of use	3
2.1.1	Intended use.	3
2.1.2	Indications	3
2.1.3	Contraindications.	3
2.2	Safety information.	3
2.2.1	Clinical user	3
2.2.2	Product	3
2.2.3	Sterility	3
2.3	Application	4
2.3.1	Selecting the right plate length for the SecureSpan plate.	4
2.3.2	Adapting the SecureSpan plate	4
2.3.3	Preparing the holes for the SecureSpan screws	4
2.3.4	Inserting the screws	4
2.3.5	Removing SecureSpan implants.	5
3.	Validated reprocessing procedure	5
3.1	General safety instructions.	5
3.2	General information	5
3.3	Single-use products.	5
3.4	Reusable products	6
3.5	Preparations at the place of use.	6
3.6	Preparing for cleaning.	6
3.7	Disassembly	6
3.7.1	Single drill guide FG868R	6
3.7.2	Depth gauge FG866R.	6
3.7.3	Screwdriver FG848R	6
3.8	Cleaning/Disinfection	7
3.8.1	Product-specific safety information on the reprocessing method	7
3.8.2	Validated cleaning and disinfection procedure.	7

3.9	Manual cleaning/disinfection	8
3.9.1	Manual cleaning with immersion disinfection.	8
3.9.2	Manual cleaning with ultrasound and immersion disinfection	9
3.10	Mechanical cleaning/disinfection	10
3.10.1	Mechanical alkaline cleaning and thermal disinfecting	10
3.11	Mechanical cleaning/disinfection with manual pre-cleaning.	11
3.11.1	Manual pre-cleaning with a brush	11
3.11.2	Manual pre-cleaning with ultrasound and brush	11
3.11.3	Mechanical alkaline cleaning and thermal disinfecting	12
3.12	Inspection	12
3.12.1	Visual inspection.	12
3.12.2	Functional test	12
3.13	Assembly	13
3.13.1	Single drill guide FG868R.	13
3.13.2	Depth gauge FG866R	13
3.13.3	Screwdriver FG848R	13
3.14	Packaging	13
3.15	Steam sterilization	13
3.16	Storage.	13
4.	Technical service.	13
5.	Disposal	13

## 1. About this document

### Note

General risk factors associated with surgical procedures are not described in these instructions for use.

### 1.1 Scope

These instructions for use apply for the products listed in the Legend above.

### Note

The applicable CE mark for the product can be seen on the label or packaging of the product.

- For article-specific instructions for use as well as information on material compatibility and lifetime see B. Braun eIFU at [eifu.bbraun.com](http://eifu.bbraun.com)

### 1.2 Safety messages

Safety messages make clear the dangers to patient, user and/or product that could arise during the use of the product. Safety messages are labeled as follows:

#### WARNING

Indicates a possible threat of danger. If not avoided, minor or moderate injury may result.

#### CAUTION

Indicates a possible threat of material damage. If not avoided, the product may be damaged.

## 2. Clinical use

### 2.1 Areas of use and limitations of use

#### 2.1.1 Intended use

The SecureSpan Laminoplasty Plating system includes the instruments listed in the Legend above.

These instruments are used to place the plates and screws to stabilize the lamina after a laminectomy procedure has been performed.

#### 2.1.2 Indications

##### Note

*The manufacturer is not responsible for any use of the product against the specified indications and/or the described applications.*

For indications, see Intended use.

#### 2.1.3 Contraindications

No known contraindications.

## 2.2 Safety information

### 2.2.1 Clinical user

#### General safety information

To prevent damage caused by improper setup or operation, and to not compromise the manufacturer warranty and liability:

- ▶ Use the product only according to these instructions for use.
- ▶ Follow the safety and maintenance instructions.
- ▶ Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge and experience.
- ▶ Store any new or unused products in a dry, clean, and safe place.
- ▶ Prior to use, check that the product is in good working order.
- ▶ Keep the instructions for use accessible for the user.

##### Note

*The user is obligated to report all severe events in connection with the product to the manufacturer and the responsible authorities of the state in which the user is located.*

#### Notes on surgical procedures

It is the user's responsibility to ensure that the surgical procedure is performed correctly.

Appropriate clinical training as well as a theoretical and practical proficiency of all the required operating techniques, including the use of this product, are prerequisites for the successful use of this product.

The user is required to obtain information from the manufacturer if there is an unclear preoperative situation regarding the use of the product.

### 2.2.2 Product

#### Product-specific safety information

The operating surgeon must be thoroughly conversant with stabilization techniques for and the biomechanical properties of the cervical spine. He must also be familiar with the theory and practice of operative techniques within the SecureSpan system.

- ▶ Always follow the instructions for use of the SecureSpan plates and screws (TA013082) and the surgical technique guide.
- ▶ To avoid damage to the working end: Carefully insert the product through the working channel (e.g. drill).

### 2.2.3 Sterility

#### Non-sterile packaged products

Art. no.	Designation
FG836R to FG840R	Trial implants
FG847R	Lamina lifter
FG848R	Locking screwdriver
FG849R	Implant holder
FG851R	Plate holder
FG861R	Drill handle
FG866R	Depth gauge
FG867R	Bone awl
FG868R	Adjustable single drill guide
LX185R	Plate bending pliers
LX159R	Plate cutters

The product is delivered in an unsterile condition.

- ▶ Clean the new product after removing its transport packaging and prior to its initial sterilization.

#### Sterile products

Art. no.	Designation
FG856SU to FG860SU FG842SU	Fixed depth drills
FG869SU	Drill for use with drill guide

The product has been sterilized by irradiation and is supplied in sterile packaging.

- ▶ Do not use products from open or damaged sterile packaging.
- ▶ Do not use the product after its use-by date.
- ▶ Do not reuse the product.

The reprocessing of the product affects its functionality. Risk of injury, illness or death due to soiling and/or impaired functionality of the product.

- ▶ Do not reprocess the product.

## 2.3 Application

### ⚠ WARNING

Risk of injury and/or malfunction!

- ▶ Prior to each use, inspect the product for loose, bent, broken, cracked, worn, or fractured components.
- ▶ Always carry out a function test prior to each use of the product.

#### 2.3.1 Selecting the right plate length for the SecureSpan plate

- ▶ Use the Trial implants **A** to select the appropriate plate and spacer size.

#### 2.3.2 Adapting the SecureSpan plate

The SecureSpan plates supplied by Aesculap are supplied as pre-bent implants. They can be bent by means of the plate benders. Further bending should be done at the existing bends.

### ⚠ CAUTION

Damage to the SecureSpan implant through excessive strain on the material!

- ▶ Always bend the SecureSpan plate in one direction only.
- ▶ Do not bend back the SecureSpan plate.
- ▶ Avoid small bending radii, bending back, notching or scratching the SecureSpan plate.

#### 2.3.3 Preparing the holes for the SecureSpan screws

The holes for the self-tapping SecureSpan screws are either created using the awl **H**, or the fixed depth drills **J**, or the drill **K** and drill guide **I**. The drills can be driven manually with drill handle **F** or with a motor system and Aesculap Intra handpiece (e.g., GD450R/GD456R).

### ⚠ CAUTION

Risk of drilling into the patient's spinal cord as a result of the cutting depth set incorrectly on the drill guide!

- ▶ Insert the drill into the drill guide and check the set cutting depth with a caliper or ruler.


Fig. 1 Drill attached to drill handle through drill guide

### Note

The maximum center punch depth with the awl **H** is 3 mm. Do not use the awl with the drill guide.

- ▶ To adjust the required depth with the depth-adjustable drill guide **I** (0 mm to 14 mm), turn guide sleeve **1** of the drill guide in direction "+" or "-". Note that the thread of the drill guide is a left-handed type. Every half-turn (= height adjustment by 0.5 mm) you will hear and feel the guide sleeve clicking into position.

#### 2.3.4 Inserting the screws

### ⚠ CAUTION

Damage to the SecureSpan plate or SecureSpan screw in the course of driving the screw into the bone through the plate, caused by incorrectly applied screwdriver!

- ▶ Fully insert the tip of the screwdriver into the cross slot of the screw.
- ▶ Lock the screwdriver to the screw through clockwise turning of the thumbwheel.

### ⚠ CAUTION

Damage to the locking thread of the SecureSpan screw, or the screwdriver tip caused by overtightening the locking thread!

- ▶ Ensure that the SecureSpan screw is firmly attached to the screwdriver.

### Note

The locking screwdriver **C** has a self-holding feature. A screw holding sleeve is not necessary.


Fig. 2

- ▶ Fully insert the tip **2** of the locking screwdriver **C** into the screw head **3**.
- ▶ Tighten the thumbwheel **4** clockwise until hand-tight.
- ▶ Apply gentle pressure on the locking screwdriver **C** while driving the SecureSpan screw into bone through the hole of the SecureSpan plate until seated.
- ▶ Detach the locking screwdriver **C** from the screw head by turning the thumbwheel **4** counterclockwise.

### 2.3.5 Removing SecureSpan implants

- ▶ Always use SecureSpan instruments to remove SecureSpan implants.

#### Note

The screwdriver **C** can also be used for screw removal.


Fig. 3 Using the depth gauge FG866R

- ▶ Using the depth gauge **G** with a SecureSpan plate **5** for the measuring example **6** of a 10 mm SecureSpan screw:
  - Loosen knurled nut **7**.
  - Slide caliper **8** of depth gauge **G** in the required direction (e.g., 10 mm, see Fig. 3).
  - Tighten knurled nut **7**.

The depth gauge is fixed at the measured depth of 10 mm.

## 3. Validated reprocessing procedure

### 3.1 General safety instructions

#### Note

Adhere to national statutory regulations, national and international standards and directives, and local, clinical hygiene instructions for sterile processing.

#### Note

For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of products.

#### Note

Mechanical reprocessing should be favored over manual cleaning as it gives better and more reliable results.

#### Note

Successful processing of this medical device can only be ensured if the processing method is first validated. The operator/sterile processing technician is responsible for this.

#### Note

If there is no final sterilization, then a virucidal disinfectant must be used.

#### Note

For up-to-date information about reprocessing and material compatibility, see B. Braun eIFU at [eifu.bbraun.com](http://eifu.bbraun.com)

The validated steam sterilization procedure was carried out in the Aesculap sterile container system.

### 3.2 General information

Dried or affixed surgical residues can make cleaning more difficult or ineffective and lead to corrosion. Therefore the time interval between application and processing should not exceed 6 h; also, neither fixating pre-cleaning temperatures >45 °C nor fixating disinfecting agents (active ingredient: aldehydes/alcohols) should be used.

Excessive measures of neutralizing agents or basic cleaners may result in a chemical attack and/or to fading and the laser marking becoming unreadable visually or by machine for stainless steel.

Residues containing chlorine or chlorides e.g. in surgical residues, medicines, saline solutions and in the service water used for cleaning, disinfection and sterilization will cause corrosion damage (pitting, stress corrosion) and result in the destruction of stainless steel products. These must be removed by rinsing thoroughly with demineralized water and then drying.

Additional drying, if necessary.

Only process chemicals that have been tested and approved (e.g. VAH or FDA approval or CE mark) and which are compatible with the product's materials according to the chemical manufacturers' recommendations may be used for processing the product. All the chemical manufacturer's application specifications must be strictly observed. Failure to do so can result in the following problems:

- Optical changes of materials, e.g. fading or discoloration of titanium or aluminum. For aluminum, the application/process solution only needs to be of pH >8 to cause visible surface changes.
- Material damage such as corrosion, cracks, fracturing, premature aging or swelling.
- ▶ Do not use metal cleaning brushes or other abrasives that would damage the product surfaces and could cause corrosion.
- ▶ Further detailed advice on hygienically safe and material-/value-preserving reprocessing can be found at [www.a-k-i.org](http://www.a-k-i.org), link to "AKI-Brochures", "Red brochure".

### 3.3 Single-use products

Art. no.	Designation
FG856SU to FG860SU FG842SU	Fixed depth drills
FG869SU	Drill for use with drill guide

- ▶ Do not reuse the product.

The reprocessing of the product affects its functionality. Risk of injury, illness or death due to soiling and/or impaired functionality of the product.

- ▶ Do not reprocess the product.

### 3.4 Reusable products

Art. no.	Designation
FG836R to FG840R	Trial implants
FG847R	Lamina lifter
FG848R	Locking screwdriver
FG849R	Implant holder
FG851R	Plate holder
FG861R	Drill handle
FG866R	Depth gauge
FG867R	Bone awl
FG868R	Adjustable single drill guide
LX185R	Plate bending pliers
LX159R	Plate cutters

Influences of the reprocessing which lead to damage to the product are not known.

A careful visual and functional inspection before the next use is the best opportunity to recognize a product that is no longer functional, see Inspection.

### 3.5 Preparations at the place of use

- ▶ If applicable, rinse non-visible surfaces preferably with deionized water, with a disposable syringe for example.
- ▶ Remove any visible surgical residues to the extent possible with a damp, lint-free cloth.
- ▶ Transport the dry product in a sealed waste container for cleaning and disinfection within 6 hours.

### 3.6 Preparing for cleaning

- ▶ Disassemble the product prior to cleaning, see Disassembly.

### 3.7 Disassembly

#### 3.7.1 Single drill guide FG868R


Fig. 4

- ▶ Remove guide sleeve 1 by turning it clockwise. Be aware that it is a left-hand thread. You will hear and feel the guide sleeve clicking into position every half turn.

#### 3.7.2 Depth gauge FG866R


Fig. 5

- ▶ Loosen knurled nut 7.
- ▶ Extract caliper 8 in the direction of the arrow.
- ▶ Unscrew and remove knurled nut 7.

#### 3.7.3 Screwdriver FG848R


Fig. 6

- ▶ Unscrew locking nut 9 of locking screwdriver C clockwise.
- ▶ Slide outer shaft 10 of the distal tip of locking screwdriver C.

### 3.8 Cleaning/Disinfection

#### 3.8.1 Product-specific safety information on the reprocessing method

Damage to or destruction of the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!

- ▶ Following the manufacturer's instructions, use cleaning and disinfecting agents
  - that are approved for high-grade steel,
  - that do not attack softeners (e.g., in silicone).

- ▶ Observe specifications regarding concentration, temperature and exposure time.
- ▶ Do not exceed the maximum allowable disinfection temperature of 95 °C.
- ▶ Use suitable cleaning/disinfecting agents if the product is disposed of in a wet condition. To prevent foaming and degradation of the efficacy of the process chemicals: prior to mechanical cleaning and disinfection, rinse the product thoroughly with running water

#### 3.8.2 Validated cleaning and disinfection procedure

Validated procedure	Specific requirements	Reference
Manual cleaning with immersion disinfection ■ FG836R to FG840R ■ FG847R to FG848R ■ FG866R to FG868R	■ Suitable cleaning brush: e.g., TA006874 and TE654202 ■ Disposable syringe 20 ml ■ Drying phase: Use a lint-free cloth or medical compressed air	Chapter Manual cleaning/disinfection and subsection: ■ Chapter Manual cleaning with immersion disinfection
Manual cleaning with ultrasound and immersion disinfection ■ FG849R ■ FG851R ■ FG861R ■ LX159R ■ LX185R	■ Suitable cleaning brush: e.g., TA006874 and TE654202 ■ Disposable syringe 20 ml ■ Drying phase: Use a lint-free cloth or medical compressed air	Chapter Manual cleaning/disinfection and subsection: ■ Chapter Manual cleaning with ultrasound and immersion disinfection
Mechanical alkaline cleaning and thermal disinfection ■ FG836R to FG840R ■ FG847R ■ FG867R	■ Place the product in a tray that is suitable for cleaning (avoiding rinsing blind spots). ■ Connect components with lumens and channels directly to the rinsing port of the injector carriage. ■ Place products in the tray with their hinges open.	Chapter Mechanical cleaning/disinfection and subsection: ■ Chapter Mechanical alkaline cleaning and thermal disinfecting
Manual pre-cleaning with brush and subsequent mechanical alkaline cleaning and thermal disinfection ■ FG849R ■ FG848R ■ FG851R ■ FG866R ■ FG868R	■ Suitable cleaning brush ■ Disposable syringe 20 ml ■ Place the product in a tray that is suitable for cleaning (avoiding rinsing blind spots). ■ Connect components with lumens and channels directly to the rinsing port of the injector carriage.	Chapter Mechanical cleaning/disinfection with manual pre-cleaning and subsection: ■ Chapter Manual pre-cleaning with a brush ■ Chapter Mechanical alkaline cleaning and thermal disinfecting
Manual pre-cleaning with ultrasound and brush, and subsequent mechanical alkaline cleaning and thermal disinfection ■ FG861R ■ LX159R ■ LX185R	■ Suitable cleaning brush ■ Disposable syringe 20 ml ■ Place the product in a tray that is suitable for cleaning (avoiding rinsing blind spots). ■ Connect components with lumens and channels directly to the rinsing port of the injector carriage.	Chapter Mechanical cleaning/disinfection with manual pre-cleaning and subsection: ■ Chapter Manual pre-cleaning with ultrasound and brush ■ Chapter Mechanical alkaline cleaning and thermal disinfecting

### 3.9 Manual cleaning/disinfection

- ▶ Prior to manual disinfecting, allow water to drip off for a sufficient length of time to prevent dilution of the disinfecting solution.
- ▶ After manual cleaning/disinfection, check visible surfaces visually for residues.
- ▶ Repeat the cleaning/disinfection process if necessary.

#### 3.9.1 Manual cleaning with immersion disinfection

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	<b>Disinfecting cleaning</b>	RT (cold)	>15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
II	<b>Intermediate rinse</b>	RT (cold)	1	-	D-W	-
III	<b>Disinfection</b>	RT (cold)	5	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
IV	<b>Final rinse</b>	RT (cold)	1	-	FD-W	-
V	<b>Drying</b>	RT	-	-	-	-

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

RT: Room temperature

\*Recommended: BBraun Stabimed fresh

- ▶ Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

#### Phase I

- ▶ Fully immerse the product in the cleaning/disinfectant for at least 15 min. Ensure that all accessible surfaces are moistened.
- ▶ Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- ▶ If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
- ▶ Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- ▶ Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

#### Phase II

- ▶ Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- ▶ Drain any remaining water fully.

#### Phase III

- ▶ Fully immerse the product in the disinfectant solution.
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- ▶ Rinse lumens at least 5 times at the beginning of the exposure time using an appropriate disposable syringe. Ensure that all accessible surfaces are moistened.

#### Phase IV

- ▶ Rinse/flush the product thoroughly (all accessible surfaces).
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during final rinse.
- ▶ Rinse lumens with an appropriate disposable syringe at least five times.
- ▶ Drain any remaining water fully.

#### Phase V

- ▶ Dry the product in the drying phase with suitable equipment (e.g. cloth, compressed air), see Validated cleaning and disinfection procedure.

### 3.9.2 Manual cleaning with ultrasound and immersion disinfection

Phase	Step	D [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	<b>Ultrasonic cleaning</b>	RT (cold)	>15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
II	<b>Intermediate rinse</b>	RT (cold)	1	-	D-W	-
III	<b>Disinfection</b>	RT (cold)	5	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
IV	<b>Final rinse</b>	RT (cold)	1	-	FD-W	-
V	<b>Drying</b>	RT	-	-	-	-

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

RT: Room temperature

\*Recommended: BBraun Stabimed fresh

- ▶ Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

#### Phase I

- ▶ Clean the product in an ultrasonic cleaning bath (frequency 35 kHz) for at least 15 min. Ensure that all accessible surfaces are immersed and acoustic shadows are avoided.
- ▶ Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- ▶ If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
- ▶ Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- ▶ Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

#### Phase II

- ▶ Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- ▶ Drain any remaining water fully.

#### Phase III

- ▶ Fully immerse the product in the disinfectant solution.
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- ▶ Rinse lumens at least five times at the beginning of the exposure time with an appropriate disposable syringe. Ensure that all accessible surfaces are moistened.

#### Phase IV

- ▶ Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during final rinse.
- ▶ Rinse lumens with an appropriate disposable syringe at least five times.
- ▶ Drain any remaining water fully.

#### Phase V

- ▶ Dry the product in the drying phase with suitable equipment (e.g. cloth, compressed air), see Validated cleaning and disinfection procedure.

### 3.10 Mechanical cleaning/disinfection

#### Note

The cleaning and disinfection device must be of tested and approved effectiveness (e.g. FDA approval or CE mark according to DIN EN ISO 15883).

#### Note

The cleaning and disinfection device used for processing must be serviced and checked at regular intervals.

#### 3.10.1 Mechanical alkaline cleaning and thermal disinfecting

Machine type: single-chamber cleaning/disinfection device without ultrasound

Phase	Step	T [°C/°F]	t [min]	Water quality	Chemistry/Note
I	Pre-rinse	<25/77	3	D-W	-
II	Cleaning	55/131	10	FD-W	<ul style="list-style-type: none"> <li>■ Concentrate, alkaline:               <ul style="list-style-type: none"> <li>- pH = 13</li> <li>- &lt;5 % anionic surfactant</li> </ul> </li> <li>■ 0.5 % working solution               <ul style="list-style-type: none"> <li>- pH ~ 11*</li> </ul> </li> </ul>
III	Intermediate rinse	>10/50	1	FD-W	-
IV	Thermal disinfection	90/194	5	FD-W	-
V	Drying	-	-	-	According to the program for cleaning and disinfection device

DW: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

\*Recommended: BBraun Helimatic Cleaner alkaline

- ▶ Check visible surfaces for residues after mechanical cleaning/disinfecting.

### 3.11 Mechanical cleaning/disinfection with manual pre-cleaning

#### Note

The cleaning and disinfection device must be of tested and approved effectiveness (e.g. FDA approval or CE mark according to DIN EN ISO 15883).

#### Note

The cleaning and disinfection device used for processing must be serviced and checked at regular intervals.

#### 3.11.1 Manual pre-cleaning with a brush

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Disinfecting cleaning	RT (cold)	>15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
II	Rinsing	RT (cold)	1	-	D-W	-

D-W: Drinking water

RT: Room temperature

\*Recommended: BBraun Stabimed fresh

- ▶ Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

#### Phase I

- ▶ Fully immerse the product in the cleaning/disinfectant for at least 15 min. Ensure that all accessible surfaces are moistened.
- ▶ Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- ▶ If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.

- ▶ Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- ▶ Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

#### Phase II

- ▶ Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.

#### 3.11.2 Manual pre-cleaning with ultrasound and brush

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Ultrasonic cleaning	RT (cold)	>15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
II	Rinsing	RT (cold)	1	-	D-W	-

D-W: Drinking water

RT: Room temperature

\*Recommended: BBraun Stabimed fresh

- ▶ Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

#### Phase I

- ▶ Clean the product in an ultrasonic cleaning bath (frequency 35 kHz) for at least 15 min. Ensure that all accessible surfaces are immersed and acoustic shadows are avoided.
- ▶ Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- ▶ If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
- ▶ Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- ▶ Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

#### Phase II

- ▶ Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.

### 3.11.3 Mechanical alkaline cleaning and thermal disinfecting

Machine type: single-chamber cleaning/disinfecting machine without ultrasound

Phase	Step	T [°C/°F]	t [min]	Water quality	Chemical
I	Prerinse	<25/77	3	D-W	-
II	Cleaning	55/131	10	FD-W	<ul style="list-style-type: none"> <li>■ Concentrate, alkaline:               <ul style="list-style-type: none"> <li>- pH ~ 13</li> <li>- &lt;5 % anionic surfactant</li> </ul> </li> <li>■ 0.5 % working solution               <ul style="list-style-type: none"> <li>- pH = 11*</li> </ul> </li> </ul>
III	Intermediate rinse	>10/50	1	FD-W	-
IV	Thermal disinfection	90/194	5	FD-W	-
V	Drying	-	-	-	According to the program for cleaning and disinfection device

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

\*Recommended: BBraun Helimatic Cleaner alkaline

- ▶ Check visible surfaces for residues after mechanical cleaning/disinfecting.

## 3.12 Inspection

- ▶ Allow the product to cool down to room temperature.
- ▶ Dry the product if it is wet or damp.

### 3.12.1 Visual inspection

- ▶ Ensure that all soiling has been removed. In particular, pay attention to mating surfaces, hinges, shafts, recessed areas, drill grooves and the sides of the teeth on rasps.
- ▶ If the product is dirty: repeat the cleaning and disinfection process.
- ▶ Check the product for damage, e.g. insulation or corroded, loose, bent, broken, cracked, worn or severely scratched and fractured components.
- ▶ Check the product for missing or faded labels.
- ▶ Check the products with long, slim shapes (in particular rotating instruments) for deformities.
- ▶ Check the cutting edges for continuity, sharpness, nicks and other damage.
- ▶ Check the surfaces for rough spots.
- ▶ Check the product for burrs that could damage tissue or surgical gloves.
- ▶ Check the product for loose or missing parts.
- ▶ Immediately put aside damaged or inoperative products and send them to Aesculap Technical Service, see Technical service.

### 3.12.2 Functional test

#### CAUTION

**Damage (metal cold welding/friction corrosion) to the product caused by insufficient lubrication!**

- ▶ **Prior to function checks, lubricate moving parts (e.g. joints, pusher components and threaded rods) with maintenance oil suitable for the respective sterilization process (e.g. for steam sterilization: STERILIT® I oil spray JG600 or STERILIT® I drip lubricator JG598).**
- ▶ Assemble disassembled products, see Assembly.
- ▶ Check that the product functions correctly.
- ▶ Check that all moving parts are working properly (e.g. hinges, locks/latches, sliding parts etc.).
- ▶ Check rotating products (e.g. reusable drills and cutters) for bends and deformities. To do this, roll the product on an even surface.
- ▶ Check for compatibility with associated products.
- ▶ Immediately put aside inoperative products and send them to Aesculap Technical Service, see Technical service.

### 3.13 Assembly

#### 3.13.1 Single drill guide FG868R

- ▶ Insert guide sleeve **1** by turning it counterclockwise. Be aware that it is a left-hand thread.  
You will hear and feel the guide sleeve clicking into position every half turn.

#### 3.13.2 Depth gauge FG866R

- ▶ Slide knurled nut **7** over the tip of caliper **8**.
- ▶ Place caliper **8** and knurled nut **7** assembly in outer shaft **6**.
- ▶ Tighten knurled nut **7**.

#### 3.13.3 Screwdriver FG848R

- ▶ Slide the outer shaft **10** over the distal tip of the inner shaft.
- ▶ Tighten locking nut **9**.

### 3.14 Packaging

- ▶ Appropriately protect products with fine working tips.
- ▶ Place the product in its holder or on a suitable tray. Ensure that sharp edges are covered.
- ▶ Package trays appropriately for the sterilization process (e.g. in Aesculap sterile containers).
- ▶ Ensure that the packaging provides sufficient protection against contamination of the product during storage.

### 3.15 Steam sterilization

#### Note

*The product can be sterilized either in disassembled or in assembled condition.*

- ▶ Check to ensure that the sterilizing agent will come into contact with all external and internal surfaces (e.g., by opening any valves and faucets).
- ▶ Validated sterilization process
  - Steam sterilization using fractional vacuum process
  - Steam sterilizer according to DIN EN 285 and validated according to DIN EN ISO 17665
  - Sterilization using fractionated vacuum process at 134 °C/holding time 5 min
- ▶ If several devices are sterilized at the same time in the same steam sterilizer: Ensure that the maximum permitted load according to the manufacturers' specifications is not exceeded.

### 3.16 Storage

- ▶ Store sterile products in germ-proof packaging, protected from dust, in a dry, dark, temperature-controlled area.
- ▶ Store sterile single-use products in germ-proof packaging in a dust-protected, dry, dark and temperature-controlled room.

## 4. Technical service

### ⚠ CAUTION

**Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.**

- ▶ **Do not modify the product.**
- ▶ **For service and repairs, please contact your national B. Braun/Aesculap agency.**

#### Service addresses

Aesculap Technischer Service

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78532 Tuttlingen / Germany

Phone: +49 7461 95-1601

Fax: +49 7461 16-2887

E-Mail: [ats@aesculap.de](mailto:ats@aesculap.de)

Other service addresses can be obtained from the address indicated above.

## 5. Disposal

### ⚠ WARNING

**Risk of infection due to contaminated products!**

- ▶ **Adhere to national regulations when disposing of or recycling the product, its components and its packaging.**

### ⚠ WARNING

**Risk of injury due to sharp-edged and/or pointed products!**

- ▶ **When disposing of or recycling the product, ensure that the packaging prevents injury by the product.**

#### Note

*The user institution is obliged to reprocess the product before its disposal, see Validated reprocessing procedure.*

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