

Aesculap Spine


- en **Instructions for use/Technical description**
activC® handle, reamer guide, chisel guide, trial implants
- USA **Note for U.S. users**
This Instructions for Use is NOT intended for United States users. Please discard. The Instructions for Use for United States users can be obtained by visiting our website at www.aesculapimplantsystems.com. If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-866-229-3002. A paper copy will be provided to you upon request at no additional cost.
- de **Gebrauchsanweisung/Technische Beschreibung**
activC® Handgriff, Fräserführung, Meißelführung, Probeimplantate
- fr **Mode d'emploi/Description technique**
Manche, guide pour fraise, guide pour burin, implants d'essai activC®
- es **Instrucciones de manejo/Descripción técnica**
Mango, guía de escariador, guía de cincel e implantes de prueba activC®
- it **Istruzioni per l'uso/Descrizione tecnica**
Impugnatura activC®, guida per la fresa, guida per lo scalpello, impianti di prova
- pt **Instruções de utilização/Descrição técnica**
Punho, guia de fresa, guia de cincel, implantes de ensaio activC®
- nl **Gebruiksaanwijzing/Technische beschrijving**
activC® handgreep, freesgeleider, beitelgeleider, testimplantaten
- da **Brugsanvisning/Teknisk beskrivelse**
activC® håndtag, rømmejernsstyr, mejselstyr, prøveimplantater
- sv **Bruksanvisning/Teknisk beskrivning**
activC® Handtag, fräsguide, mejselguide, provimplantat
- fi **Käyttöohje/Tekninen kuvaus**
activC® -kahva, -kalvinohjain, -taltaohjain, -koeimplantit
- lv **Lietošanas instrukcijas/tehniskais apraksts**
activC® rokturis, rīvurbja vadotne, kalnu vadotne, izmēģinājuma implanta
- lt **Naudojimo instrukcija/techninis aprašas**
activC® rankena, plėtiklio kreiptuvas, kalto kreiptuvas, bandomieji implantai
- ru **Инструкция по применению/Техническое описание**
Пробные имплантаты, направляющая для долота, направляющая для фрезы, рукоятка activC®
- cs **Návod k použití/Technický popis**
activC® – rukojeť, vodítko výstružníku, vodítko dláta, zkušební implantáty
- pl **Instrukcja użytkowania/Opis techniczny**
Uchwyt, prowadnica rozwiertaka, prowadnica dłuta, implanty próbne activC®
- sk **Návod na použitie/Technický opis**
activC® – rukoväť, vodiaci prvok na výstružník, vodiaci prvok na dláto, skúšobné implantáty
- hu **Használati útmutató/Műszaki leírás**
activC® nyél, lyuktágító-vezető, vésővezető, próbaimplantátumok
- sl **Navodila za uporabo/Tehnični opis**
Ročaj, vodilo za povrtalo, poskusni vsadki activC®
- hr **Upute za uporabu/Tehnički opis**
Ručka, vodilica razvrtača, vodilica dljeta, probni implantati activC®
- ro **Manual de utilizare/Descriere tehnică**
Mâner, ghidaj pentru alezor, ghidaj pentru daltă, implanturi de probă activC®
- bg **Упътване за употреба/Техническо описание**
Ръкохватка, водач на разширител, водач на длето, пробни импланти activC®
- tr **Kullanım Kılavuzu/Teknik açıklama**
activC® tutamağı, rayba kılavuzu, keski kılavuzu, deneme implantları
- el **Οδηγίες χρήσης/Τεχνική περιγραφή**
Λαβή activC®, οδηγός γλυφάνου, οδηγός σμίλης, δοκιμαστικά εμφυτεύματα

FW647R

FW645R

Legend

- 1 Trial implant
- 2 Handle
- 3 Reamer guide
- 4 Mounting block (for reamer guide)
- 5 Reamer
- 6 Reamer handpiece
- 7 HiLAN pneumatic motor
- 8 activC chisel insert
- 9 Chisel holder
- 10 Chisel guide
- 11 Locking pipes
- 12 Adjusting wheel
- 13 Locking clip

1. About this document

Note

General risk factors associated with surgical procedures are not described in these instructions for use.

1.1 Scope

These instructions for use apply for activC handle, reamer guide, chisel guide and trial implants.

- ▶ For article-specific instructions for use as well as information on material compatibility and lifetime see B. Braun eIFU at eifu.bbraun.com

1.2 Safety messages

Safety messages make clear the dangers to patient, user and/or product that could arise during the use of the product. Safety messages are labeled as follows:

⚠ DANGER

Indicates a possible threat of danger. If not avoided, death or serious injury may result.

⚠ WARNING

Indicates a possible threat of danger. If not avoided, minor or moderate injury may result.

⚠ CAUTION

Indicates a possible threat of material damage. If not avoided, the product may be damaged.

2. Clinical use

2.1 Available sizes

Art. no.	Designation
FW789R ¹	Trial implant size XS, H5
FW790R ¹	Trial implant size XS, H6
FW791R ¹	Trial implant size S, H5
FW792R ¹	Trial implant size S, H6
FW793R ¹	Trial implant size M, H5
FW794R ¹	Trial implant size M, H6
FW795R ¹	Trial implant size L, H5
FW796R ¹	Trial implant size L, H6
FW797R ¹	Trial implant size XL, H5
FW798R ¹	Trial implant size XL, H6
FW799R ¹	Trial implant size XXL, H5
FW800R ¹	Trial implant size XXL, H6
FW8701 ¹	Handle
FW871R ²	Reamer guide (with mounting block)
GB771R ²	Hi-Line XS handpiece angled XL-I
GE700SU ⁴	Hi-Line XS activC drill XL-I D1, 5
FW784R ³	activC chisel guide
FW786R ³	activC chisel holder
FW787R ⁴	activC chisels

¹ Trays FW645R and FW647R

² only tray FW647R

³ only tray FW645R

⁴ single-use products

2.2 Areas of use and limitations of use

2.2.1 Intended use

The trial implants are used in conjunction with the handpiece to verify the size, height and shape of the activC intervertebral disk prosthesis to be implanted, and in conjunction with the reamer guide or the chisel guide to prepare the fin bed for the inferior implant plate of the activC intervertebral disk prosthesis.

2.2.2 Indications

Note

The manufacturer is not responsible for any use of the product against the specified indications and/or the described applications.

For indications, see Intended use.

2.2.3 Contraindications

No known contraindications.

2.3 Safety information

2.3.1 Clinical user

General safety information

To prevent damage caused by improper setup or operation, and to not compromise the manufacturer warranty and liability:

- ▶ Use the product only according to these instructions for use.
- ▶ Follow the safety and maintenance instructions.
- ▶ Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge and experience.
- ▶ Store any new or unused products in a dry, clean, and safe place.
- ▶ Prior to use, check that the product is in good working order.
- ▶ Keep the instructions for use accessible for the user.

Note

The user is obligated to report all severe events in connection with the product to the manufacturer and the responsible authorities of the state in which the user is located.

Notes on surgical procedures

It is the user's responsibility to ensure that the surgical procedure is performed correctly.

Appropriate clinical training as well as a theoretical and practical proficiency of all the required operating techniques, including the use of this product, are prerequisites for the successful use of this product.

The user is required to obtain information from the manufacturer if there is an unclear preoperative situation regarding the use of the product.

2.3.2 Product

Product-specific safety information

Risk of injury caused by incorrect operation of the product!

- ▶ Attend appropriate product training before using the product.
- ▶ For information about product training, please contact your national B. Braun/Aesculap agency.
- ▶ Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge, or experience.
- ▶ Only combine Aesculap products with each other.
- ▶ Always adhere to applicable standards.
- ▶ To prevent damage at the working end: carefully insert the product through the working channel (e.g. trocar).
- ▶ Only use the product under vision.

2.3.3 Sterility

Non-sterile packaged products

The product is delivered in an unsterile condition.

- ▶ Clean the new product after removing its transport packaging and prior to its initial sterilization.

Sterile products

Art. no.	Designation
GE700SU	Hi-Line XS activC drill XL-I D1, 5
FW787R	activC Chisels

The product has been gamma-sterilized and is supplied in sterile packaging.

- ▶ Do not use products from open or damaged sterile packaging.
- ▶ Prior to each use, inspect the product for loose, bent, broken, cracked or fractured components.
- ▶ Do not use the product if it is damaged or defective. Set aside the product if it is damaged.
- ▶ Do not use the product after its use-by date.

Single-use products

The reprocessing of the product affects its functionality. Risk of injury, illness or death due to soiling and/or impaired functionality of the product.

- ▶ Do not reprocess the product.

2.4 Application

⚠ WARNING

Risk of injury and/or malfunction!

- ▶ Prior to each use, inspect the product for loose, bent, broken, cracked, worn, or fractured components.
- ▶ Always carry out a function test prior to each use of the product.

⚠ WARNING

Risk of injury when using the product beyond the field of view!

- ▶ Apply the product only under visual control.

Note

Follow OP-Manual no. O31302 and instructions for use "activC intervertebral disk prosthesis", see TA011995.

2.4.1 Mounting the trial implant of the selected size on the handle

- ▶ Turn locking pipes 11 of handle 2 to "open" position  see Fig. 1.
- ▶ Mount the trial implant with the selected size, height and shape onto the handle 2. Note the "CRANIAL" and "CAUDAL" markings on the handle and on the trial implant.
- ▶ Turn locking pipes 11 of handle 2 to "closed" position .
- ▶ Check the trial implant for firm seating.

⚠ WARNING

Compression of the spinal canal and of other posterior elements caused by trial implants inserted too deeply!

- ▶ Prior to insertion of trial implants, adjust the depth stop to minimum insertion depth.
- ▶ Insert trial implants in the intervertebral space under X-ray control.
- ▶ Use adjusting wheel 12 to set the depth stop of the trial implant to minimum insertion depth.
- ▶ Insert the trial implant in the intervertebral space under X-ray control. Ensure correct cranial and caudal alignment of the trial implant.
- ▶ Align the markings of the trial implant center to the midline marking.
- ▶ If necessary, correct the position of the depth stop with adjusting wheel 12 until the intended position is reached. Ensure that depth stop abuts to the anterior edge of the vertebral body.

⚠ WARNING

Increased risk of migration if the selected implant size is too small!

- ▶ Select the trial implant or implant size, respectively, to achieve maximum coverage of the vertebral body endplates.
- ▶ On the X-ray, in AP and lateral view, check to ensure that the size, position and alignment of the trial implant are correct.

⚠ WARNING

Irritation of the facet joints and strain on spinal cord and nerve roots due to excessive or insufficient distraction!

- ▶ Select the correct trial implant/ implant height.
- ▶ Check for correct height of the trial implant in a lateral X-ray.
- ▶ On the intraoperative X-ray image, in AP and lateral view, check the size, height and position of the trial implant.
- ▶ With the trial implant properly positioned, turn locking pipes 11 to "open" position  and pull off handle 2 carefully in axial direction from the trial implant.
The trial implant remains in the intervertebral space.

Note


The fin bed can either be prepared by reaming or chiseling. This can be carried out either with the reamer guide with handpiece or the chisel guide with chisel holder, depending on which set is used.

2.4.2 Fin bed reaming

Assembling the reamer guide

Note

The handle and the reamer guide are fitted with identical locking mechanisms.

- ▶ Turn locking pipes 11 of reamer guide 3 to "open" position .
- ▶ Mount the reamer guide 3 onto the trial implant in the intervertebral space. Note the "CRANIAL" and "CAUDAL" markings on the reamer guide and on the trial implant.
- ▶ Turn locking pipes 11 of reamer guide 3 to "closed" position .
Reamer guide 3 is fixed to the trial implant.
- ▶ Check again, by X-ray imaging, the position of the trial implant.

Assembling the reamer handpiece

⚠ WARNING

Damage to nerves, spinal cord and other tissue because of fin bed reamed too deeply due to using too long a reamer handpiece or reamer!

- ▶ Only use reamer handpiece Hi-Line XS XL-1 (GB771R).
- ▶ Only use activC reamer (GE700SU).
- ▶ Prior to mounting, verify that the safety stop of the pneumatic motor is in position "0".
- ▶ Assemble the compressed-air tube, the HiLAN pneumatic motor and reamer handpiece 6 as described and illustrated in the respective instructions for use, see Fig. 2.
- ▶ Hold locking clip 13 of mounting block 4 pressed down and push the mounting block over the reamer handpiece 6, see Fig. 3.
- ▶ Slide on mounting block 4 as far as it will go and release locking clip 13, see Fig. 4.
- ▶ Pull at the mounting block 4 to check whether it is engaged properly. The arrow marking on the locking clip 13 must point to the "closed"  pictogram.
- ▶ Mount reamer 5 as described in the instructions for use of reamer handpiece 6.

⚠ WARNING

Risk of injury to blood vessels, nerves, spinal cord or other tissue if wrong reamer is used!

- ▶ Only use the activC reamer GE700SU.
- ▶ Check for technical faults prior to use.
- ▶ Do not use bent reamers.
- ▶ Do not reuse reamers after first use.

⚠ DANGER

Infection hazard for patients and/or users and impairment of product functionality due to reuse. Risk of injury, illness or death due to contamination and/or impaired functionality of the product!

- ▶ Do not reprocess the product.

Note

The reaming of the fin bed determines the position and orientation of the disk prosthesis in the plane of the intervertebral disk. The fin bed must only be prepared for the inferior implant plate.

⚠ WARNING

Malpositioning of the intervertebral disk prosthesis due to failure to observe the midline!

- ▶ Check the central position of the trial implant by means of the midline marking and on the AP-view X-ray image.
- ▶ Make certain that the position of the trial implant does not change during reaming.

⚠ WARNING

Risk of injury caused by sharp chisel!

- ▶ To avoid injuries to critical blood vessels, nerves and other tissue, as well as injuries to the operating surgeon, always work with appropriate care when a reamer is clamped in the machine.
- ▶ Do not activate the motor before the handpiece with the coupled reamer is properly positioned in the guide hole and the guiding slot of the reamer guide.
- ▶ A: Insert the reamer handpiece 6 with the reamer clamped in into the front guide hole of the reamer guide 3, see Fig. 5.
- ▶ B: Insert the laterally protruding pins of the mounting block 4 into the rear guide slots of the reamer guide 3, see Fig. 5.
- ▶ With the hand control not actuated, push the safety stop to .
The system is ready for operation.

Note

The reaming of the fin bed determines the position and orientation of the disk prosthesis in the plane of the intervertebral disk.

⚠ WARNING

Malpositioning of the intervertebral disk prosthesis due to failure to observe the midline!

- ▶ Check the central position of the trial implant by means of the midline marking and on the AP-view X-ray image.
- ▶ Make certain that the position of the trial implant does not change during reaming.

⚠ WARNING

Risk of injury to nerves and spinal cord caused by too deep insertion of the reamer!

- ▶ Use activC instruments only.
- ▶ Use the instruments according to their intended use.
- ▶ For reaming the holes, apply the reamer only in combination with a correctly selected and positioned trial implant and with the reamer guide coupled.
- ▶ Take care that the depth stop of the trial implant abuts the anterior edge of the vertebral body.
- ▶ Always maintain lateral X-ray control while reaming the holes.


- ▶ To ream the first hole under lateral X-ray control: With the motor running, push reamer handpiece 6 at an even speed and with mild forward pressure into the first guiding slot down to the stop.
- ▶ Pull back reamer handpiece 6 and, after moving the pins of mounting block 4 in the second guiding slot, repeat the reaming procedure under lateral X-ray control.
The reaming depth will be limited by the stop on the reamer guide 3.
- ▶ With the hand control not actuated, shift the safety stop on the pneumatic motor to position "0".
- ▶ Cautiously pull out reamer handpiece 6 first from the rear and then from the front reamer guide 3.
- ▶ Remove reamer guide 3 and the trial implant from the vertebral disk space. Clean the fin bed with nerve hooks from the discectomy set.

2.4.3 Chiseling the fin bed

Assembling the chisel guide

Note

The handle and the chisel guide are fitted with identical locking mechanisms.

- ▶ Turn locking pipes 11 of chisel guide 10 to "open" position .
- ▶ Mount the chisel guide 10 onto the trial implant in the intervertebral space. Note the "CRANIAL" and "CAUDAL" markings on the handle and on the trial implant.
- ▶ Turn locking pipes 11 of chisel guide 10 to "closed" position .
Chisel guide 10 is fixed to the trial implant.
- ▶ Check again, by X-ray imaging, the position of the trial implant.

⚠ WARNING

Damage to nerves, spinal cord and other tissue because of fin bed reamed too deeply due to using too long a chisel.

- ▶ Only use chisel holder FW786R.
- ▶ Only use chisel FW787R.
- ▶ Mount the chisel onto the chisel holder 9, see Fig. 6.
- ▶ To do this mount the chisel sideways onto the chisel holder and push it backwards until it engages. Pull at the chisel, to check that it is engaged properly.

⚠ WARNING

Risk of injury to blood vessels, nerves, spinal cord or other tissue if wrong chisel is used!

- ▶ Only use chisel FW787R.
- ▶ Do not use bent chisels. Do not reuse chisels after first use.

Note

The chiseling of the fin bed determines the position and orientation of the disk prosthesis in the plane of the intervertebral disk. The fin bed must only be prepared for the inferior implant plate.

⚠ WARNING

Malpositioning of the intervertebral disk prosthesis due to failure to observe the midline!

- ▶ Check the central position of the trial implant by means of the midline marking and on the AP-view X-ray image.
- ▶ Make certain that the position of the trial implant does not change during reaming.

⚠ WARNING

Risk of injury caused by sharp chisel!

- ▶ To avoid injuries to critical blood vessels, nerves and other tissue, as well as injuries to the operating surgeon, always work with appropriate care when a chisel is clamped in the chisel holder.
- ▶ A: Carefully push the chisel holder 9 with the chisel through the rear chisel guide, see Fig. 7.
- ▶ B: Carefully push the chisel holder 9 with the chisel through the front chisel guide, see Fig. 7.

Note

The chiseling of the fin bed determines the position and orientation of the disk prosthesis in the plane of the intervertebral disk.

⚠ WARNING

Malpositioning of the intervertebral disk prosthesis due to failure to observe the midline!

- ▶ Check the central position of the trial implant by means of the midline marking and on the AP-view X-ray image.
- ▶ Make certain that the position of the trial implant does not change during chiseling.

⚠ WARNING

Risk of injury to nerves and spinal cord caused by too deep insertion of the reamer!

- ▶ Use activC instruments only.
- ▶ Use the instruments according to their intended use.
- ▶ Only use the fin in combination with a correctly selected and positioned trial implant and with the chisel guide coupled.
- ▶ Take care that the depth stop of the trial implant abuts the anterior edge of the vertebral body.
- ▶ Always maintain lateral X-ray control while chiseling the fin.
- ▶ Carefully mount the chisel guide with chisel on the vertebral body.
- ▶ With slight pressure and at an even speed, chisel under lateral X-ray control. The chiseling depth is limited by the stop on the chisel guide 10.
- ▶ Carefully pull the chisel backwards out of the chisel guide.
- ▶ Remove chisel guide 10 and the trial implant from the vertebral disk space.
- ▶ Clean the fin bed with nerve hooks from the discectomy set.

3. Validated reprocessing procedure

3.1 Single-use products

Art. no.	Designation
GE700SU	Hi-Line XS activC drill XL-I D1, 5
FW787R	activC chisels

- ▶ Do not reuse the product.
The reprocessing of the product affects its functionality. Risk of injury, illness or death due to soiling and/or impaired functionality of the product.
- ▶ Do not reprocess the product.

3.2 Reusable products

Influences of the reprocessing which lead to damage to the product are not known.

A careful visual and functional inspection before the next use is the best opportunity to recognize a product that is no longer functional, see Inspection.

3.3 Preparations at the place of use

- ▶ If applicable, rinse non-visible surfaces preferably with deionized water, with a disposable syringe for example.
- ▶ Remove any visible surgical residues to the extent possible with a damp, lint-free cloth.
- ▶ Transport the dry product in a sealed waste container for cleaning and disinfection within 6 hours.

3.4 Preparing for cleaning

- ▶ Disassemble the product immediately after use, as described in the respective instructions for use.
- ▶ Disassemble the product prior to cleaning, see Disassembly.

3.5 Disassembly

3.5.1 Uncoupling the reamer guide/handle from the trial implant

- ▶ Turn locking pipes 11 of reamer guide 3 or of handle 2 to "open" position .
- ▶ Pull off the trial implant in axial direction.
- ▶ Decouple reamer 5 from reamer handpiece 6 as described in instructions for use TA011267.

3.5.2 Detaching the mounting block from the reamer handpiece

- ▶ Hold locking clip 13 of mounting block 4 pressed down (arrow marking points to "open" position ) and pull off the mounting block from reamer handpiece 6.

3.5.3 Uncoupling the chisel guide/handle from the trial implant

- ▶ Turn locking pipes 11 of chisel guide 10 or of handle 2 to "open" position .
- ▶ Pull off the trial implant in axial direction.
- ▶ Raise the chisel by the tip on the instrument and push out forwards.

3.6 Cleaning/Disinfection

3.6.1 Product-specific safety information on the reprocessing method

Damage to or destruction of the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!

- ▶ Following the manufacturer's instructions, use cleaning and disinfecting agents
 - that are approved for example aluminum, plastic materials, and high-grade steel,
 - that do not attack softeners (e.g., in silicone).
- ▶ Observe specifications regarding concentration, temperature and exposure time.
- ▶ Do not exceed the maximum allowable disinfection temperature of 95 °C.
- ▶ Use suitable cleaning/disinfecting agents if the product is disposed of in a wet condition. To prevent foaming and degradation of the efficacy of the process chemicals: prior to mechanical cleaning and disinfection, rinse the product thoroughly with running water

3.6.2 Validated cleaning and disinfection procedure

Validated procedure	Specific requirements	Reference
Manual cleaning with immersion disinfection ■ FW786R	<ul style="list-style-type: none"> ■ Cleaning brush: TA011944 ■ Disposable syringe 20 ml ■ When cleaning instruments with movable hinges, ensure that these are in an open position and, if applicable, move the hinge while cleaning. ■ Drying phase: Use a lint-free cloth or medical compressed air 	Chapter Manual cleaning/disinfection and subsection: ■ Chapter Manual cleaning with immersion disinfection
Manual cleaning with ultrasound and immersion disinfection ■ FW789R-FW800R ■ FW870 ■ FW871R ■ FW784R	<ul style="list-style-type: none"> ■ Cleaning brush: TA011944 ■ Disposable syringe 20 ml ■ When cleaning instruments with movable hinges, ensure that these are in an open position and, if applicable, move the hinge while cleaning. ■ FW789R-FW800R: Clean reamer guide 3, chisel guide 10 and handle 2 with locking pipes 11 in closed position  and do not move them during cleaning. ■ Drying phase: Use a lint-free cloth or medical compressed air 	Chapter Manual cleaning/disinfection and subsection: ■ Chapter Manual cleaning with ultrasound and immersion disinfection
Mechanical alkaline cleaning and thermal disinfection ■ FW786R	<ul style="list-style-type: none"> ■ Place the product on a tray that is suitable for cleaning (avoid rinsing blind spots). 	Chapter Mechanical cleaning/disinfection and subsection: ■ Chapter Mechanical alkaline cleaning and thermal disinfection
Manual pre-cleaning with ultrasound and brush, and subsequent mechanical alkaline cleaning and thermal disinfection ■ FW789R-FW800R ■ FW870 ■ FW871R ■ FW784R	<ul style="list-style-type: none"> ■ Cleaning brush: TA011944 ■ Disposable syringe 20 ml ■ Place the product on a tray that is suitable for cleaning (avoid rinsing blind spots). ■ Connect components with lumens and channels directly to the rinsing port of the injector carriage. ■ To flush the product: Use a flushing nozzle or flushing sleeve. ■ Keep working ends open for cleaning. ■ Place the product on the tray with all product links and joints open. ■ Clean reamer guide 3, chisel guide 10 and handle 2 with locking pipes 11 in closed position  and do not move them during cleaning. 	Chapter Mechanical cleaning/disinfection with manual pre-cleaning and subsection: ■ Chapter Manual pre-cleaning with ultrasound and brush ■ Chapter Mechanical alkaline cleaning and thermal disinfection

3.7 Manual cleaning/disinfection

- ▶ Prior to manual disinfecting, allow water to drip off for a sufficient length of time to prevent dilution of the disinfecting solution.
- ▶ After manual cleaning/disinfection, check visible surfaces visually for residues.
- ▶ Repeat the cleaning/disinfection process if necessary.

3.7.1 Manual cleaning with immersion disinfection

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Disinfecting cleaning	RT (cold)	>15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
II	Intermediate rinse	RT (cold)	1	-	D-W	-
III	Disinfection	RT (cold)	5	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
IV	Final rinse	RT (cold)	1	-	FD-W	-
V	Drying	RT	-	-	-	-

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, microbiological, at least of drinking water quality)

RT: Room temperature

*Recommended: BBraun Stabimed fresh

- ▶ Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

Phase I

- ▶ Fully immerse the product in the cleaning/disinfectant for at least 15 min. Ensure that all accessible surfaces are moistened.
- ▶ Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- ▶ If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
- ▶ Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- ▶ Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

Phase II

- ▶ Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- ▶ Drain any remaining water fully.

Phase III

- ▶ Fully immerse the product in the disinfectant solution.
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- ▶ Rinse lumens at least 5 times at the beginning of the exposure time using an appropriate disposable syringe. Ensure that all accessible surfaces are moistened.

Phase IV

- ▶ Rinse/flush the product thoroughly (all accessible surfaces).
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during final rinse.
- ▶ Rinse lumens with an appropriate disposable syringe at least five times.
- ▶ Drain any remaining water fully.

Phase V

- ▶ Dry the product in the drying phase with suitable equipment (e.g. cloth, compressed air), see Validated cleaning and disinfection procedure.

3.7.2 Manual cleaning with ultrasound and immersion disinfection

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Ultrasonic cleaning	RT (cold)	>15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
II	Intermediate rinse	RT (cold)	1	-	D-W	-
III	Disinfection	RT (cold)	5	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
IV	Final rinse	RT (cold)	1	-	FD-W	-
V	Drying	RT	-	-	-	-

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

RT: Room temperature

*Recommended: BBraun Stabimed fresh

- ▶ Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

Phase I

- ▶ Clean the product in an ultrasonic cleaning bath (frequency 35 kHz) for at least 15 min. Ensure that all accessible surfaces are immersed and acoustic shadows are avoided.
- ▶ Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- ▶ If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
- ▶ Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- ▶ Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

Phase II

- ▶ Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- ▶ Drain any remaining water fully.

Phase III

- ▶ Fully immerse the product in the disinfectant solution.
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- ▶ Rinse lumens at least five times at the beginning of the exposure time with an appropriate disposable syringe. Ensure that all accessible surfaces are moistened.

Phase IV

- ▶ Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during final rinse.
- ▶ Rinse lumens with an appropriate disposable syringe at least five times.
- ▶ Drain any remaining water fully.

Phase V

- ▶ Dry the product in the drying phase with suitable equipment (e.g. cloth, compressed air), see Validated cleaning and disinfection procedure.

3.8 Mechanical cleaning/disinfection

Note

The cleaning and disinfection device must be of tested and approved effectiveness (e.g. FDA approval or CE mark according to DIN EN ISO 15883).

Note

The cleaning and disinfection device used for processing must be serviced and checked at regular intervals.

3.8.1 Mechanical alkaline cleaning and thermal disinfecting

Machine type: single-chamber cleaning/disinfecting machine without ultrasound

Phase	Step	D [°C/°F]	t [min]	Water quality	Chemical/Note
I	Prerinse	<25/77	3	D-W	-
II	Cleaning	55/131	10	FD-W	<ul style="list-style-type: none"> ■ Concentrate, alkaline: <ul style="list-style-type: none"> - pH = 13 - <5 % anionic surfactant ■ working solution 0.5% <ul style="list-style-type: none"> - pH = 11*
III	Intermediate rinse	>10/50	1	FD-W	-
IV	Thermal disinfecting	90/194	5	FD-W	-
V	Drying	-	-	-	According to the program for cleaning and disinfection device

DW: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

*Recommended: BBraun Helimatic Cleaner alkaline

- ▶ Check visible surfaces for residues after mechanical cleaning/disinfecting.

3.9 Mechanical cleaning/disinfection with manual pre-cleaning

Note

The cleaning and disinfection device must be of tested and approved effectiveness (e.g. FDA approval or CE mark according to DIN EN ISO 15883).

Note

The cleaning and disinfection device used for processing must be serviced and checked at regular intervals.

3.9.1 Manual pre-cleaning with ultrasound and brush

Phase	Step	D [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Ultrasonic cleaning	RT (cold)	>15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
II	Rinsing	RT (cold)	1	-	D-W	-

D-W: Drinking water

RT: Room temperature

*Recommended: BBraun Stabimed fresh

- ▶ Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

Phase I

- ▶ Clean the product in an ultrasonic cleaning bath (frequency 35 kHz) for at least 15 min. Ensure that all accessible surfaces are immersed and acoustic shadows are avoided.
- ▶ Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- ▶ If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
- ▶ Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- ▶ Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

Phase II

- ▶ Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.

3.9.2 Mechanical alkaline cleaning and thermal disinfecting

Machine type: single-chamber cleaning/disinfection device without ultrasound

Phase	Step	T [°C/°F]	t [min]	Water quality	Chemical
I	Prerinse	<25/77	3	D-W	-
II	Cleaning	55/131	10	FD-W	<ul style="list-style-type: none"> ■ Concentrate, alkaline: <ul style="list-style-type: none"> - pH = 13 - <5 % anionic surfactant ■ working solution 0.5% <ul style="list-style-type: none"> - pH = 11*
III	Intermediate rinse	>10/50	1	FD-W	-
IV	Thermal disinfecting	90/194	5	FD-W	-
V	Drying	-	-	-	According to the program for cleaning and disinfection device

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

*Recommended: BBraun Helimatic Cleaner alkaline

- ▶ Check visible surfaces for residues after mechanical cleaning/disinfecting.

3.10 Inspection

- ▶ Allow the product to cool down to room temperature.
- ▶ Dry the product if it is wet or damp.

3.10.1 Visual inspection

- ▶ Ensure that all soiling has been removed. In particular, pay attention to mating surfaces, hinges, shafts, recessed areas, drill grooves and the sides of the teeth on rasps.
- ▶ If the product is dirty: repeat the cleaning and disinfection process.
- ▶ Check the product for damage, e.g. insulation or corroded, loose, bent, broken, cracked, worn or severely scratched and fractured components.
- ▶ Check the product for missing or faded labels.
- ▶ Check the products with long, slim shapes (in particular rotating instruments) for deformities.
- ▶ Check the product for damage to the spiral element.
- ▶ Check the cutting edges for continuity, sharpness, nicks and other damage.
- ▶ Check the surfaces for rough spots.
- ▶ Check the product for burrs that could damage tissue or surgical gloves.
- ▶ Check the product for loose or missing parts.
- ▶ Immediately put aside damaged or inoperative products and send them to Aesculap Technical Service, see Technical service.

3.10.2 Functional test

⚠ CAUTION

Damage (metal cold welding/friction corrosion) to the product caused by insufficient lubrication!

- ▶ Prior to function checks, lubricate moving parts (e.g. joints, pusher components and threaded rods) with maintenance oil suitable for the respective sterilization process (e.g. for steam sterilization: STERILIT® I oil spray JG600 or STERILIT® I drip lubricator JG598).

- ▶ Assemble disassembled products.
- ▶ Check that the product functions correctly.
- ▶ Check that all moving parts are working property (e.g. hinges, locks/latches, sliding parts etc.).
- ▶ Check for compatibility with associated products.
- ▶ Immediately put aside inoperative products and send them to Aesculap Technical Service, see Technical service.

3.11 Packaging

- ▶ Appropriately protect products with fine working tips.
- ▶ Store products with ratchet locks fully opened or locked no further than in the first notch.
- ▶ Place the product in its holder or on a suitable tray. Ensure that sharp edges are covered.
- ▶ Package trays appropriately for the sterilization process (e.g. in Aesculap sterile containers).
- ▶ Ensure that the packaging provides sufficient protection against contamination of the product during storage.

3.12 Steam sterilization

Note

The product may only be sterilized in disassembled condition.

Note

To avoid breakage due to stress crack corrosion, sterilize the instruments with the lock fully open or locked no further than on the first ratchet tooth.

- ▶ Check to ensure that the sterilizing agent will come into contact with all external and internal surfaces (e.g., by opening any valves and faucets).
- ▶ Validated sterilization process
 - Disassemble the product
 - Steam sterilization using fractional vacuum process
 - Steam sterilizer according to DIN EN 285 and validated according to DIN EN ISO 17665
 - Sterilization using fractionated vacuum process at 134 °C/holding time 5 min
- ▶ If several devices are sterilized at the same time in the same steam sterilizer: Ensure that the maximum permitted load according to the manufacturers' specifications is not exceeded.

3.13 Storage

- ▶ Store sterile products in germ-proof packaging, protected from dust, in a dry, dark, temperature-controlled area.
- ▶ Store sterile single-use products in germ-proof packaging in a dust-protected, dry, dark and temperature-controlled room.

3.13.1 Ambient conditions

The following environmental conditions apply to the transport and storage of the product:

Temperature	-10 °C to 50 °C
Relative humidity	10 % to 90 %
Atmospheric pressure	500 hPa to 1 060 hPa

4. Technical service

⚠ CAUTION

Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

- ▶ Do not modify the product.
- ▶ For service and repairs, please contact your national B. Braun/Aesculap agency.

Service addresses

Aesculap Technischer Service

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78532 Tuttlingen / Germany

Phone: +49 7461 95-1601

Fax: +49 7461 16-2887

E-Mail: ats@aesculap.de

Other service addresses can be obtained from the address indicated above.

5. Disposal

⚠ WARNING

Risk of infection due to contaminated products!

- ▶ Adhere to national regulations when disposing of or recycling the product, its components and its packaging.

⚠ WARNING

Risk of injury due to sharp-edged and/or pointed products!

- ▶ When disposing of or recycling the product, ensure that the packaging prevents injury by the product.

Note

The user institution is obliged to reprocess the product before its disposal, see Validated reprocessing procedure.

6. Symbols on product and packaging

	"Closed" position
	"Open" position
	Ready

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