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







1. About this document

Note

General risk factors associated with surgical procedures are not described in these instructions for use.

1.1 Scope

These instructions for use apply for activL revision instruments and distraction forceps.

Note

The applicable CE mark for the product can be seen on the label or packaging of the product.

Note

Instructions for use and further information about B. Braun / AESCULAP products can be found on the B. Braun eIFU website at eifu.bbraun.com

1.2 Safety messages

Safety messages make clear the dangers to patient, user and/or product that could arise during the use of the product. Safety messages are labeled as follows:

⚠ DANGER

Indicates a possible threat of danger. If not avoided, death or serious injury may result.

⚠ WARNING

Indicates a possible threat of danger. If not avoided, minor or moderate injury may result.

⚠ CAUTION

Indicates a possible threat of material damage. If not avoided, the product may be damaged.

2. Clinical use

2.1 Areas of use and limitations of use

2.1.1 Intended use

The activL revision instruments and distraction forceps are used for the distraction/mobilization of two adjacent lumbar vertebral bodies or prostheses plates, as well as for the repositioning and revision of activL intervertebral disk prostheses.

activL angled distraction forceps FW960R

The activL angled distraction forceps is used for the distraction/mobilization of two adjacent vertebral bodies.

activL angled revision distraction forceps FW965R

The activL angled revision distraction forceps is used for distracting the two activL metal plates of an implanted activL prosthesis in order to remove the PE inlay in corrective surgery.

activL revision instrument for endplates FW966R/FW967 and handle FW998R

The activL revision instrument for endplates and the corresponding handle are used for the repositioning and revision of activL intervertebral disk prostheses.

activL revision instrument for PE inlay FW968R

The activL revision instrument for PE inlay is used for the revision of the activL PE inlay.

2.1.2 Indications

Note

The manufacturer is not responsible for any use of the product against the specified indications and/or the described applications.

For indications, see Intended use.

2.1.3 Contraindications

No known contraindications.

2.2 Safety information

Use of the activL revision instruments and distraction forceps requires precise knowledge about spine surgery and the stabilization and biomechanical situation at the spine.

Surgical application of the activL revision instruments and distraction forceps is described in detail in the associated operating manual.

2.2.1 Clinical user

General safety information

To prevent damage caused by improper setup or operation, and to not compromise the manufacturer warranty and liability:

- ▶ Use the product only according to these instructions for use.
- ▶ Follow the safety and maintenance instructions.
- ▶ Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge and experience.
- ▶ Store any new or unused products in a dry, clean, and safe place.
- ▶ Prior to use, check that the product is in good working order.
- ▶ Keep the instructions for use accessible for the user.

Note

The user is obligated to report all severe events in connection with the product to the manufacturer and the responsible authorities of the state in which the user is located.

Notes on surgical procedures

It is the user's responsibility to ensure that the surgical procedure is performed correctly.

Appropriate clinical training as well as a theoretical and practical proficiency of all the required operating techniques, including the use of this product, are prerequisites for the successful use of this product.

The user is required to obtain information from the manufacturer if there is an unclear preoperative situation regarding the use of the product.

2.2.2 Sterility

The product is delivered in an unsterile condition.

- ▶ Clean the new product after removing its transport packaging and prior to its initial sterilization.

2.3 Application

⚠ WARNING

Risk of injury and/or malfunction!

- ▶ Prior to each use, inspect the product for loose, bent, broken, cracked, worn, or fractured components.
- ▶ Always carry out a function test prior to each use of the product.

2.3.1 Distracting/mobilizing vertebral bodies with the activL angled distraction forceps FW960R

⚠ DANGER

Risk of fracture of the vertebral body endplates and/or injury to surrounding structures due to overdistraction!

- ▶ Apply proper care when distracting.
- ▶ Avoid excessive distraction.

⚠ CAUTION

Damage to the product caused by actuation of the locking lever!

- ▶ To release lock 2, gently press the handle parts 1 together and actuate locking lever 3.

Note

Mobilization of the intervertebral disk space is a critical process and should be performed under X-ray control.

Note

The intervertebral space distraction is often wedge-shaped. Sufficient load relief in the posterior intervertebral disk space is required so that parallel distraction can be performed over the entire intervertebral space.

- ▶ Insert the working end 4 into the intervertebral space, see Fig. 2. When doing so, ensure that the instrument is pushed as dorsally as possible into the intervertebral disk space.
- ▶ Distract/mobilize the distraction forceps FW960R by pressing the handle parts 1 together as required.

Note

Lock 2 engages automatically when handle parts 1 are released.

2.3.2 Removing the activL angled distraction forceps FW960R

- ▶ Gently press the handle parts 1 together and actuate locking lever 3, see Fig. 1. This releases the lock 2 of the distraction forceps FW960R.

2.3.3 Applying the activL angled revision distraction forceps FW965R

⚠ DANGER

Injury to surrounding blood vessels and soft tissue!

Damage to the product or to the activL disk prosthesis!

- ▶ To protect the blood vessels, always use a retractor in conjunction with this product.
- ▶ Always insert the product into the two holes at the right-hand side of the prosthesis.

⚠ CAUTION

Damage to the product caused by actuation of the locking lever!

- ▶ To release the lock, gently press the handle parts together and actuate locking lever 3.

- ▶ Insert pins 5 of the distraction forceps FW965R into the holes on the right-hand side of the implant, see Fig. 4.
- ▶ Distract the distraction forceps FW965R by pressing the handle parts 1 together as required.

Note

Lock 2 engages automatically when handle parts 1 are released.

2.3.4 Removing the activL angled revision distraction forceps FW965R

- ▶ Gently press the handle parts 1 together and actuate locking lever 3, see Fig. 3. This releases the lock 2 of the distraction forceps FW965R.

2.3.5 Removing the PE inlay with activL revision instrument for PE inlay FW968R from the intervertebral disk prosthesis

- ▶ Insert the jaw part 6 of the revision instrument FW968R between the two endplates 7 of the intervertebral disk prosthesis, see Fig. 5. Ensure the following:
 - The orientation of the revision instrument FW968R in the caudal direction corresponds with the "caudal" marking on the revision instrument.
 - The edges of the jaw part 6 are pointing in the direction of the inferior endplates.
- ▶ Push the revision instrument along the lower endplate, until the jaw part 6 engages on the PE inlay.
- ▶ Turn the adjustment knob 8 clockwise to fix the PE inlay in place on jaw part 6, see Fig. 6.
- ▶ Lift the revision instrument FW968R with fixated PE inlay over the edge of the inferior endplate and remove from the endplate by pulling in the anterior direction.
- ▶ To remove the PE inlay from the revision instrument FW968R:
 - Rotate the adjustment knob 8 in a counterclockwise direction.
 - Pull the PE inlay from the jaw part 6.

2.3.6 Reposition/remove endplates with activL revision instrument for endplates FW966R/FW967R

⚠ DANGER

Risk of injury to nerves and surrounding structures due to uncontrolled posterior pressure of the hooks of the anchoring rods!

- ▶ Apply proper care when manipulating the anchoring rods.
- ▶ Do not apply excessive pressure.

- ▶ Screw handle FW998R onto handpiece of revision instrument FW966R/FW967R, see Fig. 7.

Note

The pins 13 on the screw connection of the anchoring rods of the revision instrument will point in the same direction as the corresponding hooks on the front part.

Note

In case of an incomplete revision of the implant, the endplates are removed one after another. The slotted hammer is also used for the revision procedure.

Note

The endplates can grow into and adhere strongly to the vertebral bodies. This can considerably complicate the revision procedure. If the revision is of a prosthesis implanted a number of years ago, the endplates may need to be removed from the intervertebral body with a chisel or elevator prior to using the revision instrument.

- ▶ Turn the anchoring rods 9 so that both hooks 10 are pointing outwards, see Fig. 8.
- ▶ Push the revision instrument FW966R/FW967R between the endplates. When doing so, ensure that the hooks reach over the posterior edges of the plates.
- ▶ Turn both anchoring rods 9 90° in the direction of the plate to be removed.
- ▶ Pull back the anchoring rods 9. When doing so, ensure that the anchoring rods engage securely in the recesses 11 on the posterior side of the plate, see Fig. 9.
- ▶ Fix the anchoring rods 9 in place by turning the lock nut 12 counterclockwise. When doing so, hold the anchoring rods so that they cannot turn anymore.
- ▶ Carefully reposition or withdraw the plate.
- ▶ To remove the endplate completely, hit it out with the slotted hammer FW579R.
- ▶ Undo the lock nut 12 by turning it clockwise.
- ▶ Turn the anchoring rods 9 90° in the anterior direction and remove the endplate.
- ▶ If applicable, remove other endplates in the same way.

3. Validated reprocessing procedure

3.1 General safety instructions

Note

Adhere to national statutory regulations, national and international standards and directives, and local, clinical hygiene instructions for sterile processing.

Note

For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of products.

Note

Mechanical reprocessing should be favored over manual cleaning as it gives better and more reliable results.

Note

Successful processing of this medical device can only be ensured if the processing method is first validated. The operator/sterile processing technician is responsible for this.

Note

If there is no final sterilization, then a virucidal disinfectant must be used.

Note

For up-to-date information about reprocessing and material compatibility, see B. Braun eFU at eifu.bbraun.com

The validated steam sterilization procedure was carried out in the Aesculap sterile container system.

3.2 General information

Dried or affixed surgical residues can make cleaning more difficult or ineffective and lead to corrosion. Therefore the time interval between application and processing should not exceed 6 h; also, neither fixating pre-cleaning temperatures >45 °C nor fixating disinfecting agents (active ingredient: aldehydes/alcohols) should be used.

Excessive measures of neutralizing agents or basic cleaners may result in a chemical attack and/or to fading and the laser marking becoming unreadable visually or by machine for stainless steel.

Residues containing chlorine or chlorides e.g. in surgical residues, medicines, saline solutions and in the service water used for cleaning, disinfection and sterilization will cause corrosion damage (pitting, stress corrosion) and result in the destruction of stainless steel products. These must be removed by rinsing thoroughly with demineralized water and then drying.

Additional drying, if necessary.

Only process chemicals that have been tested and approved (e.g. VAH or FDA approval or CE mark) and which are compatible with the product's materials according to the chemical manufacturers' recommendations may be used for processing the product. All the chemical manufacturer's application specifications must be strictly observed. Failure to do so can result in the following problems:

- Optical changes of materials, e.g. fading or discoloration of titanium or aluminum. For aluminum, the application/process solution only needs to be of pH >8 to cause visible surface changes.
- Material damage such as corrosion, cracks, fracturing, premature aging or swelling.
- ▶ Do not use metal cleaning brushes or other abrasives that would damage the product surfaces and could cause corrosion.
- ▶ Further detailed advice on hygienically safe and material-/value-preserving reprocessing can be found at www.a-k-i.org, link to "AKI-Brochures", "Red brochure".

3.3 Reusable products

Influences of the reprocessing which lead to damage to the product are not known.

A careful visual and functional inspection before the next use is the best opportunity to recognize a product that is no longer functional, see Inspection.

3.4 Preparations at the place of use

- ▶ If applicable, rinse non-visible surfaces preferably with deionized water, with a disposable syringe for example.
- ▶ Remove any visible surgical residues to the extent possible with a damp, lint-free cloth.
- ▶ Transport the dry product in a sealed waste container for cleaning and disinfection within 6 hours.

3.5 Preparing for cleaning

- ▶ Open up product with hinges.

FW966R, FW967R, FW968R:

- ▶ Disassemble the product prior to cleaning, see Disassembly.

3.6 Disassembly

3.6.1 activL revision instrument for endplates FW966R/FW967R

- ▶ Unscrew handle FW998R of revision instrument FW966R/FW967R.
- ▶ Turn the anchoring rods 9 so that the hooks 10 are positioned in the direction of the slots 14.
- ▶ Slide the anchoring rods 9 behind the handpiece until the hooks 10 are positioned on the middle section of the revision instrument FW966R/FW967R. When doing so, slide the hooks 10 through the slots 14, see Fig. 10.
- ▶ Turn the lock nuts 12 counterclockwise until they are fully undone.
- ▶ Slide the lock nuts 12 along the anchoring rods behind the handpiece.

3.6.2 activL revision instrument for PE inlay FW968R

- ▶ Turn the adjustment knob 8 counterclockwise and gently pull on the jaw part 6 until the jaw part comes out, see Fig. 11.

3.7 Cleaning/Disinfection

3.7.1 Product-specific safety information on the reprocessing method

Damage to or destruction of the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!

- ▶ Following the manufacturer's instructions, use cleaning and disinfecting agents that are approved for, e.g., aluminum, plastics, and high-grade steel, according to the manufacturer's instructions.
- ▶ Observe specifications regarding concentration, temperature and exposure time.
- ▶ Do not exceed the maximum allowable disinfection temperature of 95 °C.
- ▶ Use suitable cleaning/disinfecting agents if the product is disposed of in a wet condition. To prevent foaming and degradation of the efficacy of the process chemicals: prior to mechanical cleaning and disinfection, rinse the product thoroughly with running water
- ▶ If the microsurgical products can be securely fixed in machines or storage devices in such a way that they will be cleaned thoroughly, clean and disinfect them mechanically.

3.7.2 Validated cleaning and disinfection procedure

Validated procedure	Specific requirements	Reference
Manual cleaning with immersion disinfection ■ FW968R ■ FW998R	<ul style="list-style-type: none"> ■ Cleaning brush: TAO11944 ■ Disposable syringe 20 ml ■ Keep working ends open for cleaning. ■ When cleaning instruments with movable hinges, ensure that these are in an open position and, if applicable, move the hinge while cleaning. ■ Drying phase: Use a lint-free cloth or medical compressed air 	<p>Chapter Manual cleaning/disinfection and subsection:</p> <ul style="list-style-type: none"> ■ Chapter Manual cleaning with immersion disinfection
Manual cleaning with ultrasound and immersion disinfection ■ FW960R ■ Fw965R ■ FW966R ■ FW967R	<ul style="list-style-type: none"> ■ Cleaning brush: TAO11944 ■ Disposable syringe 20 ml ■ Keep working ends open for cleaning. ■ When cleaning instruments with movable hinges, ensure that these are in an open position and, if applicable, move the hinge while cleaning. ■ Drying phase: Use a lint-free cloth or medical compressed air 	<p>Chapter Manual cleaning/disinfection and subsection:</p> <ul style="list-style-type: none"> ■ Chapter Manual cleaning with ultrasound and immersion disinfection
Manual pre-cleaning with brush and subsequent mechanical alkaline cleaning and thermal disinfection ■ FW968R ■ FW998R	<ul style="list-style-type: none"> ■ Cleaning brush: TAO11944 ■ Disposable syringe 20 ml ■ Place the product on a tray that is suitable for cleaning (avoid rinsing blind spots). ■ Keep working ends open for cleaning. ■ Place the product on the tray with all product links and joints open. 	<p>Chapter Mechanical cleaning/disinfection with manual pre-cleaning and subsection:</p> <ul style="list-style-type: none"> ■ Chapter Manual pre-cleaning with a brush ■ Chapter Mechanical alkaline cleaning and thermal disinfecting
Manual pre-cleaning with ultrasound and brush, and subsequent mechanical alkaline cleaning and thermal disinfection ■ FW966R ■ FW967R	<ul style="list-style-type: none"> ■ Cleaning brush: TAO11944 ■ Disposable syringe 20 ml ■ Place the product on a tray that is suitable for cleaning (avoid rinsing blind spots). ■ Keep working ends open for cleaning. ■ Place the product on the tray with all product links and joints open. 	<p>Chapter Mechanical cleaning/disinfection with manual pre-cleaning and subsection:</p> <ul style="list-style-type: none"> ■ Chapter Manual pre-cleaning with ultrasound and brush ■ Chapter Mechanical alkaline cleaning and thermal disinfecting
Manual pre-cleaning with ultrasound and brush, and subsequent mechanical alkaline cleaning and thermal disinfection ■ FW960R ■ Fw965R	<p>Deviating manual pre-cleaning:</p> <ul style="list-style-type: none"> ■ Phase I: Fully immerse the product in cold water for at least 5 min and mobilize non-rigid components at least two times. Ensure that all accessible surfaces are moistened. ■ Phase II: Clean the product with a suitable cleaning brush under running water until all discernible residues have been removed from the surface. Mobilize non-rigid components, such as set screws, links, etc. during cleaning. ■ Phase III: Clean the product in an ultrasonic cleaning bath (frequency 35 kHz) for at least 15 min. Ensure that all accessible surfaces are immersed and acoustic shadows are avoided. <p>Mechanical cleaning</p> <ul style="list-style-type: none"> ■ Place the product on a tray that is suitable for cleaning (avoid rinsing blind spots). ■ Connect components with lumens and channels directly to the rinsing port of the injector carriage. ■ Keep working ends open for cleaning. ■ Place the product on the tray with all product links and joints open. 	<p>Chapter Mechanical cleaning/disinfection with manual pre-cleaning and subsection:</p> <ul style="list-style-type: none"> ■ Chapter Manual pre-cleaning with ultrasound and brush ■ Chapter Mechanical alkaline cleaning and thermal disinfecting

3.8 Manual cleaning/disinfection

- ▶ Prior to manual disinfecting, allow water to drip off for a sufficient length of time to prevent dilution of the disinfecting solution.
- ▶ After manual cleaning/disinfection, check visible surfaces visually for residues.
- ▶ Repeat the cleaning/disinfection process if necessary.

3.8.1 Manual cleaning with immersion disinfection

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Disinfecting cleaning	RT (cold)	>15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
II	Intermediate rinse	RT (cold)	1	-	D-W	-
III	Disinfection	RT (cold)	5	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
IV	Final rinse	RT (cold)	1	-	FD-W	-
V	Drying	RT	-	-	-	-

D-W:

Drinking water

FD-W: Fully desalinated water (demineralized, microbiological, at least of drinking water quality)

RT:

Room temperature

*Recommended: B Braun Stabimed fresh

- ▶ Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

Phase I

- ▶ Fully immerse the product in the cleaning/disinfectant for at least 15 min. Ensure that all accessible surfaces are moistened.
- ▶ Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- ▶ If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
- ▶ Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- ▶ Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

Phase II

- ▶ Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- ▶ Drain any remaining water fully.

Phase III

- ▶ Fully immerse the product in the disinfectant solution.
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- ▶ Rinse lumens at least 5 times at the beginning of the exposure time using an appropriate disposable syringe. Ensure that all accessible surfaces are moistened.

Phase IV

- ▶ Rinse/flush the product thoroughly (all accessible surfaces).
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during final rinse.
- ▶ Rinse lumens with an appropriate disposable syringe at least five times.
- ▶ Drain any remaining water fully.

Phase V

- ▶ Dry the product in the drying phase with suitable equipment (e.g. cloth, compressed air), see Validated cleaning and disinfection procedure.

3.8.2 Manual cleaning with ultrasound and immersion disinfection

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Ultrasonic cleaning	RT (cold)	>15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
II	Intermediate rinse	RT (cold)	1	-	D-W	-
III	Disinfection	RT (cold)	5	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
IV	Final rinse	RT (cold)	1	-	FD-W	-
V	Drying	RT	-	-	-	-

D-W: Drinking water
FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

RT: Room temperature

*Recommended: BBraun Stabimed fresh

- ▶ Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

Phase I

- ▶ Clean the product in an ultrasonic cleaning bath (frequency 35 kHz) for at least 15 min. Ensure that all accessible surfaces are immersed and acoustic shadows are avoided.
- ▶ Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- ▶ If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
- ▶ Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- ▶ Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

Phase II

- ▶ Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- ▶ Drain any remaining water fully.

Phase III

- ▶ Fully immerse the product in the disinfectant solution.
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- ▶ Rinse lumens at least five times at the beginning of the exposure time with an appropriate disposable syringe. Ensure that all accessible surfaces are moistened.

Phase IV

- ▶ Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during final rinse.
- ▶ Rinse lumens with an appropriate disposable syringe at least five times.
- ▶ Drain any remaining water fully.

Phase V

- ▶ Dry the product in the drying phase with suitable equipment (e.g. cloth, compressed air), see Validated cleaning and disinfection procedure.

3.9 Mechanical cleaning/disinfection with manual pre-cleaning

Note

The cleaning and disinfection device must be of tested and approved effectiveness (e.g. FDA approval or CE mark according to DIN EN ISO 15883).

Note

The cleaning and disinfection device used for processing must be serviced and checked at regular intervals.

3.9.1 Manual pre-cleaning with a brush

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Disinfecting cleaning	RT (cold)	>15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
II	Rinsing	RT (cold)	1	-	D-W	-

D-W: Drinking water

RT: Room temperature

*Recommended: BBraun Stabimed fresh

- ▶ Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

Phase I

- ▶ Fully immerse the product in the cleaning/disinfectant for at least 15 min. Ensure that all accessible surfaces are moistened.
- ▶ Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- ▶ If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
- ▶ Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- ▶ Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

Phase II

- ▶ Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.

3.9.2 Manual pre-cleaning with ultrasound and brush

Phase	Step	D [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Ultrasonic cleaning	RT (cold)	>15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
II	Rinsing	RT (cold)	1	-	D-W	-

D-W: Drinking water

RT: Room temperature

*Recommended: BBraun Stabimed fresh

- ▶ Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

Phase I

- ▶ Clean the product in an ultrasonic cleaning bath (frequency 35 kHz) for at least 15 min. Ensure that all accessible surfaces are immersed and acoustic shadows are avoided.
- ▶ Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- ▶ If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
- ▶ Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- ▶ Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

Phase II

- ▶ Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.

3.9.3 Mechanical alkaline cleaning and thermal disinfecting

Machine type: single-chamber cleaning/disinfection device without ultrasound

Phase	Step	T [°C/°F]	t [min]	Water quality	Chemical
I	Prerinse	<25/77	3	D-W	-
II	Cleaning	55/131	10	FD-W	<ul style="list-style-type: none"> ■ Concentrate, alkaline: <ul style="list-style-type: none"> - pH = 13 - <5 % anionic surfactant ■ working solution 0.5% <ul style="list-style-type: none"> - pH = 11*
III	Intermediate rinse	>10/50	1	FD-W	-
IV	Thermal disinfecting	90/194	5	FD-W	-
V	Drying	-	-	-	According to the program for cleaning and disinfection device

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

*Recommended: BBraun Helimatic Cleaner alkaline

- ▶ Check visible surfaces for residues after mechanical cleaning/disinfecting.

3.10 Inspection

- ▶ Allow the product to cool down to room temperature.
- ▶ Dry the product if it is wet or damp.

3.10.1 Visual inspection

- ▶ Ensure that all soiling has been removed. In particular, pay attention to mating surfaces, hinges, shafts, recessed areas, drill grooves and the sides of the teeth on rasps.
- ▶ If the product is dirty: repeat the cleaning and disinfection process.
- ▶ Check the product for damage, e.g. insulation or corroded, loose, bent, broken, cracked, worn or severely scratched and fractured components.
- ▶ Check the product for missing or faded labels.
- ▶ Check the products with long, slim shapes (in particular rotating instruments) for deformities.
- ▶ Check the product for damage to the spiral element.
- ▶ Check the cutting edges for continuity, sharpness, nicks and other damage.
- ▶ Check the surfaces for rough spots.
- ▶ Check the product for burrs that could damage tissue or surgical gloves.
- ▶ Check the product for loose or missing parts.
- ▶ Immediately put aside damaged or inoperative products and send them to Aesculap Technical Service, see Technical service.

3.10.2 Functional test

⚠ CAUTION

Damage (metal cold welding/friction corrosion) to the product caused by insufficient lubrication!

▶ **Prior to function checks, lubricate moving parts (e.g. joints, pusher components and threaded rods) with maintenance oil suitable for the respective sterilization process (e.g. for steam sterilization: STERILIT® I oil spray JG600 or STERILIT® I drip lubricator JG598).**

- ▶ Assemble disassembled products, see Assembly.
- ▶ Check that the product functions correctly.
- ▶ Check that all moving parts are working properly (e.g. hinges, locks/latches, sliding parts etc.).
- ▶ Check for compatibility with associated products.
- ▶ Immediately put aside inoperative products and send them to Aesculap Technical Service, see Technical service.

3.11 Assembly

3.11.1 activL revision instrument for endplates FW966R/FW967R

- ▶ Slide the lock nuts 12 along the anchoring rods 9 as far as the threaded rods, see Fig. 10.
- ▶ Turn the lock nuts 12 clockwise onto the threaded rods.
- ▶ Turn the anchoring rods 9 so that the hooks 10 are positioned in the direction of the slots 14.
- ▶ Slide the hooks 10 through the slots 14.

3.11.2 activL revision instrument for PE inlay FW968R

- ▶ Insert the jaw part 6 into the tubular shaft 15 of the revision instrument FW968R until it engages, see Fig. 11. When doing so, ensure that the "caudal" marking on the tubular shaft is aligned with the edge of the jaw part 6.
- ▶ Push the jaw part 6 in the direction of the adjustment knob 8 and turn the adjustment knob clockwise until the jaw part is positioned just in front of the tubular shaft. When doing so, ensure that the jaw part is not compressed by the tubular shaft.

3.12 Packaging

- ▶ Appropriately protect products with fine working tips.
- ▶ Store products with ratchet locks fully opened or locked no further than in the first notch.
- ▶ Place the product in its holder or on a suitable tray. Ensure that sharp edges are covered.
- ▶ Package trays appropriately for the sterilization process (e.g. in Aesculap sterile containers).
- ▶ Ensure that the packaging provides sufficient protection against contamination of the product during storage.

3.13 Steam sterilization

Note

To avoid breakage due to stress crack corrosion, sterilize the instruments with the lock fully open or locked no further than on the first ratchet tooth.

- ▶ Check to ensure that the sterilizing agent will come into contact with all external and internal surfaces (e.g., by opening any valves and faucets).
- ▶ Validated sterilization process
 - Steam sterilization using fractional vacuum process
 - Steam sterilizer according to DIN EN 285 and validated according to DIN EN ISO 17665
 - Sterilization using fractionated vacuum process at 134 °C/holding time 5 min
- ▶ If several devices are sterilized at the same time in the same steam sterilizer: Ensure that the maximum permitted load according to the manufacturers' specifications is not exceeded.

3.14 Storage

- ▶ Store sterile products in germ-proof packaging, protected from dust, in a dry, dark, temperature-controlled area.

4. Technical service

⚠ CAUTION

Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

- ▶ **Do not modify the product.**
- ▶ **For service and repairs, please contact your national B. Braun/Aesculap agency.**

Service addresses

Aesculap Technischer Service

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78532 Tuttlingen / Germany

Phone: +49 7461 95-1601

Fax: +49 7461 16-2887

E-Mail: ats@aesculap.de

Other service addresses can be obtained from the address indicated above.

5. Disposal

⚠ WARNING

Risk of infection due to contaminated products!

- ▶ **Adhere to national regulations when disposing of or recycling the product, its components and its packaging.**

⚠ WARNING

Risk of injury due to sharp-edged and/or pointed products!

- ▶ **When disposing of or recycling the product, ensure that the packaging prevents injury by the product.**

Note

The user institution is obliged to reprocess the product before its disposal, see Validated reprocessing procedure.

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