| | | |
|----|---|------------|
| DE | BASIX; BASIX MIS, MOSS 100, MOSS VRS; MOSS VRS MIS Gebrauchsanweisung | Seite 2 |
| EN | BASIX; BASIX MIS, MOSS 100, MOSS VRS; MOSS VRS MIS Instruction for Use | Page 22 |
| FR | BASIX; BASIX MIS, MOSS 100 MOSS VRS; MOSS VRS MIS Mode d'emploi | Page 41 |
| IT | BASIX; BASIX MIS, MOSS 100, MOSS VRS; MOSS VRS MIS Istruzioni per l'uso | Pagina 61 |
| ES | BASIX; BASIX MIS, MOSS 100, MOSS VRS; MOSS VRS MIS Instrucciones de uso | Página 81 |
| PT | BASIX; BASIX MIS, MOSS 100, MOSS VRS; MOSS VRS MIS Instrucoes de utilizacao | Página 101 |
| EL | BASIX; BASIX MIS, MOSS 100, MOSS VRS; MOSS VRS MIS Οδηγίες χρήσης | Σελίδα 121 |
| SK | BASIX; BASIX MIS, MOSS 100, MOSS VRS; MOSS VRS MIS Názov systému | Strana 141 |
| CZ | BASIX; BASIX MIS, MOSS 100, MOSS VRS; MOSS VRS MIS Návod k použití | Strana 160 |
| NL | BASIX; BASIX MIS, MOSS 100, MOSS VRS; MOSS VRS MIS Gebruiksaanwijzing | Pagina 179 |

DEUTSCH

1. ZU DIESEM DOKUMENT

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ACHTUNG

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⚠ VORSICHT

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1.1 Depictions in the document

ATTENTION

Warns of the possibility of material damage.

CAUTION

Warns of the possibility of minor personal injury.

WARNING

Warns of the possibility of serious or fatal injury.

DANGER

Warns of the imminent danger of serious or fatal injury.

CAUTION

This is a warning symbol which warns of dangers that could arise throughout the entire subsequent process.

- In this case, the warning symbol indicates possible consequences.
- This is a precautionary measure against the danger.

1.2 Accompanying documentation

www.biedermann.com/eifu

For the implant systems that are taken into consideration in these instructions for use, there are the following surgical techniques that can be used to support the application:

- | | |
|-----------------------------------|-----------|
| - BASIX Surgical Technique | GI-0153 |
| - BASIX MIS Surgical Technique | GI-0162 |
| - MOSS 100 Surgical Technique | PD-038079 |
| - MOSS VRS Surgical Technique | GI-038994 |
| - MOSS VRS MIS Surgical Technique | GI-041021 |

1.3 Instruments

Under normal conditions, the instruments are intended for an uninterrupted use for a period of less than 60 minutes. The period of application is therefore defined as temporary in accordance with paragraph 1.1, Annex IX of Directive 93/42/EEC and Paragraph 1.1 of Chapter I of Annex VIII of the European Medical Devices Regulation ((EU) 2017(745)).

The instruments are to be used for implantation or explantation.

2. INTENDED PURPOSE

The BASIX, BASIX MIS, MOSS 100, MOSS VRS und MOSS VRS MIS implant system is a comprehensive thoracolumbosacral spinal system that offers posterior clinical solutions.

The BASIX, BASIX MIS, MOSS 100, MOSS VRS und MOSS VRS MIS implant system is a universal set of instruments and implants that are indicated for the treatment of significant mechanical instability or deformity of the spine which requires fusion with instrumentation.

DANGER

Danger for patients, users or third parties.

- The implant system must only be used in line with the intended purpose and the indications.
- The relative contraindications must be taken into account in deciding whether to use the implant system.
- The implant system must not be used if absolute contraindications apply.

3. INDICATIONS

The BASIX, BASIX MIS, MOSS 100, MOSS VRS and MOSS VRS MIS Spinal System is intended to provide immobilization and stabilization of spinal segments as an adjunct to fusion in the treatment of acute and chronic instabilities or deformities of the thoracic, lumbar and sacral spine.

In case of skeletally mature patients, the following indications must be given:

- degenerative disc disease (DDD is defined as back pain of discogenic origin with degeneration of the disc confirmed by history and radiographic studies.)
- spondylolisthesis
- trauma (i. e. fracture or dislocation)
- spinal stenosis
- deformities or curvatures (i. e. scoliosis, kyphosis and/or lordosis)
- tumor
- pseudoarthrosis
- failed previous fusion

When used for posterior non-cervical pedicle screw fixation in pediatric patients, the MOSS 100 and MOSS VRS System is indicated as an adjunct to fusion to treat adolescent idiopathic scoliosis. The MOSS 100 and MOSS VRS System is intended to be used with autograft and/or allograft. Pediatric pedicle screw fixation is limited to a posterior approach.

WARNING

Failure / breaking of pedicle screws with injection channel because of mechanical weakness due to the injection channel.

- Use of pedicle screws with injection channel only in patients with osteoporosis
- Pedicle screws with injection channel must be cemented.

4. CONTRAINDICATIONS

The systems must not be used in the following cases:

- Overall inoperability
- Advanced osteoporosis
- Active infection
- The system is not indicated for use in patients on whom it is not possible to carry out imaging procedures intraoperatively (fluoroscopy or computer navigation).

Relative contraindications

In the following cases, there is an increased risk of complications or implant failure. The systems may only be used once the surgeon has weighed these risks (benefit-risk ratio).

- Morbid obesity
- Comorbidities or cotreatments which can affect bone fusion (e.g., chemotherapy/rheumatoid arthritis)
- Failure on the patient's part to comply with the prescribed postoperative precautionary measures due to a lack of reliability (e.g, due to addictions or mental illnesses).

5. COMPLICATIONS AND SIDE EFFECTS

The following complications have been identified in the scope of a clinical evaluation of the available medical literature:

- Superficial infections
- Deep infections
- Damage to blood vessels
- Blood loss
- Haematomas
- Deep venous thrombosis
- Injuries to the pedicle cortex
- Incorrect placement of screws
- Loss of correction
- Pseudarthroses
- Nerve root irritation
- Neurological deficits
- (Partial) paralysis
- Cerebrospinal fluid leak / dural sac lesions
- Implant breakage
- Implant loosening
- Implant migration
- Infections of the urogenital tract
- Incisional hernia
- Retrograde ejaculation
- Metallosis

Additional Adverse Events for Pediatric Patients:

- Inability to use pedicle screw fixation due to anatomic limitations (pedicle dimensions, distorted anatomy).
- Pedicle screw malpositioning, with or without neurological or vascular injury.
- Proximal or distal junctional kyphosis.
- Pancreatitis.

6. WARNINGS

The following warnings are the results of the risk management assessment performed in accordance with DIN EN ISO 14971 and refer to causes which can lead to risk situations and patient injury in the form of complications or side effects.

General

A satisfactory result is not achieved with every surgical intervention.

It is possible that the pre-surgical symptoms may not be alleviated. In exceptional cases, a worsening of the situation is also possible.

Patient selection

The following patient selection factors can have a decisive influence on the long-term success of treatment:

- Patients should be selected whilst taking the purpose, indications, relative contraindications and the contraindications into account.
- Patient's age:
 - The implant system must mainly be used on patients whose skeletal growth is concluded.
 - The use of pedicle screw fixation in the pediatric population may present additional risks when patients are of smaller stature and skeletally immature. Pediatric patients may have smaller spinal structures (pedicle diameter or length) that may preclude the use of pedicle screws or increase the risk of pedicle screw malpositioning and neurological or vascular injury. Patients not skeletally mature that undergo spinal fusion procedures may have longitudinal spinal growth, or may be at risk for rotational spinal deformities (the "crankshaft phenomenon") due to continued differential growth of the anterior spine.
 - The safety and effectiveness of these devices has not been established for use as part of a pediatric growing rod construct. These devices are only intended to be used when definitive fusion is being performed at all instrumented levels. Other adverse events related to pedicle screw fixation, such as screw or rod bending, breakage, or loosening, may also occur in pediatric patients. Pediatric patients may be at increased risk for device related injury because of their smaller stature.
- Patient's body weight:
 - Obesity or adiposity can lead to excessive loading of the implant system.
- Patient's occupation or activity:
 - Physically demanding activities, e.g. heavy lifting, excessive bending, etc., can place excessive loads on the implant system.

WARNINGS

- Patient's mental state:
 - Psychological illnesses can contribute to a patient ignoring necessary precautionary measures, which can lead to excessive loading of the implant system.
- Alcohol or drug abuse:
 - Alcohol or drug abuse can contribute to a patient ignoring necessary precautionary measures which can lead to excessive loading of the implant system.
- Bone quality:
 - The user must check the bone quality of the vertebral segments to be treated. In the case of a poor condition of the bones or the possibility of progressive osteoporotic weakening of the bone substance, there is the risk of loosening implants or migration of the implant into affected segments. If the bones are in poor condition (osteoporosis), application of bone cement should be considered to secure bone anchorage.

Qualifications

Users may only use the systems if they meet the following qualifications:

- Surgeons specialized in orthopedics, neurology or spinal surgery.
- Familiarity with the state of the art in instrumentation of the spine.
- Trained in the use of the respective system.

Stability of the anterior spine

Biomechanical studies have shown that the majority of the axial load is borne by the anterior column of the spine. If there is no stable anterior support, this can lead to excess strain on the bone anchor or the system.

- If there is insufficient load transfer via the anterior column of the spine, a surgical procedure to support stability by way of, e.g. intercorporeal fusion, vertebral replacement, etc. should be considered in addition to internal spinal fixation.

Material intolerance

The implants are made of titanium alloy. In addition, CoCr rods are available.

In patients who are intolerant to one of the materials used in the respective system, the site of the operation may become infected or anaphylactic shock may be triggered.

- If it is suspected that the patient has a material intolerance, appropriate tests must be performed prior to the operation.

Patient briefing

Failure to brief the patient sufficiently can lead to excessive loading of the implant.

- Inform the patient before the surgery of any possible complications during and after the operation.
- Inform the patient of appropriate physical exertion and of any restrictions after the operation.
- Warn the patient to avoid full weight bearing until bone healing has been proven.

As-delivered condition

All of the components of the system are delivered in a non-sterile condition.

- The user must reprocess the components before use.

Compatibility of system and materials

The implant system is not compatible with systems from other manufacturers.

- The use of the implants is performed with the appropriate, related instruments of the system.
- A combination with instruments or implant systems of alternative manufacturers is not provided.

Manipulation of implant components

Manipulation of the implant components can result in damage which may reduce the resilience of the implant system.

- Do not change or readjust implants mechanically unless this is expressly intended.
- Only rods may be readjusted. When doing so, only use the intended instrument set and keep manipulation to a minimum. Never bend rods back into shape.
- Avoid defects (notches , scratches) during beding, as this causes surface defects which can lead to internal concentrations of stress.
- Rods with instrument interface for connection to a rod inserter must not be bent in this area.

Damaged implants or instruments

Damaged implants (i.e. notches, scratches, deformities) or instruments (instruments whose function is no longer present) can not guarantee proper functioning.

- Dispose damaged implants.
- Replace damaged instruments.

Implant selection and positioning

The size and shape of the human bone determine the size and stability of the implant. Incorrectly selected implants and incorrectly positioned implants can impede the function of the implant and damage anatomical structures.

- Check the correct implant size using the appropriate equipment or using imaging techniques (pre- and intraoperative x-ray, computer tomography).
- Check the correct position of the implant intraoperatively using imaging techniques.

Avoiding contamination

Contaminated components of the implant system can cause infections.

- All components must be cleaned, disinfected and sterilized before each use.
- Only use validated cleaning, disinfection and sterilization procedures specific to the equipment and product.
- Only use equipment which has been maintained and checked.
- Comply with the validated parameters in each cycle.
- Prior to the operation, the components must be removed from the container under aseptic conditions.
- Ensure that the implants are not contaminated or damaged during storage, transport and removal from the caddy or container.

Reuse of implants

Reuse of implants which have already come into contact with body tissue or bodily fluids is not permitted for the following reasons:

- Used implants can have defects which are not visible
- Biological cross-contamination
- Risk of infection for the patient

PRECAUTIONARY MEASURES

- Death of the patient

Dispose of implants which have come into contact with body tissue or bodily fluids in line with standard hospital practice

Aftercare

Delayed or the absence of spondylodesis can put excessive loads on the implant.

- Evaluate the progression of spondylodesis. In the case of a delay in or absence of spondylodesis, adapt the subsequent treatment strategy.
- Check the position of the implant, e.g. x-ray, CAT scan.

Revision and lifetime

The implants of the systems can be explanted following successful spondylodesis. In this case, the ideal lifespan of the implant is the period of time up to spondylodesis.

Magnetic resonance imaging

The implant systems have not been evaluated for safety and compatibility in the MR environment. They have not been tested for heating, migration, or image artifact in the MR environment. The safety of the implant systems in the MR environment is unknown. Scanning a patient who has one of these devices may result in patient injury.

Radiation exposure

Radiation exposure is increased due to the minimally invasive procedure.

- To keep the exposure to radiation as low as possible, a risk assessment should be performed to evaluate the advantages and disadvantages compared with an open procedure.

Final tightening

Incorrect final tightening can damage the implant and endanger the patient.

- For final tightening, always use the torque limiting T-handle.
- Always use a counter-holder during final tightening to accept the forces and torques that occur.

Pedicle Screws with injection channel

Improper use of pedicle screws with injection channel can overload the implant system or cause injury of the patient.

Therefore observe the following when using pedicle screws with injection channel:

- Evaluate the use of augmentation during the preoperative planning phase or intraoperatively, consider the bone quality.
- Perform bilateral augmentation.
- An imaging techniques must be used for control during augmentation process.
- Always use pedicle screws with injection channel with vertebroplasty cement.
- Apply the vertebroplasty in accordance with the cement manufacturer's specifications.
- Observe the Instruction for use of the cement manufacturer.

7. PRECAUTIONARY MEASURES

The implantation of pedicle screw spinal systems should be performed only by experienced spinal surgeons with specific training in the use of this pedicle screw spinal system because this is a technically demanding procedure presenting a risk of serious injury to the patient.

8. REPROCESSING

Please also observe the valid legal provisions in your country and the hygiene regulations of the medical practice or hospital. This applies in particular to the various requirements as regards effective prion inactivation (not relevant for USA).

WARNING

Procedures for prion inactivation lacking.

- With intended use, contact with patient areas relevant in terms of prions cannot be entirely ruled out in all cases.
 - High risk to the patient due to insufficient prion inactivation.
- Observe the applicable national requirements. These requirements may also include procedures which go beyond the present reprocessing instructions.

8.1 Implants – Cleaning and Disinfection

WARNING

Impermissible reuse of damaged and contaminated implants.

- Patients can become infected following contact with contaminated implants.
 - Damaged implants may fail.
- Sort out damaged implants and implants which come into contact with the patient or his/her bodily secretions and dispose of them according to the requirements of the hospital.

8.1.1 Pre-treatment directly after use

WARNING

Cross-contamination of implants through pre-treatment directly after use.

- Patients can become infected following contact with contaminated implants.
- Do not pre-treat implants directly after use.

8.1.2 Pre-treatment in the Central Sterile Services Department

WARNING

Cross-contamination of implants through pre-treatment.

- Patients can become infected following contact with contaminated implants.
- Do not pre-treat implants.

8.1.3 Manual cleaning and disinfection

WARNING

Residual contamination due to inadequate manual cleaning and manual disinfection.

- Patients can become infected following contact with inadequately cleaned and disinfected implants.
- Do not clean or disinfect instruments manually.

8.1.4 Machine cleaning and disinfection

WARNING

Residual contamination on contact surfaces.

- Patients can become infected following contact with contaminated implants.
 - Residual cleaning agents can damage implants.
- Implants and caddies must be machine cleaned and disinfected separately from each other.

ATTENTION

Damage and contamination of implants due to the use of an unsuitable disinfection method.

- Do not use chemical disinfection.

Machine cleaning and disinfection process

Selection criteria for the operating resources for the machine cleaning and disinfection process, see chapter 8.7.2 „Operating resources for machine cleaning and disinfection of implants and instruments“.

- Remove the implants from the caddy and place in a cleaning tray (for details about positioning see table in chapter 8.1.8 „Instructions for the implant groups“).
- Place the cleaning tray in the cleaning and disinfecting device.
- Check that the implants are not in contact with each other.
- Start the program.
- Remove the implants after the end of the program.

8.1.5 Inspection

Visually inspect the implants for the following defects:

- Corrosion
- Surface damage
- Chipping
- Soiling

WARNING

Impermissible reuse of damaged and soiled implants.

- Patients can become infected following contact with contaminated implants.
 - Damaged implants may fail.
- Sort out damaged and soiled implants and dispose of them according to the requirements of the hospital.

8.1.6 Care of the implants

WARNING

Damage and contamination of implants through the use of care products.

- Do not apply oil to implants.

8.1.7 Packaging

- If available, sort implants into the associated caddy.
- Place the products / caddy in the sterilization tray.
- Place the sterilization tray in the sterilization container.
- When using sterilization containers, comply with national regulations; they include EN ISO / ANSI/AAMI/ISO 11607 / EN 868-8.
- For containers without integrated filter system:
The sterilization container must be wrapped in two layers of 1-ply polypropylene wrap using sequential envelope folding techniques.
or
Place sterilization container into rigid sterilization container with integrated filter system and close lid (container in container system).

For requirements on sterilization containers, see chapter 8.1.4 „Machine cleaning and disinfection“.

8.1.8 Instructions for the implant groups

| Implant groups (for classification, see system-specific surgical instructions) | Machine cleaning and disinfection |
|--|--|
| Group Y | Standard procedure |
| Group Z | Positioning in fully extended state |

8.2 Instruments - Cleaning and Disinfection

8.2.1 Pre-treatment directly after use

WARNING

Residual contamination on contact surfaces.

- Patients can become infected following contact with contaminated instruments.
 - Residual cleaning agents can damage instruments.
- Disassemble instruments for pre-treatment as far as possible.
- Make sure that no soiling is allowed to dry onto instruments.

8.2.2 Pre-treatment in the Central Sterile Services Department

WARNING

Residual contamination on contact surfaces.

- Patients can become infected following contact with contaminated instruments.
 - Residual cleaning agents can damage instruments.
- Disassemble instruments for pre-treatment as far as possible.
- Make sure that no soiling is allowed to dry onto instruments.

Pre-treatment procedure

Selection criteria for the operating resources for the pre-treatment procedure, see chapter 8.7.1 „Operating resources for pre-treatment of instruments“.

1. Disassemble the instruments as far as possible.
2. Make up a fresh cleaning solution in the ultrasonic bath.
When doing so, observe the requirements of the cleaning agent manufacturer, e.g. temperature, concentration, exposure time and final rinsing.

3. Take instruments out of the transport trough and manually remove coarse soiling from the instruments within max. 2 hours. When doing so, use a fresh single-use cloth soaked in cleaning solution.
4. Place the instruments in the ultrasonic bath so that they are covered with cleaning solution. The instruments should not come into contact with each other. Treat the instruments in the ultrasonic bath as follows:
 - Carefully brush the instruments.
 - Rinse out the instruments. Rinse out the instruments according to the relevant instructions in the chapter 8.2.8 „Instructions for instrument groups“. For group allocation of the instruments, see the surgical instructions relating to the implant system.
 - At the same time, manipulate the moving parts and joints at least 3 times.
 - In the case of instruments with catch mechanisms, disengage where possible.
5. Activate the ultrasound for the exposure time specified by the manufacturer.
6. Rinse off the instruments under running water at least 3 times for at least 1 minute, turning several times.
7. Deposit the instruments in a cleaning tray.

8.2.3 Manual cleaning and disinfection

WARNING

Residual contamination due to inadequate cleaning and disinfection.

- Patients can become infected following contact with inadequately cleaned and disinfected instruments.
- Do not use manual cleaning and manual disinfection.

8.2.4 Machine cleaning and disinfection

WARNING

Residual contamination on contact surfaces.

- Patients can become infected following contact with contaminated instruments.
- Residual cleaning agents can damage instruments.
- Disassemble instruments as far as possible.
- Trials and caddies must be machine-cleaned and disinfected separately from each other.

ATTENTION

Damage and contamination of instruments due to the use of an unsuitable disinfection method.

- Do not perform chemical disinfection.

Machine cleaning and disinfection process

Selection criteria for the operating resources for the machine cleaning and disinfection process, see chapter 8.7.2 „Operating resources for machine cleaning and disinfection of implants and instruments“.

- Place instruments in the cleaning tray with the opening facing down diagonally (for further details see table in chapter 8.2.8 „Instructions for instrument groups“).
- Place the cleaning tray in the cleaning and disinfecting device.
- Check that the instruments are not in contact with each other.

- Start the program.
- Remove the instruments after the end of the program.

8.2.5 Inspection

Visually inspect the instruments for the following defects:

- Corrosion
- Surface damage
- Chipping
- Soiling: Clean and disinfect soiled instruments again.

WARNING

Impermissible reuse of damaged and soiled instruments.

- Patients can become infected following contact with contaminated instruments.
 - Damaged instruments may fail.
- Send back damaged instruments or dispose of them.
- Clean and disinfect soiled instruments again or send them back or dispose of them.

8.2.6 Care of the instruments

For the European market:

Selection criteria for the operating resources for the oiling, see chapter 8.7.3 „Λειτουργικοί πόροι για τη φροντίδα των οργάνων“. Grease all threaded surfaces and joints of the instruments sufficiently with instrument oil (based on white paraffin oil) to maintain their functionality.

For the US market:

Optional lubrication of the instruments can ensure their long-term functionality. This is at the discretion of the user. Selection criteria for the operating resources see chapter 8.7.3 „Operating resources for the care of instruments“.

8.2.7 Packaging

- Assemble instruments which need to be assembled for sterilization (see associated set images and surgical instructions).
- Carry out a functional check for the instruments.
- Place the products in the sterilization tray.
- Place the sterilization tray in the sterilization container.
- When using sterilization containers, comply with national regulations; they include EN ISO / ANSI/AAMI/ISO 11607 / EN 868-8.
- For containers without integrated filter system:
The sterilization container must be wrapped in two layer of 1-ply polypropylene wrap using sequential envelope folding techniques.
or
Place sterilization container into rigid sterilization container with integrated filter system and close lid (container in container system).

For requirements on sterilization containers, see chapter 8.7.4 „Operating resources for sterilization“.

8.2.8 Instructions for instrument groups

| Instrument groups (for classification, see system-specific surgical instructions) | Pre-treatment | Machine cleaning and disinfection |
|---|--|---|
| Standard | see Chapter 8.2.2 | Standard procedure |
| Group A | Thoroughly clean the lumen / cannulation with a brush at least 3 times | Place above a suitable rinsing mandrel |
| Group B | Brush carefully at least 3 times around the joint | Place above a suitable rinsing mandrel |
| Group C | Rinse out cannulation with syringe min. 5 ml and cannula with Ø 1.2 mm at least 3 times | Standard procedure |
| Group D | Brush out cannulation at least 3 times and rinse out with syringe min. 5 ml at least 3 times | Place above a suitable rinsing mandrel |
| Group E | Brush out cannulation at least 3 times and rinse out with syringe min. 40 ml at least 3 times | Place above a suitable rinsing mandrel |
| Group F | Rinse out blind lumen actively by immersing at least 3 times | Standard procedure |
| Group G | Rinse out actively at least 3 times by immersing Rinse blind lumen / variable lumen with syringe min. 5 ml at least 3 times | Standard procedure |
| Group H | Brush variable lumen at least 3 times Rinse out actively at least 3 times by immersing | Standard procedure |
| Group I | Standard procedure Remove lid | Standard procedure |
| Group S1 | Rinse out lumen / blind lumen with syringe min. 5 ml at least 3 times | Standard procedure |
| Group S2 | Rinse out blind lumen with syringe min. 2 ml and cannula with Ø 1.2 mm at least 3 times | Standard procedure |
| Group S3 | Brush out cannulation at least 3 times and rinse out with syringe min. 10 ml at least 3 times | Place above a suitable rinsing mandrel |
| Group S4 | Rinse out cannulation with syringe min. 5 ml and cannula with Ø 0.9 mm at least 3 times | Standard procedure |
| Group S5 | Rinse out annular gap with a syringe min. 2 ml at least 3 times | Standard procedure |
| Group S6 | No brush Rinse actively by immersing at least 3 times. When rinsing, bend open the tight gaps at least 3 times and also rinse inner lumen | Standard procedure |
| Group S7 | No brush Manipulate the moving parts backwards and forwards at least 5 times Rinse out actively at least 3 times by immersing | Positioning in fully extended state Place above a suitable rinsing mandrel |

8.3 Sterilization

WARNING

Use of impermissible sterilization procedures.

- Damage to implants and instruments.
- Do not use the following sterilization procedures:

- Flash sterilization procedures
- Gravity-displacement method
- Hot air sterilization
- Radiation sterilization
- Formaldehyde sterilization
- Ethylene oxide sterilization
- Plasma sterilization

WARNING

If maximum weight for sterilizing unit is exceeded, as per DIN EN 285 defined (Weight indication without sterilization container):

- Implants and instruments will not be sufficiently sterile.
- Do not exceed the following weights for container **with** integrated filter systems:
- 1/1 sterilizing unit max. 10 kg
 - 3/4 sterilizing unit max. 7 kg
 - 1/2 sterilizing unit max. 5 kg

Note: For the container in container systems, the inner container (container without integrated filter system) is included in the sterilizing unit.

- Loading of sterilization containers **without** integrated filter system according to the given marking.
- If no marking is available, the products must be stored in one layer without mutual contact.

Sterilization procedure

Selection criteria for the operating resources for the sterilization procedure and requirements on the sterilizer, see chapter 8.7.4 „Operating resources for sterilization“.

- Place the sterilization container in the sterilizer.
- Steam sterilization as fractionated vacuum method with at least three vacuum cycles and adequate product drying.

Recommended validated sterilization parameters:

| Time | Temperature | Drying time |
|-------------|-------------------------|---------------------|
| 4 minutes | At least 132 ° (270 °F) | At least 40 minutes |
| 3 minutes | At least 134 ° (273 °F) | At least 40 minutes |

8.4 Storage

Store implants and instruments in sterile packaging under the following ambient conditions:

- dry
- dust-free
- protected from damage

8.5 Reusability

- Implants can be reprocessed, as long as they are undamaged and not contaminated, and due appropriate care is taken.
- Reuse of the instruments is envisaged. To maintain the relevant functions of the instruments, they must be treated with care and handled according to the instructions.

8.6 Material resilience

WARNING

Damage and contamination of implants and instruments due to unsuitable cleaning and disinfecting agents.

- Residues can result in risks for the patient.
 - Damaged implants and instruments may fail.
- Use neutral, enzymatic or slightly alkaline cleaning agents
- Cleaning and disinfecting agents which contain the following substances must NOT be used:
- Organic, mineral or oxidising acids with a minimum permissible pH value of 5.5
 - Strong alkalis, a pH value of max. 11 is permissible
 - Oxidisers (e.g. hydrogen peroxide)
 - Halogens (e.g. chlorine, iodine, bromine)
 - Aromatic or halogenated hydrocarbons
 - Corrosion inhibitors
- Do not use rinse aids or neutralisers.

WARNING

The use of incorrect temperatures may damage implants, instruments, caddies, transport racks and sterilization containers.

- Damaged implants and instruments may fail.
- Implants, instruments, caddies, transport racks and sterilization containers should only be exposed to a maximum temperature of 142 °C (288 °F).

8.7 Technical data

ATTENTION

Damage to implants and instruments through the use of incorrect operating resources.

- As regards operating resources, the requirements given in chapter 8.2.6 „Φροντίδα των οργάνων“ must also be observed.

8.7.1 Operating resources for pre-treatment of instruments

(according to chapter 8.2.2 section „Instruments“)

The suitability of the instruments for pre-treatment has been verified by an independent, accredited and officially recognized testing laboratory using the following cleaning agents:

| Medium | Name | Manufacturer |
|----------------|---------------------------|------------------------------------|
| Cleaning agent | Neodisher mediclean forte | Dr. Weigert GmbH & Co. KG, Hamburg |

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| Properties of operating resources | |
|-----------------------------------|--|
| Cleaning agent ¹ | <ul style="list-style-type: none"> - Approved for cleaning instruments made of metal and plastic - Compatible with the materials of the instruments - Suitable for use with ultrasound - Free of alcohols and simple aldehydes, otherwise protein fixation |
| Water | <ul style="list-style-type: none"> - To max. 45 °C (113 °F) - Sterile or low-germ with max. 10 CFU/ml - Low in endotoxins at max. 0.25 EU/ml → e.g. purified water or highly purified water |
| Additional aids | <ul style="list-style-type: none"> - Absorbent and low-linting single-use cloths - Soft brushes with varying diameters - Syringes with varying volumes (for details see chapter 8.2.8) - Cannulas of varying sizes (for details see chapter 8.2.8) - For the group allocation of instruments, observe the surgical instructions relevant to the implant system - Ultrasonic bath |

8.7.2 Operating resources for machine cleaning and disinfection of implants and instruments

(according to chapter 8.1.4 section „Implants“ and chapter 8.2.4 section „Instruments“)

The suitability of the implants and instruments for machine cleaning and disinfection has been verified by an independent, accredited and officially recognized test lab using the following devices and cleaning agents:

| Device | Method | Manufacturer |
|---------------|----------------------|------------------------------------|
| RDG G 7836 CD | Thermal disinfection | Miele & Cie. GmbH & Co., Gütersloh |

| Medium | Name | Manufacturer |
|----------------|---------------------------|------------------------------------|
| Cleaning agent | Neodisher mediclean forte | Dr. Weigert GmbH & Co. KG, Hamburg |

The procedure followed corresponds to the procedure described in chapter 8.1.4 „Machine cleaning and disinfection“ and chapter 8.2.4 „Machine cleaning and disinfection“.

| Properties of operating resources | |
|-----------------------------------|---|
| Cleaning agent system | <ul style="list-style-type: none"> - Approved for cleaning implants and instruments made of metal and plastic - Compatible with the materials of the implants and instruments |

¹ If you are using a cleaning and disinfection agent, for example for health and safety reasons, please bear in mind that it should be aldehyde-free (otherwise blood contaminants will become fixed), should offer proven efficacy (e.g. VAH/DGFM or FDA/EPA approval/clearance/registration or CE certification), is suitable for disinfection of the instruments and is compatible with them (see chapter 8.2.6 „Φροντίδα των οργάνων“). Please note that the disinfectants used in pre-treatment only serve to protect staff and cannot replace the subsequent – post-cleaning – disinfection step.

Properties of operating resources

| | |
|----------------------------------|--|
| Cleaning and disinfecting device | <ul style="list-style-type: none"> - With proven efficacy (e.g. VAH/DGHM approval or FDA approval or CE certification as per EN ISO 15883) - Compliance with national regulations - Sufficient number of rinsing cycles - Use water that is sterile or low-germ and low in endotoxins <ul style="list-style-type: none"> - Germs: Max. 10 CFU/ml - Endotoxins: max. 0.25 EU/ml → e.g. purified water or highly purified water - Air for drying: filtered, oil-free, low-germ, low-particle - Suitable for thermal disinfection with A0 value of ≥ 3000 - Suitable for implants and instruments made of metal and plastic - Always perform thermal disinfection |
|----------------------------------|--|

8.7.3 Operating resources for the care of instruments

(according to chapter 8.2.6 „Care of the instruments“)

The suitability of the procedure for the care of instruments has been verified by an independent, accredited and officially recognized test laboratory using the following instrument oil:

| Medium | Name | Manufacturer |
|----------------|-----------------|-------------------------------|
| Instrument oil | Spezial Ölspray | Dr. Schumacher GmbH, Malsfeld |

For the care of the instruments an oil for surgical instruments can be used to lubricate all moving parts of the instruments. This is at the discretion of the user.

Properties of operating resources

| | |
|----------------|-------------------------------|
| Instrument oil | – based on white paraffin oil |
|----------------|-------------------------------|

⚠ WARNING

Damage and contamination of instruments due to use of unsuitable instrument oils.

- Residues can result in risks for the patient.
 - Damaged instruments may fail.
- Do not use oils with other additives such as silicones or corrosion inhibitors.
- Use instrument oil based on white paraffin oil.

8.7.4 Operating resources for sterilization

(according to chapter 8.3 „Sterilization“)

The suitability of the implants and instruments for sterilization has been verified by an independent, accredited and officially recognised test laboratory.

Properties of operating resources

| | |
|------------------|---|
| Steam sterilizer | <ul style="list-style-type: none"> - Design as per EN 13060 or EN 285 or ANSI/AAMI ST79 for USA: FDA cleared - Validation as per EN ISO 17665 - Maximum sterilization temperature plus tolerance acc. to DIN EN ISO 17665: 138 °C (280 °F) |
|------------------|---|

| Properties of operating resources | |
|---|---|
| Sterilization container | <ul style="list-style-type: none">- Compliance with national regulations, including EN ISO / ANSI/AAMI/ISO 11607 / EN 868-8 for USA: FDA cleared- Suitable for steam sterilization: To withstand temperatures up to at least 142 °C (288 °F), sufficiently steam permeable- Protection of products and sterilization trays from mechanical damage |
| Wrapping material for containers without integrated filter system | <ul style="list-style-type: none">- With proven efficacy (FDA approval or CE certification)- With compliance with national regulations |

8.7.5 Implant properties

| Feature | Property |
|-----------------------------|--|
| Material | The implants may be made of the following materials: <ul style="list-style-type: none">- Ti-6Al-4V in accordance with ASTM F 136- CoCr in accordance with ASTM F 1537 |
| As-delivered condition | Non-sterile |
| Maximum ambient temperature | 142 °C (288 °F) |

9. NOTICE FOR THE SCOPE OF EU REGULATION (MDR) 2017/745, ANNEX I, CHAP. III z)

All serious incidents that have occurred in connection with the product must be reported to the manufacturer and the competent authorities of the member state of the European Community in which the user and/or the patient is established, in accordance with EU Regulation (MDR) 2017/745.

Biedermann Motech has drawn up an SSCP, which can be obtained by the product-specific Basic UDI-DI: 42508696000010009E when the EUDAMED database becomes fully functional (<https://ec.europa.eu/tools/eudamed>).

10. SYMBOLS

The following table defines the symbols used for the product, labelling or the instruction for use of the system according to ISO 15223-1.

| Symbol | Title | Description of symbol (acc. ISO 15223-1) |
|--------|---|---|
| | Manufacturer | Indicates the medical device manufacturer. |
| | Date of manufacture | Indicates the date on which the medical device was manufactured. |
| | Medical device | Indicates that the item in question is a medical device. |
| | Catalogue number | Indicates the manufacturer's catalogue number so that the medical device can be identified. |
| | Batch code | Indicates the manufacturer's batch code so that the batch or lot can be identified. |
| | Caution | Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself. |
| | Consult electronic instructions for use www.biedermann.com/eifu | Indicates the need for the user to consult the instructions for use. www.biedermann.com/eifu |
| | Do not re-use | Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure. |
| | Non-sterile | Indicates a medical device that has not been subjected to a sterilization process. |

The following symbol is used according to the FDA „Guidance for Industry – Alternative to Certain Prescription Device Labeling Requirements“.

| Symbol | Title | Description |
|--------|---------|--|
| | Rx only | Caution: Federal law restricts this device to sale in the USA by or in the order of a physician. |

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| | |
|-----------------------------------|-----------|
| - BASIX Surgical Technique | GI-0153 |
| - BASIX MIS Surgical Technique | GI-0162 |
| - MOSS 100 Surgical Technique | PD-038079 |
| - MOSS VRS Surgical Technique | GI-038994 |
| - MOSS VRS MIS Surgical Technique | GI-041021 |

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Biedermann Motech GmbH & Co. KG