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VERTEBROPLASTY NEEDLE INSTRUCTIONS FOR USE

Produced by:

ADRIA SRL
Via Modena 46
40017 San Giovanni in Persiceto (BO)
Italy

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TEKNIMED S.A.S
8, rue du Corps Franc Pommiers
65500 VIC EN BIGORRE - FRANCE
Tél. (33)5 62 96 88 38
Fax (33)5 62 96 28 72
www.teknimed.com

GENERAL INFORMATION

Teknimed products must be implanted by qualified operators with in-depth knowledge and great surgical experience, and a full understanding of operating techniques specific to using Teknimed products.

The surgeon is responsible for complications or harmful consequences that could result from an incorrect indication or operating technique, improper use of the equipment and a failure to observe safety recommendations given in the instructions for use. These complications may be attributed neither to the manufacturer nor to the relevant Teknimed representative.

DESCRIPTION
Confort & Confort Biopsy trocars are single-use medical devices used for surgical procedures involving vertebral bone. The medical device is sterile and ready to use.

COMPOSITION
The Trocar Confort comprises:

- a 20G breakable-base chiba needle with its mandrin
- a trocar with cannulated mandrin
- a cement pump
- a trocar with mandrin canulated
- a 5-point bone biopsy cannula
- a sampling extraction mandrin

INDICATIONS
The trocar is designed for vertebroplasty techniques. Its use is proposed for bone diseases.

PRECAUTIONS FOR USE
The surgeon must have been trained in the vertebroplasty technique and must adhere closely to it. This is packaged and sterilised, strictly for single-use. Do not reuse, retreat or resterilise. Reusing, retreatring or resterilising the device may compromise its structural integrity and/or lead to the device malfunctioning, which can itself cause lesions, infections or death of the patient. Retreatring or resterilising single-use devices can also create a contamination risk and/or cause infections or infections cross-contaminated from one patient to another. Contamination of the device can lead to lesions, infections or death of the patient.

CONTRAINDICATIONS
Procedures other than those stated in the INDICATIONS section.
Proven infection or the risk of infection.
Haemorrhagic diathesis.
Patients susceptible to allergic reactions to the

materials making up the product.

WARNING
If using an orthopaedic hammer, apply it only to the anvil section. The trocar must be fitted with its mandrin when introducing it or advancing into the vertebral bone. The trocar must not be pushed further into the vertebral bone after the first injection of cement. Once the trocar is in place, do not apply force that could move it. Immediate mobilisation must be avoided to permit biological healing of bone and soft tissue.

POSITIONING
Trocar with chiba needle (Trocar Confort & Trocar Confort Biopsy)
The patient is placed in an appropriate position for the indication to be treated.
The 20G chiba needle contained in the kit is used to administer local anaesthetic.

- Put this needle in place under X-ray guidance until it comes into contact with the bone. Withdraw the internal mandrin and inject the anaesthetic locally.
- After injecting the anaesthetic, withdraw the breakable base; the needle will be used to guide the trocar fitted with its cannulated mandrin.
- The Trocar Confort assembled with its cannulated mandrin is advanced using the guide needle until it is positioned in the vertebra.
- The guide needle is withdrawn.
- To inject cement, the cannulated mandrin must be injected during the pasty phase of polymerisation to avoid venous leaks. Injection is carried out under fluoroscopic guidance.

Injection of cement (Trocar Confort & Trocar Confort Biopsy)
When the trocar is correctly positioned in the vertebral bone, the bone cement is injected under pressure directly into the fractured vertebra to stabilise it.
The injection is carried out using injectors that are connected to the Trocar Confort.
Once injection is completed, the cement pump is inserted into the trocar to inject the residual volume of cement.
The trocar is then withdrawn from the bone.

Biopsy needle (Trocart Confort Biopsy)
The bone biopsy cannula is inserted through the trocar.
The contact with the area to be sampled, the biopsy cannula is introduced into the bone and a biopsy is taken using a rotating movement.
The advance and depth of the biopsy are monitored under X-ray guidance.
Once the biopsy has been taken, the bone biopsy cannula is removed and the bone sample withdrawn using an extraction mandrin.
This bone biopsy needle is intended for use on bone tissue that readily lends itself to conventional sample collection in a conventional manner.
In the presence of a particularly hard bone (condensing bone lesions), it is advisable to use a material more suitable for this type of circumstance.

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INFORMATION FOR THE PATIENT
The doctor should make the patient aware of the potential consequences of the factors mentioned in the section on Contraindications. Specifically, those that could hinder success of the operation, as well as possible complications that could arise. The patient must also be informed about measures to be taken in order to reduce the potential consequences of these factors.

PACKAGING – STORAGE
Each product is packaged in sealed double bags and sterilised using ethylene oxide. Sterility is guaranteed until the expiry date for use. If the packaging has not been opened or damaged, Resterilisation of the product is strictly forbidden. SINGLE USE only.
Each product must be stored unopened in its original packaging in a clean, dry place at room temperature.

DISCARDING
The device or its components must be disposed of in accordance with current local waste regulations.

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