SPAZIATORE INTERSOMATICO LOMBARE IN TITANIO / LUMBAR TITANIUM INTERBODY FUSION CAGE / TITAN LUMBALER (LENDEN) FUSIONS-CAGE / ESPACIADOR INTERSOMÁTICO LUMBAR DE TITANIO / ESPACEUR INTERSOMATIQUE LOMBAIRE EN TITANE / ESPAÇADOR INTERSOMATICO LOMBAR DE TITÂNIO

ISTRUZIONI PER L'USO / INSTRUCTION FOR USE / GEBRAUCHSANWEISUNG / INSTRUCCIONES DE USO / MODE D'EMPLOI / INSTRUÇÃO DE USO

Rev. 01 – 03/09/2021

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## LUMBAR TITANIUM INTERBODY FUSION CAGE

#### IMPORTANT INFORMATION ABOUT THE LUMBAR TITANIUM INTERBODY FUSION CAGE SYSTEM

The LUMBAR TITANIUM INTERBODY FUSION CAGE is intended to recreate and maintain distance between vertebrae to support biological fusion in the thoracic, lumbar or lumbar-sacral spine zone. As complementary device, it should NOT be used as stand-alone. The

Tsunami Medical S.r.l.  
Sede Legale: Via E. Giorgi, 27 - 41124 Modena (MO) - Italy  
Sede Produttiva: Via XXV Aprile, 22 - 41037 Mirandola (MO) - Italy  
Tel.: +39 0535 38397, Fax: +39 0535 38399  
e-mail: [info@tsunamimed.com](mailto:info@tsunamimed.com) - [www.tsunamimedical.com](http://www.tsunamimedical.com)

**SPAZIATORE INTERSOMATICO LOMBARE IN TITANIO / LUMBAR TITANIUM INTERBODY FUSION CAGE / TITAN LUMBALER (LENDEN) FUSIONS-CAGE / ESPACIADOR INTERSOMÁTICO LUMBAR DE TITANIO / ESPACEUR INTERSOMATIQUE LOMBAIRE EN TITANE / ESPAÇADOR INTERSOMATICO LOMBAR DE TITÂNIO**

**ISTRUZIONI PER L'USO / INSTRUCTION FOR USE / GEBRAUCHSANWEISUNG / INSTRUCCIONES DE USO / MODE D'EMPLOI / INSTRUÇÃO DE USO**

Rev. 01 – 03/09/2021

device must always be within the limits of the vertebral cortices. If this indication is not followed, the device can be crushed because it protrudes from the spinal column and may not withstand the forces involved in such conditions

## DESCRIPTION

The LUMBAR TITANIUM INTERBODY FUSION CAGE consists of cages in variable sizes and shapes. The dimensions of implants are designed based on anatomical conditions and decisions hereto are made by physicians.

LUMBAR TITANIUM INTERBODY FUSION CAGE components should be used in combination with other spinal systems or other fixation systems in order to obtain stabilization.

All components of LUMBAR TITANIUM INTERBODY FUSION CAGE cannot be re-used under any circumstances.

LUMBAR TITANIUM INTERBODY FUSION CAGE is designed to be applied for posterior, posterior-lateral, anterior and lateral approach.

In particular, this instruction for use is applicable for codes:

ACT Cage for transforaminal arthrodesis "TLIF"

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ACTH Expandable lumbar cage

ACXH Expandable extra-lateral cage and variable lordosis for lumbar arthrodesis

ACPH Variable lordosis expandable cage for posterolateral arthrodesis

MM Expandable anterior cage for lumbar arthrodesis

CT Expandable anterolateral cage for lumbar arthrodesis

MMJ Self-locking expandable anterior cage for lumbar arthrodesis

ACXJ Self-locking extra-lateral cage for lumbar arthrodesis

ACAJ Self-locking anterior cage for lumbar arthrodesis

ACJ Self-locking anterior cage for lumbar arthrodesis

ACTZ Lumbar Cage

## MATERIALS

The entire system is made of medical grade titanium described by ISO 5832-3 or 10993-5 or ASTM F2026 or ASTM F136 standards. Tsunami Medical warrants that all devices are manufactured from one of the foregoing material specifications.

## INDICATIONS

**LUMBAR TITANIUM INTERBODY FUSION CAGE** is intended for lumbar interbody fixation for the following indications:

1. Degenerative disc disease.
2. Spondylolisthesis.
3. Spinal stenosis.
4. Trauma.
5. Tumour.
6. Pseudo-arthritis.
7. Instability of motion segments.

## CONTRAINDICATIONS

Contraindications include, but are not limited to:

<p><b>Tsunami Medical S.r.l.</b>  <b>Sede Legale: Via E. Giorgi, 27 - 41124 Modena (MO) - Italy</b>  <b>Sede Produttiva: Via XXV Aprile, 22 - 41037 Mirandola (MO) - Italy</b>  <b>Tel.: +39 0535 38397, Fax: +39 0535 38399</b>  <b>e-mail: <a href="mailto:info@tsunamimed.com">info@tsunamimed.com</a> - <a href="http://www.tsunamimedical.com">www.tsunamimedical.com</a></b></p>	<p>5</p>
--	----------

**SPAZIATORE INTERSOMATICO LOMBAR IN TITANIO / LUMBAR TITANIUM INTERBODY FUSION CAGE  
/TITAN LUMBALER (LENDEN) FUSIONS-CAGE / ESPACIADOR INTERSOMÁTICO LUMBAR DE TITANIO /  
ESPACEUR INTERSOMATIQUE LOMBAIRE EN TITANE / ESPAÇADOR INTERSOMATICO LOMBAR DE TITÂNIO**

**ISTRUZIONI PER L'USO / INSTRUCTION FOR USE / GEBRAUCHSANWEISUNG /  
INSTRUCCIONES DE USO / MODE D'EMPLOI / INSTRUÇÃO DE USO**

Rev. 01 – 03/09/2021

1. Risk of infection or infection in progress, or fever or inflammation.
2. Obesity.
3. Pregnancy.
4. Mental illness.
5. Allergy on any system components.
6. Any anatomical, medical or surgical conditions which may preclude potential or intentional benefits of spinal implants application.
7. Bone, joints or ligaments conditions such but not limited as: osteopenia, bone absorption, osteomalacia. Osteoporosis is relative contraindications a must be carefully evaluated prior surgery.
8. Implants size, shape or anchorage functionality might be not enough to achieve expected clinical results.
9. Combination with implants from other manufacturers.
10. Potential risk of unexpected patient's anatomy destruction, interference with neurological, functional or other deficits.
11. Any risk of patient's unwillingness to follow postoperative instructions.
12. Any other not described in indications.

Caution: in case of reuse there is a danger of cross contamination ; any reuse is therefore not permitted.

**POTENTIAL ADVERSE EVENTS**

Possible adverse events which might occur after spinal surgery with or without instrumentation including, but are not limited to:

1. Disassembly, bending, and/or breakage of any or all of the system's components.
2. Migration of any of the system's components.
3. Pressure on the skin from component parts in patients with inadequate tissue coverage.
4. Tissue or nerve damage caused by improper positioning and placement of implants or improper use of instruments.
5. Dura leakage, distortion or damage.
6. Neurologic dysfunctions and/or physiological dysfunctions like paresthesia, radiculopathy, paralysis, hypertension, or any others related to surgery in general associated to anesthesia.
7. Infections.
8. Loss of urinary functions.
9. Permanent or temporary or developing sexual dysfunctions.
10. Postoperative change of body curvature, change of physiological range of movement.
11. Pseudo-arthrosis or non-fusion or delayed fusion.
12. Bone loss or overgrowth, or any other bone malformations.
13. Permanent or temporary limitation or inability to perform daily activities.
14. Change in mental behavior.
15. Permanent or temporary or developing respiratory problems.
16. Permanent or temporary or developing cardiovascular deteriorations or dysfunctions.
17. Death.

In some cases, additional surgery or surgeries might be necessary to correct or change potential adverse events.

**WARNINGS**

The effectiveness and safety of interbody fixation is only applicable for certain conditions with which require the fusion supported by medical device. The device might be supportive for such mechanical instability like deformity, fracture, listhesis, dislocation, tumour, pseudo-arthrosis. The effectiveness and safety for any other conditions are unknown.

**PRECAUTIONS**

**SPAZIATORE INTERSOMATICO LOMBARRE IN TITANIO / LUMBAR TITANIUM INTERBODY FUSION CAGE  
/TITAN LUMBALER (LENDEN) FUSIONS-CAGE / ESPACIADOR INTERSOMÁTICO LUMBAR DE TITANIO /  
ESPACEUR INTERSOMATIQUE LOMBAIRE EN TITANE / ESPAÇADOR INTERSOMATICO LOMBAR DE TITÂNIO**

**ISTRUZIONI PER L'USO / INSTRUCTION FOR USE / GEBRAUCHSANWEISUNG /  
INSTRUCCIONES DE USO / MODE D'EMPLOI / INSTRUÇÃO DE USO**

Rev. 01 – 03/09/2021

The LUMBAR TITANIUM INTERBODY FUSION CAGE are complementary implants to posterior fixation systems. The applications of pedicle screw and/or interbody cages should be performed by experienced surgeons with specific training in use of LUMBAR TITANIUM INTERBODY FUSION CAGE. The spinal screw fixation system and/or interbody cage system should not be considered as sole spinal support. No implants can withstand body loads without bone support. Therefore bends, breakages, loosening, disassembling may occur over the time. A successful result is not always achievable. The factors as proper preoperative and operative procedure, comprehensive knowledge of surgical techniques, proper selection of implant's size and type are considerably important in treatment process. Patients with obesity, smokers and alcoholics are at risk for non-fusion. Also, patients with weak muscle or bone conditions and/or nervous system dysfunctions are poor candidates for spinal fusion. Prior or during or after the surgery, in order to evaluate or check the positioning of the implants or patient's anatomy or any other patient or implant correction, X-ray or CT or any other diagnostic examinations may be necessary.

The proper, patient's individual implant selection in terms of type, size, shape or design is vital to successful surgery performance. Proper implants and instruments' handling is crucial. Extensive bending or contouring should be avoided. Sharp cutting edges, reversed bending, scratching or notching may generate internal stress, which may weak the implants or construct.

**IMPORTANT:** All necessary information about surgery, potential risks, benefits and adverse effects should be conveyed to the patient prior to surgery.

#### **PRE-OPERATIVE**

1. Patients that meet the criteria described in the indications should only be selected.
2. Patients' conditions should be checked prior to surgery; any required diagnostics should be performed.
3. The efficient and adequate implants and instruments inventory must be secured and be available during the surgery.
4. All implants, instruments and any other components should be cleaned and sterilized before use. Any implants, instruments or components delivered in sterile packaging must be checked on sterility and expiration date of sterility prior surgery.
5. Implants and instruments should be stored in certain conditions to warrant the sterility and protection from any contamination or corrosive environment.
6. It's highly recommended that all personnel interacting with any mechanical components of the spinal system should be familiar with all components before use.

#### **INTRA-OPERATIVE**

1. Extreme caution should be taken when working close to or around the spinal cord and nerve roots.
2. Whenever possible or required, intra-operative diagnostic systems should be used to facilitate surgery.
3. Breakage, bends, scratch, slippage, part loosening or improper use of any implant or instrument during the surgery may cause injury to OR personal or patient.
4. It's very important to follow the surgical technique carefully. Proper application of any instrument or implant may facilitate surgery.
5. Before closing of soft tissue, double check of all implants' positioning, geometrical relations, and fixing, tightening or mounting maneuvers for all screws, nuts or other fixing parts have been performed. Image diagnostics is highly recommended at this stage.

#### **POST-OPERATIVE**

The physician's post-operative directions and warnings to the patient, and the corresponding patient's compliance, are extremely important.

1. Detailed instructions about the use and limitations of the device should be communicated with the patient.
2. The patient should be warned to avoid falls or sudden jolts in spinal position.
3. The patient should be warned for this possibility and instructed to limit and restrict physical activities, especially lifting and twisting movements and any type of sports participation. The patient should be advised not to smoke, or to consume alcohol or non-steroids or use anti-inflammatory medications such as aspirin during the bone graft healing.
4. As a precaution, before patients with implants receive any subsequent surgery (such as surgical dental procedures), prophylactic antibiotics may be considered, especially for high-risk patients.

**SPAZIATORE INTERSOMATICO LOMBARRE IN TITANIO / LUMBAR TITANIUM INTERBODY FUSION CAGE / TITAN LUMBALER (LENDEN) FUSIONS-CAGE / ESPACIADOR INTERSOMÁTICO LUMBAR DE TITANIO / ESPACEUR INTERSOMATIQUE LOMBAIRE EN TITANE / ESPAÇADOR INTERSOMATICO LOMBAR DE TITÂNIO**

**ISTRUZIONI PER L'USO / INSTRUCTION FOR USE / GEBRAUCHSANWEISUNG / INSTRUCCIONES DE USO / MODE D'EMPLOI / INSTRUÇÃO DE USO**

Rev. 01 – 03/09/2021

- Any retrieved implants should be treated in such a manner that reuse in another surgical procedure is not possible. As with all orthopaedic implants, the LUMBAR TITANIUM INTERBODY FUSION CAGE components should never be reused under any circumstances.

#### PACKAGING

LUMBAR TITANIUM INTERBODY FUSION CAGE is sterile packaged; in order to control the correct sterilization, make sure that the package arrives properly closed. Do not use if the blister is open or damaged.

#### STORAGE

The components of the LUMBAR TITANIUM INTERBODY FUSION CAGE must be handled with care to prevent damage. Store in designated trays and in areas which provide protection from dust, insects, chemical vapours and extreme changes in temperature and humidity. Sterile parts must be stored in original packages as a precaution for prevented to any damages.

#### WARRANTY AND PRODUCT COMPLAINTS

Every spinal product from Tsunami Medical is guaranteed to be free of defects in handling and materials when used properly for its intended purpose. Any implant or instrument delivered from Tsunami Medical proving to be defective after proper use following the instructions provided and used for its intended purpose will be replaced or repaired, at Tsunami discretion, without charge.

Any customer or user, who has any complaint or who has experienced any dissatisfaction in the product quality, durability, safety or reliability should inform Tsunami Medical or its dedicated Distributor in the customer's territory. When filing a complaint, please provide the component(s) name, part number, lot number(s), your name and address and a detailed description of the complaint. Please provide the report directly to Tsunami Medical or to its dedicated Distributor.

#### CUSTOMER SERVICE

Additional information about the LUMBAR TITANIUM INTERBODY FUSION CAGE or any other Tsunami Medical's Spinal System is available from Tsunami Medical. You may also want to contact your local Tsunami Medical Sales Representative.

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