obelisc™ vertebral body replacement

Gebrauchsanweisung	de	Deutsch	2
Instructions for use	en	English	7
Instructions for use	en	USA	11
Mode d'emploi	fr	Français	12
Istruzioni per l'uso	it	Italiano	17
Instrucciones de uso	es	Español	22
Instruções de utilização	pt	Português	27
Brugsanvisning	da	Dansk	32
Bruksanvisning	sv	Svenska	36
Käyttöohje	fi	Suomi	40
Kullanım Kılavuzu	tr	Türkçe	44
Οδηγίες χρήσης	el	Ελληνικά	48
Gebruiksaanwijzing	nl	Nederlands	53
Bruksanvisning	no	Norsk	58
Инструкция по применению	ru	Русский	62
Instrucțiunile de utilizare	ro	Română	67
Инструкции за употреба	bg	Български	72
Návod k použití	cs	Čeština	77
Instrukcja obsługi	lv	Latviski	81
Lietošanas instrukcija	pl	Polski	85
Návod na použitie	sl	Slovenščina	90
Navodila za uporabo	sk	Slovenčina	94
사용 설명서	ko	한국어	98
使用说明	zh	中文	102

SYSTEM:	obelisc™ vertebral body replacement	
LANGUAGE:	en	English The only version applicable to users in the U.S. is the version intended especially for the United States.
	Read all instructions carefully. Failure to follow instructions, warnings and precautionary measures may have severe consequences or result in patient injury.	
	Take into careful consideration all product materials before clinical use. The required product materials prepared as overall documentation of the system are: system related instructions for use, surgical technique and, if applicable, supplements, assembly and disassembly instructions as well as "Processing manual implants and instruments" UH 1100.	

1 Intended use

The obelisc vertebral body replacement is used for surgical reconstruction of substance defects of the anterior, thoracic and lumbar spine in humans. Additional stabilizing instrumentation is necessary, for example, with a pedicle screw-rod-system.

2 Product description

The obelisc vertebral body replacement is intended for surgical reconstruction of substance defects of the anterior thoracic and lumbar spine. Defect expansion can be performed in situ with the obelisc implant. The obelisc center pieces are continuously expandable, so enabling a defect of up to 132 mm to be bridged. Thus one or more vertebral bodies can be replaced. The required height of the obelisc implant is adjusted precisely in situ via a bevel gear drive in the inserter. The implant maintains the set height. This can only be changed by turning the bevel gear further. Final fixation is performed by means of a locking screw. Thanks to its space-saving expansion mechanism, the obelisc implant can be used in a minimally invasive manner. Obelisc consists of a center piece with corresponding end pieces, and is available in various sizes. The different obelisc end pieces enable a maximum supporting surface. A free choice of approach is possible thanks to the variable assembly of the end pieces. The teeth on the end pieces facilitate an improved grip in the vertebral body end plates. The implant only absorbs the axial compression force and always needs to be combined with a posterior or ventral stabilization system.

3 Indications

Complete or incomplete corpectomy due to destruction of a vertebral body by e.g. tumor, fracture or inflammation.

4 Contraindications

- Patients with acute infection, whether superficial or deep
- Patients with fever or leukocytosis
- Patients with obesity (according to WHO - World Health Organization)
- Patients with a history of material allergy or who tend to react to foreign bodies
- The physician must consider carefully before treating patients who are in a generally unfavorable medical or psychological state and who could be made worse by the procedure
- Patients with inadequate bone quality or quantity (e.g. severe osteoporosis, osteopenia, osteomyelitis)
- Pregnancy

5 Possible side effects

Possible complications associated with the implant are:

- Implant loosening, dislocation and/or implant failure
- Local or systemic reactions due to material intolerance
- Failure of bony fusion
- Adjacent segment degeneration

These possible complications may lead to further surgical interventions (implant removal or stabilization renewal).

Other complications which may be caused by the surgical procedure itself and independent of the implant are:

- General surgical risks and complications

- Anesthesia and blood transfusion risks as well as positioning damage
- Lung complications
- Infection
- Wound healing disorders
- Cardiovascular complications, such as blood loss, thrombosis, embolism, coagulopathy
- Gastrointestinal complications, such as gastritis, ileus, ulcer
- Neurological complications, such as spinal cord or neural root lesions with temporary or permanent sensory and/or motor disturbances (bladder and rectal disorders, sexual dysfunction)
- Intraoperative injury to blood vessels, massive hemorrhage, stroke, cerebral hemorrhage with potential life-threatening consequences
- Injury to organs adjacent to the spine, such as cervical, thoracic and abdominal organs, depending on the region operated on

6 Warnings and precautions

6.1 General

- The user has to ensure that the most current versions of the complete product materials provided as overall documentation of the system are on hand and considered. These are also available at: www.ifu.ulrichmedical.com. A printed copy may be delivered within 7 days, upon request.
- This product must be used only by an experienced spine surgeon, under sterile conditions and in an operating room where C-arm is available.
- The product is delivered non-sterile and must be sterilized before its first use. Instructions for processing can be found in the "Processing manual implants and instruments" UH 1100, prepared as part of the overall documentation of the system.
- The treating physician is responsible for making the proper selection of patients, and for acquiring the training with the system and experience required for implant selection and placement. The treating physician is also responsible for determining whether to leave the implant following surgery or to remove it.
- The product must be carefully handled and stored. Implants with any kind of damage or scratches must not be used, as product stability and resistance to fatigue may be affected.
- Implants of the corresponding ulrich medical systems must be used exclusively with the system-specific instruments intended for this purpose, unless otherwise indicated.
- Connecting ulrich medical implants to implant components of other manufacturers is not permitted. A combination with other implants from ulrich medical is not permitted.
- Given the presence of passive layers, corrosion of metal implants is very minimal. Nevertheless, there may be accelerated material fatigue and the possibility of a break, as well as increased release of metal components in the body. Corrosion is promoted by the components made of different metals coming into contact and by damage to the implant's surface. The ability to combine the various metal implants used in the system has been assured. Direct contact with metal implants from other manufacturers is not allowed.

6.2 Preoperative

- The obelisc system is only intended to absorb axial compression force. Additional stabilizing instrumentation is required to absorb the other forces (e.g. transverse forces with angled end pieces) and torques.
- Implantation should only be considered when all other treatment options have been carefully weighed and ruled out. A successful implant is also inferior to the healthy moving element(s) of the spine. On the other hand, an implant may be an efficient replacement of one or more severely deformed or degenerated moving elements in order to eliminate pain and achieve good load bearing capacity.
- The user must carefully weigh the benefits against the risks when using it in a non-mature spine.
- The risks of the procedure and the use of the implant, including any revisions which may be necessary, should be explained to the patient in detail.
- The treating physician should conduct with the patient a detailed discussion of the surgical results that may be anticipated, particularly with regard to potential physical limitations of the implant. The degree of postoperative activity affects the life time of the implant and the stability of the implant in the bone. Thus, the patient must be informed of the limitations and risks inherent in daily activities and made aware of special guidelines for movement. Smoking can inhibit fusion and may predispose to an increased rate of pseudoarthrosis. Patients should be informed about the possible negative effects of smoking on the success of a fusion procedure. The treating physician must assess whether the patient can understand and follow such instructions. Special attention should be given to postoperative discussion and the need for regular medical monitoring.
- Ensure that suitable pre-, intra- and postoperative diagnostic procedures are performed to evaluate and ensure proper implant selection and placement. For example, direct contact between the implant and delicate tissues (vessel walls, dura etc.) should be avoided in order to reduce the risk of injury in the event that revision surgery is necessary.
- Errors in implant selection could result in premature clinical implant failure. The number of segments to be treated must be accurately determined. The shape and composition of human bones limits the size and durability of the implant.

6.3 Intraoperative

- Vertebral body resection includes the neighboring disks.
- After the implant has been in contact with a patient, has been soiled or used, it must not be reused and must be discarded. Even if the exterior of the implant appears unchanged, prior stress may have caused damage that can lead to implant failure.
- In order to achieve fusion, bone or bone graft material should be grafted to the implant.

6.4 Postoperative

- The patient should be instructed to inform his or her treating physician immediately regarding any unusual changes in the area where the surgery was performed.
- The patient should be monitored for any changes to the implant area. The treating physician should assess the potential for clinical implant failure and discuss with the patient the necessary measures to promote further healing.
- The inserted implant serves to reconstruct substance defects of the surgical site over a maximum two-year healing process. After the surgical site is fused, the implant is firmly anchored in the bone. The implant is therefore not intended to be removed unless there are complications, implant failure, or a delayed healing phase (no fusion within 2 years) that require implant removal. The decision to do so should only be made after a meticulous risk to benefit assessment by a medical specialist.
- The following complications may occur due to the presence of the implant in the body:
 1. Corrosion, with local tissue inflammation or pain
 2. Implant migration, possibly leading to injuries
 3. Risk of additional injuries, caused by post-surgical trauma
 4. Distortion, loosening or breakage
 5. Pain, symptoms or unusual sensitivity due to the presence of the implant
 6. Risk of infection or inflammation
 7. Bone loss due to stress shielding
 8. Potential, unknown or unexpected long-term effects
- A prolonged healing phase, unsuccessful bony fusion or subsequent bone resorption or trauma can place undue stress on the implant, which, in turn, could lead to implant failure.
- Implant failure is possible even after successful fusion.
- The implant removal must be followed up accordingly.
- Explanted implants must not be reused.
- The patient should be informed about the necessary behavior after the implantation procedure, including measures to be taken in case of an event involving the implant. An implant card should be given to the patient.

7 Recommended application

Further description of the obelisc application is included in the surgical technique, prepared as part of the overall documentation of the system. It is also available at www.ifu.ulrichmedical.com.

8 Material information

The obelisc implants are manufactured from a titanium alloy according to ISO 5832-3 and ASTM F136.

The devices are biocompatible, corrosion-resistant and non-toxic in the intended use as per EN ISO 10993-1.

9 MRI information

Ulrich medical obelisc implants are classified as "MR Conditional" according to ASTM standard F2503.

The "MR Conditional" components are tested according to the ASTM standards: F2052; F2182; F2213 and F2119. A patient with an obelisc implant can undergo an MRI examination under the following circumstances:

- field strengths of 1.5 T and 3.0 T
- highest field gradient of 30 T/m (3000 G/cm) or less
- maximum specific absorption rate (SAR) of 2 W/kg (normal operating mode) for a maximum scan duration of 15 minutes per individual measurement, or a total scan duration of 30 minutes.

In the non-clinical test, heating of a maximum of 5.6°C was measured at 3.0 T and a 15-minute scan duration with 2 W/kg according to the calorimetric measurement. At 1.5 T and a 15 minute-scan duration, heating of a maximum of 5.3°C was measured with 2 W/kg according to the calorimetric measurement.

These scan conditions allow low-risk examination of patients. To minimize heating, the scan duration should be kept as short as possible and the SAR as low as possible. Heating can be reduced by having breaks between the individual scans or by lowering the SAR.

Artifacts: MR imaging in the area of the implants can be impaired by artifacts. Under experimental conditions, artifacts up to 30 mm were detected radially around the implant.

The scans were performed with the following parameters:

This document is valid only on the date printed. If unsure of the print date, re-print to ensure use of the latest version of the IFU. The onus resides with the user to ensure that the most up-to-date IFU is used.

- FFE sequence: TR 100 ms, TE 15 ms, flip angle 30°

- SE sequence: TR 500 ms, TE 20 ms, flip angle 70°

The SE sequence showed reduced artifacts (≤ 11 mm).

The attending physician should conduct a careful risk/benefit assessment.

10 Packaging and storage

The product is delivered non-sterile and must be sterilized before its first use. Before a system is used, all components should be meticulously checked for completeness, damages and defects. Damaged components must not be used.

Corrosion may occur when instruments are stored under unfavorable conditions. In order to avoid this, they should be stored in a dry, dust-free area. Significant temperature fluctuations are to be avoided, so that no moisture (condensation) accumulates on the instruments.

When directly exposed to metal, chemical substances may destroy this metal or release corrosive fumes. Therefore, instruments must not be stored together with chemical substances.

The system-specific trays are to be used to store the instruments.

11 Identification and traceability

Implants are labeled with the catalogue number and batch code on the packaging label and, if technically feasible, these are marked on the implant itself.

To ensure traceability, these data must be correspondingly documented in the hospital.

12 Disposal of used products

Follow country-specific regulations regarding disposal of hospital waste when disposing of used medical devices.


13 Cleaning, disinfection and sterilization

The product is delivered non-sterile and must be sterilized before its first use.

For cleaning, disinfection and sterilization of non-sterile implants and instruments the "Processing manual implants and instruments" UH 1100 that was provided in the overall documentation of the system must be used. For individual instruments which are correspondingly marked in the surgical technique, the "Assembly and disassembly instructions with special cleaning instructions" that was provided in the overall documentation of the system have to be considered.

These are also available at: www.ifu.ulrichmedical.com.

14 Labeling and symbols

 Catalogue number	 Caution	Qty. Quantity
 Batch code	 Do not re-use	 Non-sterile
 Manufacturer	 Date of manufacture	 Do not use if package is damaged
 Consult electronic instructions for use www.ifu.ulrichmedical.com	 Consult instructions for use	
 Caution: U.S. federal law restricts this device to sale by or on the order of a physician	 MR Conditional – A product which does not demonstrate any deliberate hazards in the specific MRI environment	
 CE mark with identification number of the notified body		
 Consult assembly and disassembly instructions with special cleaning instructions		

This document is valid only on the date printed. If unsure of the print date, re-print to ensure use of the latest version of the IFU. The onus resides with the user to ensure that the most up-to-date IFU is used.

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